

Withdrawn Paper

Warning Notice

This paper has been withdrawn. It has been superseded by the document identified below.

Withdrawal Date January 14, 2025


Original Release Date May 21, 2024

Superseding Document

Title Reference Measurement Procedure for the Absolute Quantification of Albumin in Urine Using Isotope Dilution-Liquid Chromatography-Tandem Mass Spectrometry (ID-LC-MS/MS) Update 1

Publication Date January 14, 2025

DOI/URL <https://doi.org/10.6028/NIST.SP.1200-31-upd1>



NIST Special Publication 1200
NIST SP 1200-31

Reference Measurement Procedure for the Absolute Quantification of Albumin in Urine Using Isotope Dilution-Liquid Chromatography-Tandem Mass Spectrometry (ID-LC-MS/MS)

Ashley Beasley-Green
N. Alan Heckert

This publication is available free of charge from:
<https://doi.org/10.6028/NIST.SP.1200-31>

NIST Special Publication 1200
NIST SP 1200-31

Reference Measurement Procedure for the Absolute Quantification of Albumin in Urine Using Isotope Dilution-Liquid Chromatography-Tandem Mass Spectrometry (ID-LC-MS/MS)

Ashley Beasley-Green
Biomolecular Measurement Division
Material Measurement Division

N. Alan Heckert
Statistical Engineering Division
Information Technology Laboratory

This publication is available free of charge from:
<https://doi.org/10.6028/NIST.SP.1200-31>

May 2024



U.S. Department of Commerce
Gina M. Raimondo, Secretary

National Institute of Standards and Technology
Laurie E. Locascio, NIST Director and Under Secretary of Commerce for Standards and Technology

NIST SP 1200-31
May 2024

Certain equipment, instruments, software, or materials, commercial or non-commercial, are identified in this paper in order to specify the experimental procedure adequately. Such identification does not imply recommendation or endorsement of any product or service by NIST, nor does it imply that the materials or equipment identified are necessarily the best available for the purpose.

Publications in the SP 1200 subseries include written procedural methods in the design and implementation of experiments that ensure successful replication of results by others. Publications may include detailed procedures, lists of required equipment and instruments, information on safety precautions, the calculation of results and reporting standards.

NIST Technical Series Policies

[Copyright, Use, and Licensing Statements](#)

[NIST Technical Series Publication Identifier Syntax](#)

Publication History

Approved by the NIST Editorial Review Board on 2024-05-21

How to Cite this NIST Technical Series Publication

Beasley-Green A, Heckert NA (2024) Reference Measurement Procedure for the Absolute Quantification of Albumin in Urine Using Isotope Dilution-Liquid Chromatography-Tandem Mass Spectrometry (ID-LC-MS/MS). (National Institute of Standards and Technology, Gaithersburg, MD), NIST Special Publication (SP) NIST SP-1200-31. <https://doi.org/10.6028/NIST.SP.1200-31>.

Author ORCID iDs

Ashley Beasley-Green: 0000-0002-2065-4218

N. Alan Heckert: 0000-0002-8430-6757

Contact Information

Ashley Beasley Green, Ph.D.

ashley.beasley@nist.gov

Abstract

Urinary excretion of albumin is a major diagnostic and prognostic marker of kidney disease; therefore, accurate measurement of urine albumin is essential to clinical diagnosis. To support the accuracy and comparability of clinical urine albumin measurements, the National Institute of Standards and Technology (NIST) has partnered with the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK) Laboratory Working Group (LWG) and the International Federation of Clinical Chemistry (IFCC) Working Group for the Standardization of Albumin Assays in Urine (WG-SAU) to develop a reference measurement system for urine albumin. NIST has developed the foundational components of the reference measurement system for urine albumin, a series of higher-order reference materials and a reference measurement procedure (RMP), to establish a traceability scheme that will link clinical urine albumin to the International System of Units (SI). The NIST RMP for urine albumin is based on the detection and quantification of signature proteotypic (trypsin) albumin peptides via isotope dilution-liquid chromatography-tandem mass spectrometry (ID-LC-MS/MS) and multiple reaction monitoring (MRM). The RMP incorporates an isotopically labeled (^{15}N) full-length recombinant human serum albumin (^{15}N -rHSA) material as the internal standard, which permits the absolute quantification of albumin in urine. Out of the 11 tryptic peptides selected to monitor the albumin in urine, 5 albumin peptides were used for the absolute quantification of albumin in urine. This publication outlines the NIST RMP for the absolute quantification of albumin in urine by ID-LC-MS/MS.

Keywords

Albumin, isotope dilution-liquid chromatography-tandem mass spectrometry (ID-LC-MS/MS), measurement uncertainty (MU), metrological traceability, multiple reaction monitoring (MRM), recombinant human serum albumin, reference measurement procedure (RMP), Standard Reference Material® (SRM), urine.

TABLE OF CONTENTS

Executive Summary	1
1. Introduction	2
2. Warning & Safety	3
2.1. Chemical and Instrumentation	3
2.2. Human Urine Material	3
3. Scope	3
4. Measurement Principle & Measurement Method	4
5. Materials	4
6. Solutions	4
7. Sample Preparation Procedure	5
7.1. Preparation of Process Samples.....	5
7.2. In-Solution Trypsin Digestion.....	5
8. LC-MS System	6
9. Quantitative Method	8
9.1. Integration of MRM Transitions.....	8
9.2. Calibration Curve.....	8
9.3. Calculation of Albumin Content.....	9
9.4. Calculation of Combined Uncertainty.....	10
9.5. Calculation of Expanded Uncertainty.....	10
10. Summary	11
References	11
Appendix A. List of Symbols, Abbreviations, and Acronyms	13
Appendix B. Data Report Template	14

List of Tables

Table 1. List of materials and instrumentation for NIST RMP.	4
Table 2. LC gradient for NIST RMP	6
Table 3. MRM segments for tandem MS analysis.	7
Table 4. MRM transition list for the NIST RMP.	7
Table 5. Description of variables for Equation 1.	9
Table 6. Description of variables for Equation 2.	10
Table 7. Description of variables for Equation 3.	10
Table 8. Description of variables for Equation 4.	10
Table 9. Description of variables for Equation 5.	10

List of Figures

Fig. 1. Proposed reference measurement system (RMS) for traceability of clinical urine albumin results using NIST SRMs and the NIST RMP. Adapted from ISO 17511 [9].....	2
Fig. 2. Detailed measurement procedure for the NIST RMP for urine albumin [10,11].	5
Fig. 3. Chromatographic profile of unlabeled (NIST SRM 2925) and labeled (IS) MRM Transitions.....	8
Fig. 4. Two-level calibration system for qt-MRM Transition 18.	9

Executive Summary

Decisions made by healthcare practitioners regarding disease diagnosis and management are heavily influenced by the validity of clinical laboratory results. It is crucial for healthcare practitioners to have access to accurate and consistent patient care information, which can only be achieved through standardization of clinical laboratory results. When clinical results are standardized, the clinical value is precise, equivalent, and not dependent on the method or laboratory used. This can be achieved by establishing metrological traceability of the results to higher-order reference materials and measurement procedures. However, when clinical results are not standardized, different values may be obtained for the same clinical sample using different methods or clinical laboratories. This can have a significant impact on patient care, from the delivery of erroneous medical decisions to inflated healthcare costs. Metrological traceability is applied to establish a traceability framework that underpins the confidence and global comparability of clinical results used in the diagnosis and management of disease.

1. Introduction

Kidney disease is a major global public health issue, with chronic kidney disease (CKD) representing one of the most prominent causes of death worldwide (12th leading cause of death in 2017). It was estimated in 2017 that more than 10 % of the general global population was diagnosed with CKD, which totals to greater than 800 million individuals [1]. Based on data from the National Health and Nutrition Examination Survey (NHANES; 2017 to 2020), an estimated 15 % of adults in the United States (37 million individuals) have been diagnosed with CKD [2,3]. The public health and economic impact of CKD and other renal diseases has led to the need for the accurate detection of kidney disease biomarkers, such as urine albumin, for early diagnosis, evaluation of treatment efficacy, and disease management.

Urinary excretion of albumin is a major diagnostic and prognostic marker of kidney disease; therefore, accurate clinical urine albumin results are essential for patient care decisions. To support the accuracy and comparability of clinical urine albumin measurements, NIST has partnered with the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK) Laboratory Working Group (LWG) and the International Federation of Clinical Chemistry (IFCC) Working Group for the Standardization of Albumin Assays in Urine (WG-SAU) to develop a reference measurement system for urine albumin. NIST has developed a series of higher-order reference materials and a reference measurement procedure (RMP) to establish a metrological traceability framework that will link clinical urine albumin to the International System of Units (SI). The traceability scheme in Figure 1, illustrates how the NIST RMP is used to support metrological traceability of clinical urine albumin results.




Fig. 1. Proposed reference measurement system (RMS) for traceability of clinical urine albumin results using NIST SRMs and the NIST RMP. Adapted from ISO 17511 [9].

The NIST RMP (Fig. 1, p.3) is intended for use in the value-assignment of albumin in NIST Standard Reference Material® (SRM) 3666 (Fig. 1, m.2), the matrix-based (human urine) secondary reference material [4,5]. Calibration of the RMP is SI traceable via the use of NIST SRM 2925 as the primary

calibrator [6]. (currently listed in JCTLM database – #C18RM1; m.2) NIST SRM 2925 serves as the primary reference material (Fig. 1, m.1) for the urine albumin reference measurement system (Fig. 1) and is currently listed in the Joint Committee for Traceability in Laboratory Medicine (JCTLM) database for reference materials [7,8]. NIST SRM 3666 is a four (4)-level (Level I to Level IV) human urine material intended for use as a secondary reference material (Fig. 1, m.2) to support the accuracy and comparability of clinical urine albumin and urine creatinine results used in clinical decisions for kidney disease [4,5].

The NIST RMP for urine albumin incorporates isotope dilution-liquid chromatography-tandem mass spectrometry (ID-LC-MS/MS) using the multiple reaction monitoring (MRM) scan mode to detect and measure proteolytic albumin peptides in urine with high selectivity and sensitivity [10, 11]. To target full-length albumin and to reduce bias introduced by sample preparation and enzymatic digestion, an isotopically labeled (^{15}N) full-length recombinant human serum albumin (^{15}N -rHSA) that is chemically equivalent to the native albumin is used as the internal standard (IS). With the incorporation of the full-length ^{15}N -rHSA IS, 11 tryptic peptides that span the sequence of albumin are measured and 5 peptides are used for the absolute quantification of albumin in urine. Using a series of measurement equations, the consensus mean mass fraction and mass concentration values are determined for albumin in urine [11].

2. Warning & Safety

2.1. Chemical and Instrumentation

Chemical hazards will vary depending on the experiment and samples being analyzed. Read the SDS for any compound being used as a solvent prior to conducting the RMP. Take all the necessary precautions for these compounds. Care should be taken with operating a LC-MS system, general hazards and controls associated with the system include: high temperatures, electrical hazards, and high voltage. Personal protective equipment (PPE) is required for the RMP.

2.2. Human Urine Material

Handle urine samples as biohazardous material capable of transmitting infectious disease. For human-derived urine products where the presence of an infectious agent may be unknown, handle human urine-based products at Biosafety Level 2 as recommended by the Centers for Disease Control and Prevention's Biosafety in Microbiological and Biomedical Laboratories (6th edition) [12].

3. Scope

The objective of the RMP is to digest albumin and use ID-LC-MS/MS to detect and quantify the proteolytic products of albumin to determine the mass fraction ($\bar{X}_{mg/g}$) and mass concentration ($\bar{X}_{mg/L}$) values and associated uncertainties of albumin in urine. No detectable effects of the human urine matrix on measurement repeatability and precision compared to a buffer system (50 mmol/L ammonium bicarbonate) were previously observed [10]. Trypsin is used to generate the proteolytic albumin peptides (unlabeled and IS) and the mass ratio is determined from the peak area ratio of the unlabeled to labeled (IS) peak area for each MRM transition. The measurement interval for the RMP aligns with the clinical ranges for albumin in urine (0 mg/L to 500 mg/L).

4. Measurement Principle & Measurement Method

The RMP is a targeted multiplexed ID-LC-MS/MS assay that incorporates a full-length ^{15}N -rHSA IS for the absolute quantification of albumin in the urine. The procedure couples ID-LC-MS/MS with the multiple reaction monitoring (MRM) MS scan mode to selectively target signature tryptic albumin peptides. The ratio of unlabeled analyte to ^{15}N -labeled IS is used to generate a calibration curve for the determination of albumin in a urine sample.

5. Materials

A list of the materials and instrumentation used in the NIST RMP are listed in Table 1.

Table 1. List of materials and instrumentation for NIST RMP.

Item	Vendor
NIST SRM 2925 Recombinant Human Serum Albumin Solution (<i>calibrator material</i>) (listed in JCTLM database – C18RM1) [6,8]	NIST
NIST SRM 3666 Albumin and Creatinine in Frozen Human Urine (<i>control material</i>) [4,5]	NIST
Internal Standard: Full-length ^{15}N - rHSA (at least 99% label incorporation)	Albumin Biosciences
Trypsin-Gold MS-grade	Promega
Dithiothreitol (DTT) (No-Weigh Format, 7.7 mg DTT per microtube)	Pierce
Iodoacetamide	Pierce
Water with 0.1 % (<i>volume fraction</i>) formic acid (high-purity LC-MS grade)	Honeywell Burdick and Jackson
Acetonitrile (ACN) with 0.1 % (<i>volume fraction</i>) formic acid polypropylene	Honeywell Burdick and Jackson
Autosampler vials (amber glass, screw cap)	Agilent
Caps for Autosampler Vials	Agilent
1.0 mL Tubes (polypropylene, sterile)	Eppendorf
6460 Triple Quadrupole Mass Spectrometer equipped with an Agilent 1290 Series LC system	Agilent
Zorbax 300 SB-C18 column (2.1 mm × 150 mm, 3.5 μm)	Agilent

6. Solutions

- a. 1 mol/L Ammonium Bicarbonate (AMBIC) (total volume – 50 mL)
 - a) Add 3.953 g of ammonium bicarbonate to 50 mL of water
- b. 50 mmol/L AMBIC (total volume – 50 mL)
 - a) Add 2.5 mL of 1 mol/L AMBIC to 47.5 mL of water
- c. 200 mmol/L AMBIC (total volume – 15 mL)
 - a) Add 3.0 mL of 1 mol/L AMBIC to 12 mL of water
- d. 50 mmol/L Acetic Acid (total volume – 50 mL)
 - a) Add 2.5 mL of 1 mol/L acetic acid to 47.5 mL of water
- e. 500 mmol/L DTT (total volume – 100 μL) (*prepare fresh during digestion procedure*)
 - a) Add 100 μL of 50 mmol/L AMBIC to 7.7 mg DTT

- f. 375 mmol/L Iodoacetamide (total volume – 132 μ L) (*prepare fresh during digestion procedure*)
 - a) Add 132 μ L of 200 mmol/L AMBIC to Iodoacetamide
- g. 1 μ g/ μ L Trypsin (total volume – 100 μ L) (*prepare fresh during digestion procedure*)
 - a) Add 100 μ L of 50 mmol/L acetic acid to 100 μ g of Trypsin

7. Sample Preparation Procedure

An illustration of the sample preparation procedure is shown in Figure 2.




Fig. 2. Detailed measurement procedure for the NIST RMP for urine albumin [10,11].

7.1. Preparation of Process Samples

1. Calibration Material: The NIST SRM 2925 is used to prepare the calibration solutions for the NIST RMP [6,8].
2. Calibration Solutions: Gravimetrically prepare calibration solutions using NIST SRM 2925 and 50 mmol/L AMBIC within analytical range (5 mg/L to 500 mg/L) of RMP. Add the volume of NIST SRM 2925 and 50 mmol/L AMBIC to the tube and record the mass of each volume (total volume of 150 μ L).
 - a. Two-curve or Three-curve calibration system that is concentration dependent
 - i. Curve #1: approx. 5 mg/L to 50 mg/L
 - ii. Curve #2: approx. 50 mg/L to 150 mg/L
 - iii. Curve #3: approx. 150 mg/L to 500 mg/L
3. Urine Specimen and Quality Control Material (NIST SRM 3666 or in-house material): Gravimetrically prepare the urine specimen. Add 150 μ L of urine specimen to the tube and record the mass.
4. Blank Sample: Gravimetrically prepare blank sample by adding 150 μ L of 50 mmol/L AMBIC to tube and record mass.
5. Internal Standard: Gravimetrically add the same amount (approx. 5 μ g to 8 μ g of protein) of internal standard to each calibration solution, quality control material, blank sample, and urine specimen.
6. Solubilization: Store Process Samples (calibration solutions, blank sample, quality control material, urine specimens) at 5 °C overnight protected from light.
7. After overnight incubation, equilibrate Process Samples to room temperature (approx. 25 °C) for 1 h prior to trypsin digestion.

7.2. In-Solution Trypsin Digestion

1. Record the start date and time for the trypsin digestion procedure
2. Incubate Process Samples at 97 °C for 10 min
3. Cool Process Samples to room temperature (approx. 25 °C) for 2 min
4. Prepare 500 mmol/L DTT stock

5. Add 500 mmol/L DTT stock to a final concentration of 5 mmol/L DTT to each Process Sample
6. Incubate Process Samples at 60 °C for 30 min
7. Prepare 375 mmol/L Iodoacetamide stock
8. Add 375 mmol/L Iodoacetamide stock to a final concentration of 15 mmol/L Iodoacetamide to each Process Sample
9. Incubate Process Samples at room temperature (approx. 25 °C) in dark for 30 min
10. Prepare 1 µg/µL Trypsin Solution
11. An approximate 1:30 mass ratio of trypsin-to-total protein is used for digestion
12. Incubate Process Samples at 37 °C for 24 h (record digestion start and end time)
13. Add 60 µL of 0.1 % (*volume fraction*) formic acid in water to quench the digestion reaction
14. Incubate Analytical Samples (digested Process Samples) at 37 °C in dark for 45 min
15. Concentrate Analytical Samples in SpeedVac (no heat) to dryness
16. Resuspend Analytical Samples in 100 µL of 0.1 % (*volume fraction*) formic acid in water

8. LC-MS System

Analysis of the digested Analytical Samples is performed on an Agilent 6460 triple quadrupole mass spectrometer in positive ion mode equipped with an Agilent 1290 Series LC system utilizing an Agilent Zorbax 300 SB-C18 column (2.1 mm × 150 mm, 3.5 µm). The column temperature is maintained at 45 °C, and the peptides are loaded onto the column at a flow rate of 200 µL/min in 97 % (*volume fraction*) mobile phase A (water with 1 mL/L formic acid) and 3 % (*volume fraction*) mobile phase B (ACN with 1 mL/L formic acid). The LC gradient is outlined in Table 2 [10,11].

Table 2. LC gradient for NIST RMP.

Time (min)	Mobile Phase A (%)	Mobile Phase B (%)	Flowrate (mL/min)
0.0	97 %	3 %	0.2
2.0	97 %	3 %	0.2
4.0	90 %	10 %	0.2
14.0	80 %	20 %	0.2
24.0	70 %	30 %	0.2
26.0	30 %	70 %	0.2
27.0	3 %	97 %	0.2
32.0	3 %	97 %	0.2
34.0	97 %	3 %	0.2
45.0	97 %	3 %	0.2

General mass spectrometric conditions: gas temperature of 300 °C; gas flow of 7 L/min; nebulizer of 20 psi (1.4×10^5 Pa); sheath gas temperature of 300 °C; sheath gas flow of 6 L/min; capillary voltage of 4000 V; and a nozzle voltage of 1500 V. The MRM segment for the MS method is outlined in Table 3 [10,11]. The segment time for the MS analysis is dependent upon the analytical column type and the LC type and configuration.

Table 3. MRM segments for tandem MS analysis.

MRM Segment Number	Segment Time (min)	Scan Mode	Divert Valve Status	Store Data (Y/N)
1	0.0	MS2 Scan	To Waste	N
2	4.0	MRM	To MS	Y
3	9.5	MRM	To MS	Y
4	12.5	MRM	To MS	Y
5	15.7	MRM	To MS	Y
6	18.5	MRM	To MS	Y
7	21.0	MRM	To MS	Y
8	29.0	MS2 Scan	To Waste	N

Analysis of the Analytical Samples is conducted in a randomized sequence with replicate measurements (4 replicates) to reduce the influence of systematic bias on the output measurements. The MRM transitions (11 peptides with 2 or 3 MRM transitions per peptide) are listed in Table 4 [10,11].

Table 4. MRM transition list for the NIST RMP.

Peptide Type – Qt or Qa ^a	MRM Transition #	Peptide	Precursor m/z	Product m/z	IS ^b Precursor m/z	IS ^b Product m/z	MS MRM Segment
Qt	Transition 3	TYETTLEK	492.75	720.30	497.23	727.40	2
Qt	Transition 7	VFDEFKPLVEEPQNLIK	682.37	900.00	689.35	909.50	7
Qt	Transition 9	FQNALLVR	480.78	685.40	487.27	695.40	4
Qt	Transition 18	LCTVATLR	467.26	660.40	472.75	669.40	3
Qt	Transition 20	YLYEIAR	464.25	651.30	469.24	659.30	4
Qa	Transition 1	DLGEENFK	476.22	723.30	481.21	731.30	3
Qa	Transition 2	DLGEENFK	476.22	229.07	481.21	231.07	3
Qa	Transition 4	TYETTLEK	492.75	265.10	497.23	279.20	2
Qa	Transition 5	VFDEFKPLVEEPQNLIK	682.37	970.50	689.35	981.50	7
Qa	Transition 6	VFDEFKPLVEEPQNLIK	682.37	712.40	689.35	721.37	7
Qa	Transition 8	FQNALLVR	480.78	276.09	487.27	279.09	4
Qa	Transition 10	QTALVELVK	500.81	488.27	506.29	493.25	5
Qa	Transition 11	QTALVELVK	500.81	587.30	506.29	593.40	5
Qa	Transition 12	RPCFSALEVDETYVPK	637.65	961.50	644.3	972.40	5
Qa	Transition 13	RPCFSALEVDETYVPK	637.65	244.17	644.3	247.16	5
Qa	Transition 14	LVAASQAALGL	507.30	189.08	513.29	191.08	6
Qa	Transition 15	LVAASQAALGL	507.30	712.40	513.29	721.40	6
Qa	Transition 16	LVNEVTEFAK	575.31	937.40	581.29	947.40	4
Qa	Transition 17	LVNEVTEFAK	575.31	694.40	581.29	701.40	4
Qa	Transition 19	LCTVATLR	467.26	274.12	472.75	276.12	3
Qa	Transition 21	YLYEIAR	464.25	277.20	469.24	279.10	4
Qa	Transition 22	AEFAEVSK	440.72	201.05	445.21	203.04	2
Qa	Transition 23	AEFAEVSK	440.72	680.32	445.21	687.30	2

^a The "Peptide Type" differentiates the Quantitative (Qt) and Qualitative (Qa) MRM Transitions. The five Qt-MRM Transitions are in bold

^b IS represents the internal standard (¹⁵N-labeled rHSA)

Figure 3 is a representative chromatographic profile of the 23 MRM transitions for both the unlabeled analyte (NIST SRM 2925) and IS (^{15}N -labeled rHSA) using the conditions and parameters listed in this document. The retention time (RT) for the unlabeled (NIST SRM 2925) MRM transitions overlaps with the RT for the corresponding IS (^{15}N -labeled rHSA) MRM transition (Fig. 3, qt-MRM Transitions labeled in bold with an asterisk*).




Fig. 3. Chromatographic profile of unlabeled (NIST SRM 2925) and labeled (IS) MRM Transitions

9. Quantitative Method

9.1. Integration of MRM Transitions

The Agilent MassHunter Quantitative Analysis software (Version B.10.00 or higher) is used to automatically integrate the MRM chromatographic peaks for both the unlabeled albumin and IS MRM transitions. Following automatic integration, manually confirm all automatic peak integrations. If manual correction of the automatic integration is required, record the corrected peak integration. A total of 46 measurements (11 peptides with 2 or 3 MRM transitions per peptide for both unlabeled and ^{15}N -labeled albumin) are collected for the 23 MRM transitions in each sample. Following integration, the instrument data (peak area, peak intensity, retention time) is imported from MassHunter into Microsoft Excel for manual quantitative assessment of the raw data. For the quantitative evaluation, the unlabeled peak area and IS peak area are used to calculate the Peak Area Ratio (PAR) for each transition and sample. The peak intensity and retention time data are used to monitor the repeatability of the MRM measurements and the performance of the LC-MS system.

9.2. Calibration Curve

Linear regression analysis (concentration ratio vs. PAR) of the calibration solutions is used to generate the calibration curve (2-level or 3-level calibration system) for each MRM transition (qt-MRM and qa-MRM transitions). An example 2-level calibration system is presented in Fig. 4. The linear equation

and coefficient of determination (R^2 value) for each calibration curve are determined using Microsoft Excel.

2-Level Calibration System for qt-MRM Transition 18




Fig. 4. Two-level calibration system for qt-MRM Transition 18.

9.3. Calculation of Albumin Content

A full description of the equations used to calculate the albumin content is presented in Ref. 11. The mass fraction value ($X_{t_{1mg/g}}$) for each MRM transition is derived from [11]:

$$X_{t_{1mg/g}} = \frac{PAR_t - B_{0,t}}{B_{1,t}} \times \left(\frac{m_U}{\frac{m_{IS}}{d_{IS}} \times C_{IS}} \times 100 \right). \quad (1)$$

Table 5. Description of variables for Equation 1.

Variable	Description
$X_{t_{1mg/g}}$	Mass fraction value of an MRM transition
PAR_t	t_1 -MRM transition PAR (unlabeled t_1 -MRM transition in urine sample to IS t_1 -MRM transition, ^{15}N -labeled IS)
$B_{0,t}$	y-intercept of the linear calibration curve for the t_1 -MRM transition
$B_{1,t}$	Slope of the linear calibration curve for the t_1 -MRM transition
m_U	Mass of the unknown urine sample
m_{IS}	Mass of the IS
d_{IS}	Density of the IS
C_{IS}	Mass concentration of the IS
m_U	Mass of the unknown urine sample

The consensus mass concentration value ($\bar{X}_{mg/L}$) is derived from [11]:

$$\bar{X}_{mg/L} = \bar{X}_{mg/g} \times (d_{urine} \times 1000). \quad (2)$$

Table 6. Description of variables for Equation 2.

Variable	Description
$\bar{X}_{mg/L}$	Consensus mass concentration value
$\bar{X}_{mg/g}$	Consensus mass fraction value
d_{IS}	Density of the IS

9.4. Calculation of Combined Uncertainty

A full description of the equations used to calculate the combined uncertainty for $\bar{X}_{mg/L}$ and $\bar{X}_{mg/g}$ is outlined in Ref. 11. The combined uncertainty for the $\bar{X}_{mg/g}$ is derived from [11]:

$$u_{mg/g} = \sqrt{u_{Type A}^2 + u_{Type B}^2}. \quad (3)$$

Table 7. Description of variables for Equation 3.

Variable	Description
$u_{Type A}$	Type A uncertainty of $\bar{X}_{mg/g}$ determined via the DSL-HHD or DSL-bootstrap methods
$u_{Type B}$	Combined Type B uncertainties [11]

The combined uncertainty for the $\bar{X}_{mg/L}$ is derived from [11]:

$$u_{mg/L} = \sqrt{d^2 u_{mg/g}^2 + \bar{X}_{mg/g} u_d^2}. \quad (4)$$

Table 8. Description of variables for Equation 4.

Variable	Description
$u_{mg/L}$	$\bar{X}_{mg/L}$ combined uncertainty
d_{urine}	Density of the urine material
u_d	Density uncertainty

9.5. Calculation of Expanded Uncertainty

The $u_{mg/g}$ and $u_{mg/L}$ values are expressed as expanded uncertainties (U), which are calculated using [11]:

$$U_{mg/g} = k \times u_{mg/g} \quad \text{or} \quad U_{mg/L} = k \times u_{mg/L}. \quad (5)$$

Table 9. Description of variables for Equation 5.

Variable	Description
$U_{mg/g}$	Expanded uncertainty for $u_{mg/g}$
$u_{mg/g}$	$\bar{X}_{mg/g}$ combined uncertainty
$U_{mg/L}$	Expanded uncertainty for $u_{mg/L}$
$u_{mg/L}$	$\bar{X}_{mg/L}$ combined uncertainty
k	Coverage factor ($k = 2$)

10. Summary

The NIST RMP is a targeted, multiplexed procedure that applies quantitative proteomics techniques for the absolute quantification of albumin in urine. The RMP incorporates a full-length IS (¹⁵N-labeled rHSA) that allows the detection of multiple signature peptides (11 peptides) that span the complete amino acid sequence of albumin. By including multiple peptides, both quantitative and qualitative information about albumin in the urine sample is observed. Applying the MS-based MRM technique offers a high degree of analytical specificity and sensitivity for the detection of albumin in urine. The triple quadrupole mass spectrometer selectively separates the precursor/product ion pair of each target peptide from the bulk digestion products; therefore, the urine sample can be analyzed with minimal sample preparation. As a result of the selective nature of the MRM method, the sensitivity of the assay is enhanced, which enables the detection of urine albumin at the lower normoalbuminuria range (lower limit of quantitation of approx. 5.0 mg/L) with high precision and accuracy. The dynamic range of the multiplexed assay traverses all three albuminuria stages (approx. 5 mg/L to 500 mg/L), and this supports the usage of the assay for early detection of microalbuminuria (30 mg/L to 300 mg/L) and the observation of renal dysfunction progression from normoalbuminuria (0 mg/L to 30 mg/L) to macroalbuminuria (> 300 mg/L).

References

- [1] Kovesdy CP Epidemiology of Chronic Kidney Disease: An Update 2022. *Kidney Int Suppl.* 2011; <https://doi.org/10.1016/j.kisu.2021.11.003>.
- [2] Centers for Disease Control and Prevention. Chronic Kidney Disease in the United States, 2021. Atlanta, GA: US Department of Health and Human Services, Centers for Disease Control and Prevention. <https://www.cdc.gov/kidneydisease/pdf/Chronic-Kidney-Disease-in-the-US-2021-h.pdf>.
- [3] Johansen KL, Chertow GM, Foley RN, et al. US Renal Data System 2020 Annual Data Report: Epidemiology of Kidney Disease in the United States. *Am J Kidney Dis.* 2021; <https://doi.org/10.1053/j.ajkd.2021.01.002>.
- [4] NIST Certificate of Analysis, SRM 3666 Albumin and Creatinine in Frozen Human Urine. 2024. <https://tsapps.nist.gov/srmext/certificates/3666.pdf>.
- [5] Beasley-Green A, Camara J, Heckert NA. Certification of Standard Reference Material 3666 Albumin and Creatinine in Frozen Human Urine. NIST Special Publication 260-238-upd1. 2023. <https://doi.org/10.6028/NIST.SP.260-238-upd1>.
- [6] Beasley-Green A, Bunk DM, Alejo W, Zhang NF. Certification of Standard Reference Material 2925 Recombinant Human Serum Albumin Solution (Primary Reference Calibrator for Urine Albumin) (Frozen). NIST Special Publication 260-199. 2020. <https://doi.org/10.6028/NIST.SP.260-199>.
- [7] The Joint Committee for Traceability in Laboratory Medicine. Available at <https://www.jctlm.org/> (accessed May 2024).
- [8] Joint Committee for Traceability in Laboratory Medicine (JCTLM) Database for Reference Materials; NIST SRM 2925 Recombinant Human Serum Albumin in Frozen Aqueous Solution, [C18RM1](#).
- [9] International Organization for Standardization (2008) ISO 17511:2020 – In vitro diagnostic medical devices — Requirements for establishing metrological traceability of values assigned to

calibrators, trueness control materials and human samples (International Organization for Standardization, Geneva, CH).

- [10] Beasley-Green A, Burris N, Bunk DM, Phinney KW. Multiplexed LC-MS/MS assay for urine albumin. *J. Proteome Res.* 2014; <https://doi.org/10.1021/pr500204c>.
- [11] Beasley-Green A, Heckert A. "Estimation of Measurement Uncertainty for the Quantification of Protein by LC-MS/MS". *Anal. Bioanal. Chem.* 2023; <https://doi.org/10.1007/s00216-023-04705-8>.
- [12] CDC/NIH; Biosafety in Microbiological and Biomedical Laboratories, 6th ed.; Meehan, P.J.; Potts, J., Eds.; U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention and National Institutes of Health; National Institutes of Health; HHS Publication (CDC) 300859 (2020); available at <https://www.cdc.gov/labs/BMBL.html> (accessed Jan 2024).

Appendix A. List of Symbols, Abbreviations, and Acronyms

ACN	Acetonitrile
AMBIC	Ammonium Bicarbonate
CKD	Chronic Kidney Disease
DTT	Dithiothreitol
HSA	Human Serum Albumin
ID-LC-MS/MS	Isotope Dilution-Liquid Chromatography-Tandem Mass Spectrometry
IFCC	International Federation of Clinical Chemistry
IFCC WG-SAU	IFCC Working Group for the Standardization of Albumin Assays in Urine
IS	Internal Standard
JCTLM	Joint Committee for Traceability in Laboratory Medicine
MRM	Multiple Reaction Monitoring
MU	Measurement Uncertainty
NHANES	National Health and Nutrition Examination Survey
NIDDK	National Institute of Diabetes and Digestive and Kidney Diseases
NIDDK LWG	NIDDK Laboratory Working Group
NIST	National Institute of Standards and Technology
PAR	Peak Area Ratio
PPE	Personal Protective Equipment
Qt	Quantitative
Qa	Qualitative
RM	Reference Material
RMP	Reference Measurement System
SRM	Standard Reference Material®

Appendix B. Data Report Template

Reference Measurement Procedure (RMP)	NIST RMP for the Absolute Quantification of Albumin in Urine Using Isotope Dilution-Liquid Chromatography-Tandem Mass Spectrometry (ID-LC-MS/MS)
Date(s) of Sample Collection (MM.DD.YY)	
Date(s) of Sample Processing (MM.DD.YY)	
Date(s) of Analysis (MM.DD.YY)	
Analyst Name	
Calibration System	<input type="checkbox"/> Two(2)-Level Calibration System <input type="checkbox"/> Three(3)-Level Calibration System

Sample Results: Mass Fraction ($\bar{X}_{mg/g}$) (units: mg/g)

Sample Name	Mass Fraction Value (mg/g), $(\bar{X}_{mg/g})$	Combined Uncertainty - Mass Fraction (mg/g), $u_{mg/g}$	Coverage factor (k)	Expanded Uncertainty - Mass Fraction (mg/g), $U_{mg/g}$

Sample Results: Mass Concentration ($\bar{X}_{mg/L}$) (units: mg/L)

Sample Name	Mass Concentration Value (mg/L), $(\bar{X}_{mg/L})$	Combined Uncertainty - Mass Concentration (mg/L), $u_{mg/L}$	Coverage factor (k)	Expanded Uncertainty - Mass Concentration (mg/L), $U_{mg/L}$