

NIST

Standard Reference Materials Catalog 1990-91

NIST Special Publication 260

U.S. DEPARTMENT OF COMMERCE

National Institute of Standards and Technology

TO ORDER

Phone: (301) 975-OSRM (6776)

Fax: (301) 948-3730

NEW NIST PUBLICATION

December 1989





UNITED STATES DEPARTMENT OF COMMERCE
National Institute of Standards and Technology
[formerly National Bureau of Standards]

Gaithersburg, Maryland 20899

Dear Colleague:

Welcome to the very first catalog of Standard Reference Materials (SRM's) to be issued by the National Institute of Standards and Technology. Of course, the materials listed herein are an extension of 84 years of service to the U.S. technical community provided by the National Bureau of Standards since 1906. Our new name is indicative of the effort we are making to extend the tradition of SRM service into new fields and into a new era of service to technology.

You will find over two hundred changes in SRM availability listed in this catalog — new additions, renewals, and revised certificates. Please carefully review the sections which relate to your work to see the new (and previously issued) SRM's that are available. If we can help you with regard to any listing please give us a call.

Many of you have told us, at SRM exhibits, or by letter, that you are happy that SRM's continue to be available through the new NIST organization. In fact, we are seeing some modest growth in our ability to supply the materials that you need. If you have a requirement that we do not currently meet please let us know about it. In the front pages of this catalog you will find a section entitled *Guide for Requesting Development of Standard Reference Materials*. You can register your need with us with just a few minutes of effort. Then, as funding permits, we can try to be responsive to your needs.

Thank you for your continued interest in Standard Reference Materials from NIST. It is our privilege to provide them to you and thereby enhance the quality of measurements and analyses made throughout the technical community.

Sincerely,

A handwritten signature in black ink that reads "William P. Reed".

William P. Reed
Acting Chief
Office of Standard Reference Materials

NIST Standard Reference Materials Catalog 1990 – 1991

R.L. McKenzie, *Editor*

Office of Standard Reference Materials
National Institute of Standards and Technology
Gaithersburg, MD 20899

CAUTION: The values shown in the catalog are nominal values only. Users should consult the certificate issued with an SRM for the certified values.



U.S. Department of Commerce
Robert A. Mosbacher, *Secretary*

National Institute of Standards and Technology
Ramond G. Kammer, *Acting Director*

Issued January 1990

TO ORDER
PHONE 301/975-6776
Telex TRT 197674NBS UT
See page 6 for
ORDERING
INSTRUCTIONS

National Institute of Standards and Technology
Special Publication 260
Supersedes NIST Spec. Publ. 260, 1988-89
163 pages (January 1990)
CODEN: XNBSAV

U.S. GOVERNMENT PRINTING OFFICE
WASHINGTON: 1990

For sale by the Superintendent of Documents
U.S. Government Printing Office, Washington, DC 20402

Contents

ABSTRACT AND KEY WORDS, 1

PROGRAM INFORMATION, 3

Definitions, 3

SRM Catalog, 4

Preparation and Availability of Standard Reference Materials, 4

Guide for Requesting Development of Standard Reference Materials, 5

ORDERING STANDARD REFERENCE MATERIALS, 6

General, 6

Terms, 7

Late Charges, 8

Proforma Invoice (Price Quotations), 8

Domestic Shipments, 8

Foreign Shipments, 9

Documentation, 9

CERTIFIED REFERENCE MATERIALS

FROM OTHER SOURCES, 10

Special Nuclear Materials, 10

International CRM's, 10

OTHER SERVICES OF THE NATIONAL INSTITUTE OF STANDARDS AND TECHNOLOGY, 11

Calibration and Related Measurement Services, 11

Office of Weights and Measures, 11

Proficiency Sample Programs, 12

Accreditation of Testing Laboratories, 12

National Center for Standards and Certification Information, 12

National Standard Reference Data System, 12

CHEMICAL COMPOSITION, 15

Ferrous Alloys, 15

Steels (Chip Form), 15

Plain Carbon Steels, 15

Low Alloy Steels, 16

Special Low Alloy Steels, 17

High Alloy Steels, 18

Stainless Steels, 19

Tool Steels, 19

Steels (Solid Form), 20

Low Alloy Steels, 20

Stainless Steels, 24

Specialty Steels, 24

High Temperature Alloys, 25

Steelmaking Alloys, 26

Cast Irons (Chip Form), 27

Cast Steels, White Cast Irons, Ductile Irons, and Blast Furnace Irons (Solid Form), 28

Nonferrous Alloys, 30

Aluminum-Base Alloys, 30

Copper-Base Alloys (Chip Form), 31

Copper-Base Alloys (Solid Form), 32

Copper "Benchmark", 33

Lead-Base Alloys, 34

Lead-Base Material, 34

Nickel-Base Alloys, 35

Carbon-Modified Silicon, 35

Trace Elements in Nickel-Base Superalloys (Chip Form), 36

Nickel Oxides (Powder Form), 36

Titanium-Base Alloys (Chip Form), 37

Zinc-Base Alloys, 38

Zirconium-Base Alloys, 38

Gases in Metals, 39

High Purity Metals, 40

Microanalytical, 41

Metals for Microanalysis, 41

Mineral Glasses for Microanalysis, 42

Glasses for Microchemical Analysis, 42

Thin Film for X-Ray Spectrometry, 43

Glass Fibers for Microanalysis, 43

Primary, Working, and Secondary Chemicals, 44

Microchemicals, 44

Spectrometric Solutions, 45

Anion Ion Chromatographic Solutions, 48

Clinical Laboratory, 49

Serum Reference Materials, 50

Biological Materials, 51

Food and Beverage, 51

Ethanol Solutions, 52

Agriculture, 53

Environmental Materials, 54

Analyzed Gases, 54

Permeation Devices, 57

Analyzed Liquids and Solids, 59

Simulated Rainwaters, 59

Alcohols in Reference Fuels, 60

Sulfur in Fossil Fuels, 61

Trace Elements, 62

Organic Constituents, 63

Industrial Hygiene, 68

Freeze-Dried Urine, 68

Thin Films for X-ray Fluorescence, 69

Materials on Filter Media, 69

Blank Filters, 69

Respirable Quartz, 70

Asbestos, 71

| | |
|--|--|
| Lubricating Materials, 72 | Melting Point, 102 |
| Metallo-Organic Compounds, 72 | Laboratory Thermometer, 103 |
| Lubricating Base Oils, 72 | Thermocouple Material, 103 |
| Catalyst Package for Lubricant Oxidation, 73 | Vapor Pressure, 103 |
| Wear-Metals in Oil, 73 | Thermal Conductivity, 104 |
| Fertilizers, 74 | Thermal Expansion, 104 |
| Ores, 75 | Thermal Resistance, 104 |
| Rocks, Minerals, and Refractories, 79 | Magnetic, 105 |
| Clays, 78 | Magnetic Susceptibility, 105 |
| Mercury in Soil, 79 | Magnetic Moment, 105 |
| Carbides, 81 | Optical, 106 |
| Glasses, 83 | Spectrophotometric, 106 |
| Cements, 84 | Reflectance, 107 |
| Portland Cement Clinkers, 85 | Specular Spectral Reflectance, 107 |
| Trace Elements, 86 | Infrared Reflectance, 107 |
| Nuclear Materials, 87 | Directional-Hemispherical Reflectance, 108 |
| Radiation Dosimetry, 87 | Refractive Index, 108 |
| Fission Track Glasses, 87 | Optical Rotation, 108 |
| Stable Isotopic Materials, 88 | Radioactivity, 109 |
| Alpha-particle, Beta-particle, Gamma-ray, and Electron-capture Solutions, 110 | Alpha-Particle Point-Sources, 111 |
| PHYSICAL PROPERTIES, 89 | Radiocarbon Dating and Ground Water Studies, 111 |
| Ion Activity, 89 | Gaseous Materials, 111 |
| pH, 89 | Gamma-ray and X-ray Point-Sources, 112 |
| pD, 89 | Low-Energy Photon Point-Sources, 112 |
| Ion-Selective Electrodes, 90 | Radium-226 Solutions, 112 |
| Metrology, 90 | Radon Analysis, 112 |
| Scanning Electron Microscope, 90 | Gamma-ray Solutions, 113 |
| Optical Microscope Linewidth-Measurement, 91 | Environmental Natural Matrix Materials for Quality Assurance Testing, 113 |
| Depth Profiling, 91 | Radiopharmaceuticals, 114 |
| Coating Thickness, 92 | Metallurgical, 115 |
| Nonmagnetic Coating on Magnetic Substrate, 92 | Abrasive Wear, 115 |
| Magnetic Coating on Magnetic Substrate, 92 | Corrosion, 115 |
| Solder Thickness, 92 | Electrochemical Potential and Thickness, 115 |
| Coating Weight, 93 | Pitting or Crevice Corrosion, 116 |
| Gold Coating on Nickel, 93 | X-ray Fluorescent Emission Target, 116 |
| Ellipsometry, 93 | X-ray Diffraction, 116 |
| Glass, 94 | Gas Transmission, 117 |
| Chemical Resistance (Durability) of Glass, 94 | Reference Fuels, 118 |
| Electrical Properties of Glass, 94 | Electrical Resistivity and Conductivity, 118 |
| Viscosity, 94 | Metals, 118 |
| Viscosity Fixpoints, 95 | Silicon, 118 |
| Relative Stress Optical Coefficient, 95 | Residual Resistivity Ratio, 119 |
| Glass Liquidus Temperature, 95 | Eddy Current, 119 |
| Density, 96 | Electrolytic Conductance, 119 |
| Microhardness, 96 | Superconducting Critical Current, 119 |
| Ultrasonics, 97 | ENGINEERING MATERIALS, 121 |
| Polymers, 97 | Standard Rubbers and Rubber-Compounding Materials, 121 |
| Molecular Weight (Melt Flow), 97 | Sizing, 122 |
| Rheology, 98 | Particle Size, 122 |
| Heat, 99 | Cement Turbidimetric and Fineness, 122 |
| Calorimetric, 99 | Surface Area of Powders, 123 |
| Combustion, 99 | Performance Standards, 123 |
| Solution, 99 | Socketed Ball Bar, 123 |
| Heat Source, 99 | Dye Penetrant Test Blocks, 123 |
| Enthalpy and Heat Capacity, 100 | Surface Roughness, 124 |
| Differential Scanning Calorimetry, 100 | Charpy V-Notch Test Blocks, 124 |
| Differential Thermal Analysis, 101 | Artificial Flaw for Eddy Current NDE, 124 |
| Superconductive Thermometric Fixed Point Devices, 101 | |
| Freezing Point Materials, 102 | |
| Defining Fixed Points, 102 | |
| Secondary Reference Points, 102 | |

Color, 125
X-ray and Photographic, 125
Magnetic Computer Storage Media, 125
Centerline Drawings for Optical Character
Recognition, 126
NIST Time Software, 127
Fire Research, 128
 Surface Flammability, 128
 Smoke Density Chamber, 128
 Flooring Radiant Panel, 128
 Tape Adhesion Testing, 128

ADDITIONAL INFORMATION, 129
 NIST Special Publications in the 260 Series, 129
 Calibration Service Contacts, 133

INDICES, 135
 Numerical Index of Standard Reference Materials
 (Name and Certificate Date), 135
 Alphabetical Index by Standard Reference
 Material (Name and Category), 145

**Program and
Sales Information** 3

Chemical Composition 15

Physical Properties 89

Engineering Materials 121

Indices 135

Abstract and Key Words

National Institute of Standards and Technology Standard Reference Materials 1990-91 Catalog

This catalog describes the Standard Reference Materials (SRM's) currently available from the National Institute of Standards and Technology (NIST), lists those in preparation, and provides ordering information. The descriptions provide *nominal* values for these SRM's. Certified values are provided in the certificates that accompany each SRM. Price Lists for SRM's are issued as separate supplements to this catalog and include new SRM's as they are issued.

Key Words: analysis, calibration, characterization, composition, concentration, materials, measurement, property, quality assurance, quality control, reference materials, Standard Reference Materials, standardization.



Look for us at the next meeting you attend! Our new exhibit booth, attended by representatives of our sales and technical staff, is scheduled to be at about a dozen professional meetings and expositions each year.

The Engineering Mechanics Building is now the home of the Standard Reference Materials Program. We look forward to continued and improved service to the technical community in our new offices and materials preparation and packaging facilities.



Program Information

The National Institute of Standards and Technology (NIST) offers for sale over 1,000 different materials through its Office of Standard Reference Materials. These materials are primarily Standard Reference Materials (SRM's) certified for their chemical composition, chemical property, or physical property, but include other reference materials. All materials bear distinguishing names and numbers by which they are permanently identified. Thus, each material bearing a given description is identical (within the specified limits) to every other sample bearing the same designation—with the exception of individually certified items, which are further identified by serial number.

Definitions

From "Terms and definitions used in connection with reference materials," ISO Guide 30-1981 (E):

1. "Reference Material (RM): A material or substance one or more properties of which are sufficiently well established to be used for the calibration of an apparatus, the assessment of a measurement method, or for assigning values to materials."
2. "Certified Reference Material (CRM): A reference material one or more of whose property values are certified by a technically valid procedure, accompanied by or traceable to a certificate or other documentation which is issued by a certifying body."

NIST Standard Reference Materials (SRM's): Certified reference materials issued by NIST. These are well-characterized materials produced in quantity to improve measurement science. SRM's are certified for specific chemical or physical properties, and are issued by NIST with certificates that report the results of the characterization and indicate the intended use of the material. They are prepared and used for three main purposes:

- (1) To help develop accurate methods of analysis (reference methods);
- (2) To calibrate measurement systems used to:
 - (a) facilitate exchange of goods,
 - (b) institute quality control,
 - (c) determine performance characteristics, or
 - (d) measure a property at the state-of-the-art limit; and
- (3) To assure the long-term adequacy and integrity of measurement quality assurance programs.

NIST certified values are obtained by one of three routes of measurement:

- (1) A previously validated reference method,
- (2) Two or more independent, reliable measurement methods, or
- (3) A network of cooperating laboratories, technically competent and thoroughly knowledgeable with the material being tested.

These measurement routes are described in "The Role of Standard Reference Materials in Measurement Systems," NBS Monograph 148, 54 pages (Jan 1975).

Reference Materials (RM's) listed in this catalog are sold by, but not certified by, NIST. They meet the ISO definition for RM's, and many meet the definition for CRM's. The documentation issued with these materials is either a:

(1) "Report of Investigation," the sole authority of which is the author of the report. RM's are intended to further scientific or technical research on that particular material. The principal consideration in issuing an RM is to provide a homogeneous material so that investigators in different laboratories are assured that they are investigating the same material.

(2) "Certificate," issued by the certifying agency (other than NIST), e.g., other national laboratories, other government agencies, other standardizing bodies, or other non-profit organizations. When deemed to be in the public interest and when alternate means of national distribution do not exist, NIST acts as the distributor for such materials. This service is available to organizations that qualify and have the reference materials that would help meet national measurement needs.



Stan Rasberry, the new Acting Director of Measurement Services, displays one of the many SRM's developed in the program of the Office of Standard Reference Materials (OSRM). The enthusiasm and decisiveness Stan exercised while Chief of OSRM will continue to be an asset to the SRM program as well as to the additional programs to which he will now provide leadership — Standard Reference Data, Physical Measurement Services and the National Voluntary Laboratory Accreditation Programs.

SRM Catalog

New catalogs of NIST Standard Reference Materials are published approximately every 2 years, listing materials available and materials in preparation, and deleting discontinued materials. Catalog supplements (Price Lists) are issued simultaneously with new catalogs and approximately every 6 months to keep the catalog current between editions. These supplements list current prices, and reflect any changes in material availability—listing new and renewed materials and dropping discontinued ones.

The numerical values given in this catalog to describe the materials' properties are NOMINAL values only and are to be used only as guides in selecting SRM's. They are NOT TO BE USED in place of the values given on the certificate issued with the materials.

Two indices are provided for user convenience. The first is a numerical index that lists the numbers, names, and certificate dates of the materials in the catalog. The second is an alphabetical index that lists categories of materials, elements, and names of materials.

Preparation and Availability of Standard Reference Materials

New and renewal SRM's are being prepared continually. These SRM's are included in the next edition of the catalog and its supplements. Prospective users that have requested that their names be added to the SRM mail list are notified as these new items become available. To have your name placed on this mail list, please write to the address given below.

Renewal SRM's are intended to be completed before the supply of an existing SRM is exhausted. This is not always possible and an SRM may be out-of-stock for a time. When this occurs, those ordering the material are so notified and possible substitutes (if any) are suggested. When a renewal is issued, customers who have ordered the previous lot are promptly notified of the price and availability of the renewal. If little demand exists or if an alternate source of supply becomes available, production of an SRM may be discontinued permanently.

Renewal SRM's are not identical to the preceding lot; however, they meet the same specifications and can be used for the same purpose. For example, the first 0.1 percent carbon Bessemer steel was prepared in 1909 (Standard Sample No. 8). Since then a number of renewals, 8a, 8b, 8c, etc., were prepared. The current SRM 8j, Bessemer Steel (Simulated), 0.1% C, represents the eleventh lot of the material. Each lot differs somewhat in detailed analysis, thus the use of the specific certificate for that lot is essential.

Guide for Requesting Development of Standard Reference Materials

The National Institute of Standards and Technology has the function to develop, produce, and distribute Standard Reference Materials (SRM's) that provide a basis for comparison of measurements on materials, and that aid in the control of production processes. To perform this function, the Office of Standard Reference Materials evaluates the requirements of science, industry, and government for carefully characterized reference materials, and directs their production and distribution.

NIST currently has over 1,000 SRM's available, about 100 new ones in preparation, and requests for the production of many others.

To be an SRM, a candidate material must meet one or more of these criteria:

1. It would permit users to attain more accurate measurements.
2. Its production elsewhere would not be economically or technically feasible.
3. It would be an industry-wide standard for commerce from a neutral source not otherwise available.
4. Its production by NIST would provide continued availability of a well-characterized material important to science, industry, or government.

NIST has recognized and responded to requests to enlarge the scope of the SRM program to include all types of well-characterized materials for use in calibrating measurement systems, or for producing scientific data that can be referred to a common base. However, the requests for new SRM's greatly exceed the Institute's capacity to produce and certify such materials. Consequently, requests for new SRM's of limited use, or for which the need is not very great, are deferred in favor of requests that clearly show a critical need. To determine which requests receive top priority, NIST needs and uses information supplied by industry and such interested organizations as the American National Standards Institute, American Nuclear Society, American Petroleum Institute, American Society for Testing and Materials, etc.

Accordingly, while NIST welcomes all requests for developing new SRM's, both NIST and industry would be helped if such requests provide information that permit objective assessment of the urgency and importance of the proposed new reference materials.

Requests for the development of new Standard Reference Materials should provide information such as listed below.

1. Short title of the proposed SRM.
2. Purpose for which the SRM would be used.
3. Reasons why the SRM is needed.
4. Special characteristics and requirements for the material. Include additional requirements and reasons if more than one SRM is necessary for standardization in this area.
5. An estimate of the probable present and future (6-10 year) demand for such an SRM in your operations and elsewhere. (National and international estimates are useful.)
6. Whether such an SRM, or a similar one, could be produced or obtained from a source other than NIST; and if so, justify its preparation by NIST.
7. Miscellaneous pertinent information to aid justification for the SRM, such as: (a) an estimate of the potential range of application, monetary significance of the measurement affected, scientific and technological significance including, when feasible, estimates of the impact upon industrial productivity, growth, quality assurance or control, and (b) supporting letters from industry leaders, trade organizations, interested committees, and others.

All such requests should be addressed to:

Office of Standard Reference Materials

ATTN: SRM Development

Room 205, Building 202

National Institute of Standards and Technology

Gaithersburg, MD 20899

Ordering Standard Reference Materials

General

Purchase orders for all SRM's should be addressed to:

Office of Standard Reference Materials
Room 205, Building 202
National Institute of Standards and Technology
Gaithersburg, MD 20899
Telephone: (301) 975-OSRM [6776]
FTS: 879-OSRM [6776]
Fax number: (301) 948-3730
Telex: TRT197674NIST UT

All orders should give the number of units, catalog number, and name of the material requested. For example: "1 each, SRM 79a, Fluorspar (Customs Grade)." The materials described in this catalog are sold only in the units listed or multiples thereof.

Acceptance of an order does not imply acceptance of any provisions set forth in the order contrary to the policy, practice, or regulations of the National Institute of Standards and Technology or the U.S. Government.

In general, orders received for "out-of-stock" material will be filled with the renewal material, if available; otherwise they will be canceled. Customers are notified when an order is canceled; and their names are placed on a notification list. This list is used when a renewal material is issued to notify customers of the price and availability of the item. Customers so notified are requested to submit a new order if they still want the item.

For some individually certified SRM's, production lots are small and may entail frequent stock outages. In these cases, the notification list is used to fill orders on a "first come, first served" basis. NOTE: For such SRM's, customers are notified that the SRM is again available and are requested to confirm their original purchase orders.



Lee Best, Manager of Sales and Marketing, has developed an excellent sales and marketing staff through the thoroughness and initiative of her leadership. We look forward to meeting many of you at one of the many exhibits Lee has scheduled for the Office of Standard Reference Materials.



Rosemary Blasingame, who displays enthusiasm and endurance in every responsibility she receives, serves as SRM assistant, performing many tasks in sales and marketing, including foreign sales.

Terms

Prices quoted are in U.S. dollars (\$), and are published in the catalog supplements (price lists). When price lists are issued, they are sent to persons or organizations on the SRM mail list. These prices are subject to change without notice and orders will be billed for the prices in effect at the time of shipment. No discounts are given on purchases of SRM's or RM's.

Remittances of the purchase price need not accompany the purchase order. Payment of invoices is expected within 30 days of the receipt of the invoice. Payment on foreign orders may be made by any of the following:

- a. Banker's draft against U.S.A. bank,
- b. Bank to bank transfer to U.S.A. bank,
- c. Cash against documents,
- d. Sight draft,
- e. International money order, or
- f. UNESCO coupons.

Letters of credit: If a letter of credit or any method of payment other than those listed above is to be used, the services of an agent in the United States must be secured to act in your behalf. Your agent would purchase the material and our invoice would indicate that the agent is the purchaser. The material would be shipped to your agent, who would tranship in accordance with your instructions.

NOTE: A customer identification number; i.e., social security number (EIN) for consumer customers; tax identification number (TIN) for commercial customers; or agency location code (ALC) for government customers, is required with each order.



Phyllis Wagner, an SRM assistant, provides consistency and dependability to the delivery of many SRM's through her expertise in regulations and procedures for the shipment of hazardous materials.



Cindy Leonard provides customer sales and services with a joyfulness and friendliness that is appreciated by our customers.

Late Charges

Unless otherwise notified, payment is due within 30 days of shipment of the order to the customer. U.S. Treasury regulations require that late charges be assessed for each 30-day period, or portion thereof, that the payment is overdue.

Proforma Invoice (Price Quotation)

Proforma service will be provided only to those requiring such service.

Domestic Shipments

Shipments of material (except for certain restricted categories and refrigerated items) intended for the United States and Canada are normally shipped prepaid, providing the parcel does not exceed the weight limitations prescribed by postal laws and regulations. Refrigerated items are shipped prepaid air express with shipping costs added to the invoice.

Donna Clark, systems analyst, displays persistence and attention to details in the development and maintenance of the OSRM computer programs for sales and inventory control.





Jody Hines (left) and Sandy Barber combine their characteristic qualities of diligence, determination, flexibility and cooperativeness to the task of updating information for the preparation and publication of the new 1990-91 SRM Catalog.

Foreign Shipments

The regulations of various nations covering the importation of SRM's differ widely; any attempt to list all possible variations would be impractical. Therefore, where shipping practices outlined below do not apply, purchasers will be informed of the best method of shipment for their countries.

Most foreign orders will be shipped by prepaid International Air Parcel Post. Exceptions are those items in restricted categories, those items requiring refrigeration, and shipments exceeding parcel post weight limits. These exceptions will be shipped FOB Gaithersburg, MD, unless an agent (shipping or brokerage firm) located in the United States is used. When an agent is required, the purchaser will be notified and will be requested to obtain the services of one and inform us of the agent's name and address. In such cases, the material will be packed for overseas shipment and will be forwarded to the agent FOB Gaithersburg, MD.

Documentation

The documents we furnish are:

- a. Two commercial invoices,
- b. Two sight drafts,
- c. Two packing slips, and
- d. An air waybill for air shipments.

(All documents are printed in English.)

If documents other than those listed above are required, the services of an agent in the United States will be needed to purchase and ship the material.

NOTE: Orders and inquiries submitted in English will be processed more rapidly than those requiring translation.

Certified Reference Materials From Other Sources

Special Nuclear Materials

On October 1, 1987, the New Brunswick Laboratory began issuing special nuclear reference materials as NBL Certified Reference Materials (CRM's). These CRM's include the plutonium and uranium assay and isotopic materials previously issued by the National Institute of Standards and Technology. All orders or inquiries should be addressed to:

U.S. Department of Energy
New Brunswick Laboratory
Attn: Reference Materials Sales
9800 S. Cass Avenue, Bldg. 350
Argonne, IL 60439
(312) 972-2767

International CRM's

Certified reference materials (CRM's) are available from many sources. The International Organization for Standardization (ISO), through its Council Committee on Reference Materials (REMCO), has prepared an international Directory of Certified Reference Materials. Inquiries may be directed to:

Dr. M. Parkany
Secretary for REMCO
International Organization for Standardization
1, Rue de Varembe
Case Postale 56
1211 Geneva 20
Switzerland

The International Union of Pure and Applied Chemistry (IUPAC), through its Commission on Physicochemical Measurements and Standards, issues a catalog of CRM's that are useful for the realization of physicochemical properties. It also has prepared a number of related documents. The current IUPAC edition is: "Physicochemical Measurements: Catalogue of Reference Materials from National Laboratories," Revised 1976, Pure & Appl. Chem., 48, 503-515 (1976).

Julie Frum makes the Division Office a pleasant and productive place by her dependability and warm spirit of hospitality.



Other Services of the National Institute of Standards and Technology

Calibration and Related Measurement Services

The measurement services of NIST include the calibration of standards, test of instruments, and certain interlaboratory testing programs. These services are described in NIST Special Publication 250, National Institute of Standards and Technology Calibration Services Users Guide, 1989 ed. [Available from the Superintendent of Documents, U.S. Government Printing Office, Washington, DC 20402.]

An abbreviated list of the services offered through this program appears under Additional Information. These services are performed at either the NIST Washington laboratories (Gaithersburg, MD) or those in Boulder, CO. For additional information on available measurement services, consult Special Publication 250 or write to:

Office of Physical Measurement Services
Room B362 Physics Building
National Institute of Standards and Technology
Gaithersburg, MD 20899

Telephone: (301) 975-2002

Requests for measurement services available in Boulder should be addressed to:

Measurement Services Clerk
National Institute of Standards and Technology
Boulder, CO 80303

Telephone: (303) 497-3753

Office of Weights and Measures

The NIST Office of Weights and Measures operates a Type Evaluation Program which provides for an evaluation of (1) prototype weighing and measuring devices to determine compliance with the requirements of NBS Handbook 44, "Specifications, Tolerances, and Other Technical Requirements for Commercial Weighing and Measuring Devices," (2) standards to determine compliance with the requirements of NBS Handbook 105-1, 105-2, 105-3, "Specifications and Tolerances for Reference Standard and Field Standard Weights and Measures." This program may be used by manufacturers and weights and measures officials in determining the acceptability of devices for commercial use or the suitability of reference and field standards. For information on programs of NIST and the States, write or telephone:

Office of Weights and Measures
Room A617 Administration Building
National Institute of Standards and Technology
Gaithersburg, MD 20899

Telephone: (301) 975-4004

Proficiency Sample Programs

General information on the Proficiency Sample Programs may be obtained from:

Materials Reference Laboratories
National Institute of Standards and Technology
Gaithersburg, MD 20899

Telephone: (301) 975-6704

Information is available on the following programs:

Proficiency Sample Programs for Hydraulic Cements and Portland Cement Concrete
Proficiency Sample Programs for Soils, Aggregates, and Bituminous Materials
Inspection of Cement and Concrete Testing Laboratories
Inspection of Soils and Bituminous Testing Laboratories

Accreditation of Testing Laboratories

General information about the National Voluntary Laboratory Accreditation Program (NVLAP) or application packages may be obtained from:

Manager, Laboratory Accreditation
Room A531 Administration Building
National Institute of Standards and Technology
Gaithersburg, MD 20899

Telephone: (301) 975-4016

National Center for Standards and Certification Information

The National Center for Standards and Certification Information (NCSCI) contains title information or full texts for more than 240,000 engineering or related standards issued by U.S. technical societies, professional organizations, and trade associations; State purchasing offices; U.S. Federal Government agencies; and major foreign national and international standardizing bodies. NCSCI publishes general and specific indices of standards. Information services which are free consist of searching Key-Word-In-Context (KWIC) Indices to determine whether any published standards, specifications, codes, test methods, or recommended practices exist for a given item or product. Inquiries should be directed to:

National Center for Standards and Certification Information
Room A629 Administration Building
National Institute of Standards and Technology
Gaithersburg, MD 20899

Telephone: (301) 975-4040

National Standard Reference Data System

The National Standard Reference Data System (NSRDS) is a nationwide program established to compile and critically evaluate quantitative physical science data and assure its availability to the technical community. The program publishes compilations of critically evaluated data, critical reviews of experimental techniques, and bibliographies. A complete list of NSRDS publications is available from the Office of Standard Reference Data (OSRD). OSRD responds to queries within the scope of the program by providing references, referrals, documentation, or data, as available. Inquiries or requests for information should be directed to:

Office of Standard Reference Data
Room A323 Physics Building
National Institute of Standards and Technology
Gaithersburg, MD 20899

Telephone: (301) 975-2208

A sampling of Environmental SRM's which illustrate the diversity of SRM matrices and the creativity that must be exercised in the packaging of SRM's to maintain the integrity of the sample and its certified properties.



Roger Brown (left) contributes generous, industrious leadership to the shipping activities. Carlton Fisher (right front), whose qualities of orderliness and neatness combine with the responsibility and availability provided by Gary Proulx (center back) and Jim Fort (right back) to create a productive, efficient shipping department that sent out about 45,000 units this past year.



A diamond saw wafering machine, one of a wide variety of equipment used in the preparation of SRM materials, is being used to cut glass discs from glass rods.



Dan Swearingen, shown here using an automated band saw to prepare samples for one of the metal SRM's, is appreciated for the dependability and diligence he applies to the activities of sample preparation.

Chemical Composition

Ferrous Alloys

Steels (Chip Form)

These SRM's are for checking chemical methods of analysis. They consist of steel alloys selected to provide a wide range of analytical values for elements. They are furnished in 150-gram units (unless otherwise noted) as chips usually sized between 0.4 to 1.2 mm, prepared from selected portions of commercial ingots.

Plain Carbon Steels

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | | | | | |
|------|-----------------------------------|---|-------|-------|--------|-------|-------------------|-------|-------|-------------------|-------|
| | | C | Mn | P | S | Si | | | | | |
| 8j | Bessemer (simulated), 0.1C | 0.081 | 0.505 | 0.095 | Grav | 0.077 | 0.058 | | | | |
| 11h | BOH, 0.2C | 0.200 | 0.510 | 0.010 | | 0.026 | 0.21 ₁ | | | | |
| 12h | BOH, 0.4C | 0.407 | 0.842 | 0.018 | | 0.027 | 0.235 | | | | |
| 13g | Carbon, 0.6C | 0.613 | 0.853 | 0.006 | | 0.031 | 0.35 ₅ | | | | |
| 14f | BOH, 0.8C | 0.753 | 0.410 | 0.009 | | 0.039 | 0.172 | | | | |
| 15h | BOH, 0.1C | 0.076 | 0.373 | 0.005 | | 0.019 | 0.008 | | | | |
| 16f | BOH, 1.1C | 0.97 | 0.404 | 0.014 | | 0.026 | 0.214 | | | | |
| 19h | BES, 0.2C | 0.215 | 0.393 | 0.016 | | 0.022 | 0.211 | | | | |
| 20g | AISI 1045 | 0.462 | 0.665 | 0.012 | | 0.028 | 0.305 | | | | |
| 105 | High-Sulfur (Carbon Only) | 0.193 | | | | | | | | | |
| 152a | BOH, 0.5C (Tin bearing) | 0.486 | 0.717 | 0.012 | | 0.030 | 0.202 | | | | |
| 178 | Basic Oxygen Furnace 0.4C | 0.395 | 0.824 | 0.012 | | 0.014 | 0.163 | | | | |
| 335 | BOH, 0.1C (Carbon only) 300 g | 0.092 | | | | | | | | | |
| 337a | BOH, 1.1C (Carbon & Sulfur) 300 g | 0.969 | | | | 0.024 | | | | | |
| 368 | AISI 1211 | 0.089 | 0.82 | 0.084 | | 0.132 | 0.007 | | | | |
| SRM | | Cu | Ni | Cr | V | Mo | Co | Ti | Sn | Al (total) | N |
| 8j | | 0.020 | 0.113 | 0.047 | 0.015 | 0.038 | | | | | |
| 11h | | 0.061 | 0.028 | 0.025 | 0.001 | | | 0.004 | | | |
| 12h | | 0.073 | 0.032 | 0.074 | 0.003 | 0.006 | | | | (0.038) | 0.006 |
| 13g | | 0.066 | 0.061 | 0.050 | 0.001 | | | | | 0.04 ₈ | |
| 14f | | 0.072 | 0.053 | 0.070 | 0.002 | 0.013 | | | | 0.060 | |
| 15g | | 0.013 | 0.017 | 0.018 | <0.001 | 0.009 | | | | 0.061 | |
| 16f | | 0.006 | 0.008 | 0.020 | 0.002 | 0.003 | 0.003 | | | | |
| 19h | | 0.466 | 0.248 | 0.173 | 0.003 | 0.038 | | | | 0.002 | |
| 20g | | 0.034 | 0.034 | 0.036 | 0.002 | 0.008 | | | | 0.040 | |
| 152a | | 0.023 | 0.056 | 0.046 | 0.001 | 0.036 | | | 0.032 | | |
| 178 | | 0.032 | 0.010 | 0.016 | 0.001 | 0.003 | | | | | |
| 368 | | 0.010 | 0.008 | 0.030 | 0.001 | 0.003 | | | | 0.010 | |

Values in parentheses are not certified, but are given for information only.

Low Alloy Steels

| SRM | Type | (Other Forms) | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|---------------------------|-------------------|---|-------------------|---------|------------|-------------------|-----------|
| | | | C | Mn | P | S | Si | Cu |
| 30f | Cr-V (SAE 6150) | | 0.490 | 0.79 | 0.011 | 0.009 | 0.283 | 0.074 |
| 32e | Ni-Cr (SAE 3140) | | 0.409 | 0.798 | 0.008 | 0.022 | 0.021 | 0.278 |
| 33e | Ni-Mo (SAE 4820) | | 0.186 | 0.525 | 0.005 | 0.009 | 0.262 | 0.070 |
| 36b | Cr-Mo | | 0.114 | 0.404 | 0.007 | 0.019 | 0.258 | 0.179 |
| 72g | AISI 4130 | | 0.278 | 0.492 | 0.009 | 0.014 | 0.223 | 0.011 |
| 100b | Manganese (SAE T1340) | | 0.397 | 1.89 | 0.023 | 0.029 | 0.210 | 0.064 |
| 106b | Cr-Mo-Al (Nitrally G) | | 0.326 | 0.506 | 0.008 | 0.016 | 0.017 | 0.274 |
| 125b | High-Silicon | 1134 | 0.028 | 0.278 | 0.029 | 0.008 | 2.89 | 0.071 |
| 129c | High-Sulfur (SAE 112) | | 0.125 | 0.769 | 0.076 | 0.245 | 0.020 | 0.013 |
| 131d | Low Carbon-Silicon (100g) | | 1218 | 0.0035 | | 0.0011 | | |
| 139b | Cr-Ni-Mo (AISI 8640) | 1222 | 0.403 | 0.778 | 0.013 | 0.019 | 0.242 | 0.097 |
| 155 | Cr-W | | 0.905 | 1.24 | 0.015 | 0.010 | 0.322 | 0.083 |
| 163 | Cr (100g) | | 0.933 | 0.897 | 0.007 | 0.027 | 0.488 | 0.087 |
| 291 | Cr-Mo (ASTM A213) | | 0.177 | 0.55 _o | 0.008 | 0.020 | 0.23 _o | 0.047 |
| 293 | Cr-Ni-Mo (AISI 8620) | | 0.222 | 0.96 _o | 0.018 | 0.022 | 0.30 _o | 0.032 |
| SRM | Ni | Cr | V | Mo | Sn | Al (total) | N | Other |
| 30f | 0.070 | 0.945 | 0.182 | | | | 0.010 | |
| 32e | 1.19 | 0.678 | 0.002 | 0.023 | (0.011) | | 0.009 | |
| 33e | 3.36 | 0.068 | (0.001) | 0.224 | (0.002) | 0.030 | | |
| 36b | 0.203 | 2.18 | 0.004 | 0.996 | | | | |
| 72g | 0.016 | 0.905 | 0.003 | 0.170 | | (0.041) | (0.008) | |
| 100b | 0.030 | 0.063 | 0.003 | 0.237 | | | 0.004 | |
| 106b | 0.217 | 1.18 | 0.003 | 0.199 | | | 1.07 | |
| 125b | 0.038 | 0.019 | | 0.008 | 0.003 | 0.329 | | Ca 0.0051 |
| 129c | 0.251 | 0.014 | 0.012 | 0.002 | | | | |
| 139b | 0.510 | 0.488 | 0.004 | 0.182 | | | 0.007 | |
| 155 | 0.100 | 0.485 | 0.014 | 0.039 | | | | W 0.517 |
| 163 | 0.081 | 0.982 | | 0.029 | | | 0.007 | |
| 291 | 0.065 | 1.33 | | 0.53 _s | | 0.002 | | |
| 293 | 0.48 _o | 0.51 _o | 0.004 | 0.20 ₄ | | 0.039 | | |

Values in parentheses are not certified, but are given for information only.

Special Low Alloy Steels

| SRM | Type | (Other forms) | Chemical Composition (Nominal Weight Percent) | | | | | | | | | | |
|------|----------------------|----------------------|---|--------------------|-------------------|-------------------|-------------------|---------------------|-------------------|---------------------|----------|-----------|---------------------|
| | | | C | Mn | P | S | Si | Cu | Ni | Cr | | | |
| 361 | AISI 4340 | 661,1095,1261a | 0.383 | 0.66 | 0.014 | 0.0143 | 0.222 | 0.042 | 2.00 | 0.69 ₄ | | | |
| 362 | AISI 94B17 (Mod) | 662,1096 | 0.160 | 1.04 | 0.041 | 0.0360 | 0.39 | 0.50 | 0.59 | 0.30 | | | |
| 363 | Cr-V (Mod) | 663,1097,1263a | 0.62 | 1.50 | 0.02 ₉ | 0.0068 | 0.74 | 0.10 | 0.30 | 1.31 | | | |
| 364 | High Carbon (Mod) | 664,1098,1264a | 0.87 | 0.25 ₅ | 0.01 | 0.0250 | 0.06 ₅ | 0.24 ₉ | 0.14 ₄ | 0.06 ₃ | | | |
| 365 | Iron, Electrolytic | 665,1099,1265a | (Supplies exhausted, see other forms.) | | | | | | | | | | |
| 2161 | A IN PREP | 1761 | (1.03) | (0.68) | (0.043) | (0.033) | (0.19) | (0.30) | (1.99) | (0.21) | | | |
| 2162 | B IN PREP | 1762 | (0.34) | (2.03) | (0.036) | (0.03) | (0.36) | (0.12) | (1.15) | (0.92) | | | |
| 2163 | C IN PREP | 1763 | (0.20) | (1.59) | (0.012) | (0.022) | (0.65) | (0.045) | (0.49) | (0.51) | | | |
| 2164 | D IN PREP | 1764 | (0.59) | (1.22) | (0.023) | (0.012) | (0.06) | (0.51) | (0.20) | (1.50) | | | |
| 2165 | E | 1765 | (0.006) | (0.14) | (0.007) | (0.004) | (0.005) | (0.002) | (0.15) | (0.05) | | | |
| 2166 | F | 1766 | (0.015) | (0.06) | (0.004) | (0.002) | (0.01) | (0.014) | (0.02) | (0.02) | | | |
| 2167 | G | 1767 | (0.051) | (0.02) | (0.005) | (0.009) | (0.02) | (0.002) | (0.001) | (0.001) | | | |
| SRM | V | Mo | W | Co | Ti | As | Sn | Al (total) | Nb | Ta | Zr | N | Ca |
| 361 | 0.011 | 0.19 | 0.017 | 0.032 | 0.020 | 0.017 | 0.010 | 0.02 ₁ | 0.022 | 0.020 | 0.009 | (0.0037) | 0.0001 ₀ |
| 362 | 0.040 | 0.068 | 0.20 | 0.30 | 0.084 | 0.09 ₂ | 0.016 | 0.09 ₃ | 0.29 | 0.20 | 0.19 | (0.00404) | 0.0002 ₁ |
| 363 | 0.31 | 0.028 | 0.046 | 0.048 | 0.050 | 0.010 | 0.10 ₄ | 0.24 | 0.049 | (0.053) | 0.049 | (0.0041) | 0.0002 ₂ |
| 364 | 0.10 ₅ | 0.49 | 0.10 | 0.15 | 0.24 | 0.05 ₂ | 0.008 | (0.008) | 0.15 ₇ | 0.11 | 0.068 | (0.0032) | 0.00003 |
| 365 | 0.0006 | 0.0050 | | 0.007 ₀ | 0.0006 | (0.0002) | (0.0002) | (0.0007) | | | | 0.0013 | |
| 2161 | (0.05) | (0.10) | | (0.03) | (0.17) | (0.01) | (0.04) | (0.05) | (0.02) | (0.05) | (0.01) | | |
| 2162 | (0.20) | (0.36) | | (0.06) | (0.1) | (0.02) | (0.04) | (0.07) | (0.07) | (0.02) | (0.03) | | |
| 2163 | (0.31) | (0.49) | | (0.09) | (0.31) | (0.05) | (0.008) | (0.05) | (0.10) | (0.01) | (0.04) | | |
| 2164 | (0.11) | (0.20) | | (0.01) | (0.03) | (0.009) | (0.015) | (0.01) | (0.04) | (0.03) | (<0.001) | | |
| 2165 | (0.004) | (0.005) | | (0.002) | (0.005) | (0.006) | (0.002) | (0.006) | (<0.002) | (0.005) | (0.001) | | |
| 2166 | (0.009) | (0.004) | | (0.003) | (0.001) | (<0.002) | (0.001) | (0.01) | (0.003) | (0.01) | (0.001) | | |
| 2167 | (0.03) | (0.02) | | (0.006) | (0.01) | (<0.002) | (0.007) | (0.004) | (0.007) | (<0.005) | (0.004) | | |
| SRM | B | Pb | Sb | Bi | Ag | Se | Te | Ce | La | Nd | Fe | | |
| 361 | 0.0003 ₇ | 0.00002 ₅ | 0.0042 | (0.0004) | 0.0004 | (0.004) | (0.0006) | 0.0040 | (0.001) | 0.0007 ₅ | (95.6) | | |
| 362 | 0.0025 | 0.0004 ₈ | 0.013 | (0.002) | 0.0011 | (0.0012) | (0.0011) | 0.0019 | (0.001) | 0.0007 ₅ | (95.3) | | |
| 363 | 0.0007 ₈ | 0.0018 ₆ | 0.002 | (0.0008) | 0.0037 | (0.00016) | (0.0009) | 0.0030 | (0.002) | 0.0012 | (94.4) | | |
| 364 | 0.0106 | 0.023 ₀ | 0.034 | (0.0009) | (0.00002) | (0.00021) | (0.0002) | 0.0005 ₇ | (0.0002) | 0.0001 ₈ | (96.7) | | |
| 365 | 0.00012 | 0.00001 ₉ | | | | | | | | | 99.90 | | |
| SRM | Mg | Zn | Pr | Ge | O | H | Au | Hf | Sr | | | | |
| 361 | 0.0002 ₆ | (0.0001) | (0.0003) | [0.006] | (0.0009) | (<0.0005) | (<0.00005) | (0.0002) | | | | | |
| 362 | 0.0006 ₈ | (0.0005) | (0.0003) | [0.002] | (0.00107) | (<0.0005) | (<0.00005) | (0.0003) | | | | | |
| 363 | 0.0006 ₂ | (0.0004) | (0.0004) | [0.010] | (0.00066) | (<0.0005) | 0.0005 | (0.0005) | | | | | |
| 364 | 0.00016 | [0.001] | (0.0001) | [0.003] | (0.0010) | (<0.0005) | 0.0001 | (0.0013) | | | (0.001) | | |

Values in parentheses are not certified, but are given for information only.

Values in brackets are approximate values from heat analysis and are given for information only.

High Alloy Steels

| SRM | Type | (Other Forms) | Chemical Composition (Nominal Weight Percent) | | | | | | |
|------|--|---------------|---|-------|--------|-------|--------|-------|-------|
| | | | C | Mn | P | S | Si | Cu | |
| 126c | High-Nickel (36% Ni) | 1158 | 0.025 | 0.468 | 0.004 | Grav | Comb | 0.194 | 0.040 |
| 344 | Cr-Ni (Mo Precipitation Hardening) | | 0.069 | 0.57 | 0.018 | | 0.005 | 0.395 | 0.106 |
| 345 | Cr-Ni (Cu Precipitation Hardening) | | 0.048 | 0.224 | 0.018 | 0.012 | 0.012 | 0.610 | 3.44 |
| 346a | Valve Steel | 1233 | 0.502 | 9.16 | 0.031 | | 0.002 | 0.219 | 0.375 |
| 348a | High Temperature Alloy (A286) Ni-Cr | 1230 | 0.044 | 0.64 | 0.023 | | 0.0007 | 0.43 | 0.14 |
| 868 | High Temperature Alloy Fe-Ni-Co | 1250 | 0.022 | 0.052 | <0.003 | | 0.0025 | 0.097 | 0.022 |

| SRM | Ni | Cr | V | Mo | Co | Ti | Al (total) | Nb | Ta | B | Fe |
|------|-------|-------|-------|-------|--------|----------|------------|--------|------------|----------|---------|
| 126c | 36.05 | 0.062 | 0.001 | 0.011 | 0.008 | | | | | | |
| 344 | 7.28 | 14.95 | 0.040 | 2.40 | | 0.076 | 1.16 | | | | |
| 345 | 4.24 | 16.04 | 0.041 | 0.122 | 0.089 | | | | | | |
| 346a | 3.43 | 21.08 | 0.096 | 0.237 | (0.05) | (<0.001) | (0.001) | 0.231 | 0.002 | | |
| 348a | 24.2 | 14.8 | 0.23 | 1.18 | 0.15 | 2.12 | 0.24 | (0.01) | Sn (0.008) | (<0.001) | N 0.415 |
| 868 | 37.78 | 0.077 | 0.077 | 0.014 | 16.1 | 1.48 | 0.99 | (0.07) | W (0.07) | 0.0055 | (55.2) |
| | | | | | | | | 2.99 | 0.003 | 0.0078 | 40.5 |

Values in parentheses are not certified, but are given for information only.



Curt Fales, a valuable member of the materials preparation staff because of his resourcefulness and attention to detail, is shown here using the diamond saw wafering machine.

Stainless Steels

| SRM | Type | (Other Forms) | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|---|-------------------|---|------------------|--------|----------|--------|-----------|
| | | | C | Mn | P | S | Si | Cu |
| 73c | Cr (SAE 420) | | 0.310 | 0.330 | 0.018 | 0.036 | 0.181 | 0.080 |
| 101g | Stainless (AISI 304 L) (100g) | | 0.0136 | 0.085 | 0.007 | 0.0078 | 1.08 | 0.029 |
| 121d | Cr-Ni-Ti (AISI 321) | 1171 | 0.067 | 1.80 | 0.019 | 0.013 | 0.54 | 0.121 |
| 123c | Cr-Ni-Nb (AISI 348) | 1172 | 0.056 | 1.7 _s | 0.024 | 0.014 | 0.59 | 0.103 |
| 133b | Cr-Mo | | 0.128 | 1.07 | 0.018 | 0.328 | 0.327 | 0.080 |
| 160b | Cr-Ni-Mo (AISI 316) | 1155 | 0.044 | 1.64 | 0.020 | 0.016 | 0.509 | 0.172 |
| 166c | Low Carbon (AISI 3162) Carbon Only (100g) | | 0.0078 | | | | | |
| 339 | Cr-Ni-Se (SAE 303Se) | | 0.052 | 0.738 | 0.129 | 0.013 | 0.654 | 0.199 |
| 343a | Cr-Ni (AISI 431) | 1219 | 0.149 | 0.42 | 0.026 | 0.001 | 0.545 | 0.162 |
| 367 | Cr-Ni (AISI 446) | 1267 | 0.093 | 0.315 | 0.018 | 0.016 | 0.58 | |
| SRM | Ni | Cr | V | Mo | Co | Ti | Nb | Ta |
| 73c | 0.246 | 12.82 | 0.030 | 0.091 | | | | |
| 101g | 10.00 | 18.46 | 0.041 | 0.004 | 0.09 | | | |
| 121d | 11.17 | 17.4 _s | | 0.165 | 0.10 | 0.342 | | |
| 123c | 11.3 _s | 17.4 _s | | 0.22 | 0.12 | | 0.65 | <0.001 |
| 133b | 0.230 | 12.63 | 0.071 | 0.052 | | | | |
| 160b | 12.26 | 18.45 | 0.047 | 2.38 | 0.101 | | | 0.001 |
| 339 | 8.89 | 17.42 | 0.058 | 0.248 | 0.096 | | | 0.247 |
| 343a | 2.16 | 15.64 | 0.056 | 0.164 | (0.04) | (<0.001) | (0.01) | (<0.0001) |
| 367 | 0.29 | 24.19 | 0.08 | | | | | 0.078 |
| | | | | | | | | 0.168 |

Values in parentheses are not certified, but are given for information only.

Tool Steels

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|----------------------|---|-------|-------|-------|-------|-------|
| | | C | Mn | P | S | Si | Cu |
| 50c | W-Cr-V | 0.719 | 0.342 | 0.022 | Grav | Comb | |
| 132b | Tool Steel (AISI M2) | 0.864 | 0.341 | 0.012 | 0.010 | 0.009 | 0.311 |
| 134a | Mo-W-Cr-V | 0.808 | 0.218 | 0.018 | 0.007 | 0.004 | 0.079 |
| 153a | Co-Mo-W-Cr-V | 0.902 | 0.192 | 0.023 | 0.007 | 0.007 | 0.185 |
| | | | | | 0.007 | 0.270 | 0.088 |
| SRM | Ni | Cr | V | Mo | W | Co | Sn |
| 50c | 0.069 | 4.13 | 1.16 | 0.082 | 18.44 | | 0.018 |
| 132b | 0.230 | 4.38 | 1.83 | 4.90 | 6.28 | 0.029 | 0.022 |
| 134a | 0.088 | 3.67 | 1.25 | 8.35 | 2.00 | | 0.012 |
| 153a | 0.168 | 3.72 | 2.06 | 8.85 | 1.76 | 8.47 | 0.024 |

Steels (Solid Form)

These SRM's are furnished in various forms. The 600 series is for microchemical methods of analysis such as electron probe microanalysis, spark source mass spectrometric analysis, and laser probe analysis. The 1100, 1200, and 1700 series are for optical emission and x-ray spectroscopic methods of analysis. These materials have been prepared to ensure high homogeneity.

NOTE: Values in parentheses are not certified, but are given for additional information on the chemical composition.

Nominal Sizes for Solid Steel SRM's:

600 Series: 3.2 mm (1/8 in) diameter, 51 mm (2 in) long.

1100, 1200, and 1700 Series: 31 mm (1 1/4 in) diameter, 19 mm (3/4 in) thick.

C indicates a chill cast sample: 31 mm (1 1/4 in) diameter, 19 mm (3/4 in) thick.

Low-Alloy Steels

| SRM | Type | (Other Forms) | Chemical Composition (Nominal Weight Percent) | | | | |
|-------|---|---------------|---|-------------------|--------------------|---------|--------------------|
| | | | C | Mn | P | S | Si |
| 1134 | High-Silicon | 125b | 0.026 | 0.277 | 0.028 | 0.009 | 2.89 |
| 1135 | High-Silicon | 179 | 0.027 | 0.094 | 0.006 | 0.026 | 3.19 |
| 1217 | Nickel | 33e | 0.186 | 0.525 | 0.005 | 0.009 | 0.262 |
| 1218 | Low Carbon and Sulfur Silicon | 131d | 0.0029 | 0.014 | (0.002) | 0.0011 | (3.2) |
| C1221 | Resulfurized/Rephosphorized, (AISI 1211) | | 0.020 | 0.102 | 0.090 | 0.112 | 0.876 |
| 1222 | Cr-Ni-Mo (AISI 8640) | 139b | 0.43 | 0.78 | 0.013 | 0.022 | 0.24 |
| 1224 | Carbon (AISI 4130) | | 0.75 | 0.41 | 0.009 | 0.039 | 0.173 |
| 1225 | Low Alloy (AISI 4130) | | 0.274 | 0.48 | 0.007 | 0.014 | 0.221 |
| 1226 | Low Alloy | | 0.085 | 0.274 | 0.0022 | 0.0044 | 0.231 |
| 1227 | Basic Open Hearth, 1% C | | 0.97 | 0.402 | 0.014 | 0.026 | 0.215 |
| 1228 | Basic Open Hearth, 0.1% C | | 0.072 | 0.365 | 0.004 | 0.018 | 0.007 |
| 1254 | Ca in Low Alloy (Si) | | (0.03) | (0.28) | (0.03) | (0.008) | (2.9) |
| *661 | AISI 4340 | | 0.39 ₂ | 0.66 | 0.015 | 0.015 | 0.223 |
| *662 | AISI 94B17 (Modified) | | 0.163 | 1.05 | 0.044 | 0.037 | 0.40 |
| *663 | Cr-V (Modified) | | 0.57 | 1.50 | 0.02 ₉ | 0.0055 | 0.74 |
| *664 | High Carbon (Modified) | | 0.87 ₁ | 0.25 ₈ | 0.010 | 0.025 | 0.066 |
| *665 | Electrolytic Iron | | 0.008 | 0.0057 | 0.002 ₅ | 0.0059 | 0.008 ₀ |
| 1269 | Line Pipe Steel (AISI 1526, Mod) | | 0.298 | 1.35 | 0.012 | 0.0061 | 0.189 |
| 1270 | Cr-Mo Low Alloy, A336 (F-22) | | 0.077 | 0.626 | 0.0065 | 0.0065 | 0.247 |
| C1285 | Low Alloy (A242 Mod) | | 0.058 | 0.332 | 0.072 | 0.020 | 0.36 |
| 1286 | Low Alloy (Hy 80) | | 0.196 | 0.152 | 0.008 | 0.017 | 0.130 |
| 1761 | Low Alloy Steel | 2161 | 1.03 | 0.678 | 0.042 | 0.035 | 0.18 |
| 1762 | Low Alloy Steel | 2162 | 0.337 | 2.00 | 0.034 | 0.030 | 0.35 |
| 1763 | Low Alloy Steel | 2163 | 0.203 | 1.58 | 0.012 | 0.023 | 0.63 |
| 1764 | Low Alloy Steel | 2164 | 0.592 | 1.21 | 0.021 | 0.012 | 0.057 |
| 1765 | Low Alloy Steel | 2165 | 0.006 | 0.144 | 0.006 | 0.0038 | 0.004 |
| 1766 | Low Alloy Steel | 2166 | 0.015 | 0.067 | 0.002 | 0.0024 | 0.010 |
| 1767 | Low Alloy Steel | 2167 | 0.052 | 0.021 | 0.004 | 0.0090 | 0.026 |

Values in parentheses are not certified, but are given for information only.

*SRM's 661, 662, 663, 664, and 665 are sold in a set only as SRM 668.

Low-Alloy Steels (Continued)

| SRM | Cu | Ni | Cr | V | Mo | W | Co | Ti |
|-------|-------------------|-------------------|--------------------|-------------------|---------|-------------------|--------------------|-------------|
| 1134 | 0.070 | 0.038 | 0.019 | | 0.008 | | | |
| 1135 | 0.056 | 0.050 | 0.022 | <0.01 | 0.014 | | | |
| 1217 | 0.070 | 3.36 | 0.068 | (0.001) | 0.224 | | (0.06) | (0.001) |
| 1218 | 0.003 | (0.002) | 0.006 | (<0.001) | (0.003) | | (0.002) | (0.004) |
| C1221 | 0.041 | 0.067 | 0.049 | (0.0007) | 0.038 | | (0.010) | (0.0014) |
| 1222 | 0.097 | 0.51 | 0.48 | 0.005 | 0.18 | | (0.016) | (0.002) |
| 1224 | 0.072 | 0.054 | 0.071 | 0.002 | 0.013 | | | |
| 1225 | | 0.018 | 0.91 | 0.004 | 0.166 | | | |
| 1226 | 0.125 | 5.42 | 0.467 | 0.0018 | 0.446 | (0.005) | 0.029 | 0.0021 |
| 1227 | 0.006 | 0.007 | 0.019 | 0.002 | 0.003 | | 0.003 | (0.0008) |
| 1228 | 0.012 | 0.018 | 0.016 | <0.001 | 0.009 | | | |
| 1254 | (0.07) | (0.04) | (0.02) | | (0.008) | | | |
| *661 | 0.042 | 1.99 | 0.69 | 0.011 | 0.19 | 0.01, | 0.032 | 0.020 |
| *662 | 0.51 | 0.60 | 0.30 | 0.04, | 0.07, | 0.21 | 0.30 | 0.084 |
| *663 | 0.09 ₈ | 0.32 | 1.31 | 0.31 | 0.030 | 0.04 ₆ | 0.048 | 0.050 |
| *664 | 0.25 ₀ | 0.14 ₂ | 0.06 ₆ | 0.10 ₆ | 0.49 | 0.10 ₂ | 0.15 | 0.23 |
| *665 | 0.0058 | 0.041 | 0.007 ₂ | 0.0006 | 0.0050 | <1 | 0.007 ₀ | 0.0006 |
| 1269 | 0.095 | 0.108 | 0.201 | 0.004 | 0.036 | (0.001) | (0.014) | (0.009) |
| 1270 | 0.114 | 0.174 | 2.34 | 0.013 | 0.956 | (0.003) | 0.038 | (0.003) |
| C1285 | 0.37 | 1.17 | 0.80 | 0.150 | 0.164 | (0.03) | 0.036 | Ce (0.0021) |
| 1286 | 0.043 | 2.81 | 1.53 | 0.0057 | 0.334 | (0.13) | 0.116 | 0.040 |
| 1761 | 0.30 | 1.99 | 0.220 | 0.053 | 0.103 | (0.02) | (0.028) | 0.18 |
| 1762 | 0.120 | 1.15 | 0.92 | 0.200 | 0.35 | (0.01) | 0.062 | 0.095 |
| 1763 | 0.43 | 0.51 | 0.50 | 0.30 | 0.50 | (0.03) | 0.095 | 0.31 |
| 1764 | 0.51 | 0.202 | 1.48 | 0.106 | 0.200 | (<0.01) | (0.01) | 0.028 |
| 1765 | (0.0015) | 0.154 | 0.051 | 0.0040 | 0.005 | (0.001) | (0.001) | 0.0055 |
| 1766 | 0.015 | 0.021 | 0.024 | 0.009 | 0.0035 | (0.001) | 0.0020 | 0.0005 |
| 1767 | 0.0014 | 0.002 | 0.002 | 0.033 | 0.021 | (0.003) | 0.005 | 0.011 |

Values in parentheses are not certified, but are given for information only.

*SRM's 661, 662, 663, 664, and 665 are sold in a set only as SRM 668.



Helen Tyler, shown here applying heat-shrink packaging, ensures the integrity of packaging activities through her dependability and alertness to the details of packaging many different materials.

Low-Alloy Steels (Continued)

| SRM | As | Sn | Al (total) | B | Pb | Ag | Ge | |
|-------|-------------------|--------------------|-------------------|---------------------|----------------------|------------|-------------|----------|
| 1134 | | 0.003 | 0.329 | | | | | |
| 1135 | | 0.004 | 0.0028 | | | | | |
| 1217 | | | 0.030 | | | | | |
| 1218 | | | 0.005 | | | | Zr (0.002) | |
| C1221 | | | 0.111 | | | | | |
| 1222 | | | 0.038 | | | | | |
| 1224 | | | 0.060 | | | | | |
| 1226 | | (0.003) | 0.054 | | (0.0001) | Nb (0.005) | Zr (0.010) | |
| 1227 | | | (0.028) | | | | Zr (0.0006) | |
| 1228 | | | 0.061 | | | | | |
| 1254 | | (0.003) | (0.33) | | | | Ca 0.0053 | |
| *661 | 0.017 | 0.010 ₁ | 0.02 ₁ | 0.0005 | 0.00002 ₅ | 0.0004 | [0.006] | |
| *662 | 0.09 ₂ | 0.016 | 0.09 ₅ | 0.0025+ | 0.0004 ₃ | (0.0010) | [0.002] | |
| *663 | 0.010 | (0.095) | 0.24 | 0.0009 ₁ | 0.0022 | (0.0038) | [0.010] | |
| *664 | 0.05 ₂ | (0.005) | (0.008) | 0.011 | 0.024 | (0.00002) | [0.003] | |
| *665 | [0.0002] | <5 | (0.0007) | 0.00013 | 0.00001 ₅ | <0.2 | <50 | |
| 1269 | (0.006) | (0.039) | 0.016 | (<0.0001) | 0.005 | (0.0002) | | |
| 1270 | (0.02) | (0.02) | (0.005) | (0.0033) | (0.0016) | (0.0001) | | |
| C1285 | (0.022) | 0.35 | (0.12) | | | Sb (0.04) | Zr (0.02) | |
| 1286 | 0.019 | 0.012 | 0.109 | (0.006) | (0.0002) | Nb (0.012) | Zr (0.021) | |
| 1761 | 0.011 | (0.05) | 0.06 | (0.002) | | | | |
| 1762 | 0.018 | 0.046 | 0.069 | (0.004) | | | | |
| 1763 | 0.055 | 0.011 | 0.043 | (0.005) | | | | |
| 1764 | 0.010 | (0.02) | 0.009 | (0.001) | | | | |
| 1765 | 0.0010 | 0.002 | (0.006) | (0.0009) | 0.0003 | 0.0002 | | |
| 1766 | 0.0035 | 0.0010 | 0.012 | (0.0004) | 0.003 | 0.0005 | | |
| 1767 | 0.0005 | 0.006 | 0.004 | (0.0010) | (<0.0005) | | | |
| SRM | O | N | H | Nb | Se | Ta | Sr | Zr |
| *661 | (0.0009) | (0.0037) | (<0.0005) | 0.022 | 0.004 | 0.020 | | 0.009 |
| *662 | (0.0011) | (0.0041) | (<0.0005) | 0.30 | [0.001] | 0.21 | | 0.20 |
| *663 | (0.0007) | (0.0041) | (<0.0005) | 0.049 | [0.0001] | (0.053) | | 0.050 |
| *664 | [0.0017] | [0.003] | (<0.0005) | 0.15 ₇ | [0.0003] | 0.11 | | 0.69 |
| *665 | <70 | <20 | <5 | <0.5 | <0.1 | <0.5 | | <0.1 |
| 1761 | | | | (0.02) | | (0.05) | | (0.01) |
| 1762 | | | | (0.07) | | (0.02) | | (0.03) |
| 1763 | | | | (0.10) | | (0.01) | | (0.04) |
| 1764 | | | | (0.04) | | (0.03) | | (<0.001) |
| 1765 | | | | (<0.002) | | (0.005) | | (0.001) |
| 1766 | | | | (0.003) | | (0.01) | | (0.001) |
| 1767 | | | | (0.007) | | (<0.005) | | (0.004) |

*SRM's 661, 662, 663, 664, and 665 are sold in a set only as SRM 668.

Values in parentheses are not certified, but are given for information only.

Values in brackets are approximate values from heat analysis and are given for information only.

Low-Alloy Steels (Continued)

| SRM | Sb | Bi | Ca | Mg | Te | Zn | |
|------|--------------------|--------------------|-----------|----------|-----------|-----------|--------|
| *661 | 0.004 ₂ | | (<0.0001) | (0.0001) | 0.0006 | (0.0001) | |
| *662 | 0.012 | (0.002) | (0.0002) | (0.0006) | (0.0005) | (0.0005) | |
| *663 | (0.0007) | (0.0008) | (<0.0001) | (0.0005) | (0.0022) | (0.0004) | |
| *664 | (0.035) | (0.0009) | (<0.0001) | (0.0001) | (0.0002) | [0.001] | |
| *665 | <0.5 | <0.5 | <0.1 | <0.2 | <0.1 | <3 | |
| SRM | Au | Ce | Hf | La | Nd | Pr | Fe |
| *661 | (<0.0005) | 0.001 ₃ | [0.0002] | 0.0004 | 0.0003 | (0.00014) | (95.6) |
| *662 | (<0.00005) | (0.0011) | [0.006] | 0.0004 | (0.0005) | (0.00012) | (95.3) |
| *663 | 0.0005 | (0.0016) | [0.0015] | 0.0006 | (0.0007) | (0.00018) | (94.4) |
| *664 | 0.0001 | | [0.005] | 0.00007 | (0.00012) | (0.00003) | (96.7) |
| *665 | <0.02 | <0.05 | <0.2 | <0.05 | <0.05 | <0.05 | <99.9 |

*SRM's 661, 662, 663, 664, and 665 are sold in a set only as SRM 668.

Values in parentheses are not certified, but are given for information only.

Values in brackets are approximate values from heat analysis and are given for information only.



John Norris applies the knowledge and skill of many years experience to the inspection of steel samples that have been analyzed by spark discharge.

Stainless Steels

| SRM | Type | Other Forms | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|--------|--|-------------|---|------------------|---------|-------------------|-------------------|-------------------|-------------------|-------------------|------------|
| | | | C | Mn | P | S | Si | Cu | Ni | Cr | |
| C1151 | Cr-Ni | | 0.039 | 2.50 | 0.017 | 0.038 | 0.38 | 0.418 | 7.29 | 22.70 | |
| C1152 | Cr-Ni | | 0.148 | 0.96 | 0.021 | 0.0064 | 0.80 | 0.102 | 10.88 | 17.81 | |
| C1153 | Cr-Ni | | 0.264 | 0.50 | 0.030 | 0.018 | 1.07 | 0.23 | 8.77 | 16.69 | |
| C1153a | Cr-Ni | | 0.225 | 0.544 | 0.030 | 0.019 | 1.00 | 0.226 | 8.76 | 16.70 | |
| C1154 | Cr-Ni | | 0.086 | 1.42 | 0.06 | 0.053 | 0.50 | 0.40 | 12.92 | 19.06 | |
| 1155 | Cr-Ni-Mo (AISI 316) | 160b | 0.046 | 1.63 | 0.020 | 0.018 | 0.50 ₂ | 0.169 | 12.1 _s | 18.4 _s | |
| 1171 | Cr-Ni-Ti | 121d | 0.067 | 1.8 _s | 0.018 | 0.01 ₃ | 0.54 | 0.121 | 11.2 | 17.4 | |
| 1172 | Cr-Ni-Nb | 123c | 0.056 | 1.7 _s | 0.025 | 0.01 ₄ | 0.59 | 0.10 ₅ | 11.3 _s | 17.4 _s | |
| 1219 | Cr-Ni (AISI 431) | 343a | 0.149 | 0.42 | 0.026 | 0.001 | 0.545 | 0.162 | 2.16 | 15.64 | |
| 1223 | Chromium Steel | 133b | 0.127 | 1.08 | 0.018 | 0.329 | 0.327 | 0.081 | 0.232 | 12.64 | |
| 1267 | AISI 446 | 367 | 0.093 | 0.315 | 0.018 | 0.015 | 0.58 | | 0.29 | 24.14 | |
| C1287 | High-Alloy, ACI HK (AISI 310, Mod.) | | 0.36 | 1.66 | 0.029 | 0.024 | 1.66 | 0.58 | 21.16 | 23.98 | |
| C1288 | High-Alloy, ACI CN-7M (A-743) | | 0.056 | 0.83 | 0.023 | 0.010 | 0.41 | 3.72 | 29.3 | 19.55 | |
| C1289 | High-Alloy, ACI CA-6NM (AISI 414 Mod.) | | 0.014 | 0.35 | 0.017 | 0.021 | 0.156 | 0.205 | 4.13 | 12.12 | |
| SRM | V | Mo | Co | Ti | N | Al | Nb | Ta | W | Pb | Zr |
| C1151 | 0.037 | 0.80 | 0.032 | (0.006) | (0.23) | (0.004) | (0.014) | (0.006) | | 0.0039 | (0.005) |
| C1152 | 0.030 | 0.43 | 0.22 | (0.011) | (0.055) | (0.004) | (0.16) | (0.001) | | 0.0047 | (0.004) |
| C1153 | 0.18 | 0.24 | 0.127 | (0.014) | (0.134) | (0.003) | (0.050) | (0.032) | | 0.0054 | (0.003) |
| C1153a | 0.176 | 0.24 | 0.127 | (0.013) | (0.11) | (0.004) | (0.48) | (0.03) | | 0.006 | (0.0001) |
| C1154 | 0.135 | 0.07 | 0.38 | (0.004) | (0.084) | (0.004) | (0.23) | (0.075) | | 0.0178 | (0.004) |
| 1155 | 0.047 | 2.38 | 0.10 ₁ | | | | | | | 0.001 | |
| 1171 | 0.16 _s | 0.10 | 0.34 | | | | | | | | |
| 1172 | 0.22 | 0.12 | | | | 0.65 | <0.001 | | | | |
| 1219 | 0.056 | 0.164 | (0.04) | (<0.001) | 0.078 | (0.001) | (0.01) | Sn (0.008) | (0.02) | (<0.0001) | B (<0.001) |
| 1223 | 0.068 | 0.053 | | | (0.05) | (<0.005) | | Sn (0.004) | | | (0.0001) |
| 1267 | 0.08 | | | | 0.17 | | | | | | |
| C1287 | 0.09 | 0.46 | 0.31 | 0.050 | (0.034) | (0.06) | (0.07) | O (0.017) | | 0.008 | (0.006) |
| C1288 | 0.086 | 2.83 | 0.10 | 0.012 | (0.028) | (0.0025) | (0.22) | O (0.029) | (0.2) | 0.0041 | (0.002) |
| C1289 | 0.007 | 0.82 | 0.035 | 0.005 | (0.017) | (0.0016) | (0.10) | O (0.027) | | 0.0005 | (0.001) |

Values in parentheses are not certified, but are given for information only.

Specialty Steels

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | | | | | | |
|------|---------------------|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|----------------|
| | | C | Mn | P | S | Si | Cu | Ni | Cr | V | Mo | W |
| 1157 | Tool (AISI M2) | 0.836 | 0.34 | 0.011 | 0.004 | 0.18 | 0.088 | 0.228 | 4.36 | 1.82 | 4.86 | 6.28 |
| 1158 | High-Nickel (Ni 36) | 0.025 | 0.468 | 0.004 | 0.005 | 0.194 | 0.039 | 36.03 | 0.062 | 0.001 | 0.010 | 0.008 |
| 1233 | Valve Steel | 0.502 | 9.16 | 0.031 | 0.002 | 0.219 | 0.375 | 3.43 | 21.08 | 0.096 | 0.237 | (0.01) N 0.415 |

High-Temperature Alloys

| SRM | Type | Other Forms | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|-------|-------------------------------------|-------------|---|-------|------------|--------|------------------|---------|------------------|------------------|----------|
| | | | C | Mn | P | S | Si | Cu | | | |
| 1199* | L 605 | | (0.14) | 1.42 | (0.005) | | 0.83 | | | | |
| 1200* | S 816 | | (0.40) | 1.34 | (0.015) | | 0.86 | | | | |
| 1230 | A 286 | 348a | 0.044 | 0.64 | 0.023 | 0.0007 | 0.43 | 0.14 | | | |
| 1244 | Inconel 600 | | 0.062 | 0.29 | 0.010 | 0.003 | 0.12 | 0.26 | | | |
| 1245 | Inconel 625 | | 0.036 | 0.18 | 0.011 | 0.001 | 0.40 | 0.37 | | | |
| 1246 | Incoloy 800 | | 0.082 | 0.91 | 0.018 | 0.001 | 0.18 | 0.49 | | | |
| 1247 | Incoloy 825 | | 0.021 | 0.38 | 0.018 | 0.002 | 0.32 | 1.75 | | | |
| 1250 | Fe-Ni-Co | 868 | 0.022 | 0.052 | <0.003 | 0.0025 | 0.097 | 0.022 | | | |
| C2400 | High-Alloy Steel, ACI (17/4 PH) | | 0.036 | 0.71 | 0.013 | 0.003 | 0.61 | 2.63 | | | |
| C2401 | High-Alloy Steel, (ACI-CD-4M-Cu) | | 0.062 | 1.03 | 0.025 | 0.027 | 0.74 | 3.17 | | | |
| C2402 | Hasteloy C | | 0.010 | 0.64 | 0.007 | 0.018 | 0.85 | 0.19 | | | |
| SRM | Ni | Cr | Mo | Co | Ti | Al | Nb | Ta | Fe | W | B |
| 1199 | 10.2 | 19.9 | (<0.02) | 51.6 | (<0.01) | | (<0.02) | | 0.6 ₅ | 15.4 | |
| 1200 | 20.0 | 19.9 | 4.0 ₀ | 42.0 | (0.03) | | 3.1 ₈ | 1.08 | 3.19 | 3.8 ₆ | |
| 1230 | 24.2 | 14.8 | 1.18 | 0.15 | 2.12 | 0.24 | (0.07) | V 0.23 | (55) | (0.07) | 0.0055 |
| 1244 | 73.2 | 15.7 | 0.20 | 0.058 | 0.25 | 0.26 | (0.14) | | 9.6 | | <0.05 |
| 1245 | 59.5 | 21.9 | 8.6 | 0.074 | 0.28 | 0.26 | 3.5 | <0.01 | 4.5 | | <0.001 |
| 1246 | 30.8 | 20.1 | 0.36 | 0.076 | 0.32 | 0.30 | (0.09) | | 46.2 | | <0.001 |
| 1247 | 43.5 | 23.4 | 2.73 | 0.089 | 0.75 | 0.060 | (0.46) | | 26.5 | | 0.002 |
| 1250 | 37.78 | 0.077 | 0.014 | 16.1 | 1.48 | 0.99 | 2.99 | 0.003 | 40.5 | V 0.077 | 0.0078 |
| C2400 | 4.07 | 17.06 | 0.23 | 0.10 | | | 0.15 | V 0.092 | | (0.1) | (0.0004) |
| C2401 | 5.46 | 25.1 | 2.13 | 0.19 | | | (0.002) | V 0.20 | | (0.18) | (0.0004) |
| C2402 | 51.5 | 16.15 | 17.1 | 1.50 | Sn (0.001) | | (<0.01) | V 0.22 | 7.3 | 4.29 | (0.0004) |

Values in parentheses are not certified, but are given for information only.

*SRM's 1199 and 1200 sold only in a set as S1199.

Bill Reed brings the qualities of friendliness, resourcefulness and thoroughness to his new responsibility as the Acting Chief of the Office of Standard Reference Materials.




Steelmaking Alloys

These SRM's are for checking chemical methods of analysis for major constituents and for selected minor elements. They are furnished as fine powders (usually <0.1 mm).

| SRM | Type | Wt/ Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|-----|---|------------------------|---|---------|-------|----------|-------------|-----------|----------|--------|----------|
| | | | C | Mn | P | S | Si | Cu | Ni | O | |
| 57a | Silicon Metal | 60 | 0.024 | 0.015 | 0.003 | 0.003 | 98.55 | 0.004 | 0.008 | (~0.3) | |
| 58a | Ferrosilicon (73% Si- Regular Grade) | 75 | 0.014 | 0.16 | 0.009 | <0.002 | 73.20 | 0.024 | 0.012 | (0.20) | |
| 59a | Ferrosilicon | 50 | 0.046 | 0.75 | 0.016 | 0.002 | 48.10 | 0.052 | 0.033 | | |
| 195 | Ferrosilicon (75% Si- High-Purity Grade) | 75 | 0.034 | 0.17 | 0.017 | 0.001 | 75.3 | 0.047 | 0.032 | (0.42) | |
| 64c | Ferrochromium, High-Carbon | 100 | 4.68 | 0.16 | 0.020 | 0.067 | 1.22 | 0.005 | 0.43 | | |
| 196 | Ferrochromium, Low- Carbon | 100 | 0.035 | (0.282) | 0.020 | 0.003 | 0.373 | | | | |
| 71 | Calcium Molybdate | 60 | | | | | | | | | |
| 90 | Ferrophosphorus | 75 | | | 26.2 | | | | | | |
| 340 | Ferroniobium | 100 | 0.061 | 1.70 | 0.036 | | 4.39 | | | | |
| 68c | Ferromanganese, High-Carbon | 100 | 6.72 | 80.04 | 0.19 | 0.008 | 0.225 | | Sn 0.063 | | |
| 689 | Ferrochromium Silicon | 100 | 0.043 | 0.32 | 0.026 | 0.002 | 39.5 | 0.013 | 0.20 | (0.06) | |
| SRM | Cr | V | Mo | Ti | Al | Nb | Zr | Ca | Fe | B | As |
| 57a | 0.024 | 0.013 | Pb <0.001 | 0.040 | 0.47 | | 0.002 | 0.17 | 0.50 | 0.001 | <0.001 |
| 58a | 0.020 | (0.002) | (0.01) | 0.051 | 0.95 | Co <0.01 | 0.002 | 0.30 | 25.23 | 0.0010 | (0.0020) |
| 59a | 0.080 | | | 0.35 | | | | 0.042 | 50.05 | 0.058 | |
| 195 | <0.01 | (0.001) | (0.01) | 0.037 | 0.046 | Co <0.01 | 0.011 | 0.053 | 23.6 | 0.0010 | (0.0024) |
| 64c | 68.00 | 0.15 | | 0.02 | | Co 0.051 | | N 0.045 | 24.98 | | |
| 196 | 70.83 | (0.12) | | | | | | | | | |
| 71 | | | 35.3 | 0.06 | | | | | 1.92 | | |
| 90 | | | | | | | | | | | |
| 340 | | | | 0.89 | | 57.51 | Ta 3.73 | | | | |
| 68c | 0.074 | | | | | | | | 12.3 | 0.021 | |
| 689 | 36.4 | 0.09 | Pb (0.004) | 0.40 | 0.049 | Co 0.034 | Bi (<0.003) | N (0.002) | 23.2 | 0.0017 | (0.009) |

Values in parentheses are not certified, but are given for information only.



This group of steel SRM's illustrates the variety of forms of materials required to meet the SRM needs of the technical community.

Cast Irons (Chip Form)

These SRM's are furnished in 150-g units (unless otherwise noted) for use in checking chemical methods of analysis.

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|------------------------|---|------------------|-------------------|-------------------|---------|--------|
| | | C | Mn | P | S | Si | Cu |
| | Total | Graphitic | | Grav | Comb | | |
| 3d | White (110g) | 2.54 | 0.40 | 0.02 ₅ | 0.05 ₂ | 1.31 | 0.043 |
| 4k | Cast | 3.2 ₂ | 2.6 ₆ | 0.82 ₅ | 0.149 | 0.043 | 1.33 |
| 5L | Cast | 2.60 | 1.98 | 0.68 | 0.284 | 0.124 | 1.82 |
| 6g | Cast | 2.85 | 2.01 | 1.05 | 0.557 | 0.124 | 1.05 |
| 7g | Cast (High Phosphorus) | 2.69 | 2.59 | 0.612 | 0.794 | 0.061 | 0.060 |
| 82b | Cast (Ni-Cr) | 2.85 | 2.37 | 0.745 | 0.025 | 0.007 | 2.10 |
| 107c | Cast (Ni-Cr-Mo) | 2.99 | 1.98 | 0.480 | 0.079 | 0.059 | 1.21 |
| 115a | Cast (Cu-Ni-Cr) | 2.62 | 1.96 | 1.00 | 0.086 | 0.064 | 2.13 |
| 122h | Cast (Car Wheel) | 3.52 | 2.82 | 0.543 | 0.311 | 0.072 | 0.513 |
| 334 | Gray Cast | 2.83 | | | | 0.043 | 0.028 |
| 338 | White Cast | 3.33 | | (0.76) | (0.054) | 0.015 | (1.82) |
| 341 | Ductile | 1.81 | 1.23 | 0.92 | 0.024 | 0.007 | 2.44 |
| 342a | Nodular | 1.86 | 1.38 | 0.274 | 0.019 | 0.006 | 2.73 |
| 890 | HC 250+V | 2.91 | | 0.62 | 0.025 | 0.015 | 0.67 |
| 891 | Ni-Hard, Type I | 2.71 | | 0.55 | 0.038 | 0.029 | 0.56 |
| 892 | Ni-Hard, Type IV | 3.33 | | 0.76 | 0.054 | 0.015 | 1.83 |
| SRM | Ni | Cr | V | Mo | Co | Ti | |
| 3d | 0.025 | 0.03 | (0.002) | (0.007) | | (0.003) | |
| 4k | 0.042 | 0.116 | 0.024 | 0.040 | Zn (<0.001) | (0.03) | |
| 5L | 0.086 | 0.148 | 0.034 | 0.020 | | 0.050 | |
| 6g | 0.135 | 0.370 | 0.056 | 0.035 | | 0.059 | |
| 7g | 0.120 | 0.048 | 0.010 | 0.012 | | 0.044 | |
| 82b | 1.22 | 0.333 | 0.027 | 0.002 | | 0.027 | |
| 107c | 2.20 | 0.693 | 0.015 | 0.83 | | 0.019 | |
| 115a | 14.49 | 1.98 | 0.014 | 0.050 | | 0.020 | |
| 122h | 0.078 | 0.052 | 0.041 | (0.003) | | 0.034 | |
| 338 | (5.5) | (10.2) | (0.04) | | (0.32) | | |
| 341 | 20.32 | 1.98 | 0.012 | 0.010 | | 0.018 | |
| 342a | 0.058 | 0.034 | | 0.006 | | 0.020 | |
| 890 | 0.397 | 32.4 | 0.45 | 0.018 | (0.03) | | |
| 891 | 4.48 | 2.23 | 0.039 | 0.27 | 0.19 | (0.01) | |
| 892 | 5.53 | 10.18 | 0.041 | 0.20 | 0.31 | (0.02) | |

Values in parentheses are not certified, but are given for information only.

Cast Irons (Chip Form) (Continued)

| SRM | As | Sn | Al (total) | Mg | N | Fe |
|------|---------|---------|------------|------------|----------|-----------|
| 4k | (0.03) | (0.004) | (0.004) | Sb(<0.001) | (0.0016) | Pb(0.001) |
| 5L | | | | | 0.005 | |
| 6g | 0.042 | | | | 0.005 | |
| 7g | 0.014 | | | | 0.004 | |
| 341 | | | | 0.068 | | |
| 342a | | | | 0.070 | | |
| 890 | (0.008) | | (<0.01) | | (0.089) | (61.8) |
| 891 | (0.004) | (<0.01) | (0.008) | | (0.012) | (88.5) |
| 892 | (0.006) | (0.02) | (0.009) | | (0.019) | (77.4) |

Values in parentheses are not certified, but are for information only.

Cast Steels, White Cast Irons, Ductile Irons, and Blast Furnace Irons (Solid Form)

These SRM's are for analysis of cast steels and cast irons by rapid instrumental methods.

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|--------|-------------------------------|---|-------|-------|----------|-------------------|-------------------|-------------------|------------------|--|
| | | C | Mn | P | S | Si | Cu | Ni | Cr | |
| C1137a | White Cast Iron | 2.86 | 0.52 | 0.087 | 0.017 | 1.15 | 0.192 | 2.17 | 0.643 | |
| 1138a | Cast Steel (No. 1) | 0.11 ₈ | 0.35 | 0.035 | 0.056 | 0.25 | 0.09 | 0.10 | 0.13 | |
| 1139a | Cast Steel (No. 2) | 0.79 ₀ | 0.92 | 0.012 | 0.013 | 0.80 | 0.47 | 0.98 | 2.1 ₈ | |
| 1144a | Blast Furnace Iron (2) | 4.32 | 1.23 | 0.084 | 0.083 | 0.18 ₂ | 0.09 ₁ | 0.06 ₃ | 0.029 | |
| C1145a | White Cast Iron | 2.92 | 0.187 | 0.215 | 0.191 | 0.271 | 0.46 | 0.62 | 0.63 | |
| C1146a | White Cast Iron | 1.97 | 1.60 | 0.55 | 0.016 | 3.93 | 1.48 | 3.07 | 2.56 | |
| C1150a | White Cast Iron | 3.32 | 0.77 | 0.078 | 0.065 | 1.35 | 0.112 | 0.097 | 0.155 | |
| C1173 | Cast Steel 3 | 0.453 | 0.174 | 0.031 | 0.092 | 1.38 | 0.204 | 4.04 | 2.63 | |
| 1173 | Ni-Cr-Mo-V Steel | 0.423 | 0.19 | 0.033 | 0.092 | 1.28 | 0.204 | 4.06 | 2.70 | |
| C1290 | High Alloy (HC-250+V) | 3.04 | 0.66 | 0.030 | 0.013 | 0.971 | 0.065 | 0.917 | 30.5 | |
| C1291 | High Alloy (Ni-Hard, Type I) | 2.67 | 1.14 | 0.028 | 0.032 | 1.34 | 0.26 | 4.34 | 2.78 | |
| C1292 | High Alloy (Ni-Hard, Type IV) | 3.47 | 0.55 | 0.049 | 0.016 | 0.59 | 0.36 | 5.04 | 11.4 | |
| C2423 | Ductile Iron | 3.76 | 0.98 | 0.27 | (0.0006) | 1.67 | 1.55 | 0.146 | 0.322 | |
| C2423a | Ductile Iron | 3.66 | 0.91 | 0.246 | (<0.001) | 1.59 | 1.61 | 0.147 | 0.322 | |
| C2424 | Ductile Iron | 2.68 | 0.268 | 0.041 | 0.024 | 3.37 | 0.125 | 0.061 | 0.13 | |
| C2424a | Ductile Iron | 2.76 | 0.207 | 0.034 | 0.016 | 3.30 | 0.099 | 0.045 | 0.15 | |
| C2425 | Ductile Iron | 3.26 | 0.76 | 0.191 | 0.012 | 2.50 | 0.47 | 0.55 | 0.092 | |
| C2425a | Ductile Iron | 3.30 | 0.72 | 0.188 | 0.010 | 2.38 | 0.47 | 0.57 | 0.085 | |

**Cast Steels, White Cast Irons, Ductile Irons,
and Blast Furnace Irons (Solid Form) (Continued)**

| SRM | V | Mo | Ti | As | Al | Te | Co |
|--------|-------------------|---------|----------|----------|----------|-------------------|---------|
| C1137a | 0.019 | 0.86 | (0.04) | | (0.007) | Mg 0.032 | |
| 1138a | 0.02 _o | 0.05 | (0.0012) | (<0.005) | (0.067) | Fe (98.7) | |
| 1139a | 0.26 | 0.51 | (0.004) | (<0.005) | (0.13) | Fe (93.0) | |
| 1144a | 0.02 _s | (0.007) | 0.32 | (0.004) | (<0.005) | 0.02 _z | |
| C1145a | 0.112 | 0.48 | 0.012 | (0.02) | (0.04) | | 0.058 |
| C1146a | 0.20 | 1.52 | 0.20 | (0.16) | (0.028) | Pb 0.0018 | 0.13 |
| C1150a | 0.040 | 0.086 | 0.040 | (0.017) | (0.005) | Pb 0.001 | 0.014 |
| C1173 | 0.42 | 1.46 | 0.037 | (0.02) | (0.005) | Pb (0.0006) | 0.064 |
| 1173 | 0.42 | 1.50 | (0.015) | | | Nb (0.045) | 0.076 |
| C1290 | 0.442 | (0.041) | | | | | |
| C1291 | 0.031 | 0.32 | | | | | |
| C1292 | 0.041 | 0.25 | | | | | |
| C2423 | 0.048 | 0.155 | 0.10 | | (0.09) | | (0.02) |
| C2423a | 0.043 | 0.159 | 0.10 | | (0.08) | | (0.02) |
| C2424 | 0.083 | 0.019 | 0.050 | | (<0.01) | | (0.05) |
| C2424a | 0.081 | 0.019 | 0.045 | | (<0.01) | | (0.05) |
| C2425 | 0.013 | 0.30 | 0.19 | | (0.02) | | (0.02) |
| C2425a | 0.013 | 0.29 | 0.20 | | (0.02) | | (0.03) |
| SRM | Mg | | Ce | | La | | B |
| C2423 | | 0.058 | | 0.036 | | 0.011 | (0.01) |
| C2423a | | 0.076 | | 0.031 | | 0.0042 | (0.01) |
| C2424 | | 0.006 | | 0.0046 | | 0.0011 | (0.002) |
| C2424a | | 0.014 | | 0.0053 | | 0.0010 | (0.001) |
| C2425 | | 0.040 | | 0.0062 | | 0.0015 | (0.10) |
| C2425a | | 0.047 | | 0.023 | | 0.0037 | (0.1) |

Values in parentheses are not certified, but are given for information only.

Nonferrous Alloys

Aluminum-Base Alloys

| SRM | Type | Wt/ Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | |
|-------|---------------------------------------|------------------------|---|-------|--------|--------|--------|----------|---------|
| | | | Mn | Si | Cu | Ni | Cr | V | |
| 87a | Al-Si (Chip) | 75 | 0.26 | 6.24 | 0.30 | 0.57 | 0.11 | <0.01 | |
| 853 | Alloy 3004 (Chip) | 30 | 1.26 | 0.18 | 0.15 | 0.004 | <0.001 | 0.017 | |
| 1240 | Alloy 3004 | Disk | 1.26 | 0.18 | 0.15 | 0.004 | <0.001 | 0.017 | |
| 1240a | Alloy 3004 | Disk | 1.27 | 0.18 | 0.15 | 0.004 | <0.001 | 0.017 | |
| 1240b | Alloy 3004 | Disk | 1.27 | 0.18 | 0.15 | 0.004 | <0.001 | 0.017 | |
| 854 | Alloy 5182 (Chip) | 30 | 0.38 | 0.16 | 0.050 | 0.020 | 0.032 | 0.016 | |
| 1241a | Alloy 5182 | Disk | 0.38 | 0.16 | 0.050 | 0.020 | 0.032 | 0.016 | |
| 1241b | Alloy 5182 | Disk | 0.38 | 0.16 | 0.050 | 0.020 | 0.032 | 0.016 | |
| 855 | Casting Alloy 356 (fine millings) | 30 | 0.057 | 7.17 | 0.13 | 0.015 | 0.013 | | |
| 1255a | Casting Alloy 356 | Disk | 0.053 | 7.22 | 0.12 | 0.017 | 0.012 | 0.024 | |
| 856 | Casting Alloy 380 (fine millings) | 30 | 0.35 | 9.21 | 3.51 | 0.37 | 0.055 | | |
| 1256a | Casting Alloy 380 | Disk | 0.38 | 9.18 | 3.51 | 0.41 | 0.055 | 0.018 | |
| C1257 | High Purity | Disk | <1.0 | <10 | <1.0 | <1.0 | <1.0 | <0.1 | |
| 858 | Alloy 6011 (modified) (fine millings) | 35 | 0.48 | 0.79 | 0.84 | 0.0006 | 0.0011 | 0.0030 | |
| 1258 | Alloy 6011 (35mm D×19mm thick) | Disk | 0.48 | 0.78 | 0.84 | 0.0006 | 0.0011 | | |
| 859 | Alloy 7075 (fine millings) | 35 | 0.078 | 0.17 | 1.59 | 0.063 | 0.176 | 0.0082 | |
| 1259 | Alloy 7075 (35mm D×19mm thick) | Disk | 0.079 | 0.18 | 1.60 | 0.063 | 0.173 | | |
| SRM | Ti | Sn | Ga | Fe | Pb | Mg | Zn | Zr | Be |
| 87a | 0.18 | 0.05 | 0.02 | 0.61 | 0.10 | 0.37 | 0.16 | | |
| 853 | 0.018 | | 0.018 | 0.50 | | 1.11 | 0.052 | 0.002 | |
| 1240 | 0.022 | | 0.018 | 0.50 | | 1.11 | 0.052 | 0.002 | |
| 1240a | 0.022 | | 0.018 | 0.50 | | 1.12 | 0.051 | 0.002 | |
| 1240b | 0.021 | | 0.018 | 0.50 | | 1.11 | 0.051 | 0.002 | |
| 854 | 0.030 | | 0.018 | 0.20 | | 4.54 | 0.051 | 0.002 | |
| 1241a | 0.032 | | 0.018 | 0.20 | | 4.54 | 0.052 | 0.002 | |
| 1241b | 0.034 | | 0.018 | 0.20 | | 4.54 | 0.051 | 0.002 | |
| 855 | 0.15 | 0.010 | | 0.16 | 0.015 | 0.37 | 0.083 | | |
| 1255a | 0.156 | 0.013 | | 0.14 | 0.017 | 0.36 | 0.083 | Sr 0.02 | |
| 856 | 0.068 | 0.10 | | 0.92 | 0.10 | 0.061 | 0.96 | | |
| 1256a | 0.084 | 0.10 | | 0.90 | 0.10 | 0.062 | 1.02 | Sr 0.020 | |
| C1257 | (<0.1) | (<0.1) | (<0.1) | 1.0 | (<0.1) | <1.0 | (<0.1) | (<0.1) | (<0.1) |
| 858 | 0.042 | | | 0.078 | | 1.01 | 1.04 | | <0.0001 |
| 1258 | (0.04) | | (0.010) | 0.079 | | 0.98 | 1.03 | | |
| 859 | 0.041 | | | 0.202 | | 2.45 | 5.46 | | 0.0026 |
| 1259 | (0.04) | | (0.022) | 0.205 | | 2.48 | 5.44 | | 0.0025 |

Values in parentheses are not certified, but are given for information only.

Copper-Base Alloys (Chip Form)

| SRM | Type | Wt/ Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|------|---|------------------------|--|----------|---------|----------|----------|-----------|-----------|---------|-----------|
| | | | Cu | Ni | Fe | Zn | Pb | | | | |
| 158a | Bronze, Silicon | 150 | 90.93 | 0.001 | 1.23 | 2.08 | 0.097 | | | | |
| 871 | Bronze, Phosphor (CDA 521) | 100 | 91.68 | <0.001 | 0.025 | 0.010 | | | | | |
| 872 | Bronze, Phosphor (CDA 544) | 100 | 87.36 | 0.003 | 4.0 | 4.13 | | | | | |
| 874 | Cupro-Nickel, 10% (CDA 706) "High-Purity" | 100 | 88.49 | 10.18 | 1.22 | 0.002 | <0.0005 | | | | |
| 875 | Cupro-Nickel, 10% (CDA 706) "Doped" | 100 | 87.83 | 10.42 | 1.45 | 0.11 | 0.0092 | | | | |
| 879 | Nickel Silver (CDA 762) | 100 | 57.75 | 12.11 | 0.0020 | 30.04 | 0.002 | | | | |
| 880 | Nickel Silver (CDA 770) | 100 | 54.51 | 18.13 | 0.004 | 27.3 | 0.002 | | | | |
| 1034 | *Unalloyed Copper | rod | (99.96%) | (0.6) | (2.0) | (<11) | (0.5) | | | | |
| 1035 | **Leaded-Tin Bronze Alloy | 50 | (78.5) | (0.75) | (0.001) | (0.25) | (13.5) | | | | |
| SRM | Mn | Sb | Sn | Cr | P | Ag | Si | Al | Te | Cd | Se |
| 158a | 1.11 | | 0.96 | | 0.026 | | 3.03 | 0.46 | | | |
| 871 | | | 8.14 | | 0.082 | | | | | | |
| 872 | | | 4.16 | | 0.26 | | | | | | |
| 874 | 0.0020 | <0.001 | 0.007 | | 0.002 | | (0.0006) | | | <0.0002 | 0.00015 |
| 875 | <0.0007 | <0.001 | 0.009 | | 0.0020 | | (0.0008) | | (<0.0001) | 0.0022 | 0.0004 |
| 879 | <0.001 | | | | | | | | | | |
| 880 | <0.001 | | | | | | | | | | |
| 1034 | (<0.1) | (0.2) | (<0.2) | (0.3) | | | (8.1) | (<2) | (<2) | (0.5) | (<1) |
| 1035 | | | (6.8) | | | | | | | | (3.3) |
| SRM | Bi | O | Co | C | Au | H | S | As | Mg | Ti | |
| 874 | <0.0002 | (0.06) | | (0.0023) | | (0.0016) | (0.0011) | (<0.0006) | (0.0002) | | (0.0001) |
| 875 | 0.003 | (0.14) | | (0.0035) | | (0.004) | (0.0011) | (0.0010) | (0.0010) | | (<0.0002) |
| 1034 | (0.2) | (363) | (0.2) | | (<0.05) | | 2.8 ppm | (0.2) | (<1) | | |
| 1035 | | (0.64) | | | | | 22.3 ppm | | P (0.004) | | |

Values in parentheses are not certified, but are given for information only.

*Values for SRM 1034 are ppm by weight.

**Sulfur value for SRM 1035 is ppm by weight.

Copper-Base Alloys (Solid Form)

The SRM's with "C" prefix are chill-cast blocks, 31 mm square, 19 mm thick; the others are wrought disks, 31 mm in diameter and 19 mm thick. Both forms have nearly identical chemical compositions.

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | | |
|------------|------------------------|---|-------------------|--------|-------------------|-------------------|-------------------|------------|
| | | Cu | Zn | Pb | Fe | Sn | Ni | Al |
| 1103 | Free-Cutting Brass | 59.27 | 35.72 | 3.73 | 0.26 | 0.88 | 0.15 | |
| 1104 | Free-Cutting Brass | 61.33 | 35.31 | 2.77 | 0.088 | 0.43 | 0.070 | |
| 1107 C1107 | Naval Brass B | 61.2 ₁ | 37.3 ₄ | 0.18 | 0.037 | 1.04 | 0.098 | |
| 1108 C1108 | Naval Brass C | 64.9 ₅ | 34.4 ₂ | 0.063 | 0.050 | 0.39 | 0.033 | |
| | | | | | | | | |
| | | | | | | | | |
| 1110 C1110 | Red Brass B | 84.5 ₉ | 15.2 ₀ | 0.033 | 0.033 | 0.051 | 0.053 | |
| 1111 C1111 | Red Brass C | 87.1 ₄ | 12.8 ₁ | 0.013 | 0.010 | 0.019 | 0.022 | |
| 1112 C1112 | Gilding Metal A | 93.3 ₈ | 6.3 ₀ | 0.057 | 0.07 ₀ | 0.12 | 0.10 ₀ | |
| 1113 C1113 | Gilding Metal B | 95.0 ₃ | 4.8 ₀ | 0.026 | 0.04 ₃ | 0.06 ₄ | 0.057 | |
| 1114 C1114 | Gilding Metal C | 96.4 ₅ | 3.4 ₇ | 0.012 | 0.01 ₇ | 0.02 ₇ | 0.021 | |
| | | | | | | | | |
| 1115 C1115 | Commercial Bronze A | 87.9 ₆ | 11.7 ₃ | 0.013 | 0.13 | 0.10 | 0.074 | |
| 1116 C1116 | Commercial Bronze B | 90.3 ₇ | 9.4 ₄ | 0.042 | 0.046 | 0.04 ₄ | 0.048 | |
| 1117 C1117 | Commercial Bronze C | 93.0 ₁ | 6.8 ₇ | 0.069 | 0.014 | 0.02 ₁ | 0.020 | |
| | | | | | | | | |
| | | | | | | | | |
| 1119 C1119 | Aluminum Brass B | 77.1 ₂ | 20.5 ₃ | 0.051 | 0.03 ₂ | | | 2.14 0.050 |
| 1275 | Cupro-Nickel (CDA 706) | 88.2 | 0.085 | 0.006 | 1.46 | 0.008 | 9.76 | 0.0005 |
| 1276 | Cupro-Nickel (CDA 715) | 67.8 | 0.038 | 0.004 | 0.56 | 0.023 | 30.5 | 0.0004 |
| | | | | | | | | |
| SRM | As | Be | Cd | Mn | P | Si | Ag | |
| 1103 | | | | | 0.003 | | | |
| 1104 | | | | | 0.005 | | | |
| | | | | | | | | |
| | | | | | | | | |
| 1108 C1106 | | | | 0.005 | | | | |
| 1108 C1108 | | | | 0.025 | | | | |
| 1112 C1112 | | | | | 0.009 | | | |
| | | | | | | | | |
| | | | | | | | | |
| 1113 C1113 | | | | | 0.008 | | | |
| 1114 C1114 | | | | | 0.009 | | | |
| 1115 C1115 | | | | | 0.005 | | | |
| 1116 C1116 | | | | | | 0.008 | | |
| 1117 C1117 | | | | | 0.002 | | | |
| | | | | | | | | |
| | | | | | | | | |
| 1275 C1119 | (0.001) | | 0.0003 | 0.42 | 0.070 | | | |
| 1276 | | | 0.0002 | 1.01 | 0.005 | (0.001) | (0.004) | |
| | | | | | 0.006 | | | |
| SRM | Te | Co | Cr | Se | Mg | B | S | Ti |
| 1275 | (0.0002) | 0.024 | (0.0002) | 0.0004 | 0.003 | (0.0009) | (0.008) | (0.0002) |
| 1276 | | 0.045 | | 0.0005 | 0.12 | | | |

Values in parentheses are not certified, but are given for information only.

Copper "Benchmark"

| SRM | | Type | Cu(Wt%) | Chemical Composition (Nominal Parts Per Million by Weight) | | | | | | | | |
|--------|---------|--------------------------|---------|--|-------|--------|--------|-------|--------|-------|---------|------|
| (Chip) | (Solid) | | | Sb | As | Bi | Cr | Co | Fe | Pb | Mn | |
| 393 | | Unalloyed Copper "O" | 99.998 | 0.25 | 0.41 | <0.1 | <0.5 | 0.02 | <1 | 0.039 | <0.01 | |
| 394 | 494 | Unalloyed Copper I | 99.908 | 4.5 | 2.6 | 0.35 | 2.0 | 0.5 | 147 | 26.5 | 3.7 | |
| 395 | 495 | Unalloyed Copper II | 99.944 | 8.0 | 1.6 | 0.50 | 6.0 | | 96 | 3.25 | 5.3 | |
| 396 | 496 | Unalloyed Copper III | 99.955 | <1 | <0.2 | 0.07 | 4.3 | 0.4 | 143 | 0.41 | 7.5 | |
| | 457 | Unalloyed Copper IV | 99.96 | 0.2 | 0.2 | 0.2 | (0.3) | (0.2) | 2.0 | 0.5 | <0.1 | |
| 398 | 498 | Unalloyed Copper V | 99.98 | 7.5 | 25 | 2.0 | (0.3) | 2.8 | 11.4 | 9.9 | (0.3) | |
| 399 | 499 | Unalloyed Copper VI | 99.79 | 30 | 47 | 10.5 | (0.5) | 0.5 | 20.0 | 114 | (0.3) | |
| 400 | 500 | Unalloyed Copper VII | 99.70 | 102 | 140 | 24.5 | (0.5) | 0.6 | 41 | 128 | (0.2) | |
| | C1251 | Phosphorized Copper VIII | 99.96 | 14 | 14.4 | (3) | 2.8 | 8.8 | (264) | 7.5 | (5) | |
| | C1252 | Phosphorized Copper IX | 99.89 | 42 | 115 | 21 | 7.4 | 90 | (35) | 60 | (17) | |
| | C1253 | Phosphorized Copper X | 99.42 | (140) | 432 | 70 | 216 | 495 | (330) | 244 | (380) | |
| 454 | | Unalloyed Copper XI | 99.84 | 24 | 46 | 19 | | (4) | (50) | 66 | | |
| SRM | | Ni | Se | Ag | S | Te | Sn | Zn | Al | Cd | Au | Mg |
| 393 | | 0.05 | <0.05 | 0.10 | <1 | <0.5 | <0.1 | <0.1 | <0.1 | <0.1 | <0.05 | <0.1 |
| 394 | 494 | 11.7 | 2.00 | 50.5 | 15 | 0.58 | 70 | 405 | (<2) | (0.5) | (0.07) | (<1) |
| 395 | 495 | 5.4 | 0.63 | 12.2 | 13 | 0.32 | 1.5 | 12.2 | (<2) | (0.4) | (0.13) | (<1) |
| 396 | 496 | 4.2 | 0.62 | 3.30 | 9.5 | (0.02) | 0.8 | 5.0 | (<2) | (0.6) | (<0.05) | (<1) |
| | 457 | 0.6 | 4.2 | 8.1 | (4) | 0.29 | <0.2 | <11 | (<2) | (<1) | (<0.05) | (<1) |
| 398 | 498 | 7.0 | 17.5 | 20.1 | (11) | 10.1 | 4.8 | 24 | (<2) | (22) | (0.1) | (<1) |
| 399 | 499 | 506 | 95 | 117 | (10) | 50 | (~90) | 45 | (<2) | (<1) | (4) | (<1) |
| 400 | 500 | 603 | 214 | 181 | (9) | 153 | (~200) | 114 | (<2) | (<1) | (10) | (<1) |
| | C1251 | 22 | 11.4 | 85 | (31) | 15 | (15) | 8 | (2) | 2 | 15.0 | (10) |
| | C1252 | 128 | 53.6 | 166.6 | (29) | 51 | (110) | 60 | (7) | 14 | 34.9 | (20) |
| | C1253 | (500) | 164 | 495 | 55 | 199 | (470) | 350 | (180) | 74 | 74.4 | (80) |
| 454 | | (150) | 479 | 286 | | 27 | 2.2 | 7 | | | 7.5 | |
| SRM | | Si | Be | B | Ca | Li | Pd | P | Ti | Zr | | |
| 393 | | <0.5 | <0.01 | <0.01 | <0.05 | <0.01 | <0.05 | <0.05 | <0.5 | <0.5 | <0.5 | |
| 394 | 494 | (<2) | | | | | | | | | | |
| 395 | 495 | (<2) | | | | | | | | | | |
| 396 | 496 | (<2) | | | | | | | | | | |
| 398 | 498 | (<2) | | | | | | | | | | |
| 399 | 499 | (<2) | | | | | | | | | | |
| 400 | 500 | (<2) | | | | | | | | | | |
| | C1251 | (15) | (<0.5) | | (4) | (0.04) | | | (0.4%) | | | |
| | C1252 | (13) | (<5) | | (6) | (0.03) | | | | | | |
| | C1253 | (350) | (12) | | (1) | (9) | | | 518 | | | |
| 454 | | | | | | | (0.1) | | | | | |

Values in parentheses are not certified, but are given for information only.

Lead-Base Alloys

| SRM | | Type | Chemical Composition (Nominal Weight Percent) | | | | | | | |
|------|------|-------------------------|---|-------|-------|------|-------|-------|-------|--------|
| Chip | Disk | | Cu | Ni | As | Sn | Sb | Bi | Ag | Fe |
| 1129 | | Solder 63Sn-37Pb | 0.16 | 0.010 | 0.055 | 62.7 | 0.13 | 0.13 | 0.075 | |
| 127b | 1131 | Solder 60Pb-40Sn | 0.011 | 0.012 | 0.01 | 39.3 | 0.43 | 0.06 | 0.01 | |
| 53e | 1132 | Bearing Metal(Pb-Sb-Sn) | 0.054 | 0.003 | 0.057 | 5.84 | 10.26 | 0.052 | | <0.001 |

Lead-Base Material

| SRM | C2416 | | C2417 | | C2418 | | | |
|---|-----------|-------------|-----------------|------------------|-----------|--|--|--|
| | Type | Bullet Lead | Lead-Base Alloy | High-Purity Lead | | | | |
| | Size | 50 mm | 50 mm | 50 mm | | | | |
| Chemical Composition (Nominal Weight Percent) | | | | | | | | |
| Sb | 0.79 | | 0.010 | | (<0.0001) | | | |
| As | 0.056 | | 0.011 | | (<0.0001) | | | |
| Bi | 0.10 | | 0.010 | | (<0.0005) | | | |
| Cu | 0.065 | | 0.010 | | (<0.0001) | | | |
| S | 0.0015 | | (<0.0005) | | — | | | |
| Ag | 0.0044 | | 0.010 | | 0.0001 | | | |
| Sn | 0.09 | | (<0.010) | | (<0.0005) | | | |
| Al | (<0.0001) | | (<0.0001) | | (<0.0001) | | | |
| Cd | (0.0002) | | (<0.0002) | | 0.0003 | | | |
| Ca | (<0.001) | | (<0.001) | | (<0.0005) | | | |
| Co | (<0.0002) | | (<0.0002) | | (<0.0005) | | | |
| Fe | (<0.0005) | | (<0.0003) | | (<0.0005) | | | |
| Mn | (<0.0005) | | (<0.0003) | | (<0.0005) | | | |
| Ni | (<0.0005) | | (<0.0005) | | (<0.0005) | | | |
| Te | (<0.0005) | | (<0.0005) | | (<0.0005) | | | |
| Zn | (<0.0005) | | (<0.0005) | | (<0.0005) | | | |

These lead-base materials are issued in the form of disks 50 mm in diameter and 16 mm thick. They are intended for use with optical emission spectrometric methods of analysis.

Nickel-Base Alloys

| SRM | Type | Wt/Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | |
|-------|-------------------------------------|--------------------|---|-------------------|--------|--------|-----------|--------|-------|-------|
| | | | C | Mn | P | S | Si | Cu | Ni | Cr |
| 349a | Waspaloy Ni-Co-Cr | 150 | 0.035 | 0.019 | 0.003 | 0.0024 | 0.018 | 0.007 | 58.1 | 19.3 |
| 882 | Ni-Cu-A | 100 | 0.006 | 0.0007 | | 0.0014 | 0.006 | 31.02 | 65.25 | |
| 864 | Inconel, 600 | 100 | 0.064 | 0.29 | 0.010 | 0.003 | 0.12 | 0.26 | 73.1 | 15.7 |
| 865 | Inconel, 625 | 100 | 0.037 | 0.18 | 0.012 | 0.001 | 0.41 | 0.36 | 59.5 | 21.9 |
| 866 | Incoloy, 800 | 100 | 0.082 | 0.92 | 0.017 | 0.001 | 0.17 | 0.49 | 30.8 | 20.1 |
| 867 | Incoloy, 825 | 100 | 0.021 | 0.39 | 0.018 | 0.002 | 0.32 | 1.74 | 43.5 | 23.4 |
| 1159 | Electronic and Magnetic Alloy Ni-Fe | Disk | 0.007 | 0.30 _s | 0.003 | 0.003 | 0.32 | 0.038 | 48.2 | 0.06 |
| 1160 | Electronic and Magnetic Alloy Ni-Mo | Disk | 0.019 | 0.55 _s | 0.003 | 0.001 | 0.37 | 0.021 | 80.3 | 0.05 |
| 1243 | Waspaloy | Disk | 0.024 | 0.019 | 0.003 | 0.0018 | 0.018 | 0.007 | 58.78 | 19.20 |
| C1248 | Ni-Cu | Disk | 0.266 | 0.31 | 0.002 | 0.0008 | 1.61 | 29.80 | 65.75 | 0.095 |
| SRM | Mo | Co | Ti | Al | B | Fe | Nb | | | |
| 349a | 4.25 | 12.46 | 3.06 | 1.23 | 0.005 | 1.15 | V | 0.012 | | |
| 882 | | | 0.57 | 2.85 | | 0.009 | | | | |
| 864 | 0.20 | 0.059 | 0.26 | 0.26 | <0.005 | 9.6 | | (0.14) | | |
| 865 | 8.6 | 0.072 | 0.28 | 0.21 | <0.001 | 4.5 | | 3.5 | | |
| 866 | 0.36 | 0.075 | 0.31 | 0.29 | <0.001 | 46.1 | | (0.09) | | |
| 867 | 2.73 | 0.089 | 0.75 | 0.062 | 0.002 | 26.6 | | (0.45) | | |
| 1159 | 0.01 _s | 0.022 | | | | 51.0 | | | | |
| 1160 | 4.3 _s | 0.054 | | | | 14.3 | | | | |
| 1243 | 4.25 | 12.46 | 3.06 | 1.23 | 0.005 | 0.79 | V | 0.12 | | |
| C1248 | 0.006 | Pb 3.8 µg/g | Sn 1.1 µg/g | 0.009 | | 2.10 | Zn 3 µg/g | | | |

Values in parentheses are not certified, but are given for information only.

Carbon Modified Silicon

SRM 1216—This SRM is intended for the calibration of instruments used to measure total elemental carbon. The SRM consists of three, one gram bottles of chemically modified microparticulate silica certified at the carbon levels of 0.70, 9.06, and 17.04% carbon, respectively.

| SRM | Type | Bottle | Percent Carbon |
|------|-------------------------|--------|----------------|
| 1216 | Carbon Modified Silicon | I | 0.070 |
| | | II | 9.06 |
| | | III | 17.04 |

Trace Elements in Nickel-Base Superalloys (Chip Form)

| SRM | Type | Wt/Unit (grams) | Nominal Trace Composition (Parts Per Million by Weight) | | | | |
|-----|----------------|--------------------|---|-------|------|------|-------|
| | | | Pb | Bi | Se | Te | Tl |
| 897 | "Tracealloy" A | 35 | 11.7 | (0.5) | 9.1 | 1.05 | 0.51 |
| 898 | "Tracealloy" B | 35 | 2.5 | (1.0) | 2.00 | 0.54 | 2.75 |
| 899 | "Tracealloy" C | 35 | 3.9 | (0.3) | 9.5 | 5.9 | 0.252 |

| SRM | Approximate Base Composition (Weight Percent) | | | | | | | | | | | |
|-----|---|--------|-------|-------|--------|-------|-------|-------|---------|--------|--------|-------|
| | C | Cr | Co | Ni | W | Nb | Al | Ti | B | Zr | Ta | Hf |
| 897 | (0.12) | (12.0) | (8.5) | (Bal) | (1.75) | (0.9) | (2.0) | (2.0) | (0.010) | (0.10) | (1.75) | (1.2) |
| 898 | (0.12) | (12.0) | (8.5) | (Bal) | (1.75) | (0.9) | (2.0) | (2.0) | (0.010) | (0.10) | (1.75) | (1.2) |
| 899 | (0.12) | (12.0) | (8.5) | (Bal) | (1.75) | (0.9) | (2.0) | (2.0) | (0.010) | (0.10) | (1.75) | (1.2) |

Values in parentheses are not certified, but are given for information only.

Nickel Oxides (Powder Form)

| SRM | Type | Wt/Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|-----|---------|--------------------|---|-------|-------|--------|-------|-------|-------|-------|-------|
| | | | Mn | Si | Cu | Cr | Co | Ti | Al | Fe | Mg |
| 671 | Oxide 1 | 25 | 0.13 | 0.047 | 0.20 | 0.025 | 0.31 | 0.024 | 0.009 | 0.39 | 0.030 |
| 672 | Oxide 2 | 25 | 0.095 | 0.11 | 0.018 | 0.003 | 0.55 | 0.009 | 0.004 | 0.079 | 0.020 |
| 673 | Oxide 3 | 25 | 0.0037 | 0.006 | 0.002 | 0.0003 | 0.016 | 0.003 | 0.001 | 0.029 | 0.003 |

| SRM | Nominal Trace Composition (Parts Per Million by Weight) | | | | | | | | | | | |
|-----|---|------|------|-------|--------|--------|--------|--------|--------|--------|--------|-------|
| | Pb | Se | Bi | As | Sn | Sb | Cd | Ga | Ag | Te | Tl | Zn |
| 671 | 16 | 2.0 | 0.07 | (59) | (2.7) | (0.4) | (0.7) | (0.8) | (0.5) | (<0.2) | (<0.1) | (160) |
| 672 | 38 | 0.40 | 0.3 | (74) | (4) | (0.5) | (1.7) | (0.4) | (0.3) | (<0.2) | (<0.1) | (140) |
| 673 | 3.5 | 0.2 | 0.06 | (0.4) | (<0.5) | (<0.5) | (0.05) | (<0.1) | (<0.1) | (0.4) | (<0.1) | (1.7) |

Values in parentheses are not certified, but are given for information only.

Titanium-Base Alloys

| SRM | Type | Wt/ Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|-------------------|------------------------|--|-------------------|-------|-------------------|---------|------|
| | | | C | Mn | Cr | Cu | Mo | |
| 173b | Al-V | 50 | 0.025 | | | 0.008 | 0.013 | |
| 176 | Al-Sn | 100 | 0.015 | 0.0008 | | 0.003 | 0.0003 | |
| 641 | Mn (A) | Disk | | 6.6 ₈ | | | | |
| 642 | Mn (B) | Disk | | 9.0 ₈ | | | | |
| 643 | Mn (C) | Disk | | 11.6 ₈ | | | | |
| 644 | Cr-Fe-Mo (A) | Disk | | | 1.03 | | 3.61 | |
| 646 | Cr-Fe-Mo (C) | Disk | | | 3.43 | | 1.11 | |
| 647 | Al-Mo-Sn-Zr | 50 | 0.006 | | | | 1.96 | |
| 648 | Al-Sn-Zr-Cr-Mo | 50 | 0.011 | | 3.84 | | 3.75 | |
| 650 | Unalloyed A | 30 | | 0.016 | 0.002 | 0.033 | 0.002 | |
| 651 | Unalloyed B | 30 | | 0.005 | 0.037 | 0.032 | 0.031 | |
| 652 | Unalloyed C | 30 | | 0.046 | 0.082 | 0.081 | 0.039 | |
| 654a | Al-V (B) | Disk | (<0.1) | (0.20) | | | (<0.05) | |
| 1133 | Al-Sn-Zr-Cr-Mo | Disk | 0.011 | | | | 3.75 | |
| SRM | Fe | Al | V | Sn | Si | N | W | Zr |
| 173b | 0.23 | 6.36 | 4.31 | (0.03) | 0.046 | 0.015 | | |
| 176 | 0.07 ₀ | 5.16 | | 2.47 | | 0.01 ₀ | | |
| 644 | 1.36 | | | | | | | |
| 646 | 2.14 | | | | | | | |
| 647 | 0.075 | 5.88 | (<0.02) | 2.02 | | (<0.01) | | 3.90 |
| 648 | 0.15 | 5.13 | | 1.98 | 0.027 | (0.01) | | 1.84 |
| 650 | 0.024 | <0.01 | 0.009 | 0.03 | 0.004 | | 1.55 | |
| 651 | 0.058 | <0.006 | 0.021 | 0.026 | 0.011 | | 0.39 | |
| 652 | 0.67 | 0.039 | 0.024 | 0.053 | 0.16 | | 0.5 | |
| 654a | (0.20) | 6.3 ₄ | 3.9 ₅ | | | | | |
| 1133 | 0.15 | 5.13 | | 1.98 | 0.027 | (0.01) | | 1.84 |

Values in parentheses are not certified, but are given for information only.

Zinc-Base Alloys

| SRM | Type | Wt/ Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | | |
|---|-------------------------|------------------------|---|--------------------|----------------------|--------|-----------|-----------------------|--------|--------------------|-------|
| | | | Mn | Cu | Ni | Sn | Al | Cd | Fe | Pb | Mg |
| 94c | Die Casting Alloy | 150 | 0.014 | 1.01 | 0.006 | 0.006 | 4.13 | 0.002 | 0.018 | 0.006 | 0.042 |
| Chemical Composition (Nominal Weight Percent) | | | | | | | | | | | |
| SRM | Type | | Cu | Al | Mg | Fe | Pb | Cd | Sn | Cr | |
| 625 | Zinc-base A-ASTM AG 40A | | 0.034 | 3.06 | 0.070 | 0.036 | 0.0014 | 0.0007 | 0.0006 | 0.0128 | |
| 626 | Zinc-base B-ASTM AG 40A | | 0.056 | 3.56 | 0.020 | 0.103 | 0.0022 | 0.0016 | 0.0012 | 0.039 _s | |
| 627 | Zinc-base C-ASTM AG 40A | | 0.132 | 3.88 | 0.030 | 0.023 | 0.0082 | 0.0051 | 0.0042 | 0.0038 | |
| 628 | Zinc-base D-ASTM AC 41A | | 0.611 | 4.59 | 0.0094 | 0.066 | 0.0045 | 0.0040 | 0.0017 | 0.0087 | |
| 629 | Zinc-base E-ASTM AC 41A | | 1.50 | 5.15 | 0.094 | 0.017 | 0.0135 | 0.0155 | 0.012 | 0.0008 | |
| 630 | Zinc-base F-ASTM AC 41A | | 0.976 | 4.30 | 0.030 | 0.023 | 0.0083 | 0.0048 | 0.0040 | 0.0031 | |
| 631 | Zinc spelter (modified) | | 0.001 _s | 0.50 | (<0.001) | 0.005 | (0.001) | 0.0002 | 0.0001 | 0.0001 | |
| SRM | Mn | Ni | Si | In | Ga | Ca | Ag | Ge | | | |
| 625 | 0.031 | 0.0184 | 0.017 | | | | | | | | |
| 626 | 0.048 | 0.047 | 0.042 | | | | | | | | |
| 627 | 0.014 | 0.0029 | 0.021 | | | | | | | | |
| 628 | 0.0091 | 0.030 | 0.008 | | | | | | | | |
| 629 | 0.0017 | 0.0075 | 0.078 | | | | | | | | |
| 630 | 0.0106 | 0.0027 | 0.022 | | | | | | | | |
| 631 | 0.0001 _s | (<0.0005) | (0.002) | 0.002 _s | (0.00 ₂) | <0.001 | (<0.0005) | (0.000 ₂) | | | |

Values in parentheses are not certified, but are given for information only.

Zirconium-Base Alloys

| SRM | Type | Wt/Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | | | | | |
|--|-------------|--------------------|---|--------|--------|-------|--------|------|-------|------|------|--------|-------|------|
| | | | C | Mn | Hf | Cu | Ni | Cr | Ti | Sn | Fe | N | Al | |
| 360b | Zircaloy-4 | 100 | 0.011 | 0.0010 | 0.008 | 0.002 | 0.0025 | 0.10 | 0.002 | 1.55 | 0.21 | 0.0045 | 0.004 | |
| Chemical Composition (Nominal Parts Per Million) | | | | | | | | | | | | | | |
| SRM 31 mm D × 9.5 mm thick | Type | | Hf | C | Cr | Cu | Fe | Mn | Mo | Ni | N | Si | Ti | W |
| 1234 | Zirconium A | | 46 | (80) | (55) | (<10) | (240) | (10) | (2) | (20) | (14) | (40) | (20) | (25) |
| 1237 | Zircaloy D | | 31 | (100) | (1510) | (<10) | (1650) | (10) | (<10) | (40) | (19) | (35) | (30) | (25) |

Values in parentheses are not certified, but are given for information only.

Gases in Metals

These SRM's are for determining hydrogen, oxygen, and nitrogen by vacuum fusion, inert gas fusion, and neutron activation methods. SRM's 1095 to 1099 are sold only in a set as SRM 1089.

| SRM | Type | Form | Oxygen (ppm) | Hydrogen (ppm) | Nitrogen (ppm) |
|-------|--|-----------|-----------------|-------------------|-------------------|
| 352c | Unalloyed titanium for hydrogen | Platelets | | (IN PREP) | |
| 1090 | Ingot iron | Rod | (491) | | (60) |
| 1091a | Stainless steel (AISI 431) | Rod | 132.2 | | |
| 1093 | Valve steel | Rod | 60 | | |
| 1094 | Maraging steel | Rod | 4. ₅ | | (71) |
| 1089 | Set of 5: 1095, 1096, 1097, 1098, and 1099 | Rods | | | |
| 1754 | Low-Alloy Steel, AISI 4320 | Rod | 24 | | 81 |

Values in parentheses are not certified, but are given for information only.



Dale Friend removes some material from the large cone blender for homogeneity analysis. Dale's dedication and availability as supervisor of the preparation and packaging activities is a valued asset of the program.

High-Purity Metals

These SRM's are for determining impurity elements in high-purity metals. (See also specific metals.)

| SRM | Type | Unit Size | Chemical Composition (Nominal Parts Per Million by Weight) | | | | | | | | |
|-------|-------------------------------|------------------------------------|--|---------|-----------|---------|---------|-------|-------|------------------|-----------|
| | | | Cu | Ni | Sn | Pb | Zr | | | | |
| 685W* | High-Purity Gold (Wire) | 1.4 mm D × 102 mm long | 0.1 | (<0.05) | (<0.07) | | | | | | |
| 685R* | High-Purity Gold (Rod) | 5.9 mm D × 25 mm long | 0.1 | (<0.05) | (<0.07) | | | | | | |
| 680a | High-Purity Platinum (Wire) | 0.51 mm D: L1 (10 cm); L2 (1 m) | 0.1 | <1 | | <1 | <0.1 | | | | |
| 681 | Doped-Platinum (Wire) | 0.51 mm D: L1 (10 cm); L2 (1 m) | 5.0 | 0.5 | | 12 | 11 | | | | |
| 682* | High-Purity Zinc | Semicirc 57 mm D | 0.042 | (<0.1) | (0.02) | | | | | | |
| 683* | Zinc Metal | Semicirc 57 mm D | 5.9 | | (0.02) | 11.1 | | | | | |
| 728 | Zinc, Intermediate Purity | Shot, 450 g | 5.7 | | (0.02) | 11.1 | | | | | |
| 726 | Selenium, Intermediate Purity | Shot, 450 g | <1 | <0.5 | <1 | <1 | Mn <0.3 | | | | |
| C1257 | Aluminum, High Purity | Disk | (<0.1) | (<0.1) | (<0.1) | (<0.1) | (<0.1) | | | | |
| SRM | Ag | Mg | In | Fe | O | Pd | Au | Rh | Ir | Cd | Tl |
| 685W* | [0.1] | (<0.2) | 0.00, | 0.3 | [2] | | | | | | |
| 685R* | [0.1] | (<0.2) | 0.00, | 0.2 | [<2] | | | | | | |
| 680a | <0.1 | <1 | | 1.3 | 4 | 0.2 | <1 | <0.2 | <0.01 | | |
| 681 | 2.0 | 12 | | 5 | 7 | 6 | 9 | 9 | 11 | | |
| 682* | (0.02) | (<0.1) | | (0.1) | (<0.5) | | | | | (0.1) | |
| 683* | 1.3 | | | 2.2 | | | | | | 1.1 | (0.2) |
| 728 | 1.1 | | | 2.7 | | | | | | 1.1 _s | (0.2) |
| 726 | <1 | <1 | S 12 | 1 | Cr <1 | Mo <0.3 | Te 0.3 | As <2 | Al <1 | B <1 | Ca <1 |
| C1257 | Si 2.0 | 5.0 | | 1.0 | Cr (<0.1) | | | | | (<0.1) | Ca (<0.1) |

*Certificate gives upper limits for other elements found to be present.

Values in parentheses are not certified, but are given for information only.

Values in brackets are subject to greater error since only one method of analysis was employed.

RM 1R—Ultra-Purity Aluminum Polycrystalline Rods

These rods are intended for use in research on the mechanical and physical properties of extremely pure aluminum; e.g., in the determination of resistivity as a function of strain at cryogenic temperatures to facilitate the design of cryogenic magnets, or superconductor stabilizing elements. Unit of issue: 4.2 mm in diameter and 25.4 mm long.

Microanalytical

These SRM's provide a highly homogeneous material at microscopic spatial resolution. They are intended primarily for use in calibration of quantitative electron probe, secondary ion mass spectrometry, spark source mass spectrometry, and laser probe microanalytical techniques.

| SRM | Type | Unit Size |
|------|----------------------------------|------------------------------|
| 470 | Mineral Glasses (K-411 & K-412) | 2 Rods: 1×1×15 mm |
| 480 | Tungsten-20% Mo Alloy | Rod: 1 mm D, 1 mm long |
| 481 | Au-Ag Set | 6 Wire: 0.5 mm D, 50 mm long |
| 482 | Au-Cu Set | 6 Wire: 0.5 mm D, 50 mm long |
| 1871 | Glasses (K-456, K-493, & K-523) | 3 Rods: 2×2×20 mm |
| 1872 | Glasses (K-453, K-491, & K-968) | 3 Rods: 2×2×20 mm |
| 1873 | Glasses (K-458, K-489 & K-963) | 3 Rods: 2×2×20 mm |
| 1874 | Glasses (K-495, K-490, & K-546) | 3 Rods: 2×2×20 mm |
| 1875 | Glasses (K-496, K-497, & K-1013) | 3 Rods: 2×2×20 mm |
| 2063 | Thin Film Mg-Si-Ca-Fe | 3 mm diameter film |
| 2064 | Thin Film Al-Si-Ca-Ti-Mn-Zn | 3 mm diameter film |

Metals for Microanalysis

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | |
|-----|-----------------------|---|--------------------|--------------------|------|------|
| | | Au | Cu | Ag | W | Mo |
| 480 | Tungsten-20% Mo Alloy | | | | 78.5 | 21.5 |
| 481 | Au 100 A | 100.0 ₀ | | | | |
| | Au-20% Ag B | 80.0 ₅ | | 19.9 ₆ | | |
| | Au-40% Ag C | 60.0 ₅ | | 39.9 ₂ | | |
| | Au-60% Ag D | 40.0 ₃ | | 59.9 ₃ | | |
| | Au-80% Ag E | 22.4 ₃ | | 77.5 ₈ | | |
| | Ag 100 F | | | 100.0 ₀ | | |
| 482 | Au 100 A | 100.0 ₀ | | | | |
| | Au-20% Cu B | 80.1 ₅ | 19.8 ₃ | | | |
| | Au-40% Cu C | 60.3 ₆ | 39.6 ₄ | | | |
| | Au-60% Cu D | 40.1 ₀ | 59.9 ₂ | | | |
| | Au-80% Cu E | 20.1 ₂ | 79.8 ₅ | | | |
| | Cu 100 F | | 100.0 ₀ | | | |

Mineral Glasses for Microanalysis

| SRM 470 | | Composition (Nominal Weight Percent) | | | | | |
|---------|------------------|--------------------------------------|-------|-----|-------|-------|--------------------------------|
| Glass | SiO ₂ | FeO | | MgO | | CaO | Al ₂ O ₃ |
| K-411 | 54.30 | | 14.42 | | 14.67 | 15.47 | — |
| K-412 | 45.35 | | 9.96 | | 19.33 | 15.25 | 9.27 |

Glasses for Microchemical Analysis

| SRM 1871 | | | SRM 1872 | | | SRM 1873 | | | SRM 1874 | | | SRM 1875 | | | |
|--------------------------------------|---------|---------|----------|---------|---------|----------|---------|----------|----------|---------|----------|----------|----------|----------|----------|
| Glass | | | Glass | | | Glass | | | Glass | | | Glass | | | |
| K-456 | K-493 | K-523 | K-453 | K-491 | K-968 | K-458 | K-489 | K-963 | K-495 | K-490 | K-546 | K-496 | K-497 | K-1013 | |
| Composition (Nominal Weight Percent) | | | | | | | | | | | | | | | |
| Pb | 65.67 | 63.28 | 63.10 | 54.21 | 54.69 | 54.74 | — | (1.32) | — | — | (1.47) | — | — | (0.86) | — |
| Si | 13.37 | (13.09) | (12.94) | — | (0.11) | — | 23.05 | (22.23) | (21.96) | — | (0.19) | — | — | (0.13) | — |
| Ge | — | — | (0.20) | 28.43 | 26.10 | 25.93 | — | — | (0.47) | — | — | (0.50) | — | — | (0.34) |
| Ba | — | — | (0.61) | — | — | (0.46) | 41.79 | 39.53 | 39.21 | — | — | (0.99) | — | — | (0.52) |
| Zn | — | — | — | — | — | — | 3.01 | 2.93 | 2.95 | — | — | — | — | — | — |
| P | — | — | (0.24) | — | — | (0.21) | — | — | (0.33) | — | — | (0.42) | 32.98 | 31.59 | 32.26 |
| Mg | — | — | (0.12) | — | — | (0.22) | — | — | (0.34) | — | — | (0.17) | 6.65 | 6.49 | 5.86 |
| Al | — | (0.13) | — | — | (0.10) | — | — | (0.11) | — | 10.89 | (10.2) | (10.1) | 6.47 | 5.97 | 6.08 |
| B | — | — | — | — | — | — | — | — | — | (23.0) | (21.5) | (21.6) | — | [0.05] | — |
| Zr | — | (0.38) | (0.33) | — | (0.26) | (0.48) | — | (0.40) | (0.61) | — | (0.53) | (0.52) | — | (0.32) | (0.45) |
| Ti | — | (0.20) | (0.21) | — | (0.14) | (0.16) | — | (0.27) | (0.32) | — | (0.31) | (0.39) | — | (0.22) | (0.21) |
| Ce | — | (0.53) | — | — | (0.59) | — | — | [0.80] | — | — | (1.46) | — | — | (0.94) | — |
| Ta | — | (0.64) | — | — | (0.52) | — | — | (0.95) | — | — | (1.02) | — | — | (0.71) | — |
| Fe | — | (0.25) | — | — | (0.17) | — | — | (0.35) | — | — | (0.38) | — | — | (0.26) | — |
| Li | — | — | — | — | — | — | — | — | — | (2.3) | (2.2) | (2.2) | — | [0.0005] | — |
| Ni | — | — | (0.25) | — | — | (0.20) | — | — | (0.33) | — | — | (0.39) | — | — | (0.31) |
| Eu | — | — | (0.73) | — | — | (0.64) | — | — | (0.95) | — | — | (1.21) | — | — | (0.53) |
| U | — | — | (0.23) | — | — | (0.05) | — | — | (0.16) | — | — | (0.24) | — | — | (0.15) |
| Th | — | — | (0.08) | — | — | (0.12) | — | — | (0.06) | — | — | (0.16) | — | — | (0.10) |
| Cr | — | — | (0.20) | — | — | (0.19) | — | — | (0.31) | — | — | (0.14) | — | — | (0.14) |
| O | (20.35) | (20.58) | (20.80) | (16.73) | (16.45) | (16.67) | (31.86) | (31.70) | (32.00) | (63.49) | (60.74) | (61.36) | (53.90)* | (52.46)* | (53.05)* |
| Total | (99.39) | (99.08) | (100.19) | (99.37) | (99.13) | (100.07) | (99.71) | (100.59) | (100.00) | (99.68) | (100.01) | (100.39) | (100.00) | (100.00) | (100.00) |

Values in parentheses are for information only, they are *not certified*.

Values in brackets were calculated from the weight of material added to the melt, they are *not certified*.

*Oxygen values in SRM 1875 were calculated by difference, not by the stoichiometry of the oxides as was done for the other glasses.

Thin Film for X-Ray Spectrometry

This SRM is for standardizing chemical analysis by x-ray spectrometry and energy loss spectrometry on the analytical electron microscope

| SRM | Type | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|--|---|-------|-------|-------|-------|---------|
| | | Mg | Si | Ca | Fe | Ar | O |
| 2063 | Thin Film Mg-Si-Ca-Fe | 8.04 | 23.89 | 12.89 | 12.43 | (0.8) | (41.95) |
| 2064 | Thin Film Al-Si-Ca-Ti-Mn-Zn (IN PREP) | | | | | | |

Values in parentheses are not certified, but are given for information only.

Glass Fibers for Microanalysis—RM 8531

| | K-456 | K-493 | K-453 | K-491 | K-458 | K-489 | K-495 | K-490 | K-496 | K-497 |
|--------------------------------|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| | Chemical Composition (Nominal Weight Percent) | | | | | | | | | |
| SiO ₂ | 28.77 | 27.89 | — | 0.19 | 49.38 | 46.76 | — | 0.42 | — | 0.27 |
| PbO | 71.23 | 69.08 | 58.72 | 59.35 | — | 1.28 | — | 1.55 | — | 0.99 |
| GeO ₂ | — | — | 41.28 | 37.98 | — | — | — | — | — | — |
| BaO | — | — | — | — | 46.80 | 43.88 | — | — | — | — |
| ZnO | — | — | — | — | 3.82 | 3.72 | — | — | — | — |
| P ₂ O ₅ | — | — | — | — | — | — | — | — | 79.54 | 76.03 |
| MgO | — | — | — | — | — | — | — | — | 9.03 | 8.64 |
| Al ₂ O ₃ | — | 0.20 | — | 0.16 | — | 0.29 | 20.00 | 18.68 | 11.43 | 10.92 |
| B ₂ O ₃ | — | 0.14 | — | 0.11 | — | 0.20 | 75.00 | 70.00 | — | 0.15 |
| ZrO ₂ | — | 0.49 | — | 0.40 | — | 0.70 | — | 0.85 | — | 0.54 |
| TiO ₂ | — | 0.32 | — | 0.26 | — | 0.46 | — | 0.55 | — | 0.35 |
| CeO ₂ | — | 0.68 | — | 0.56 | — | 0.98 | — | 1.19 | — | 0.76 |
| Ta ₂ O ₅ | — | 0.88 | — | 0.72 | — | 1.26 | — | 1.53 | — | 0.98 |
| Fe ₂ O ₃ | — | 0.32 | — | 0.26 | — | 0.046 | — | 0.55 | — | 0.35 |
| Li ₂ O | — | 0.001 | — | 0.001 | — | 0.002 | 5.00 | 4.67 | — | 0.001 |

Primary, Working, and Secondary Chemicals

These SRM's are high-purity chemicals defined as primary, working, and secondary standards in accordance with recommendations of the Analytical Chemistry Section of the International Union of Pure and Applied Chemistry [Ref. Analyst 90, 251 (1965)]. These definitions are as follows:

Primary Standard:

a commercially available substance of purity 100 ± 0.02 percent (Purity 99.98 + percent).

Working Standard:

a commercially available substance of purity 100 ± 0.05 percent (Purity 99.95 + percent).

Secondary Standard:

a substance of lower purity which can be standardized against a primary grade standard.

| SRM | Type | Wt/Unit (grams) | Certified Use | Purity Stoichiometric |
|------|---------------------------------|--------------------|---|-----------------------------|
| 17d | Sucrose | 60 | Polarimetric Value | (99.9) ^a |
| 40h | Sodium Oxalate | 60 | Reductometric Value | 99.972 |
| 41c | Dextrose (D-Glucose) | 70 | Reductometric Value | 99.9 |
| 83d | Arsenic Trioxide | 60 | Reductometric Value | 99.9926 |
| 84j | Potassium Hydrogen Phthalate | 60 | Acidimetric Value | 99.996 |
| 136e | Potassium Dichromate | 60 | Oxidimetric Value | 99.984 |
| 350a | Benzoic Acid | 30 | Acidimetric Value | 99.9958 |
| 723a | Tris(hydroxymethyl)aminomethane | 50 | Basimetric Value | 99.9703 |
| 951 | Boric Acid | 100 | Acidimetric and Boron Isotopic Value | 100.00 |
| 987 | Strontium Carbonate | 1 | Assay and Isotopic | 99.98 |
| 999 | Potassium Chloride | 60 | Assay Standard for: Potassium Chloride | 99.98 ₁ 99.99 |

^aSucrose=Moisture <0.02 percent, Ash <0.005 percent.

Microchemical

| SRM | Type | Wt/ Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | | | |
|------|----------------------|------------------------|---|------|-------|----|----|-------|-------|--------------------|
| | | | C | H | N | Br | Cl | F | S | CH ₃ O- |
| 141c | Acetanilide | 2 | 71.09 | 6.71 | 10.36 | | | | | |
| 142 | Anisic Acid | 2 | | | | | | | | 20.40 |
| 143c | Cystine | 2 | 29.99 | 5.03 | 11.66 | | | | | 26.69 |
| 148 | Nicotinic Acid | 2 | 58.54 | 4.09 | 11.38 | | | | | |
| 2141 | Urea | 2 | | | 46.63 | | | | | |
| 2142 | o-Bromobenzoic Acid | 2 | | | 39.80 | | | | | |
| 2143 | p-Fluorobenzoic Acid | 2 | | | | | | | 13.54 | |
| 2144 | m-Chlorobenzoic Acid | 2 | | | | | | 22.62 | | |

Spectrometric Solutions

These SRM's are intended as standard stock solutions for use in atomic absorption spectrometry, optical emission (plasma) spectrometry, or any other analytical technique that requires aqueous solutions for calibrating instruments. Each SRM is a single element solution of 50 mL with a concentration of 10 mg/mL, except where noted.

| SRM | Element | Acid Concentration |
|------|--------------|------------------------------|
| 3101 | Aluminum | HCl 10% |
| 3102 | Antimony | HCl 50% |
| 3103 | Arsenic | HCl 15% |
| 3104 | Barium | HCl 10% |
| 3105 | Beryllium | HCl 10% |
| 3106 | Bismuth | HNO ₃ 10% |
| 3107 | Boron (5.00) | Water |
| 3108 | Cadmium | HNO ₃ 10% |
| 3109 | Calcium | HCl 10% |
| 3110 | Cerium | HNO ₃ 10% |
| 3111 | Cesium | HCl 1% |
| 3112 | Chromium | HCl 10% |
| 3113 | Cobalt | HNO ₃ 10% |
| 3114 | Copper | HNO ₃ 10% |
| 3115 | Dysprosium | HCl 10% |
| 3116 | Erbium | HCl 10% |
| 3117 | Europium | HCl 10% |
| 3118 | Gadolinium | HCl 10% |
| 3119 | Gallium | HCl 10% |
| 3120 | Germanium | Oxalic Acid 10% |
| 3121 | Gold | HCl 10% |
| 3122 | Hafnium | HNO ₃ 10% + HF 2% |
| 3123 | Holmium | HCl 10% |
| 3124 | Indium | HCl 10% |
| 3125 | Iridium | IN PREP |
| 3126 | Iron | HCl 10% |
| 3127 | Lanthanum | HCl 10% |
| 3128 | Lead | HNO ₃ 10% |
| 3129 | Lithium | HCl 1% |
| 3130 | Lutetium | HCl 10% |
| 3131 | Magnesium | HCl 10% |
| 3132 | Manganese | HNO ₃ 10% |
| 3133 | Mercury | HNO ₃ 10% |
| 3134 | Molybdenum | HCl 10% |
| 3135 | Neodymium | HCl 10% |
| 3136 | Nickel | HNO ₃ 10% |
| 3137 | Niobium | 5% HNO ₃ + 2% HF |
| 3138 | Palladium | HCl 10% |
| 3139 | Phosphorus | HCl 0.05% |
| 3140 | Platinum | HCl 10% |
| 3141 | Potassium | HCl 1% |
| 3142 | Praseodymium | HCl 10% |
| 3143 | Rhenium | HNO ₃ 10% |
| 3144 | Rhodium | IN PREP |
| 3145 | Rubidium | HCl 1% |

Spectrometric Solutions (Continued)

| SRM | Element | Acid Concentration |
|------|-----------------|-------------------------------------|
| 3146 | Ruthenium | IN PREP |
| 3147 | Samarium | HCl 10% |
| 3148 | Scandium | HCl 10% |
| 3149 | Selenium | HNO ₃ 10% |
| 3150 | Silicon | Water |
| 3151 | Silver | HNO ₃ 10% |
| 3152 | Sodium | HCl 1% |
| 3153 | Strontium | HCl 10% |
| 3154 | Sulfur | H ₂ SO ₄ 0.1% |
| 3155 | Tantalum | 5% HNO ₃ + 2% HF |
| 3156 | Tellurium | HCl 10% |
| 3157 | Terbium | HCl 10% |
| 3158 | Thallium | HNO ₃ 10% |
| 3159 | Thorium | HNO ₃ 10% |
| 3160 | Thulium | HCl 10% |
| 3161 | Tin | HCl 60% |
| 3162 | Titanium | HCl 40% |
| 3163 | Tungsten | 7% HNO ₃ + 4% HF |
| 3164 | Uranium | HNO ₃ 10% |
| 3165 | Vanadium (5.00) | HNO ₃ 10% |
| 3166 | Ytterbium | HCl 10% |
| 3167 | Yttrium | HCl 10% |
| 3168 | Zinc | HCl 10% |
| 3169 | Zirconium | 10% HNO ₃ + 2% HF |



Dolly Downs helps implement the technical and editorial review process for all certificates in addition to performing many other reporting and documentation activities that help keep track of the many certification projects. Her sincerity and neatness are instrumental in keeping all these activities functioning smoothly.

SRM's now cover a wide range of materials, including organics in shale oil and elemental composition of glasses for microanalysis.



Multielement Spectrometric Solutions

| Element | Source, Purity | Concentration, $\mu\text{g/mL}$ |
|--|-------------------------------------|---------------------------------|
| SRM 3171 Multielement Mix A Standard Solution | | |
| Aluminum | Metal, (99.99+) | 100.0 ± 0.5 |
| Beryllium | Metal, (99.9) | 10.0 ± 0.1 |
| Cadmium | Metal, (99.99+) | 100.0 ± 0.5 |
| Chromium | Metal, (99.99+) | 100.0 ± 0.5 |
| Iron | Metal, (99.99+) | 100.0 ± 0.5 |
| Magnesium | Metal, (99.99) | 100.0 ± 0.5 |
| Manganese | Metal, (99.99) | 100.0 ± 0.5 |
| Nickel | Metal, (99.999) | 100.0 ± 0.5 |
| Potassium | KCl, (99.98) | 500.0 ± 2.5 |
| Sodium | NaCl, (99.98) | 100.0 ± 0.5 |
| SRM 3172 Multielement Mix B Standard Solution | | |
| Arsenic | As_2O_3 , (99.9926) | 200.0 ± 1.0 |
| Barium | BaCO_3 , (99.99) | 10.0 ± 0.1 |
| Calcium | CaCO_3 , (99.99) | 10.0 ± 0.1 |
| Cobalt | Metal, (99.99) | 100.0 ± 0.5 |
| Copper | Metal, (99.998) | 100.0 ± 0.5 |
| Lead | Metal, (99.99) | 100.0 ± 0.5 |
| Selenium | Metal, (99.99) | 500.0 ± 2.5 |
| Silver | Metal, (99.999+) | 100.0 ± 0.5 |
| Strontium | SrCO_3 , (99.99) | 10.0 ± 0.1 |
| Zinc | Metal, (99.99+) | 100.0 ± 0.5 |
| SRM 3173 Multielement Mix C Standard Solution | | |
| IN PREP | | |
| SRM 3174 Multielement Mix D Standard Solution | | |
| Aluminum | Metal, (99.99+) | 100.0 ± 0.5 |
| Beryllium | Metal, (99.9) | 100.0 ± 0.5 |
| Boron | H_3BO_3 , (99.99) | 100.0 ± 0.5 |
| Cadmium | Metal, (99.999) | 100.0 ± 0.5 |
| Gold | Metal, (99.999) | 100.0 ± 0.5 |
| Hafnium | Metal, (99.95) | 100.0 ± 0.5 |
| Iron | Metal, (99.99) | 100.0 ± 0.5 |
| Lead | Metal, (99.995) | 100.0 ± 0.5 |
| Titanium | Metal, (99.99) | 50.0 ± 0.25 |
| Zirconium | Metal, (99.84) | 100.0 ± 1.0 |

Anion Ion Chromatographic Solutions

These SRM's are single-component solutions prepared gravimetrically for use in anion ion chromatography, or any other technique that requires aqueous standard solutions for calibration on control materials.

| SRM | Anion | Wt/Unit (mL) | Concentration ($\mu\text{g/g}$) |
|------|-----------|-----------------|--------------------------------------|
| 3181 | Sulfate | 50 | 1000 |
| 3182 | Chloride | 50 | 1000 |
| 3183 | Fluoride | 50 | 1000 |
| 3184 | Bromide | IN PREP | 50 |
| 3185 | Nitrate | IN PREP | 50 |
| 3186 | Phosphate | IN PREP | 50 |



John Savoy, shown here placing labels on bottles of a new SRM, performs his responsibilities in packaging materials with the consistency and orderliness necessary to accurately package and label a large number of different materials.

Clinical Laboratory

These SRM's are for calibrating apparatus and validating analytical methods used in clinical and pathology laboratories. See also: Spectrophotometric SRM's and Temperature SRM's.

| SRM | Type | Associated Publications | Purity % | Wt/Unit |
|---------|---|-------------------------|----------------------|-------------------------|
| 900 | Antiepilepsy Drug Level Assay (phenytoin, ethosuximide, phenobarbital, and primidone) | | 4 drugs 3 levels | Set of 4 vials |
| 909 | Human Serum | | # | Set of 6 vials |
| 910 | Sodium Pyruvate | | 98.7 | 25 g |
| 911b | Cholesterol | | 99.8 | 2 g |
| 912a | Urea | | 99.9 | 25 g |
| 913 | Uric Acid | | 99.7 | 10 g |
| 914a | Creatinine | | 99.7 | 10 g |
| 915 | Calcium Carbonate | SP 260-36 | 99.9 | 20 g |
| 916a | Bilirubin | | 98.3 | 100 mg |
| 917a | D-Glucose (Dextrose) | | 99.9 | 25 g |
| 918 | Potassium Chloride | SP 260-63 | 99.9 | 30 g |
| 919a | Sodium Chloride | SP 260-60 | IN PREP | 30 g |
| 920 | D-Mannitol | | 99.8 | 50 g |
| 921 | Cortisol (Hydrocortisone) | | 98.9 | 1 g |
| 922 | Tris(hydroxymethyl) aminomethane | | 99.99 | 25 g |
| 923 | Tris(hydroxymethyl) aminomethane HC1 | | 99.69 | 35 g |
| 924 | Lithium Carbonate | SP 260-69 | 100.0 | 30 g |
| 925 | VMA (4-hydroxy-3-methoxymandelic acid) | | 99.4 | 1 g |
| 926 | Bovine Serum Albumin (Powder) | | * | 5 g |
| 927a | Bovine Serum Albumin (7% Solution) | | * | 10 vials, 2.15 mL ea. |
| 928 | Lead Nitrate | | 100.00 | 30 g |
| 929 | Magnesium Gluconate Dihydrate | | (100.1) | 5 g |
| 937 | Iron Metal | | 99.90 | 50 g |
| 938 | 4-Nitrophenol | | (99.75) | 15 g |
| 956 | Electrolytes in Serum for ISE | IN PREP | | |
| 968 | Fat Soluble Vitamins in Human Serum | | 3 levels | Set of 6 vials |
| 998 | Angiotensin I (Human) | | 94.1 | 0.5 mg |
| 1589 | Polychlorinated Biphenyls (PCB's) in Human Serum (as Aroclor 1260) | | — | Set of 3 bottles |
| 1595 | Tripalmitin | | 99.5 | 2 g |
| 1598 | Inorganic Constituents in Bovine Serum | | | Set of 2 vials |
| 1599 | Anticonvulsant Drug Level Assay (valproic acid and carbamazepine) | | 2 drugs/ 3 levels | Set of 4 vials |
| 1700a | Blood Gas: CO ₂ -10%, Bal N ₂ | | — | Cylinder, 20 cubic feet |
| 1701a | Blood Gas: CO ₂ -5%, O ₂ -12%, Bal N ₂ | | — | Cylinder, 20 cubic feet |
| 1702a | Blood Gas: CO ₂ -5%, O ₂ -20%, Bal N ₂ | | — | Cylinder, 20 cubic feet |
| 1703a | Blood Gas: CO ₂ -10%, O ₂ -7%, Bal N ₂ | | — | Cylinder, 20 cubic feet |
| 1951a | Cholesterol in Human Serum (Frozen) | IN PREP | | Set of 6 bottles |
| 1952a | Cholesterol in Human Serum (Freeze-dried) | IN PREP | | Set of 3 bottles |
| RM 8430 | Aspartate Aminotransferase (AST) Human Erythrocyte Source | | — | Set of 3 bottles |

*Conforms to NCCLS specification ACC-1.

Electrolytes, selected organics.



The SRM program has continually expanded in response to needs for certified reference materials. The SRMs shown here have primary relevance to clinical and health related areas of analysis.

Serum Reference Materials

These materials are for calibrating instrumentation and evaluating the reliability of analytical methods for the determination of major, minor, and trace constituents in blood serum, plasma, and similar biological fluids.

| Constituent | Concentrations | | | |
|-------------|-----------------------|----------|-----------------------|--------|
| | SRM 909 (Procedure A) | | SRM 909 (Procedure B) | |
| Cadmium | 1.46 | ng/mL/g | 1.24 | ng/mL |
| Calcium | 3.560 | mmo1/L/g | 3.013 | mmo1/L |
| Chloride | 128.0 | mmo1/L/g | 108.4 | mmo1/L |
| Chromium | 108 | ng/mL/g | 91.3 | ng/mL |
| Cholesterol | 4.326 | mmo1/L/g | 3.66 | mmo1/L |
| Copper | 1.29 | µg/mL/g | 1.10 | µg/mL |
| Creatinine | 0.179 | mmo1/L/g | 0.152 | mmo1/L |
| Glucose | 7.38 | mmo1/L/g | 6.25 | mmo1/L |
| Iron | 2.34 | µg/mL/g | 1.98 | µg/mL |
| Lead | 23.7 | ng/mL/g | 20.0 | ng/mL |
| Lithium | 1.945 | mmo1/L/g | 1.65 | mmo1/L |
| Magnesium | 1.425 | mmo1/L/g | 1.21 | mmo1/L |
| Potassium | 4.155 | mmo1/L/g | 3.52 | mmo1/L |
| Sodium | 158.4 | mmo1/L/g | 134.1 | mmo1/L |
| Urea | 11.39 | mmo1/L/g | 9.64 | mmo1/L |
| Uric Acid | 0.567 | mmo1/L/g | 0.480 | mmo1/L |
| Vanadium | 3.19 | ng/mL/g | 2.70 | ng/mL |

Biological Materials

These SRM's are intended for use in the calibration of apparatus and methods used in the analysis of biological materials for major, minor, and trace constituents.

Food and Beverage

| SRM | 1548 | 1549 | 1566a* | 1567a | 1568a | 1569 | 1577a | RM 50 |
|-----------------|-------------|--|---------------|--------------|--------------|---------------|---------------|----------------|
| Type | Total Diet | Non-fat Milk Powder | Oyster Tissue | Wheat Flour | Rice Flour | Brewers Yeast | Bovine* Liver | Albacore* Tuna |
| Unit Size | | 100 g | | | 80 g | 50 g | 50 g | 70 g |
| ELEMENTS | | Nominal Composition in $\mu\text{g/g}$, unless otherwise noted. | | | | | | |
| Aluminum | 202.5 | (2) | | 5.7 | 4.4 | | (2) | (4.39) |
| Antimony | (0.01) | (0.00027) | | | (0.0005) | | (0.003) | |
| Arsenic | 14.0 | (0.0019) | | (0.006) | 0.29 | | 0.047 | (3.3) |
| Bromine | | (12) | | (6) | (8) | | (9) | |
| Cadmium | 4.15 | 0.0005 | | 0.026 | 0.022 | | 0.44 | |
| Calcium | 0.196 wt. % | 1.30 wt. % | 0.17 wt. % | 0.0191 wt. % | 0.011 | | 120 | |
| Chlorine | 0.829 wt. % | 1.09 wt. % | 0.87 wt. % | (565) | (300) | | 0.28 wt. % | |
| Chromium | 1.43 | 0.0026 | | — | — | 2.12 | | |
| Cobalt | 0.57 | (0.0041) | | (0.006) | (0.018) | | 0.21 | |
| Copper | 66.3 | 0.7 | 2.6 | 2.1 | 2.4 | | 158 | |
| Fluorine | (240) | (0.20) | | — | — | | | |
| Iodine | 4.46 | 3.38 | | (0.0009) | (0.009) | | | |
| Iron | 539 | 1.78 | 32.2 | 14.1 | 7.4 | | 194 | |
| Lead | 0.371 | 0.019 | | (<0.020) | (<0.010) | | 0.135 | |
| Magnesium | 0.118 wt. % | 0.120 wt. % | 0.055 wt. % | 0.040 wt. % | 0.056 wt. % | | 600 | |
| Manganese | 12.3 | 0.26 | 5.1 | 9.4 | 20.0 | | 9.9 | |
| Mercury | 0.0642 | 0.0003 | | (0.0005) | 0.0058 | | 0.004 | (0.95) |
| Molybdenum | | (0.34) | | 0.48 | 1.46 | | 3.5 | |
| Nickel | 2.25 | | 0.41 | — | (0.16) | | | |
| Nitrogen | 6.81 wt. % | | 3.5 wt. % | — | — | | (10.7%) | |
| Phosphorus | 0.623 wt. % | 1.06 wt. % | 0.32 wt. % | 0.134 wt. % | 0.153 wt. % | | 1.11 wt. % | |
| Potassium | 0.790 wt. % | 1.69 wt. % | 0.60 wt. % | 0.133 wt. % | 0.1280 wt. % | | 0.996 wt. % | |
| Rubidium | (3) | (11) | | 0.68 | 6.14 | | 12.5 | |
| Selenium | 2.21 | 0.11 | 0.25 | 1.1 | 0.38 | | 0.71 | (3.6) |
| Silver | 1.68 | (<0.0003) | | — | — | | 0.04 | |
| Sodium | 0.417 wt. % | 0.497 wt. % | 0.63 wt. % | 6.1 | 6.6 | | 0.243% | |
| Strontium | 11.1 | | | — | — | | 0.138 | |
| Sulfur | 0.862 wt. % | 0.351 wt. % | 0.28 wt. % | 0.165 wt. % | 0.120 wt. % | | 0.78% | |
| Tellurium | | | | — | (<0.002) | | | |
| Thallium | | | | — | — | | (0.003) | |
| Thorium | (0.04) | | | — | | | | |
| Tin | (3) | (<0.02) | | (0.0033) | (0.0047) | | | |
| Uranium | 0.132 | | | (.0003) | (0.0003) | | 0.00071 | |
| Vanadium | 4.68 | | | (0.011) | (0.007) | | 0.099 | |
| Zinc | 830 | 46.1 | 30.8 | 11.6 | 19.4 | | 123 | (13.6) |

Values in parentheses are not certified, but are given for information only.

*Indicates freeze-dried.

Food and Beverage (Continued)

SRM 1845
Type Cholesterol in Whole Egg Powder

Unit Size 3 bottles, 8 g each
Cholesterol 19.0 mg/g

1563
Cholesterol and Fat Soluble Vitamins in Coco-
nut Oil
Set of 10
IN PREP



Bob Alvarez, a project manager for OSRM, friendly and outgoing, provides a rich background of experience and knowledge in the SRM program. Bob provides leadership for many of the organic, clinical, chemical and environmental SRM's.

Ethanol Solutions See also: Alcohol in Reference Fuels

| SRM | Type | Certified Constituent | Wt/Unit |
|------------|----------------------------|---|--|
| 1590 | Stabilized Wine | Ethanol: 18.57% by volume | Set of 10, 10-mL vials |
| 1828 | Ethanol-Water Solutions | Ethanol: 95.629 wt% Ethanol: 0.2992 wt% Ethanol: 0.1487 wt% | Set: 1, 15-mL vial 2, 3-mL vials 2, 3-mL vials |

Agricultural

| SRM | 1515 | 1547 | 1572 | 1573a | 1575 | 2695 | RM 8412 | RM 8413 |
|----------------|--------------|--|---------------|---------------|--------------|---------------------------|---------------|----------------|
| Type | Apple Leaves | Peach Leaves | Citrus Leaves | Tomato Leaves | Pine Needles | Fluoride in Vegetation | Corn Stalk | Corn Kernel |
| Unit Size | IN PREP | IN PREP | 70 g | IN PREP | 70 g | IN PREP | 34 g | 47 g |
| ELEMENT | | Nominal Composition in $\mu\text{g/g}$, unless otherwise noted. | | | | | | |
| Aluminum | | | 92 (0.04) | | 545 (0.2) | | | (4) |
| Antimony | | | 3.1 | | 0.21 | | | |
| Arsenic | | | 21 | | | | | |
| Barium | | | | | | | | |
| Boron | | | | | | | | |
| Bromine | | | (8.2) | | (9) | | | |
| Cadmium | | | 0.03 | | (<0.5) | | | |
| Calcium | | | 3.15% | | 0.41% | | (2160) | (42) |
| Cerium | | | (0.28) | | (0.4) | | | |
| Cesium | | | (0.098) | | | | | |
| Chlorine | | | (414) | | | | (2440) | (450) |
| Chromium | | | 0.8 | | 2.6 | | | |
| Cobalt | | | (0.02) | | (0.1) | | | |
| Copper | | | 16.5 | | 3.0 | | (8) | (3.0) |
| Europium | | | (0.01) | | (0.006) | | | |
| Fluorine | | | | | | | | (0.24) |
| Iodine | | | 1.84 | | | | (0.65) | |
| Iron | | | 90 | | 200 | | (139) | (23) |
| Lanthanum | | | (0.19) | | (0.2) | | | |
| Lead | | | 13.3 | | 10.8 | | | |
| Magnesium | | | 0.58% | | | | (1600) | (990) |
| Manganese | | | 23 | | 675 | | (15) | (4.0) |
| Mercury | | | 0.08 | | 0.15 | | | |
| Molybdenum | | | 0.17 | | | | | |
| Nickel | | | 0.6 | | (3.5) | | | |
| Nitrogen | | | (2.86%) | | (1.2%) | | (6970) | (13750) |
| Phosphorous | | | 0.13% | | 0.12% | | | |
| Potassium | | | 1.82% | | 0.37% | | (17350) | (3570) |
| Rubidium | | | 4.84 | | 11.7 | | | |
| Samarium | | | (0.052) | | | | | |
| Scandium | | | (0.01) | | (0.03) | | | |
| Selenium | | | (0.025) | | | | | |
| Sodium | | | 160 | | | | (0.016) | (0.004) |
| Strontium | | | 100 | | 4.8 | | (28) | |
| Sulfur | | | 0.407% | | | | (12) | |
| Tellurium | | | (0.02) | | | | | |
| Thallium | | | (<0.01) | | (0.05) | | | |
| Thorium | | | | | 0.037 | | | |
| Tin | | | (0.24) | | | | | |
| Uranium | | | (<0.15) | | 0.020 | | | |
| Zinc | | | 29 | | | | (32) | (15.7) |

Values in parentheses are not certified, but are given for information only.

Environmental Materials

Analyzed Gases

These SRM's are for calibrating apparatus used to measure various components of gas mixtures and atmospheric pollutants. All cylinders conform to the appropriate DOT specifications.

| SRM | Type | Certified Component | Nominal Concentration |
|-------|---|--|----------------------------------|
| 1811 | Aromatic Organic Gases in Nitrogen | | |
| | Benzene | C ₆ H ₆ | 0.25 μmole/mole |
| | Toluene | C ₆ H ₅ CH ₃ | 0.25 μmole/mole |
| | Chlorobenzene | C ₆ H ₅ Cl | 0.25 μmole/mole |
| | Bromobenzene | C ₆ H ₅ Br | 0.25 μmole/mole |
| 1812 | Aromatic Organic Gases in Nitrogen | | |
| | Benzene | C ₆ H ₆ | 10 μmole/mole |
| | Toluene | C ₆ H ₅ CH ₃ | 10 μmole/mole |
| | Chlorobenzene | C ₆ H ₅ Cl | 10 μmole/mole |
| | Bromobenzene | C ₆ H ₅ Br | 10 μmole/mole |
| 1813 | Aliphatic Organic Gases in Nitrogen | | |
| | Carbon Tetrachloride | CCl ₄ | 0.25 μmole/mole |
| | Chloroform | CHCl ₃ | 0.25 μmole/mole |
| | Tetrachloroethylene | CCl ₂ CCl ₂ | 0.2 μmole/mole |
| | Vinyl Chloride | CH ₂ CHCl | 0.25 μmole/mole |
| 1814 | Aliphatic Organic Gases in Nitrogen | | |
| | Carbon Tetrachloride | CCl ₄ | 10 μmole/mole |
| | Chloroform | CHCl ₃ | 10 μmole/mole |
| | Tetrachloroethylene | CCl ₂ CCl ₂ | 10 μmole/mole |
| | Vinyl Chloride | CH ₂ CHCl | 10 μmole/mole |
| 1805 | Benzene in Nitrogen | C ₆ H ₆ | 0.25 μmole/mole |
| 1806 | Benzene in Nitrogen | C ₆ H ₆ | 10 μmole/mole |
| 1700a | Blood Gas: CO ₂ -10%, Bal N ₂ | CO ₂ in N ₂ | 10% |
| 1701a | Blood Gas: CO ₂ -5%, O ₂ -12%, Bal N ₂ | CO ₂ & O ₂ in N ₂ | 5% Carbon Dioxide, 12% Oxygen |
| 1702a | Blood Gas: CO ₂ -5%, O ₂ -20%, Bal N ₂ | CO ₂ & O ₂ in N ₂ | 5% Dioxide, 20% Oxygen |
| 1703a | Blood Gas: CO ₂ -10%, O ₂ -7%, Bal N ₂ | CO ₂ & O ₂ in N ₂ | 10% Carbon Dioxide, 7% Oxygen |
| 1670 | Carbon Dioxide in Air | CO ₂ | 330 μmole/mole |
| 1671 | Carbon Dioxide in Air | CO ₂ | 340 μmole/mole |
| 1672 | Carbon Dioxide in Air | CO ₂ | 350 μmole/mole |
| 2607 | Carbon Dioxide and Nitrous Oxide in Air | | |
| | Carbon Dioxide | CO ₂ | 340 μmole/mole |
| | Nitrous Oxide | N ₂ O | 300 μmole/mole |
| 2608 | Carbon Dioxide and Nitrous Oxide in Air | | |
| | Carbon Dioxide | CO ₂ | 340 μmole/mole |
| | Nitrous Oxide | N ₂ O | 300 μmole/mole |
| 2609 | Carbon Dioxide and Nitrous Oxide in Air | | |
| | Carbon Dioxide | CO ₂ | 380 μmole/mole |
| | Nitrous Oxide | N ₂ O | 330 μmole/mole |
| 2610 | Carbon Dioxide and Nitrous Oxide in Air | | |
| | Carbon Dioxide | CO ₂ | 380 μmole/mole |
| | Nitrous Oxide | N ₂ O | 330 μmole/mole |

Analyzed Gases (Continued)

| SRM | Type | Certified Component | Nominal Concentration |
|-------|---|--|---|
| 1674b | Carbon Dioxide in Nitrogen | CO ₂ | 7.0 mole percent |
| 1675b | Carbon Dioxide in Nitrogen | CO ₂ | 14.0 mole percent |
| 2619a | Carbon Dioxide in Nitrogen | CO ₂ | 0.5 mole percent |
| 2620a | Carbon Dioxide in Nitrogen | CO ₂ | 1.0 mole percent |
| 2621a | Carbon Dioxide in Nitrogen | CO ₂ | 1.5 mole percent |
| 2622a | Carbon Dioxide in Nitrogen | CO ₂ | 2.0 mole percent |
| 2623a | Carbon Dioxide in Nitrogen | CO ₂ | 2.5 mole percent |
| 2624a | Carbon Dioxide in Nitrogen | CO ₂ | 3.0 mole percent |
| 2625a | Carbon Dioxide in Nitrogen | CO ₂ | 3.5 mole percent |
| 2626a | Carbon Dioxide in Nitrogen | CO ₂ | 4.0 mole percent |
| 2612a | Carbon Monoxide in Air | CO | 10 μmole/mole |
| 2613a | Carbon Monoxide in Air | CO | 20 μmole/mole |
| 2614a | Carbon Monoxide in Air | CO | 45 μmole/mole |
| 1677c | Carbon Monoxide in Nitrogen | CO | 10 μmole/mole |
| 1678c | Carbon Monoxide in Nitrogen | CO | 50 μmole/mole |
| 1679c | Carbon Monoxide in Nitrogen | CO | 100 μmole/mole |
| 1680b | Carbon Monoxide in Nitrogen | CO | 500 μmole/mole |
| 1681b | Carbon Monoxide in Nitrogen | CO | 1000 μmole/mole |
| 2635a | Carbon Monoxide in Nitrogen | CO | 25 μmole/mole |
| 2636a | Carbon Monoxide in Nitrogen | CO | 250 μmole/mole |
| 2637a | Carbon Monoxide in Nitrogen | CO | 2500 μmole/mole |
| 2638a | Carbon Monoxide in Nitrogen | CO | 5000 μmole/mole |
| 2639a | Carbon Monoxide in Nitrogen | CO | 1 mole percent |
| 2640 | Carbon Monoxide in Nitrogen | CO | 2 mole percent |
| 2641 | Carbon Monoxide in Nitrogen | CO | 4 mole percent |
| 2642a | Carbon Monoxide in Nitrogen | CO | 8 mole percent |
| 2725 | Carbon Monoxide and Propane in Nitrogen | CO C ₃ H ₈ | 1.6 mole percent 600 μmole/mole |
| 2726 | Carbon Monoxide and Propane in Nitrogen | CO C ₃ H ₈ | 8 mole percent 3000 μmole/mole |
| 2727 | Carbon Monoxide, Propane and Carbon Dioxide in Nitrogen | CO C ₃ H ₈ CO ₂ | 1.6 mole percent 600 μmole/mole 11 mole percent |
| 2728 | Carbon Monoxide, Propane and Carbon Dioxide in Nitrogen | CO C ₃ H ₈ CO ₂ | 8 mole percent 3000 μmole/mole 14 mole percent |
| 1658a | Methane in Air | CH ₄ | 1 μmole/mole |
| 1659a | Methane in Air | CH ₄ | 10 μmole/mole |
| 1660a | Methane-Propane in Air | CH ₄ C ₃ H ₈ | 4 μmole/mole 1 μmole/mole |
| 1683b | Nitric Oxide in Nitrogen | NO | 50 μmole/mole |
| 1684b | Nitric Oxide in Nitrogen | NO | 100 μmole/mole |
| 1685b | Nitric Oxide in Nitrogen | NO | 250 μmole/mole |
| 1686b | Nitric Oxide in Nitrogen | NO | 500 μmole/mole |
| 1687b | Nitric Oxide in Nitrogen | NO | 1000 μmole/mole |
| 2627a | Nitric Oxide in Nitrogen | NO | 5 μmole/mole |
| 2628a | Nitric Oxide in Nitrogen | NO | 10 μmole/mole |
| 2629a | Nitric Oxide in Nitrogen | NO | 20 μmole/mole |

Analyzed Gases (Continued)

| SRM | Type | Certified Component | Nominal Concentration | |
|-------|------------------------------------|---|-----------------------|------------------------------|
| 2630 | Nitric Oxide in Nitrogen | NO | 1500 | $\mu\text{mole}/\text{mole}$ |
| 2631 | Nitric Oxide in Nitrogen | NO | 3000 | $\mu\text{mole}/\text{mole}$ |
| 2654 | Nitrogen Dioxide in Air | NO ₂ | 500 | $\mu\text{mole}/\text{mole}$ |
| 2655 | Nitrogen Dioxide in Air | NO ₂ | 1000 | $\mu\text{mole}/\text{mole}$ |
| 2656 | Nitrogen Dioxide in Air | NO ₂ | 2500 | $\mu\text{mole}/\text{mole}$ |
| 2657a | Oxygen in Nitrogen | O ₂ | 2 | mole percent |
| 2658a | Oxygen in Nitrogen | O ₂ | 10 | mole percent |
| 2659a | Oxygen in Nitrogen | O ₂ | 21 | mole percent |
| 1665b | Propane in Air | C ₃ H ₈ | 3 | $\mu\text{mole}/\text{mole}$ |
| 1666b | Propane in Air | C ₃ H ₈ | 10 | $\mu\text{mole}/\text{mole}$ |
| 1667b | Propane in Air | C ₃ H ₈ | 50 | $\mu\text{mole}/\text{mole}$ |
| 1668b | Propane in Air | C ₃ H ₈ | 100 | $\mu\text{mole}/\text{mole}$ |
| 1669b | Propane in Air | C ₃ H ₈ | 500 | $\mu\text{mole}/\text{mole}$ |
| 2645a | Propane in Nitrogen | C ₃ H ₈ | 500 | $\mu\text{mole}/\text{mole}$ |
| 2646a | Propane in Nitrogen | C ₃ H ₈ | 1000 | $\mu\text{mole}/\text{mole}$ |
| 2647a | Propane in Nitrogen | C ₃ H ₈ | 2500 | $\mu\text{mole}/\text{mole}$ |
| 2648a | Propane in Nitrogen | C ₃ H ₈ | 5000 | $\mu\text{mole}/\text{mole}$ |
| 2649 | Propane in Nitrogen | C ₃ H ₈ | 1 | mole percent |
| 2650 | Propane in Nitrogen | C ₃ H ₈ | 2 | mole percent |
| 2651 | Propane in Nitrogen and Oxygen | C ₃ H ₈ | 0.01 | mole percent |
| 2652 | Propane in Nitrogen and Oxygen | O ₂ | 5.0 | mole percent |
| 1661a | Sulfur Dioxide in Nitrogen | C ₃ H ₈ | 0.01 | mole percent |
| 1662a | Sulfur Dioxide in Nitrogen | O ₂ | 10.0 | mole percent |
| 1663a | Sulfur Dioxide in Nitrogen | SO ₂ | 500 | $\mu\text{mole}/\text{mole}$ |
| 1664a | Sulfur Dioxide in Nitrogen | SO ₂ | 1000 | $\mu\text{mole}/\text{mole}$ |
| 1665a | Sulfur Dioxide in Nitrogen | SO ₂ | 1500 | $\mu\text{mole}/\text{mole}$ |
| 1666a | Sulfur Dioxide in Nitrogen | SO ₂ | 2500 | $\mu\text{mole}/\text{mole}$ |
| 1693a | Sulfur Dioxide in Nitrogen | SO ₂ | 50 | $\mu\text{mole}/\text{mole}$ |
| 1694a | Sulfur Dioxide in Nitrogen | SO ₂ | 100 | $\mu\text{mole}/\text{mole}$ |
| 1696 | Sulfur Dioxide in Nitrogen | SO ₂ | 3500 | $\mu\text{mole}/\text{mole}$ |
| 1804 | Ambient Toxic Organics in Nitrogen | (Fifteen components, call for details.) | | |
| 1808 | Tetrachloroethylene in Nitrogen | C ₂ Cl ₄ | 0.25 | $\mu\text{mole}/\text{mole}$ |
| 1809 | Tetrachloroethylene in Nitrogen | C ₂ Cl ₄ | 10 | $\mu\text{mole}/\text{mole}$ |
| 2730 | Hydrogen Sulfide in Nitrogen | H ₂ S | 5 | $\mu\text{mole}/\text{mole}$ |
| 2731 | Hydrogen Sulfide in Nitrogen | H ₂ S | 20 | $\mu\text{mole}/\text{mole}$ |

Permeation Devices

These SRM's are for calibrating air pollution monitoring apparatus, and may be used to verify air pollution analytical methods and procedures. Each tube is individually certified.

SRM's 1625, 1626, and 1627 are certified over the temperature range of 20 to 30 °C. SRM 1629a is calibrated at 25.0 °C only; cannot be shipped by air.

| SRM | Type | Tube Length (cm) | Permeation Rate ($\mu\text{g}/\text{min}$) at 25 °C | Typical Concentrations (ppm) Flow Rates (liters per minute) | | |
|------|--------------------------------|------------------|---|--|--------|--------|
| | | | | 1 | 5 | 10 |
| 1625 | Sulfur Dioxide Permeation Tube | 10 | 2.8 | 1.07 | 0.214 | 0.107 |
| 1626 | Sulfur Dioxide Permeation Tube | 5 | 1.4 | 0.535 | 0.107 | 0.0535 |
| 1627 | Sulfur Dioxide Permeation Tube | 2 | 0.56 | 0.214 | 0.0428 | 0.0214 |



Gerald Roderick (left) and Bill Dorko, of the Gas and Particulate Science Division, connect some gas cylinders to a sampling manifold in preparation for certification measurements.



Mark Cronise checks the certificate and supplemental documentation for a gas composition SRM against the SRM and serial numbers in preparation for shipment to a customer. Mark's initiative in implementing inventory control and processing procedures have benefited the SRM program.



Wayne Wolf, a Research Associate with OSRM, is providing initiative and resourcefulness in the development and management of food and nutritional SRM's. Here he examines some materials that are to be distributed to the technical divisions for certification analyses.

Analyzed Liquids and Solids

These SRM's are for analysis of materials for constituents of interest in health or environmental problems.
See also: Clinical SRM's and Industrial Hygiene SRM's.

| SRM | Type | Unit Size | Certified Element | | | | |
|-------|--|-----------|-------------------|----------------|--------|------------|------------|
| | | | Lead | Nickel | Sulfur | Mercury | Vanadium |
| 1579 | Powdered Lead Base Paint | 35 g | 11.87% | | | | |
| 1618 | Vanadium and Nickel in Residual Fuel Oil | 100 mL | | 75 µg/g (4.3%) | | | 423 µg/g |
| 1630 | Trace Mercury in Coal | 50 g | | | | 0.13 µg/g | |
| 1638b | Lead in Reference Fuel | 12 vials | | 767 µg/g | | | |
| 1641b | Mercury in Water (µg/mL) | 6×20 mL | | | | 1.52 µg/mL | |
| 2712 | Lead in Reference Fuel | 6×20 mL | 11.4 µg/g | | | | |
| 2713 | Lead in Reference Fuel | 6×20 mL | 19.4 µg/g | | | | |
| 2714 | Lead in Reference Fuel | 6×20 mL | 28.1 µg/g | | | | |
| 2715 | Lead in Reference Fuel | 6×20 mL | 784 µg/g | | | | |
| 8505 | Vanadium in Crude Oil | 250 mL | | | | | (390 µg/g) |

Simulated Rainwaters

These materials were developed to aid in the analysis of acidic rainwater by providing stable, homogeneous material as control standards at two levels of acidity.

NOTE: Values in parentheses are not certified.


| SRM | Type | Unit of Issue | |
|-------------------------------------|---------|--------------------------------|---------|
| Constituent Element/Parameter | IN PREP | Set 4: 2-50mL each of 2 levels | |
| | | 2694-I | 2694-II |
| pH, 25 °C | | (Target Values) | |
| Specific Conductance (µS/cm, 25 °C) | 26 | 4.3 | 3.6 |
| Acidity, meq/L | 0.05 | 26 | 130 |
| Fluoride, mg/L | 0.05 | 0.05 | 0.3 |
| Chloride, mg/L | (0.2) | (0.2) | (1.0) |
| Nitrate, mg/L | 0.5 | 0.5 | 7.0 |
| Sulfate, mg/L | 2 | 2 | 11 |
| Sodium, mg/L | 0.2 | 0.2 | 0.4 |
| Potassium, mg/L | 0.05 | 0.05 | 0.1 |
| Ammonium, mg/L | — | — | (1.0) |
| Calcium, mg/L | 0.01 | 0.01 | 0.05 |
| Magnesium, mg/L | 0.02 | 0.02 | 0.05 |

Values in parentheses are not certified, but are given for information only.

Alcohols in Reference Fuels

These SRM's are for calibrating instruments and validating methods used to determine various alcohols in gasoline. Each SRM is issued as a set of sealed 20-mL ampoules.

| SRM | Type | Nominal Concentration in Weight Percent | | | Methanol and t-Butanol |
|------|----------------------------|---|----------|---------|------------------------|
| | | Wt/Unit | Methanol | Ethanol | |
| 1829 | Alcohols in Reference Fuel | Set (6) | 0.335 | 11.39 | 10.33 + 6.63 |
| 1837 | Methanol and t-Butanol | Set (5) | | | 10.33 + 6.63 |
| 1838 | Ethanol | Set (5) | | 11.39 | |
| 1839 | Methanol | Set (5) | 0.335 | | |



Tom Gills, Program Manager for Production and Certification, has implemented many improvements in the production and certification procedures and documentation through his initiative and determination to achieve these objectives.



These SRM's represent the large number of environmental SRM's that are now available.

Sulfur in Fossil Fuels

| SRM | Type | Unit Size | Sulfur Wt.% | Furnace Ash Wt.% | HHV2 | |
|-------|-------------------------------|-----------|-------------|------------------|---------------------|-------------------------|
| | | | | | MJ•Kg ⁻¹ | (Btu•lb ⁻¹) |
| 1616 | Sulfur in Kerosine | 100 mL | 0.0152 | | | |
| 1617 | Sulfur in Kerosine | 100 mL | 0.169 | | | |
| 1619 | Sulfur in Residual Fuel Oil | 100 mL | 0.719 | | | |
| 1620b | Sulfur in Residual Fuel Oil | IN PREP | 4.3 | (Target Value) | | |
| 1621c | Sulfur in Residual Fuel Oil | 100 mL | 1.040 | | | |
| 1622c | Sulfur in Residual Fuel Oil | 100 mL | 2.012 | | | |
| 1623b | Sulfur in Residual Fuel Oil | IN PREP | 0.30 | (Target Value) | | |
| 1624b | Sulfur in Distillate Fuel Oil | IN PREP | 0.20 | (Target Value) | | |
| 2717 | Sulfur in Residual Fuel Oil | IN PREP | 3.5 | (Target Value) | | |
| 2682 | Coal (Sub-bituminous) | 50 g | 0.47 | 6.37 | | |
| 2683a | Coal (Bituminous) | 50 g | 1.89 | 6.85 | 32.40 | (13930) |
| 2684a | Coal (Bituminous) | 50 g | 3.05 | 11.09 | 29.03 | (12480) |
| 2685 | Coal (Bituminous) | 50 g | 4.62 | 16.53 | 27.38 | (11770) |
| 2692 | Sulfur in Coal, 1% | 50 g | 1.115 | 7.9 | 30.61 | (13160) |
| 2720 | Synthetic Sulfur in Oil | IN PREP | | | | |
| 8590 | High Sulfur Gas Oil Feed | IN PREP | | | | |

NOTE: The calorific values (MJ•Kg⁻¹) may decrease upon the aging or normal oxidation of the coals. NIST will continue to monitor these calorific values and report any substantive change to the purchaser.

Trace Elements

| SRM | 1632b | 1633a | 1634b | 1635 | 1643b | 1646 | 1648 | 2689 | 2690 | 2691 | 2704 |
|--------------|---|--------------------|-------------|-------------------------|-------------|----------------------------|---------------------------|--------------------|--------------------|--------------------|------------------------------|
| Type | Coal (Bituminous) | Coal Fly Ash | Fuel Oil | Coal (Subbituminous) | Water | Estua- rine Sediment | Urban Particu- late | Coal Fly Ash | Coal Fly Ash | Coal Fly Ash | Buffalo River Sediment |
| Unit Size | 55 g | 75 g | 100 mL | 75 g | 950 mL | 75 g | 2 g | 30 g | 30 g | 30 g | 50 g |
| ELEMENT | Nominal Concentrations in $\mu\text{g/g}$, unless otherwise noted. | | | | | | | | | | |
| Aluminum | 0.855% | 14.3% | (16) | (0.32%) | | 6.25% | 3.42% | 12.94% | 12.35% | 9.81% | 6.11 |
| Antimony | (0.24) | 6.8 | | (0.14) | | (0.4) | (45) | | | | 3.79 |
| Arsenic | 3.72 | 145 | 0.12 | 0.42 | (49) ng/g | 11.6 | 115 | | | | 23.4 |
| Barium | 67.5 | (0.15%) | (1.3) | | 44 ng/g | | (737) | | | | 414 |
| Beryllium | | (12) | | | 19 ng/g | (1.5) | | | | | |
| Bismuth | | | | | (11) ng/g | | | | | | |
| Bromine | (17) | | | | B (94) ng/g | | | | | | |
| Cadmium | 0.0573 | 1.00 | | 0.03 | 20 ng/g | 0.36 | (500) | | | | 3.45 |
| Calcium | 0.204% | 1.11% | (15) | | | 0.83% | 75 | | | | 2.60 |
| Carbon | 78.11% | | | | | | | 2.18% | 5.71% | 18.45% | 3.348 |
| Cerium | (9) | (180) | | (3.6) | | (80) | (55) | | | | (72) |
| Cesium | (0.44) | (11) | | | | (3.7) | (3) | | | | (6) |
| Chlorine | (1260) | | | | | | (0.45%) | | | | (<0.01) |
| Chromium | (11) | 196 | (0.7) | 2.5 | 18.6 ng/g | 76 | 403 | | | | 135 |
| Cobalt | 2.29 | (46) | 0.32 | (0.65) | 26 ng/g | 10.5 | (18) | | | | 14.0 |
| Copper | 6.28 | 118 | | 3.6 | 21.9 ng/g | 18 | 609 | | | | 98.6 |
| Europium | (0.17) | (4) | | (0.06) | | (1.5) | (0.8) | | | | (1.3) |
| Gallium | | (58) | | (1.05) | | | (1.4) | | | | (15) |
| Germanium | | | | | | | | | | | |
| Hafnium | (0.43) | (8) | | (0.29) | | | (4.4) | | | | (8) |
| Hydrogen | 5.07% | | | | | | | | | | |
| Indium | | | | | | | (1.0) | | | | |
| Iodine | | | | | | | (20) | | | | |
| Iron (Total) | 0.759% | 9.4% | 31.6 | 0.239% | 99 ng/g | 3.35% | 3.91% | 9.32% | 3.57% | 4.42% | (2) |
| Lanthanum | (5.1) | | | | | (42) | | | | | 4.11 |
| Lead | 3.67 | 72.4 | (2.8) | 1.9 | 23.7 ng/g | 28.2 | 0.655% | | | | 161 |
| Lithium | (10) | | | | | (49) | | | | | (50) |
| Magnesium | 0.0383% | 0.455% | | | | 1.09% | (0.8%) | | | | 1.20 |
| Manganese | 12.4 | 179 | 0.23 | 21.4 | 28 ng/g | 375 | (860) | 0.61% | 1.53% | 3.12% | 555 |
| Mercury | | 0.16 | (<0.001) | | | 0.063 | | (0.03%) | (0.03%) | (0.02%) | 1.44 |
| Molybdenum | (0.9) | (29) | | | 85 ng/g | (2.0) | | | | | |
| Nickel | 6.10 | 127 | 28 | 1.74 | 49 ng/g | 32 | 82 | | | | 44.1 |
| Nitrogen | 1.56% | | | | | | | | | | |
| Phosphorus | | | | | | 0.054% | (1.4%) | 0.10% | 0.52% | 0.51% | 0.0998 |
| Potassium | 0.0748% | 1.88% | | | | 1.05% | | 2.20% | 1.04% | 0.34% | 2.00 |
| Rubidium | 5.05 | 131 | | | | (87) | (52) | | | | (100) |
| Samarium | (0.87) | | | | | (4.4) | | | | | (6.7) |
| Scandium | (1.9) | (40) | | | | (7) | | | | | (12) |
| Selenium | 1.29 | 10.3 | 0.18 | (0.63) | 9.7 ng/g | (10.8) | (0.6) | 27 | | | (1.1) |
| Silicon | (1.4%) | 22.8% | | 0.9 | | (31%) | | | 24.06% | 25.85% | 29.08 |
| Silver | | | | | | | | | | | |
| Sodium | 0.0515% | 0.17% | (90) | (0.24%) | 9.8 ng/g | (2.0%) | (6) | 0.425% | 0.25% | 0.24% | 1.09% |
| Strontium | (102) | 830 | | | 227 ng/g | | | (0.07%) | (0.20%) | (0.27%) | (0.547) |
| Sulfur | 1.89% | (0.18%) | 2.80% | 0.33% | | (0.96%) | (5%) | | 0.15% | 0.83% | (0.4) |
| Tellurium | | | | | | (0.5) | | | | | |
| Thallium | | 5.7 | | | 8.0 ng/g | (10) | | | | | 1.2 |
| Thorium | 1.342 | 24.7 | | 0.62 | | (0.51%) | (7.4) | | | | (9.2) |
| Titanium | 0.0454% | (0.8%) | | (0.02%) | | | (0.40%) | | 0.75% | 0.52% | 0.90% |
| Tungsten | (0.48) | | | | | | (4.8) | | | | 0.457 |
| Uranium | 0.436 | 10.2 | | 0.24 | | | | 5.5 | | | 3.13 |
| Vanadium | (14) | 297 | 55.4 | 5.2 | 45.2 ng/g | 94 | 140 | | | | 95 |
| Zinc | 11.89 | 220 | 3.0 | 4.7 | 66 ng/g | 138 | 0.476% | | | | 438 |

Values in parentheses are not certified, but are given for information only.

Organic Constituents

| SRM | Type | Unit of Issue |
|-------|---|--------------------------|
| 1491 | Aromatic Hydrocarbons in Hexane | 5 ampoules |
| 1492 | Chlorinated Pesticides in Hexane | 5 ampoules |
| 1493 | Polychlorinated Biphenyls in 2,2,4 Trimethylpentane | 5 ampoules |
| 1580 | Shale Oil | Set of 5, 2mL/ampoules |
| 1581 | Polychlorinated Biphenyls in Oil | Set of 8, 5mL/ampoules |
| 1582 | Petroleum Crude Oil | Set of 5, 2mL/ampoules |
| 1583 | Chlorinated Pesticides in <i>Isooctane</i> | Set of 6, 2mL/ampoules |
| 1584 | Phenols in Methanol | Set of 5, 2mL/ampoules |
| 1585 | Chlorinated Biphenyls | Set of 5, 1.2mL/ampoules |
| 1586 | Isotopically Labelled Priority Pollutants | Set of 6, 2mL/ampoules |
| 1587 | Nitro PAH in Methanol | Set of 4, 1mL/ampoules |
| 1588 | Organics in Cod Liver Oil | Set of 5, 1.2mL/ampoules |
| 1589 | Polychlorinated Biphenyls in Human Serum (as Aroclor 1260) | Set of 3 |
| 1596 | Dinitropyrene Isomers and 1-Nitropyrene in Methylene Chloride | Set of 5, 1.3mL/ampoules |
| 1597 | Complex Mixture of Polycyclic Aromatic Hydrocarbons | Set of 4, 1.2mL/ampoules |
| 1614 | Dioxin (2,3,7,8 TCDD) in <i>Isooctane</i> | Set of 6, 1.2mL/ampoules |
| 1639 | Halocarbons (in Methanol) | Set of 5, 1.5mL/ampoules |
| 1647a | Priority Pollutant PAH (in Acetonitrile) | Set of 5, 1.2mL/ampoules |
| 1649 | Urban Dust/Organics | 10 grams |
| 1650 | Diesel Particulate Matter | Set of 5, 100mg/ampoules |
| 1939 | Polychlorinated Biphenyls in Sediments | IN PREP |
| 1940 | Polychlorinated Biphenyls in Sediments | IN PREP |
| 1941 | Organics in Marine Sediment | IN PREP |
| 1974 | Organics in Mussel Tissue | IN PREP |
| 1975 | Diesel Particulate Bioassay | IN PREP |



Anna Carroll, punctual and accurate, receives and processes customers orders by phone and mail.

Organic Constituents (Continued)

| SRM | 1491 | 1580 | 1582 | 1597 | | 1647 | 1649 | 1650 |
|---|---------------------|---------------------|---------------------|-------------------|----------------------|----------------------|---------------------|---------------------|
| Constituents | ($\mu\text{g/g}$) | ($\mu\text{g/g}$) | ($\mu\text{g/g}$) | (μg) | ($\mu\text{g/mL}$) | ($\mu\text{g/mL}$) | ($\mu\text{g/g}$) | ($\mu\text{g/g}$) |
| Anthracene | 11.69 | | | 101 | 87.4 | 3.29 | | |
| Benz[a]anthracene | 5.37 | | 3.0 | 98.6 | 85.3 | 5.03 | 2.6 | 6.5 |
| Benzo[e]pyrene | 10.14 | 21 | 1.1 | 95.8 | 82.9 | 5.30 | 2.9 | 1.2 |
| Benzo[e]pyrene | 8.40 | 18 | | | | | | (10) |
| Fluoranthene | 8.84 | 54 | 2.5 | | | 10.1 | 7.1 | 51 |
| o-Cresol | | 385 | | | | | | |
| Phenol | | 407 | | | | | | |
| Perylene | 10.65 | 3.4 | 31 | 26.1 | 22.6 | | | (0.13) |
| Pyrene | 8.81 | 104 | | 235 | 204 | 9.84 | | 48 |
| 2,6-Dimethylphenol | | 175 | | | | | | |
| Benzo[f]quinoline (5,6-Benzoquinoline) | | 16 | | | | | | |
| Naphthalene | 10.30 | | | 1160 | 1000 | 22.5 | | |
| Acenaphthylene | 10.40 | | | | | 19.1 | | |
| Acenaphthene | 10.89 | | | | | 21.0 | | |
| 1-Nitropyrene | | | | | | | 19 | |
| Fluorene | 10.87 | | | | | 4.92 | | |
| Phenanthrene | 10.48 | | 101 | 462 | 400 | 5.06 | | (71) |
| Chrysene | 10.50 | | | 71.7 | 62.0 | 4.68 | | (22) |
| Benzo[b]fluoranthene | 7.85 | | | | | 5.11 | | |
| Benzo[k]fluoranthene | 8.33 | | | | | 5.02 | | (2.1) |
| Benzo[ghi]perylene | 7.90 | | | 53.7 | 46.5 | 4.01 | 4.5 | 2.4 |
| Dibenz[a,h]anthracene | 7.74 | | | | | 3.68 | | |
| Indeno[1,2,3-cd]pyrene | 9.40 | | | 602 | 52.1 | 4.06 | 3.3 | (0.23) |
| Dibenzothiophene | | | 33 | | | | | |
| Fluoranthene | | | | 322 | 278 | | | |
| Tryphénylene | | | | 12.1 | 10.5 | | | |
| 1-Methylnaphthalene | 12.4 | | | | | | | |
| Biphenyl | 10.46 | | | | | | | |
| 2,6-Dimethylnaphthalene | 10.8 | | | | | | | |
| 2,3,5-Trimethylnaphthalene | 9.9 | | | | | | | |
| 1-Methylphenanthrene | 10.4 | | | | | | | |

Values in parentheses are not certified, but are given for information only.

SRM 1639—Certified Concentration of Halocarbons at 23 ± 3 °C.

| Compound | Concentration, $\text{ng}/\mu\text{L}$ |
|----------------------|--|
| Chloroform | 6235 |
| Chlorodibromomethane | 124.6 |
| Bromodichloromethane | 389.9 |
| Bromoform | 86.5 |
| Carbon Tetrachloride | 157.0 |
| Trichloroethylene | 85.8 |
| Tetrachloroethylene | 40.6 |

Organic Constituents (Continued)

SRM 1492 Chlorinated Pesticides in Hexane

| Pesticide | Concentration (ng/g) |
|--------------------|----------------------|
| Hexachlorobenzene | 308 |
| gamma-HCH | 310 |
| Heptachlor | 299 |
| Aldrin | 304 |
| Heptachlor epoxide | 307 |
| cis-Chlordane | 305 |
| trans-Nonachlor | 297 |
| Dieldrin | 307 |
| Mirex | 306 |
| 2,4'-DDE | 303 |
| 4,4'-DDE | 306 |
| 2,4'-DDD | 299 |
| 4,4'-DDD | 296 |
| 2,4'-DDT | 307 |
| 4,4'-DDT | 302 |

SRM 1581 Polychlorinated Biphenyls in Oils

| Matrix | Aroclor Type | Concentration ($\mu\text{g/g}$) |
|-----------------|--------------|-----------------------------------|
| Motor Oil | 1242 | 100 |
| Motor Oil | 1260 | 100 |
| Transformer Oil | 1242 | 100 |
| Transformer Oil | 1260 | 100 |

SRM 1583 Chlorinated Pesticides in 2,2,4-Trimethylpentane

| Pesticide | Concentrations | |
|---------------------|---------------------|-----------------------------|
| | ($\mu\text{g/g}$) | ($\mu\text{g/mL}$, 23 °C) |
| Y-BHC (Lindane) | 1.11 | 0.77 |
| d-BHC | 0.76 | 0.53 |
| Aldrin | 0.86 | 0.59 |
| Heptachlor Epoxide | (0.997) | |
| 4,4'-DDE (p,p'-DDE) | 1.23 | 0.85 |
| 4,4'-DDT (p,p'-DDT) | 1.90 | 1.31 |

SRM 1584 Priority Pollutant Phenols in Methanol

| Compound | Concentration ($\mu\text{g/mL}$, 23 °C) |
|-----------------------|---|
| 2-Chlorophenol | 64.4 |
| Phenol | 29.7 |
| 2-Nitrophenol | 25.2 |
| 2,4-Dimethylphenol | 51.6 |
| 2,4-Dichlorophenol | 35.6 |
| 4-Chloro-m-cresol | 27.4 |
| 2,4,6-Trichlorophenol | 20.4 |
| 4-Nitrophenol | 20.7 |
| 4,6-Dinitro-o-cresol | 20.1 |
| Pentachlorophenol | 15.4 |
| 2,4-Dinitrophenol | (22.4) |

Organic Constituents (Continued)

SRM 1585 Chlorinated Biphenyls in 2,2,4-Trimethylpentane (Iso octane)

| Chlorinated Biphenyl | Concentration | |
|-----------------------------------|---------------|-----------------|
| | µg/g | µg/mL at 23.0°C |
| 4-Chlorobiphenyl | 43.3 | 29.9 |
| 4,4'-Dichlorobiphenyl | 9.53 | 6.57 |
| 2,4,4'-Trichlorobiphenyl | 3.70 | 2.55 |
| 2,2',5,5'-Tetrachlorobiphenyl | 7.72 | 5.32 |
| 3,3',4,4'-Tetrachlorobiphenyl | 6.62 | 4.56 |
| 2,2',4,5,5'-Pentachlorobiphenyl | 5.24 | 3.61 |
| 2,2',3,4,4',5'-Hexachlorobiphenyl | 2.37 | 1.63 |
| 2,2',4,4',5,5'-Hexachlorobiphenyl | 3.06 | 2.11 |

SRM 1586 Isotopically Labeled and Unlabeled Priority Pollutants in Methanol

| Compound | Concentrations (µg/g) | |
|----------------------------|-----------------------|------------------|
| | 1586-1 (unlabeled) | 1586-2 (labeled) |
| Carbon tetrachloride | 128.5 | 124.4 |
| Benzene | 101.1 | 99.0 |
| Chlorobenzene | 133.0 | 144.0 |
| Phenol | 117.0 | 116.0 |
| Nitrobenzene | 126.0 | 134.5 |
| 2-Nitrophenol | 103.6 | 101.9 |
| 2,4-Dichlorophenol | 102.5 | 82.2 |
| Naphthalene | 126.5 | 126.6 |
| Bis(2-ethylhexyl)phthalate | 63.9 | 60.4 |
| Benzo[a]pyrene | 49.2 | 44.1 |

SRM 1587 Nitrated Polycyclic Aromatic Hydrocarbons in Methanol

| Compound | Concentrations | |
|--------------------------|----------------|----------------|
| | (µg/g) | (µg/mL, 23 °C) |
| 2-Nitrofluorene | 9.67 | 7.64 |
| 9-Nitroanthracene | 5.01 | 3.96 |
| 3-Nitrofluoranthene | 9.24 | 7.30 |
| 1-Nitropyrene | 8.95 | 7.07 |
| 7-Nitrobenz[a]anthracene | 9.27 | 7.32 |
| 6-Nitrochrysene | 8.13 | 6.42 |
| 6-Nitrobenzo[a]pyrene | (6.1) | (4.8) |

Organic Constituents (Continued)

SRM 1588 Organics in Cod Liver Oil

| Compound | Concentration (ng/g) |
|-------------------|----------------------|
| Hexachlorobenzene | 148 |
| alpha-HCH | 86 |
| trans-Chlordane | 50 |
| cis-Chlordane | 158 |
| trans-Nonachlor | 209 |
| Dieldrin | 150 |
| 4,4'-DDD | 277 |
| 4,4'-DDE | 641 |
| 2,4'-DDT | 156 |
| 4,4'-DDT | 529 |
| PCB 101 | 261 |
| PCB 138 | 276 |
| PCB 170 | 45 |
| PCB 180 | 107 |

SRM 1589 PCB's in Human Serum

| Compound | Concentrations | |
|--------------|----------------|---------|
| | (ng/g) | (ng/mL) |
| Aroclor 1260 | 106.0 | 107.9 |

SRM 1596 Dinitropyrene Isomers and 1-Nitropyrene in Methylene Chloride

| Compound | Concentration | |
|-------------------|---------------|-------|
| | µg/g | µg/mL |
| 1-Nitropyrene | 4.38 | 5.81 |
| 1,3-Dinitropyrene | 2.10 | 2.79 |
| 1,6-Dinitropyrene | 4.82 | 6.39 |
| 1,8-Dinitropyrene | 7.90 | 10.48 |

SRM 1614 Dioxin (2,3,7,8-TCDD in Isooctane)

| Compound | Concentrations | |
|-------------------------------|----------------|----------------|
| | (ng/g) | (ng/mL, 23 °C) |
| 2,3,7,8-TCDD | 98.3 | 67.8 |
| 2,3,7,8-TCDD- ¹³ C | 95.6 | 65.9 |

GC/MS System Performance

These SRM's are for evaluating the sensitivity of gas chromatographic/mass spectrometry (GC/MS) instrumentation. They consist of two concentrations each of methyl stearate and benzophenone.

| SRM | Type | Concentrations (ng/µL) | | Unit Size |
|------|--------------------------|------------------------|--------------|------------------|
| | | Methyl Stearate | Benzophenone | |
| 1543 | GC/MS System Performance | 0.99; 4.98 | 1.01; 5.01 | 1 Set, 4 vials |
| 8443 | GC/MS System Performance | 0.99; 4.98 | 1.01; 5.01 | 5 Sets, 20 vials |

Industrial Hygiene

These SRM's were developed for industrial hygiene analyses to provide reference materials for toxicology research and for monitoring human exposure to selected toxic elements.

Freeze-Dried Urine

SRM's 2670, 2671a, and 2672a consist of freeze-dried urine in 30 mL serum bottles. The freeze-dried urine SRM's are to be reconstituted by the addition of 20 mL of pure water to each bottle. Each unit contains a set of four bottles, two bottles each at normal and elevated levels.

| SRM | 2670 | | 2671a | | 2672a | |
|-----------|-----------|----------------|-----------|----------------|-----------|----------------|
| | Type | Toxic Metals | | Fluoride | | Mercury |
| | Low Level | Elevated Level | Low Level | Elevated Level | Low Level | Elevated Level |
| Aluminum | (0.18) | µg/mL | (0.18) | µg/mL | | |
| Arsenic | (0.06) | µg/mL | 0.48 | µg/mL | | |
| Beryllium | (<0.0005) | µg/mL | (0.033) | µg/mL | | |
| Cadmium | (0.00040) | µg/mL | 0.088 | µg/mL | | |
| Calcium | 0.105 | mg/mL | 0.105 | mg/mL | | |
| Chloride | (4.4) | mg/mL | (4.4) | mg/mL | | |
| Chromium | (0.013) | µg/mL | 0.085 | µg/mL | | |
| Copper | 0.13 | µg/mL | 0.37 | µg/mL | | |
| Gold | (0.008) | ng/mL | (0.24) | µg/mL | 0.55 mg/L | 5.7 mg/L |
| Lead | (0.01) | µg/mL | 0.109 | µg/mL | | |
| Magnesium | 0.063 | mg/mL | 0.063 | mg/mL | | |
| Manganese | (0.03) | µg/mL | (0.33) | µg/mL | | |
| Mercury | (0.002) | µg/mL | 0.105 | µg/mL | | |
| Nickel | (0.07) | µg/mL | (0.30) | µg/mL | | |
| Platinum | (0.008) | ng/mL | (0.12) | µg/mL | | |
| Potassium | (1.5) | mg/mL | (1.5) | mg/mL | | |
| Selenium | 0.030 | µg/mL | 0.46 | µg/mL | | |
| Sodium | 2.62 | mg/mL | 2.62 | mg/mL | | |
| Sulfate | (1.3) | mg/mL | (1.3) | mg/mL | | |
| Vanadium | | | (0.12) | µg/mL | | |

Values in parentheses are not certified, but are given for information only.

Freeze-Dried Urine (Continued)

| SRM | 1507 | 1508 | RM 8444 |
|-------------------------|--------------|------------------|--|
| Type | THC in Urine | Cocaine in Urine | Cotinine in Urine |
| Unit Size | 3 bottles | IN PREP | 4 vials |
| Certified Concentration | 20 ng/mL | | Blank level 0.8 ng/g Low level 54 ng/g High level 488 ng/g |

Thin Films for X-ray Fluorescence

These SRM's are for standardizing x-ray spectrometers. They may be useful in elemental analysis of particulate matter collected on filter media, and where x-ray spectrometer calibration functions are determined using thin film standards. Each SRM is individually certified and consists of a silica-base glass film (0.5 μm thick) deposited on a 47 mm diameter polycarbonate filter, mounted on an aluminum ring.

| SRM | Type | Chemical Composition (Nominal $\mu\text{g}/\text{cm}^2$) | | | | | | | | | | | |
|------|-----------------|---|----|----|----|----|----|----|----|----|----|---|----|
| | | Al | Ca | Co | Cu | Fe | Pb | K | Mn | Si | Ti | V | Zn |
| 1832 | Thin-Glass Film | 15 | 20 | 1 | 2 | | | | 5 | 36 | | 5 | |
| 1833 | Thin-Glass Film | | | | | 15 | 17 | 18 | | 35 | 14 | | 4 |

Materials on Filter Media

These SRM's consist of potentially hazardous materials deposited on filters to be used to determine the levels of these materials in industrial atmospheres.

| SRM | Type | Unit Size | Material Certified | Quantity Certified ($\mu\text{g}/\text{filter}$) | | | |
|-------|---------------------------------------|-------------|----------------------|--|---------------|---------------|------------------|
| | | | | I | II | III | IV |
| 2676c | Metals on Filter Media | Set of 8 | Cadmium | 0.954 | 2.83 | 10.09 | (<0.01) |
| | | | Lead | 7.47 | 14.92 | 29.81 | (<0.01) |
| | | | Manganese | 2.11 | 9.92 | 19.85 | (<0.01) |
| | | | Zinc | 9.99 | 49.68 | 99.28 | (<0.01) |
| 2677 | Beryllium and Arsenic on Filter Media | 2 sets of 4 | Beryllium Arsenic | 0.052 0.103 | 0.256 1.07 | 1.03 10.5 | <0.001 <0.002 |
| 2679a | Quartz on Filter Media | Set of 4 | Quartz Clay | <2 (370) | 30.8 (370) | 80.2 (370) | 202.7 (370) |

Values in parentheses are not certified, but are given for information only.

Blank Filters

"These SRM's, except for SRM's 2678-81, consist of potentially hazardous materials deposited on filter media to be used to determine the levels of these materials in industrial atmospheres. SRM's 2678-81 are blank filters that have been analyzed at NIST and are to be used in the assessment of system blank and in the evaluation of detection limits."

| SRM | Type | Diameter | Pore Size | Nominal | |
|------|----------------------------|----------|-------------|-------------|------|
| | | | | Filter wt g | |
| 2678 | Cellulose Acetate Membrane | 47 mm | 0.45 micron | | 0.09 |
| 2680 | Cellulose Acetate Membrane | 37 mm | 0.80 micron | | 0.05 |
| 2681 | Ashless | 42.5 mm | — | | 0.14 |



Ray McKenzie, a project manager for OSRM, provides responsive, creative leadership for a variety of physical properties and engineering SRM's.

Respirable Quartz

These SRM's consists of quartz powders that are in the respirable size range. They are intended for use in determining the level of quartz in an industrial atmosphere by x-ray diffraction.

| SRM | Type | Constituent Certified | Amount |
|------|--------------|------------------------------------|--------|
| 1878 | Alpha Quartz | 95.5% Crystalline α -quartz | 5 g |
| 1879 | Cristobalite | 98.0% Crystalline Cristobalite | 5 g |

Asbestos

These SRM's consist of four 3×3 mm sections of a 0.4 mm pore size polycarbonate filter containing chrysotile fibers mixed with an urban dust or with grunerite fibers. They are intended for use in evaluating the techniques used to count and identify chrysotile and grunerite asbestos fibers in filter samples by transmission electron microscopy.

| SRM | Type | Fiber Loading |
|-------|-------------------------------------|--|
| 1876a | Chrysotile Asbestos | 37 fibers/0.01 mm ² |
| 8410 | Chrysotile Asbestos Research Filter | 7.9 fibers/0.01 mm ² |
| 8411 | Mixed asbestos Research Filter | 138 fibers/0.01 mm Chrysotile 43 fibers/0.01 mm Grunerite |

SRM 1866, Bulk Asbestos-Common, is a set of the three common bulk mine-grade asbestos materials and one synthetic glass fiber sample. There are 4-5 grams of each sample in the set of materials. The three asbestos types are Chrysotile, Grunerite (Amosite), and Riebeckite (Crocidolite). The optical properties of each of these materials observed by polarized light microscopy have been characterized so that these samples may serve as primary calibration standards for the identification of asbestos types in building material. SRM 1867, Bulk Asbestos-Uncommon, will contain three materials: Anthophylite, Actinolite, and Tremolite and will be used in the same manner and for the same purpose as SRM 1866.

| SRM | Type | Asbestos Type |
|------|----------------------------------|-------------------------------------|
| 1866 | Bulk Asbestos-Common | Chrysotile, Grunerite, Riebeckite |
| 1867 | Bulk Asbestos-Uncommon (IN PREP) | Anthophylite, Actinolite, Tremolite |



SRM 1866, Bulk Asbestos (a set of four materials), our newest asbestos SRM, joins our other available asbestos SRM's — Chrysotile Filter (SRM 1876a) and Research Filters (8410 and 8411).

Lubricating Materials

Metallo-Organic Compounds

These SRM's are for preparing solutions in oils of known and reproducible concentrations of metals. Certificates give directions for preparing a solution of known concentration in lubricating oil.

| SRM | Type | Constituent Certified | | |
|-------|---|-----------------------|---------------|-----------------|
| | | Element | (Wt. percent) | Wt/Unit (grams) |
| 1075a | Aluminum 2-ethylhexanoate | Aluminum | 8.07 | 5 |
| 1051b | Barium cyclohexanebutyrate | Barium | 28.7 | 5 |
| 1053a | Cadmium cyclohexanebutyrate | Cadmium | 24.8 | 5 |
| 1078b | Tris (1-phenyl-1,3-butanediono)chromium (III) | Chromium | 9.6 | 5 |
| 1080a | Bis(1-phenyl-1,3-butanediono)copper (II) | Copper | 16.37 | 5 |
| 1079b | Tris (1-phenyl-1,3-butanediono)iron (III) | Iron | 10.45 | 5 |
| 1059c | Lead cyclohexanebutyrate | Lead | 37.5 | 5 |
| 1060a | Lithium cyclohexanebutyrate | Lithium | 4.1 | 5 |
| 1061c | Magnesium cyclohexanebutyrate | Magnesium | 6.45 | 5 |
| 1065b | Nickel cyclohexanebutyrate | Nickel | 13.89 | 5 |
| 1071b | Triphenyl phosphate | Phosphorus | 9.48 | 5 |
| 1066a | Octaphenylcyclotetrasiloxane | Silicon | 14.14 | 5 |
| 1077a | Silver 2-ethylhexanoate | Silver | 42.60 | 5 |
| 1069b | Sodium cyclohexanebutyrate | Sodium | 12.0 | 5 |
| 1070a | Strontium cyclohexanebutyrate | Strontium | 20.7 | 5 |
| 1057b | Dibutyltin bis (2-ethylhexanoate) | Tin | 22.95 | 5 |
| 1052b | Bis(1-phenyl-1,3-butanediono)oxovanadium (IV) | Vanadium | 13.01 | 5 |
| 1073b | Zinc cyclohexanebutyrate | Zinc | 16.66 | 5 |

Lubricating Base Oils

Each of these SRM's consists of a series of five concentrations (5 bottles, 20 g each) of a single element in a base oil.

| SRM | Type | Element | Concentration ($\mu\text{g/g}$) | | | | |
|------|----------------------------------|---------|-----------------------------------|------|-----------|------|-------|
| | | | I | II | III | IV | V |
| 1818 | Chlorine in Lubricating Base Oil | Cl | 29 | 63 | 78 | 231 | 558 |
| 1819 | Sulfur in Lubricating Base Oil | S | 299 | 1070 | 2865 | 6030 | 10550 |
| 1836 | Nitrogen in Lubricating Base Oil | N | | | (IN PREP) | | |

Catalyst Package for Lubricant Oxidation

SRM 1817b is intended primarily for use in evaluating the oxidation stability of lubricating oils, i.e., automotive crankcase lubricants. The SRM contains: (1) an oxidized/nitrated fuel fraction, (2) a metal naphthenate mixture, and (3) distilled water. The metal naphthenate mixture has the following weight percentages of metal naphthenates: lead-82, iron-7, copper-4, manganese-3.5, and tin-3.5. SRM 1817b is available as a kit of 5 ampoules of each of the three components. The fuel and metal catalysts are sealed under inert atmosphere to ensure their stabilities.

Wear-Metals in Oil

| SRM | 1083 | 1085 |
|---|-------------------|-------------------------------|
| Type | Base Oil (ppm) | Wear-Metals in Oil 300 ppm |
| Unit Size | 150 mL | 85 mL |
| ELEMENT (Values in $\mu\text{g/g}$) | | |
| Aluminum | (<0.5) | 296 |
| Chromium | (<0.02) | 298 |
| Copper | (<0.5) | 295 |
| Iron | (<1) | 300 |
| Lead | (<0.04) | (305) |
| Magnesium | (<0.1) | 297 |
| Molybdenum | (<0.01) | 292 |
| Nickel | (<0.4) | 303 |
| Silicon | (<1) | (308) |
| Silver | (<0.05) | (291) |
| Sulfur | (980) | (4806) |
| Tin | (<0.4) | 296 |
| Titanium | (<5) | 300 |

Values in parentheses are not certified, but are given for information only.

Catalyst Characterization Materials


| SRM | Type | Wt/Unit (grams) |
|--------|--|----------------------|
| RM8580 | Nickel on Fresh Alumina (NiO 2.92 Wt %) | 25 gm |
| RM8581 | Nickel on Fresh Alumina (NiO 14.35 Wt %) | 25 gm |
| RM8582 | Nickel on Fresh Alumina (NiO 54.97 Wt %) | 25 gm |
| RM8583 | Platinum on Fresh Alumina (Pt 0.5760 Wt %) | 50 gm |
| RM8584 | Platinum on Fresh Alumina (Pt 0.540 Wt %) | 50 gm |
| RM8586 | FCC Catalyst | 25 gm |
| RM8587 | FCC Catalyst | 10 gm |
| RM8588 | Faujasite Y | 25 gm |
| RM8589 | FCC Catalysts | Set of 6, 50 gm each |

Fertilizers

These SRM's are intended for use in the fertilizer industry as working standards for the determination of the certified constituents.

| SRM | Type | Wt/Unit (grams) | Composition (Nominal Weight Percent) | | | | | |
|------|--------------------------------|--------------------|--------------------------------------|-------|-------|-------------------------------|------------------|-------|
| | | | N | P | K | P ₂ O ₅ | K ₂ O | CaO |
| 193 | Potassium Nitrate | 90 | 13.85 | | 38.66 | | | |
| 194 | Ammonium Dihydrogen Phosphate | 90 | 12.15 | 26.92 | | | | |
| 200 | Potassium Dihydrogen Phosphate | 90 | | 22.74 | 28.76 | | | |
| 120c | Phosphate Rock (Florida) | 90 | | | | 33.34 | 0.110 | 48.02 |
| 694 | Phosphate Rock (Western) | 90 | | | | 30.2 | 0.51 | 43.6 |

| SRM | Composition (Nominal Weight Percent) | | | | | | | | | | | |
|-----|--------------------------------------|-----|--------------------------------|--------------------------------|------|-------------------|--------|------------------|--------------------------------|-------|---------|-------------------------------|
| | SiO ₂ | F | Fe ₂ O ₃ | Al ₂ O ₃ | MgO | Na ₂ O | MnO | TiO ₂ | Cr ₂ O ₃ | CdO | U | V ₂ O ₅ |
| 694 | 11.2 | 3.2 | 0.79 | 1.8 | 0.33 | 0.86 | 0.0116 | (0.11) | (0.10) | 0.015 | 0.01414 | 0.31 |



SRM's 1845, Cholesterol in Egg Powder and 1566a, Oyster Tissue are two of the important and popular food related SRM's now available.


Ores

| SRM | 79a | 180 | 181 | 182 | 183 | | | | |
|---------------------|-----------------------------|--------------------------|--------------------------------------|---------------------------|-----------------------------|--------------------|------------------|----------|-----------|
| Type | Fluorspar, Customs Grade | Fluorspar, High Grade | Lithium Ore (Spodumene) | Lithium Ore (Petalite) | Lithium Ore (Lepidolite) | | | | |
| Unit Weight | 120g | 120 g | 45 g | 45 g | 45 g | | | | |
| Constituents | | | | | | | | | |
| CaF ₂ | 97.39% | 98.80% | 6.3 ₉ % | 4.3 ₄ % | 4.1 ₂ % | | | | |
| Li ₂ O | | | | | | | | | |
| SRM | 330 | | 331 | | | | | | |
| Type | Copper, Ore Mill Heads | | Copper, Ore Mill Tails | | | | | | |
| Unit Weight | 100 g | | 100 g | | | | | | |
| Constituents | | | | | | | | | |
| Cu | 0.84% | | | 0.091% | | | | | |
| Re | 0.30 ppm | | | 0.04 ppm | | | | | |
| Mo | 0.018% | | | 0.0022% | | | | | |
| Au | (0.093 ppm) | | | (0.034 ppm) | | | | | |
| Ag | (1.51 ppm) | | | (0.243 ppm) | | | | | |
| SRM | Type | Wt/ Units (grams) | Constituent (Nominal Weight Percent) | | | | | | |
| 277 | Tungsten Concentrate | 100 | WO ₃ | Ca | Fe | Pb | Mn | | |
| 2430 | | 100 | 70.26 | (0.37) As 0.002 | (7.4) (1.0) | (0.07) Bi 0.078 | (10.0) (0.12) | | |
| SRM | Mo | Nb | O ₂ | P | Si | S | Ta | Sn | Ti |
| 277 | (0.06) | (1.00) | (21.4) | (0.03) | (0.85) | (0.25) | (0.20) | (0.54) | (2.2) |
| 2430 | 0.22 | Cu (0.01) | Al (0.4) | 0.017 | Mg (0.5) | 0.26 | (<0.01) | K (0.16) | Na (0.02) |

Values in parentheses are not certified, but are given for information only.

Ores (Continued)

| SRM | 27f | 690 | 691 | 692 | 693 |
|--|---------------------|---------------------|------------------------|-----------------------|--------------------|
| Type | Iron Ore, Sibley | Iron Ore, Canada | Iron Oxide, Reduced | Iron Ore, Labrador | Iron Ore, Nimba |
| Unit Weight | 100 g | 150 g | 100 g | 150 g | 150 g |
| Constituents (Nominal Weight Percent) | | | | | |
| Al ₂ O ₃ | 0.82 | 0.18 | 1.22 | 1.41 | 1.02 |
| CaO | 0.039 | 0.20 | 0.63 | 0.023 | 0.016 |
| Co | | | 0.030 | | |
| Cu | | | 0.032 | | |
| Total Fe | 65.97 | 66.85 | 90.8 | 59.58 | 65.11 |
| MgO | 0.019 | 0.18 | 0.52 | 0.035 | 0.013 |
| MnO | 0.011 | 0.23 | 0.043 | 0.46 | 0.091 |
| P | 0.041 | 0.011 | 0.006 | 0.039 | 0.056 |
| K ₂ O | 0.008 | 0.0030 | | 0.039 | 0.0028 |
| SiO ₂ | 4.17 | 3.71 | 3.7 | 10.14 | 3.87 |
| Na ₂ O | 0.012 | 0.003 | 0.186 | 0.008 | 0.0028 |
| S | 0.005 | 0.003 | 0.008 | 0.005 | 0.005 |
| TiO ₂ | 0.019 | 0.022 | 0.27 | 0.045 | 0.035 |



Pictured here are four of the more than 70 spectrometric standard solutions now available either as single element or multielement solutions.

Ores (Continued)

| SRM | 69b | 696 | 697 | 698 | 699 | 120c | 694 | 25d | 670 | 600 | 1835 |
|--|-------------------|------------------|--------------------|-------------------|---------------------------|-------------------------|-------------------------|---------------|------------|---------------------|------------|
| Type | Bauxite, Arkansas | Bauxite, Surinam | Bauxite, Dominican | Bauxite, Jamaican | Alumina (Reduction Grade) | Phosphate Rock, Florida | Phosphate Rock, Western | Manganese Ore | Rutile Ore | Bauxite, Australian | Borate Ore |
| Unit Weight | 60 g | 60 g | 60 g | 60 g | 60 g | 90 g | 90 g | 100 g | 90 g | 90 g | 60 g |
| Constituents (Nominal Weight Percent) | | | | | | | | | | | |
| Al ₂ O ₃ | 48.8 | 54.5 | 45.8 | 48.2 | | 1.30 | 1.8 | 5.32 | | 40.0 | 3.474 |
| BaO | (0.008) | (0.004) | (0.015) | (0.008) | | | | (0.21) | | | 0.0497 |
| CdO | | | | | | 0.0010 | 0.015 | | | | |
| CaO | 0.13 | 0.018 | 0.71 | 0.62 | 0.036 | 48.02 | 43.6 | | | 0.22 | 21.622 |
| Co | (0.0001) | (0.00009) | (0.0013) | (0.0045) | | F 3.2 | | (0.052) | | | |
| Cr ₂ O ₃ | 0.011 | 0.047 | 0.100 | 0.080 | 0.0002 | | (0.10) | | 0.23 | 0.024 | |
| Fe ₂ O ₃ | 7.14 | 8.70 | 20.0 | 19.6 | 0.013 | 1.02 | 0.79 | | 0.86 | 17.0 | 1.141 |
| MgO | 0.085 | 0.012 | 0.18 | 0.058 | 0.0006 | 0.32 | 0.33 | | | 0.05 | 3.411 |
| MnO | 0.110 | 0.004 | 0.41 | 0.38 | 0.0005 | 0.027 | 0.0116 | Mn 51.78 | | 0.013 | 0.0333 |
| P ₂ O ₅ | 0.118 | 0.050 | 0.97 | 0.37 | 0.0002 | 33.34 | 30.2 | 0.25 | | 0.039 | |
| K ₂ O | 0.068 | 0.009 | 0.062 | 0.010 | | 0.110 | 0.51 | 0.93 | 0.23 | 1.261 | |
| SiO ₂ | 13.43 | 3.79 | 6.81 | 0.69 | 0.014 | 5.5 | 11.2 | 2.52 | 0.51 | 20.3 | 18.408 |
| Na ₂ O | (0.025) | (0.007) | (0.036) | (0.015) | 0.59 | 0.52 | 0.86 | | | 0.022 | 3.484 |
| SO ₃ | 0.63 | 0.21 | 0.13 | 0.22 | | | | | | 0.19 | 1.477 |
| TiO ₂ | 1.90 | 2.64 | 2.52 | 2.38 | | 0.103 | (0.11) | 0.13 | 96.16 | 1.31 | 0.1332 |
| U | | | | | | | 0.01414 | | | | |
| V ₂ O ₅ | 0.028 | 0.072 | 0.063 | 0.064 | 0.0005 | | 0.31 | | 0.66 | 0.060 | |
| ZnO | 0.0035 | 0.0014 | 0.037 | 0.029 | 0.013 | | (0.19) | | | | 0.003 |
| ZrO ₂ | 0.29 | 0.14 | 0.065 | 0.061 | | 0.010 | | | 0.84 | 0.060 | |
| Ga ₂ O ₃ | | | | | | | | | | | |
| Li ₂ O | | | | | 0.002 | | | | 14.28 | | |
| Available Oxygen | | | | | | | | | | | |
| Moisture | | | | | | | | (0.96) | | | |
| Loss on Ignition | 27.2 | 29.9 | 22.1 | 27.3 | | | | | | | |

Values in parentheses are not certified, but are given for information only.

Clays

| SRM | 97b | 98b | 679 |
|------------------|---------------|---------------|--------------|
| Type | Flint Clay | Plastic Clay | Brick Clay |
| Unit Weight | 60 g | 60 g | 75 g |
| Element | | | |
| Al | 20.76 wt. % | 14.30 wt. % | 11.01 wt. % |
| Ba | (0.018) wt. % | (0.07) wt. % | |
| Ca | 0.0249 wt. % | 0.0759 wt. % | 0.1628 wt. % |
| Ce | | | (105) |
| Cs | (3.4) | (16.5) | (9.6) |
| Cr | 227 | 119 | 109.7 |
| Co | (3.8) | (16.3) | (26) |
| Eu | (0.84) | (1.3) | (1.9) |
| Hf | (13) | (7.2) | (4.6) |
| Fe | 0.831 wt. % | 1.18 wt. % | 9.05 wt. % |
| Li | 550 | 215 | 71.7 |
| Mg | 0.113 wt. % | 0.358 wt. % | 0.7552 wt. % |
| Mn | 47 | 116 | (1730) |
| P | (200) | (300) | (750) |
| K | 0.513 wt. % | 2.81 wt. % | 2.433 wt. % |
| Rb | (33) | (180) | (190) |
| Sc | (22) | (22) | (22.5) |
| Si | 19.81 wt. % | 26.65 wt. % | 24.34 wt. % |
| Na | 0.0492 wt. % | 0.1496 wt. % | 0.1304 wt. % |
| Sr | 84 | 189 | 73.4 |
| Th | (36) | (21) | (14) |
| Ti | 1.43 wt. % | 0.809 wt. % | 0.577 wt. % |
| Zn | (87) | (110) | (150) |
| Zr | (0.05) wt. % | (0.022) wt. % | |
| Sb | (2.2) | (1.6) | |
| Loss on Ignition | (13.3 wt. %) | (7.5 wt. %) | |

Values in parentheses are not certified, but are given for information only.

Karen Applestein provides the alertness and attention to detail that ensures that each customers order or service request is handled punctually and accurately.



Rocks, Minerals, and Refractories

| SRM | 1c | 88b | 70a | 99a | 81a | 165a | 1413 | 2709 | 2710 |
|--|---------------------------|-----------------------|-------------------|-------------------|---------------|-----------------------|---------------------------|-------------------------------------|-------------------------------------|
| Type | Lime-stone, argilla-ceous | Limestone, dolomitic | Feld-spar, potash | Feld-spar, soda | Glass sand | Glass sand (low iron) | Glass sand (high alumina) | Agricultural Soil (Baseline Levels) | Agricultural Soil (Elevated Levels) |
| Unit Weight | 50 g | 75 g | 40 g | 40 g | 75 g | 75 g | 75 g | IN PREP | IN PREP |
| Constituents (Nominal Weight Percent) | | | | | | | | | |
| Al ₂ O ₃ | 1.30 | 0.336 | 17.9 | 20.5 | 0.66 | 0.059 | 9.90 | | |
| BaO | | CO ₂ 46.37 | 0.02 | 0.26 | | | 0.12 | | |
| CaO | 50.3 | 30.12 | 0.11 | 2.14 | | | 0.74 | | |
| Cr ₂ O ₃ | | | | | 46 µg/g | (1) µg/g | | | |
| Fe ₂ O ₃ | 0.55 | 0.277 | 0.07 _s | 0.06 _s | 0.082 | 0.012 | 0.24 | | |
| MgO | 0.42 | 21.03 | | 0.02 | | | 0.06 | | |
| MnO | 0.025 | 0.0160 | | | | | | | |
| P ₂ O ₅ | 0.04 | 0.0044 | | 0.02 | | | | | |
| K ₂ O | 0.28 | 0.1030 | 11.8 | 5.2 | | | 3.94 | | |
| Rb ₂ O | | | 0.06 | | | | | | |
| SiO ₂ | 6.84 | 1.13 | 67.1 | 65.2 | | | 82.77 | | |
| Na ₂ O | 0.02 | 0.0290 | 2.5 _s | 6.2 | | | 1.75 | | |
| SrO | 0.030 | 0.0076 | | | | | | | |
| TiO ₂ | 0.07 | (0.016) | 0.01 | 0.007 | 0.12 0.034 | 0.011 0.006 | 0.11 | | |
| ZrO ₂ | | | | | | | | | |
| Loss on Ignition | 39.9 | (46.98) | 0.40 | 0.26 | | | | | |

Values in parentheses are not certified, but are given for information only.

Mercury in Soil

| SRM | RM 8406 | RM 8407 | RM 8408 |
|-------------|-----------------|-----------------|-----------------|
| Type | Mercury in Soil | Mercury in Soil | Mercury in Soil |
| Unit Weight | 110 µg/g | 50 µg/g | 0.07 µg/g |

Rocks, Minerals, and Refractories (Continued)

| SRM | 154b | 278 | 688 | 76a | 77a | 78a |
|--|------------------|---------------|-------------|--|--|--|
| Type | Titanium Dioxide | Obsidian Rock | Basalt Rock | Burnt Refractory (Al ₂ O ₃ -40%) | Burnt Refractory (Al ₂ O ₃ -60%) | Burnt Refractory (Al ₂ O ₃ -70%) |
| Unit Weight | 90 g | 35 g | 60 g | 75 g | 75 g | 75 g |
| Constituents (Nominal Weight Percent) | | | | | | |
| Al ₂ O ₃ | | 14.15 | 17.36 | 38.7 | 60.2 | 71.7 |
| CaO | (~0.01) | 0.983 | (12.17) | 0.22 | 0.05 | 0.11 |
| Cr | | | 332 µg/g | | | |
| Cu | | 5.9 µg/g | | | | |
| FeO | | 1.36 | 7.64 | | | |
| Fe ₂ O ₃ | (0.006) | 2.04 | 10.35 | 1.6 _o | 1.0 _o | 1.2 |
| Pb | | 16.4 µg/g | 3.3 µg/g | | | |
| Li ₂ O | | | | 0.042 | 0.02 ₅ | 0.12 |
| MgO | (~0.01) | (0.23) | (8.4) | 0.52 | 0.38 | 0.70 |
| MnO | | 0.052 | 0.167 | | | |
| Ni | | 3.6 µg/g | | | | |
| P ₂ O ₅ | (0.04) | 0.036 | 0.134 | 0.12 _o | 0.092 | 1.3 |
| K ₂ O | | 4.16 | 0.187 | 1.33 | 0.09 _o | 1.22 |
| Rb | | 127.5 µg/g | 1.91 µg/g | | | |
| SiO ₂ | (0.01) | 73.05 | 48.4 | 54.9 | 35.0 | 19.4 |
| Na ₂ O | | 4.84 | 2.15 | 0.07 | 0.037 | 0.078 |
| Sr | | 63.5 µg/g | 169.2 µg/g | | | |
| SrO | | | | 0.037 | 0.009 | 0.25 |
| Th | | 12.4 µg/g | 0.33 µg/g | | | |
| TiO ₂ | 99.74 | 0.245 | 1.17 | 2.0 _o | 2.6 _o | 3.2 _o |
| Tl | | 0.54 µg/g | | | | |
| U | | 4.58 µg/g | | | | |
| Loss on Ignition | | | | (0.34) | (0.22) | (0.42) |

Values in parentheses are not certified, but are given for information only.

Rocks, Minerals, and Refractories (Continued)

| SRM | 103a | 198 | 199 |
|--|-------------------|--------------|--------------|
| Type | Chrome Refractory | Silica Brick | Silica Brick |
| Unit Weight | 60 g | 45 g | 45 g |
| Constituents (Nominal Weight Percent) | | | |
| Al ₂ O ₃ | 29.96 | 0.16 | 0.48 |
| CaO | 0.69 | 2.71 | 2.41 |
| Cr ₂ O ₃ | 32.06 | | |
| FeO | 12.43 | | |
| Fe ₂ O ₃ | | 0.66 | 0.74 |
| Li ₂ O | | 0.001 | 0.002 |
| MgO | 18.54 | 0.07 | 0.13 |
| MnO | 0.11 | 0.008 | 0.007 |
| P ₂ O ₅ | 0.01 | 0.022 | 0.015 |
| K ₂ O | | 0.017 | 0.094 |
| SiO ₂ | 4.63 | | |
| Na ₂ O | | 0.012 | 0.015 |
| TiO ₂ | 0.22 | 0.02 | 0.06 |
| ZrO ₂ | 0.01 | | |
| Loss on Ignition | | 0.21 | 0.17 |

Values in parentheses are not certified, but are given for information only.

Carbides

| SRM | Type | Wt/Unit (grams) | Chemical Composition (Nominal Weight Percent) | | | | | |
|------|------------------|--------------------|---|---------|--------|------|----------------|----------------|
| | | | SiC | Total C | Free C | Fe | O ₂ | N ₂ |
| 112b | Silicon Carbide | 80 | 97.37 | 29.43 | 0.26 | 0.13 | | |
| 276a | Tungsten Carbide | 75 | | 6.11 | (0.02) | | (0.03) | (0.003) |


Values in parentheses are not certified, but are given for information only.

Cemented Carbides

SRM's 887-9 are fine powders prepared from sintered Tungsten Carbide bare materials.

| SRM | 887 | 888 | 889 |
|-----------------------------|--------------------------------|-------------------------------------|--|
| Type | Cemented Carbide (W83-Co10) | Cemented Carbide (W64- Co25-Ta5) | Cemented Carbide (W75- Co9-Ta5-Ti4) |
| Unit Weight | 100 g | 100 g | 100 g |
| Chemical Composition | | | |
| Cobalt | 10.35 | 24.7 | 9.50 |
| Tantalum | | 4.77 | 4.60 |
| Titanium | | | 4.03 |
| Carbon | (5.5) | (4.6) | (6.0) |
| Tungsten | Bal | Bal | Bal |

Values in parentheses are not certified, but are given for information only.



The portland and calcium aluminate SRM's pictured here are recent additions to very popular "630" portland cement composition SRM's. A new series—SRM's 1884-1889—are now available to replace the "630" series as those supplies are exhausted.

Glasses

| SRM | 89 | 91 | 92 | 93a | 620 | 621 | 1411 | 1412 | 1830 | 1831 | 1834 |
|---|-------------|-------------|----------------------------|--------------------------|---------------------------|----------------------------|-------------------|-----------------|---------------------------|---------------------------|-----------------|
| Type | Lead-Barium | Opal Powder | Low-Boron Soda-Lime Powder | High-Boron Borosilicate | Soda-Lime, Flat | Soda-Lime, Container | Soft Borosilicate | Multi Component | Soda-Lime, Float | Soda-Lime, Sheet | Fused Ore Glass |
| Unit Size | | | | | | | | | | | |
| | 45 g | 45 g | 45 g | Wafer 32 mm D×6 mm | 3 platelets 35×35×3 mm | 3 disks 38 mm D×5 mm | 10 platelets | 8 platelets | 3 platelets 38×38×6 mm | 3 platelets 37×37×3 mm | IN PREP |
| Constituent (Nominal Weight Percent) | | | | | | | | | | | |
| SiO ₂ | 65.35 | 67.5 | (75.0) | 80.8 | 72.08 | 71.13 | 58.04 | 42.38 | 73.07 | 73.08 | |
| PbO | 17.50 | 0.10 | | | | | | 4.40 | | | |
| Al ₂ O ₃ | 0.18 | 6.01 | | 2.2 _s | 1.80 | 2.76 | 5.68 | 7.52 | 0.12 | 1.21 | |
| Fe ₂ O ₃ | 0.049 | 0.079 | | 0.028 | 0.043 | 0.040 | 0.050 | (0.031) | 0.121 | 0.087 | |
| ZnO | | 0.08 | (0.2) | | | | 3.85 | 4.48 | | | |
| CdO | | | | | | | | 4.38 | | | |
| MnO | 0.088 | (0.008) | | | | | | | | | |
| TiO ₂ | 0.01 | 0.019 | | 0.01 _s | 0.018 | 0.014 | 0.02 | | 0.011 | 0.019 | |
| ZrO ₂ | 0.005 | 0.009 | | 0.04 _s | | 0.007 | | | | | |
| CaO | 0.21 | 10.49 | (8.3) | 0.01 | 7.11 | 10.71 | 2.18 | 4.53 | 8.56 | 8.20 | |
| BaO | 1.40 | | | | | 0.12 | 5.00 | 4.67 | | | |
| Li ₂ O | | | | | | | | (4.50) | | | |
| MgO | 0.03 | (0.008) | (0.1) | 0.00 _s | 3.69 | 0.27 | 0.33 | (4.69) | 3.90 | 3.51 | |
| K ₂ O | 8.40 | 3.24 | (0.6) | 0.01 _s | 0.41 | 2.01 | 2.97 | 4.14 | 0.04 | 0.33 | |
| Na ₂ O | 5.70 | 8.47 | (13.1) | 3.9 _s | 14.39 | 12.74 | 10.14 | 4.69 | 13.75 | 13.32 | |
| B ₂ O ₃ | | | 0.70 | 12.5 _s | | | 10.94 | 4.53 | | | |
| P ₂ O ₅ | 0.23 | 0.023 | | | | | | | | | |
| As ₂ O ₅ | 0.36 | 0.10 | | | | | | | | | |
| As ₂ O ₃ | 0.03 | 0.09 | | | 0.056 | 0.030 | | | | | |
| SO ₃ | 0.03 | | | | 0.28 | 0.13 | | | 0.26 | 0.25 | |
| Cl | 0.05 | 0.015 | | 0.06 _s | | | 0.09 | 4.55 | | | |
| SrO | | | | | | | | | | | |
| F | | 5.73 | | | | | | | | | |
| Loss on Ignition | 0.32 | | (0.42) | | | | | | | | |

Values in parentheses are not certified, but are given for information only.

Cements

These SRM's are for x-ray spectroscopic and chemical analysis of portland cements and related materials. Each unit consists of three sealed vials each containing approximately 5 g of material.

| SRM | 1880 | 1881 | 1882 | 1883 | 1884 | 1885 | 1886 |
|---|-------------------|-------------------|--------|--------|----------|-----------|-----------|
| Type | BLACK | WHITE | ORANGE | SILVER | IVORY | TURQUOISE | CRANBERRY |
| Unit Weight | 15 g | 15 g | 15 g | 15 g | 15 g | 15 g | 15 g |
| Constituent (Nominal Weight Percent) | | | | | | | |
| CaO | 63.1 ₄ | 58.6 ₈ | 37.6 | 27.8 | 64.01 | 62.14 | 67.43 |
| SiO ₂ | 19.8 ₂ | 22.2 ₅ | 3.40 | 0.35 | 23.19 | 21.24 | 22.53 |
| Al ₂ O ₃ | 5.0 ₃ | 4.1 ₉ | 38.6 | 71.2 | 3.31 | 3.68 | 3.99 |
| Fe ₂ O ₃ | 2.91 | 4.68 | 15.8 | 0.08 | 3.30 | 4.40 | 0.31 |
| SO ₃ | 3.37 | 3.6 ₅ | | | 1.67 | 2.22 | 2.04 |
| MgO | 2.6 ₉ | 2.62 | 1.25 | 0.29 | 2.32 | 4.02 | 1.60 |
| K ₂ O | 0.91 | 1.17 | 0.12 | (0.01) | 0.51 | 0.83 | 0.16 |
| TiO ₂ | 0.23 | 0.2 ₃ | 1.83 | (0.01) | 0.16 | 0.20 | 0.19 |
| Na ₂ O | 0.28 | 0.04 | (0.06) | 0.32 | 0.13 | 0.38 | 0.02 |
| SrO | 0.06 | 0.11 | | | 0.048 | 0.037 | 0.11 |
| P ₂ O ₅ | 0.29 | 0.09 | | | 0.12 | 0.10 | 0.025 |
| Mn ₂ O ₃ | 0.08 | 0.26 | | | 0.11 | 0.12 | 0.013 |
| F | 0.10 | 0.09 | | | (0.03) | (0.05) | (0.01) |
| ZnO | 0.01 | 0.01 | | | (0.02) | (0.03) | (<0.01) |
| Cr ₂ O ₃ | Cl 0.02 | | | | (<0.01) | (<0.04) | (<0.01) |
| Cl | | | | | (0) | (0.02) | (0) |
| L.O.I. (Loss-on-Ignition) | | | | | 1.17 | 0.74 | 1.73 |
| Total | | | | | (100.05) | (100.19) | (100.02) |

Values in parentheses are not certified, but are given for information only.

Cements (continued)

| SRM | 1887 | 1888 | 1889 | 8486 | 8487 | 8488 |
|---|----------|----------|----------|-------------------------|-------------------------|-------------------------|
| Type | BROWN | PURPLE | GRAY | Portland Cement Clinker | Portland Cement Clinker | Portland Cement Clinker |
| Unit Weight | 15 g | 15 g | 15 g | 3×10g | 3×10g | 3×10g |
| Constituent (Nominal Weight Percent) | | | | | | |
| CaO | 62.88 | 63.78 | 65.08 | (63.36) | (67.20) | (66.50) |
| SiO ₂ | 19.98 | 20.86 | 20.44 | (22.48) | (21.43) | (22.68) |
| Al ₂ O ₃ | 5.59 | 5.35 | 5.61 | (4.70) | (5.53) | (4.90) |
| Fe ₂ O ₃ | 2.16 | 3.18 | 2.67 | (3.60) | (1.98) | (4.07) |
| SO ₃ | 4.61 | 3.16 | 2.68 | (0.27) | (0.88) | (0.31) |
| MgO | 1.26 | 0.71 | 1.38 | (4.73) | (1.48) | (0.98) |
| K ₂ O | 1.27 | 0.56 | 0.32 | (0.42) | (0.72) | (0.35) |
| TiO ₂ | 0.27 | 0.29 | 0.21 | (0.25) | (0.27) | (0.24) |
| Na ₂ O | 0.10 | 0.14 | 0.11 | (0.10) | (0.14) | (0.11) |
| SrO | 0.07 | 0.07 | 0.20 | (0.05) | (0.11) | (0.13) |
| P ₂ O ₅ | 0.075 | 0.085 | 0.15 | (0.06) | (0.29) | (0.08) |
| Mn ₂ O ₃ | 0.072 | 0.025 | 0.24 | (0.10) | (0.04) | (0.03) |
| F | (0.11) | (0.02) | (0.04) | | | |
| ZnO | (0.01) | (0.01) | (<0.01) | | | |
| Cr ₂ O ₃ | (<0.01) | (0.01) | (0.01) | | | |
| Cl | (0.007) | (0.015) | (0.002) | | | |
| L.O.I. (Loss on-Ignition) | 1.49 | 1.79 | 0.92 | (0.16) | (0.17) | (0.21) |
| Total | (99.908) | (100.05) | (100.04) | (100.28) | (100.20) | (100.60) |

Values in parentheses are not certified, but are given for information only.

Portland Cement Clinkers

These RM's are intended primarily for use in the determination of the abundance of major phases in cement clinkers, i.e., the percentages of alite (C₃S)*, belite (C₂S)*, aluminate (C₃A)*, and ferrite ((C₂(A,F))*.

[*Note: Cement chemist's notation: C=CaO, S=SiO₂, A=Al₂O₃, F=Fe₂O₃.]

| RM | Concentration (wt. %) | | |
|----------------|-----------------------|----------|----------|
| | Phase | 8486 | 8487 |
| Alite | (58.47) | (73.39) | (64.97) |
| Belite | (23.18) | (7.75) | (18.51) |
| Aluminate | (1.15) | (12.09) | (4.34) |
| Ferrite | (13.68) | (3.27) | (12.12) |
| Free CaO | (0.18) | (2.45) | (0.00) |
| Periclase | (3.21) | (0.09) | (0.05) |
| Alkali Sulfate | (0.14) | (0.98) | (0.03) |
| Total | (100.01) | (100.02) | (100.02) |

Values in parentheses are not certified, but are given for information only.

Trace Elements

The SRM's are for trace chemical analysis, specifically for calibrating instruments and evaluating analytical techniques used to determine trace elements in inorganic matrices.

| SRM | 607 | 610–611 | 612–613 | 614–615 | 616–617 |
|----------------------|--------------------------------------|-------------------------------------|-------------------------|-------------------------|-------------------------|
| Type | Trace Elements in Potassium Feldspar | Trace Elements in Glass | Trace Elements in Glass | Trace Elements in Glass | Trace Elements in Glass |
| Concentration | ($\mu\text{g/g}$) | (500 ppm) | (50 ppm) | (1 ppm) | (0.02 ppm) |
| Wafer Thickness | | 610 3 mm 611 1 mm | 612 3 mm 613 1 mm | 614 3 mm 615 1 mm | 616 3 mm 617 1 mm |
| Unit of Issue | 5 g | 6 Wafers | 6 Wafers | 6 Wafers | 6 Wafers |
| Element | | Nominal Concentrations (ppm) | | | |
| Antimony | | | | (1.06) | (0.078) |
| Barium | | | (41) | | |
| Boron | | (351) | (32) | (1.30) | (0.20) |
| Cadmium | | | | (0.55) | |
| Cerium | | | (39) | | |
| Cobalt | | (390) | (35.5) | (0.73) | |
| Copper | | (444) | (37.7) | 1.37 | (0.80) |
| Dysprosium | | | (35) | | |
| Erbium | | | (39) | | |
| Europium | | | (36) | (0.99) | |
| Gadolinium | | | (39) | | |
| Gallium | | | | (1.3) | (0.23) |
| Gold | | (25) | (5) | (0.5) | (0.18) |
| Iron | | 458 | 51 | (13.3) | (11) |
| Lanthanum | | | (36) | (0.83) | (0.034) |
| Lead | | 426 | 38.57 | 2.32 | 1.85 |
| Manganese | | 485 | (39.6) | | |
| Neodymium | | | (36) | | |
| Nickel | | 458.7 | 38.8 | (0.95) | |
| Potassium | | (461) | (64) | 30 | 29 |
| Rubidium | 523.90 | 425.7 | 31.4 | 0.855 | (0.100) |
| Samarium | | | (39) | | |
| Scandium | | | | (0.59) | (0.026) |
| Silver | | (254) | 22.0 | 0.42 | |
| Strontium | 65.485 | 515.5 | 78.4 | 45.8 | 41.72 |
| Thallium | | (61.8) | (15.7) | (0.269) | (0.0082) |
| Thorium | | 457.2 | 37.79 | 0.748 | 0.0252 |
| Titanium | | (437) | (50.1) | (3.1) | (2.5) |
| Uranium | | 461.5 | 37.38 | 0.823 | 0.0721 |
| Ytterbium | | | (42) | | |
| Zinc | | (433) | | | |

In addition to the elements listed above, the glass SRM's contain the following 25 elements: As, Be, Bi, Cs, Cl, F, Ge, Hf, Hg, Li, Lu, Mg, Nb, P, Pr, Se, S, Te, Tb, Tm, Sn, W, V, Y, and Zr.

NOTE: Glass—Nominal Composition: 72% SiO₂, 12% CaO, 14% Na₂O, and 2% Al₂O₃.

Values in parentheses are not certified, but are given for information only.

Nuclear Radiation Monitoring

Radiation Dosimetry

This SRM is a cobalt-in-aluminum alloy wire 0.5 mm in diameter and 1 meter long for use as a neutron density monitor standard.

| SRM | Identification (Batch Name) | Cobalt Content (Weight Percent) |
|-----|---|---------------------------------|
| 953 | Neutron density monitor wire (Co in Al) | 0.116 |

Fission Track Glass

This SRM contains uranium at three concentration levels, and will aid laboratories performing fission track analyses in interlaboratory comparisons of data and in monitoring neutron fluences. The material was irradiated in the NIST 10 Megawatt Research Reactor, at two different neutron energies.

The SRM unit contains four unirradiated glass wafers and two irradiated wafers.

| SRM | Uranium Content (ng/g) | U (Atom Percent) | Reactor Position | Cu Foil | Au Foil |
|------|------------------------|------------------|------------------|--------------|--------------|
| 963a | 0.823 | 0.2792 | RT-3: RT-4: | 41.2 39.5 | 45.8 43.0 |

Stable Isotopic Materials

The isotopic composition of these SRM's has been determined by mass spectrometry.

| SRM | Isotopic Reference Standards | Element Certified | Wt/Unit (grams) |
|------|--|-------------------|-----------------------------|
| 951 | Boric Acid | Boron | 100 |
| 952 | Boric Acid, 95% Enriched ^{10}B | Boron | 0.25 |
| 975 | Sodium Chloride | Chlorine | 0.25 |
| 976 | Copper Metal | Copper | 0.25 |
| 977 | Sodium Bromide | Bromine | 0.25 |
| 978a | Silver Nitrate | Silver | 0.25 |
| 979 | Chromium Nitrate | Chromium | 0.25 |
| 980 | Magnesium Metal | Magnesium | 0.25 |
| *981 | Lead Metal, Natural | Lead | 1.0 |
| *982 | Lead Metal, Equal Atom (206/208) | Lead | 1.0 |
| *983 | Lead Metal, Radiogenic (92%-206) | Lead | 1.0 |
| 984 | Rubidium Chloride, assay and isotopic | Rubidium | 0.25 |
| 985 | Potassium Chloride, assay and isotopic | Potassium | 1.0 |
| 986 | Nickel | Nickel | 1.0 |
| 987 | Strontium Carbonate, assay and isotopic | Strontium | 1.0 |
| 989 | Rhenium, assay and isotopic | Rhenium | pkg. (50) |
| 990 | Silicon, assay and isotopic | Silicon | wafer, 3 cm \times 0.2 cm |
| 991 | Lead-206 Spike, assay and isotopic | Lead | 15 |
| 994 | Gallium Metal, isotopic | Gallium | 0.25 |
| 997 | Thallium Metal, isotopic | Thallium | 0.25 |

*Sold as a set containing SRM 981, 982, and 983.



John Moody, a research chemist in the Inorganic Analytical Research Division, prepares some botanical materials in the clean room for trace analysis.

Physical Properties

ION ACTIVITY

pH

These SRM's are used to prepare solutions of known hydrogen ion concentration to calibrate commercial pH instruments. SRM's 186Id and 186IId, 191a and 192a, and 922 and 923 are certified for use as admixtures only. SRM's 186Id and 186IId may be used to prepare a solution with a pH of 6.863 at 25 °C, or a physiological buffer solution with a pH of 7.41 at 25 °C.

| SRM | Type | pH(S) Values (at 25 °C) | Wt/Unit (grams) |
|--------|---|----------------------------|--------------------|
| 185f | Potassium hydrogen phthalate | 4.006 | 60 |
| 186Id | Potassium dihydrogen phosphate | (see above) | 30 |
| 186IId | Disodium hydrogen phosphate | | 30 |
| 187c | Sodium tetraborate decahydrate (Borax) | 9.180 | 30 |
| 188 | Potassium hydrogen tartrate | 3.557 | 60 |
| 189a | Potassium tetroxalate | 1.681 | 65 |
| 191a | Sodium bicarbonate | 10.011 | 25 |
| 192a | Sodium carbonate | | 30 |
| 922 | Tris(hydroxymethyl)aminomethane | 7.699 | 25 |
| 923 | Tris(hydroxymethyl)aminomethane hydrochloride | | 35 |

pD

These SRM's are for the preparation of solutions of known deuterium-ion concentration to calibrate pH indicating equipment to indicate pD data. SRM's 2186I and 2186II, and 2191a and 2192a are certified for use as admixtures only.

| SRM | Type | pD(S) Values (at 25 °C) | Wt/Unit (grams) |
|--------|--------------------------------|-------------------------|--------------------|
| 2181 | Hepes/pH Buffer | 7.503 | 60 |
| 2182 | Hepes/pH Buffer | 7.516 | 60 |
| 2185 | Potassium hydrogen phthalate | 4.518 | 60 |
| 2186I | Potassium dihydrogen phosphate | 7.428 | 30 |
| 2186II | Disodium hydrogen phosphate | | 30 |
| 2191a | Sodium bicarbonate | 10.732 | 30 |
| 2192a | Sodium carbonate | | 30 |

Ion-Selective Electrodes

These SRM's are certified for the calibration of ion-selective electrodes and have conventional ionic activities based on the Stokes-Robinson hydration theory for ionic strengths greater than 0.1 mole per liter.

| SRM | Type | Certified Property | Wt/Unit (grams) |
|------|--------------------|--------------------|--------------------|
| 2201 | Sodium Chloride | pNa, pCl | 125 |
| 2202 | Potassium Chloride | pK, pCl | 160 |
| 2203 | Potassium Fluoride | pF | 125 |

METROLOGY

Scanning Electron Microscope (SEM)

These SRM's are for calibrating the magnification scale and evaluating the performance of Scanning Electron Microscopes. SRM 484f has spacings of 0.5, 1, 2, 5, 10, 30, and 50 μm and can be used to calibrate the magnification scale of an SEM from 1000 to 20,000 X to an accuracy of 5 percent or better. SRM 2069a consists of graphitized natural fibers with smooth and uniform edges on an SEM specimen mount.

| SRM | Type | Size |
|-------|----------------------------|-----------------------|
| 484f | SEM Magnification Standard | 11 mm D, 0.65 mm high |
| 2069a | SEM Performance Standard | 12 mm D, 3 mm peg |



Robbin Frazer creates a warm, friendly atmosphere in the certification and production office where she helps perform the documentation and records keeping for SRM certification projects.

Optical Microscope Linewidth-Measurement

These SRM's are for use in calibrating optical microscopes used to measure the widths of opaque lines and clear spaces on integrated-circuit photomasks. They can also be used to calibrate line spacings and line-to-space ratios. The accuracy of a measured linewidth or line spacing is $\pm 0.05 \mu\text{m}$ or better. They are not for use with partially transmitting materials, in reflected light with opaque materials, or in a scanning electron microscope. SRM 475 is made with anti-reflective chromium on a borosilicate glass substrate. SRM 476 is made with bright chromium.

| SRM | Type | Spacings | Size |
|-----|---|-------------------------|-------------------------------------|
| 473 | Optical Linewidth Standard IN PREP | | |
| 475 | Linenwidth Measurement Standard | 0.5 to 12 μm | 6.35 \times 6.35 \times 0.15 cm |
| 476 | Linenwidth Measurement Standard IN PREP | 0.5 to 12 μm | 6.35 \times 6.35 \times 0.15 cm |

Depth Profiling

This SRM is for calibrating equipment used to measure sputtered depth and erosion rates in surface analysis. SRM 2135c consists of nine alternating metal thin-film layers—five layers of pure chromium and four of pure nickel—on a polished silicon (100) substrate. It is certified for total chromium and total nickel thickness, for individual layer uniformity, for Ni/Cr bi-layer uniformity, and for individual layer thickness. The nominal thicknesses for Cr and Ni are 53 and 66 nm, respectively.

| SRM | Type | Unit/Size |
|-------|---|----------------------------------|
| 2135c | Ni-Cr Thin-Film Depth Profile Standard | IN PREP |
| 2136 | Cr/CrO Thin-Film Depth Profile Standard | 1 \times 2.54 \times 0.04 cm |
| 2137 | Boron Implant in Silicon Depth Profile | IN PREP |

Lorna Sniegoski, research chemist in the Organic Analytical Research Division, uses a microbalance to weigh an analyte to prepare a solution for analysis of an SRM.



COATING THICKNESS

These magnetic type thickness SRM's are 30×30 mm for calibrating coating thickness gages used to measure the thickness of nonmagnetic coatings on steel, or nickel on steel. The steel substrates have the properties of AISI 1010 steel and the nickel coatings have the properties of an annealed Watts nickel electrodeposited free of cobalt and iron.

These SRM's may be used to measure the thickness of paint and other organic coatings on steel, as well as zinc (galvanized) and other nonmagnetic metallic coatings.

Nonmagnetic Coating on Magnetic Substrate (Cu and Cr on Steel)

| SRM | Unit Size | Nominal Coating Thickness | |
|-------|-----------|---------------------------|--------------------|
| | | micrometer | milliinch (mil) |
| 1321 | Set of 4 | 34, 37, 42, 47 | 7.4, 1.5, 1.7, 1.9 |
| 1322 | Set of 4 | 53, 60, 69, 80 | 2.1, 2.4, 2.8, 3.2 |
| 1323 | Set of 4 | 4, 112, 135, 167 | 3.8, 4.5, 5.4, 6.7 |
| 1357 | Set of 3 | 6, 20, 48 | 0.24, 0.8, 1.9 |
| 1358 | Set of 3 | 80, 225, 1000 | 3.1, 10, 39 |
| 1359 | Set of 4 | 48, 140, 505, 800 | 2.0, 5.5, 20, 32 |
| 1360 | Set of 4 | 2.5, 6, 12, 20 | 0.1, 0.2, 0.5, 0.8 |
| 1361a | Set of 4 | 6, 12, 25, 48 | 0.2, 0.5, 1.0, 2.0 |
| 1362a | Set of 4 | 40, 80, 140, 205 | 1.6, 3.1, 5.5, 7.9 |
| 1363a | Set of 4 | 255, 385, 505, 635 | 9.8, 16, 20, 26 |
| 1364a | Set of 4 | 800, 1000, 1525, 1935 | 32, 39, 59, 79 |

Magnetic Coating on Magnetic Substrate (Nickel on Steel)

| SRM | Unit Size | Nominal Coating Thickness | |
|--------|-----------|---------------------------|--------------------|
| | | micrometer | milliinch (mil) |
| 1365a* | Set of 4 | 3, 9, 15, 20 | 0.1, 0.4, 0.6, 0.8 |
| 1366a* | Set of 4 | 25, 34.5, 40, 50 | 1.0, 1.4, 1.6, 2.0 |

*SRM's 1365a and 1366a will be reissued as SRM's 1331-1334 and 1335-1338, respectively.

Solder Thickness Standards

| SRM | Type | Unit Size |
|------|------------------|-----------|
| 2321 | Solder Thickness | IN PREP |
| 2322 | Solder Thickness | IN PREP |

COATING WEIGHT

The gold coating SRM's are 15×15 mm and were measured by beta-ray backscatter and x-ray fluorescence techniques relative to NIST gold coating materials for which the average weights per unit area were determined by weight and area measurements. These SRM's are for calibrating equipment used to measure weight per unit area of gold coating of equivalent purity.

Gold Coating on Nickel

| SRM | Unit Size | Nominal Coating Weight (mg/cm ²) | Nominal Coating Thickness | |
|-------|-----------|--|---------------------------|------------------|
| | | | micrometer | microinch |
| 1379 | 1 each | 0.35 | 0.175 | 7 |
| 1380 | 1 each | 0.55 | 0.275 | 11 |
| 1387 | 1 each | 2.2 | 1.4 | 45 |
| 1399b | Set of 4 | 1.5, 3.0, 6.0, 14.0 | 0.8, 1.5, 3, 7 | 30, 60, 120, 280 |

ELLISSOMETRY

Each of these SRM's is certified for the ellipsometric parameters of delta (Δ) and psi (Ψ) and the derived thickness and refractive index of the silicon dioxide layer on the silicon wafer.

| SRM | Type | Unit Size | Nominal Thickness |
|------|---------------------|-----------------|-------------------|
| 2531 | Thin Film Thickness | IN PREP | 50 nm |
| 2532 | Thin Film Thickness | 76 mm Dia wafer | 100 nm |
| 2533 | Thin Film Thickness | IN PREP | 200 nm |
| 2534 | Thin Film Thickness | IN PREP | 25 nm |

Carmelita Davis exhibits a warm, expressive personality and dependability as she handles the Materials Safety Date Sheets along with other responsibilities.



Glass

Chemical Resistance (Durability) of Glass

These SRM's are for checking test methods and calibrating equipment used to determine the resistance of glass containers to chemical attack. The values below represent the volume of fiftieth-normal sulfuric acid used to titrate to the methyl-red end point the alkaline extract from a crushed sample of glass after exposure to high-purity water at 121 °C.

| SRM | Type | Unit of Issue | mL of N/50 H ₂ SO ₄ |
|-----|------------------|---------------|---|
| 622 | Soda-lime-silica | 2.2 kg | 7.67 |
| 623 | Borosilicate | 2.2 kg | 0.34 |

Electrical Properties of Glass

SRM 624 is for checking test methods and for calibrating equipment used to determine the dc volume resistivity of glass per ASTM C657. SRM 774 is for checking methods used to determine dielectric constant and ac loss characteristics of insulating materials per ASTM D150.

| SRM | Type | Unit of Issue | Approximate Value |
|-----|--------------------------------------|---------------|---|
| 624 | Lead-silica, for dc resistivity | 200 g | $\log_{10}\rho \sim 9.9 \Omega\text{-cm}$ |
| 774 | Lead-silica, for dielectric constant | 5×5×2.5 cm | K ~ 7.47 |

Viscosity

SRM's 710a, 711, and 717 are rectangular bars for checking the performance of high-temperature viscosity equipment (rotating cylinders) and low-temperature viscosity equipment (fiber elongation, beam-bending, parallel-plates, etc.).

| SRM | (IN PREP) | Temperature (°C) at Viscosity (poises) | | | | | | | | | | |
|------|-----------|--|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|------------------|------------------|------------------|
| | | 10 ² | 10 ³ | 10 ⁴ | 10 ⁵ | 10 ⁶ | 10 ⁷ | 10 ⁸ | 10 ⁹ | 10 ¹⁰ | 10 ¹¹ | 10 ¹² |
| 710a | (IN PREP) | | | | | | | | | | | |
| 711 | | 1327.1 | 1072.8 | 909.0 | 794.7 | 710.4 | 645.6 | 594.3 | 552.7 | 518.2 | 489.2 | 464.5 |
| 717 | | 1545.1 | 1248.8 | 1059.4 | 927.9 | 831.2 | 757.1 | 698.6 | 651.1 | 611.9 | 579.0 | 550.9 |

Viscosity Fixpoints

| SRM | Type of Glass | Unit of Issue | Softening Point °C | Annealing Point °C | Strain Point °C |
|------|---|---------------|--------------------|--------------------|-----------------|
| 709 | Extra Dense Lead | | | | |
| 710a | Soda Lime-Silica, type 523/586 IN PREP | 500 g | 384 | 328 | 311 |
| 711 | Lead-Silica, type 617/366 | 1.3 kg | 602 | 432 | 392 |
| 712 | Mixed Alkali Lead Silicate $\frac{1}{4}$ in patties (6 pcs.) | 225 g | 528 | 386 | 352 |
| 713 | Dense Barium Crown 620/603 $1\frac{1}{8}$ in diam $\times \frac{5}{8}$ in thick gobs (4 pcs.) | 225 g | 738 | 631 | 599 |
| 714 | Alkaline Earth Alumina Silicate $\frac{1}{4}$ in diam cane (16 pcs.—6 in long) | 225 g | 908 | 710 | 662 |
| 715 | Alkali-Free Aluminosilicate $\frac{1}{4}$ in diam cane (13 pcs.—6 in long) | 200 g | 961 | 764 | 714 |
| 716 | Neutral, $\frac{1}{2}$ in diam cane (6 pcs.—6 in long) | 250 g | 794 | 574 | 530 |
| 717 | Borosilicate, 4.2 cm \times 4.2 cm \times 12.5 cm bar | 450 g | 720 | 516 | 471 |

Relative Stress Optical Coefficient

These glasses are for calibrating instruments used to measure relative stress optical coefficient per ASTM C770. They are rectangular bars.

| SRM | Type of Glass | Unit of Issue | Relative Stress Optical Coefficient at $\lambda = 546.1$ nm | | |
|-----|------------------|---------------|---|---|--|
| 708 | Lead-Silica, A | 625 g | Glass A | $C = 2.857$ Brewsters, $10^{-12}m^2/N$ | |
| | Borosilicate, B | 275 g | Glass B | $C = 3.652$ Brewsters, $10^{-12}m^2/N$ | |
| 709 | Extra dense Lead | 500 g | | $C = -1.359$ Brewsters, $10^{-12}m^2/N$ | |

Glass Liquidus Temperature

This SRM is for checking test methods and for calibrating equipment used to determine the liquidus temperature of glass by the gradient furnace methods per ASTM C829.

| SRM | Type | Unit of Issue | Temperature, °C |
|-----|--|---------------|-----------------|
| 773 | Soda-lime-silica, for liquidus temperature 2.5 \times 2.5 \times 0.6 cm | 60 g | 990 |

Density

SRM's 211c, 2211, and 2213 are certified for density (air saturated at 1 atm) at 20, 25, and 30 °C, and may be used to calibrate pycnometers and density balances.

SRM's 1840 and 1841a are certified for density at 20 °C and may be used to determine the density of solids and liquids by means of hydrostatic weighing.

| SRM | Type | Density 20 °C (g/cm ³) | Amount |
|-------|---|------------------------------------|--------|
| 211c | Toluene | 0.86686 | 5 mL |
| 2211 | Toluene | 0.86686 | 8 mL |
| 2213 | 2,2,4 Trimethylpentane (<i>Isooctane</i>) | 0.691929 | 25 mL |
| 1841a | Silicon | 2.329 | 200 g |
| 1825 | Glass Density | 2.20185 | 22 g |
| 1826 | Glass Density | 2.54938 | 21 g |
| 1827 | Glass Density | 3.0495 | 20 g |

Microhardness

These SRM's are for use in calibrating and checking the performance of microhardness testers. These test blocks were made by electroforming the test metal on a steel substrate. The hardness numbers for 1893 through 1896 are each certified at loads of 25, 50, and 100-gram force, while 1905, 1906, and 1907 are certified for 300, 500, and 1000 gram-force, respectively.

| SRM | Type | Load | Hardness | Size |
|------|-------------------------|-----------------|------------------------|----------------|
| 1893 | Bright Copper (Knoop) | 25, 50, 100 g-f | 125 kg/mm ² | 12.5 mm square |
| 1894 | Bright Copper (Vickers) | 25, 50, 100 g-f | 125 kg/mm ² | 12.5 mm square |
| 1895 | Bright Nickel (Knoop) | 25, 50, 100 g-f | 600 kg/mm ² | 12.5 mm square |
| 1896 | Bright Nickel (Vickers) | 25, 50, 100 g-f | 600 kg/mm ² | 12.5 mm square |
| 1905 | Bright Nickel (Knoop) | 300 g-f | 600 kg/mm ² | 12.5 mm square |
| 1906 | Bright Nickel (Knoop) | 500 g-f | 600 kg/mm ² | 12.5 mm square |
| 1907 | Bright Nickel (Knoop) | 1000 g-f | 600 kg/mm ² | 12.5 mm square |



Gina Montgomery, a friendly, outgoing member of the customer sales and service staff, maintains and updates the SRM customer mail lists with valued dependability.

Ultrasonics

SRM 1855 is for point-by-point calibration of apparatus used to measure ultrasonic power. SRM 1856 is a displacement-measuring transducer to be used to determine the size and character of surface vibrations.

| SRM | Type | Frequency Range | Unit |
|------|------------------------------|-----------------|------|
| 1855 | Ultrasonic Power Transducer | 1.6 to 21.6 MHz | Each |
| 1856 | Acoustic Emission Transducer | 0.1 to 1 MHz | Each |

Polymers

Molecular Weight (Melt Flow)

| SRM | Type | Wt/Unit (grams) |
|------|--|--------------------|
| 705 | Polystyrene, narrow molecular weight distribution, $M_w \approx 179,300$, $M_w/M_n \approx 1.07$ | 5 |
| 706 | Polystyrene, broad molecular weight distribution, $M_w \approx 257,800$, $M_w/M_n \approx 2.1$ | 18 |
| 1474 | Polyethylene Melt Flow | IN PREP |
| 1475 | Polyethylene, linear, $M_w \approx 52,000$, $M_w/M_n \approx 2.9$ | 50 |
| 1476 | Polyethylene, branched | 50 |
| 1478 | Polystyrene, narrow molecular weight distribution, $M_w \approx 37,400$ ($M_w/M_n \approx 1.04$) | 2 |
| 1479 | Polystyrene, narrow molecular weight distribution, $M_w \approx 1,050,000$ | 2 |
| 1480 | Polyurethanes Low MW | IN PREP |
| 1481 | Polyurethanes HI MW | IN PREP |
| 1482 | Polyethylene, linear, $M_w \approx 13,600$ ($M_w/M_n \approx 1.19$) | 1 |
| 1483 | Polyethylene, linear, $M_w \approx 32,100$ ($M_w/M_n \approx 1.11$) | 1 |
| 1484 | Polyethylene, linear, $M_w \approx 119,600$ ($M_w/M_n \approx 1.19$) | 1 |
| 1487 | Poly(methylmethacrylate), $M_n \approx 6,000$ | 2 |
| 1488 | Poly(methylmethacrylate), $M_n \approx 29,000$ | 2 |
| 1489 | Poly(methylmethacrylate), $M_n \approx 115,000$ | 2 |
| 1496 | Polyethylene Gas Pipe Resin (Unpigmented) | 2 lbs. |
| 1497 | Polyethylene Gas Pipe Resin (Pigmented) | 20 lbs. |
| 1923 | Polystyrene Sulfonates | IN PREP |
| 1924 | Polystyrene Sulfonates | IN PREP |
| 1925 | Polystyrene Sulfonates | IN PREP |
| 8450 | Polyethylene Piping $\frac{1}{2}$ in | 2 Pieces |
| 8451 | Polyethylene Piping 2 in | 2 Pieces |
| 8452 | Polyethylene Piping 4 in | 2 Pieces |
| 8453 | Polyethylene Piping Socket T | 3 Pieces |
| 8454 | Polyethylene Piping Butt T | 3 Pieces |

These materials are certified for the properties indicated in the table.

| Property | Method | 705 | 706 | 1475 | 1476 | 1478 | 1479 | 1482 | 1483 | 1484 | 1487 | 1488 | 1489 |
|---|--|-----|-----|------|------|------|------|------|------|------|------|------|------|
| Molecular Weight: Weight Average | (Light Scattering) (Sed. Equili.) (Gel Permeation Chromatography- GPC) | X | X | X | | X | X | X | X | X | X | X | X |
| Number Average | (Osmometry) (GPC) | X | | X | | X | X | X | X | | | X | |
| Molecular Weight Distribution | (GPC) | | | X | | | | | | | | | |
| Limiting Viscosity No. Benzene 25 °C Benzene 35 °C Cyclohexane 35 °C 1-Chloronaphthalene 130 °C 1,2,4-trichlorobenzene 130 °C Decahydronaphthalene 130 °C Tetrahydrofuran 25 °C Toluene 25 °C | (Capillary Viscometer) | X | X | | | X | | | | | | | X |
| Melt Flow | (ASTM) | | | X | X | | | | | | | | |
| Density | (ASTM) | | | X | X | | | | | | | | |
| Heat Capacity | (Adiabatic) | X | | X | | | | | | | | | |

Rheology

This SRM is for calibrating instruments used in polymer technology and science to determine rheological properties of polymer melts or solutions. It is certified for Rate of Shear, Viscosity, and First Normal Stress Difference at 25 °C.

| SRM | Type | Unit size |
|------|------------------------------------|-----------|
| 1490 | Polyisobutylene Solution in Cetane | 250 mL |

Heat

Calorimetric

These SRM's are intended to relate the gain or loss of energy and work experienced during a chemical reaction or by change of temperature to the units of energy and work as defined by the International System of Units (SI). The unit for energy and work under this system is the joule, which is related to the calorie by the equation: 4.184 joule=1 calorie.

Combustion Calorimetric

| SRM | Type | Approximate Heat of Combustion (MJ/kg) | Unit Amount |
|-------|---|--|-------------|
| 39i | Benzoic Acid | 26.434 | 30 g |
| 2213 | 2,2,4-Trimethylpentane (<i>Isooctane</i>) | 47.712 | 25 mL |
| 1656 | Thianthrene | 33.480 | 30 g |
| 1657 | Synthetic Refuse Derived Fuel | 13.87 | 100 g |
| 2151 | Nicotinic Acid | 22.184 | 25 g |
| 2152 | Urea | 10.536 | 25 g |
| 2683a | Coal, Bituminous: %S= 1.85; %Ash= 6.85 | 32.30 (13887 Btu/lb) | 50 g |
| 2684a | Coal, Bituminous: %S= 3.00; %Ash= 11.09 | 28.92 (12433 Btu/lb) | 50 g |
| 2685 | Coal, Bituminous: %S= 4.62; %Ash= 16.53 | 27.38 (11770 Btu/lb) | 50 g |

NOTE: The calorific values (MJ/kg) may decrease upon the aging or normal oxidation of the coals. NIST will continue to monitor these calorific values and report any substantive change to the purchaser.

Solution Calorimetric

| SRM | Type | Heat of Solution (MJ/kg) | Wt/Unit (grams) |
|------|--|------------------------------------|-----------------|
| 724a | Tris(hydroxymethyl)aminomethane (Hydrochloric Acid and Sodium Hydroxide Solution Calorimetry) | Evolved 0.24576 Absorbed 0.1418 | 50 |
| 1655 | Potassium Chloride (Water Solution Calorimetry) | Absorbed (0.235) | 30 |

Heat Source Calorimetric

| SRM | Type | Heat of Evolution (MJ/kg) | Wt/Unit (grams) |
|------|--|---------------------------|-----------------|
| 1651 | Zirconium-barium chromate heat source powder | 1.46 | 50 |
| 1652 | Zirconium-barium chromate heat source powder | 1.632 | 50 |
| 1653 | Zirconium-barium chromate heat source powder | 1.762 | 50 |

Enthalpy and Heat Capacity

| SRM | Type | Temperature Range (K) | Unit Size |
|------|----------------------|-----------------------|-----------|
| 705 | Polystyrene, powder | 10–350 | 5 g |
| 1475 | Polyethylene, powder | 5–360 | 50 g |

Differential Scanning Calorimetry

These SRM's are for calibrating differential scanning calorimeters, differential thermal analyzers, and similar instruments.

| SRM | Type | Melting Temperature | Enthalpy of Fusion | Unit of Issue (mm) |
|-------|--------------------|---------------------|----------------------------|--------------------|
| 2220 | Tin (99.9995%) | 505.10 K | 7147 J/mol (107.46 J/g) | 25×25×0.127 |
| 2221a | Zinc (99.999%) | 692.74 K | 7026 J/mol (60.22 J/g) | 25×25×0.0508 |
| 2222 | Biphenyl (99.984%) | 342.41 K | 120.41 J/g | 1 g |
| 2223 | Potassium Nitrate | IN PREP | | |
| 2225 | Mercury | 234.30 | 11.469 J/g | 2.5g |

This SRM is for evaluating methods of determining purity by differential scanning calorimetry. It consists of phenacetin and phenacatin doped with p-aminobenzoic acid.

| SRM | Type | Dopant Level (p-ABA, mol%) | Unit |
|------|-------------------------|----------------------------|----------------------|
| 1514 | Thermal Analysis Purity | 0, 0.7, 2, 5 | Set of 4, 0.5 g/vial |

Differential Thermal Analysis

GM's 754, 757, 758, 759, 760, and 761 have been issued by NIST in cooperation with the International Confederation of Thermal Analysis as standards for calibrating differential thermal analysis, differential scanning calorimetry, and thermogravimetry equipment under operating conditions.

| GM | Material | | Peak Temp. | Unit |
|-----|-----------------------------|------------------------------------|------------------------|---------------------|
| 754 | Polystyrene | (glass transition) | 105 °C | 10 g |
| 757 | 1,2-Dichloroethane | (melting point) | -32 °C | 4 mL |
| | Clycohexane | (transition point) | -83 °C | 4 mL |
| | Phenyl Ether o-Terphenyl | (melting point) (melting point) | 7 °C 30 °C 58 °C | 4 mL 4 mL 5 g |
| 758 | Potassium Nitrate | (transition point) | 128 °C | 10 g |
| | Indium | (melting point) | 157 °C | 3 g |
| | Tin | (melting point) | 232 °C | 3 g |
| | Potassium Perchlorate | (transition point) | 300 °C | 10 g |
| | Silver Sulfate | (transition point) | 430 °C | 3 g |
| 759 | Potassium Perchlorate | (transition point) | 300 °C | 10 g |
| | Silver Sulfate | (transition point) | 430 °C | 3 g |
| | Quartz | (transition point) | 573 °C | 3 g |
| | Potassium Sulfate | (transition point) | 583 °C | 10 g |
| | Potassium Chromate | (transition point) | 665 °C | 10 g |
| 760 | Quartz | (transition point) | 573 °C | 3 g |
| | Potassium Sulfate | (transition point) | 583 °C | 10 g |
| | Potassium Chromate | (transition point) | 665 °C | 10 g |
| | Barium Carbonate | (transition point) | 810 °C | 10 g |
| | Strontium Carbonate | (transition point) | 925 °C | 10 g |
| 761 | Permanorm 3 | (magnetic transition) | 259 °C | 1 g |
| | Nickel | (magnetic transition) | 353 °C | 1 g |
| | Mumetal | (magnetic transition) | 381 °C | 1 g |
| | Permanorm 5 | (magnetic transition) | 454 °C | 1 g |
| | Trafoperm | (magnetic transition) | 750 °C | 1 g |

Superconductive Thermometric Fixed Point Device

Each device is composed of small cylinders of high purity material mounted in a threaded copper stud and enclosed by a mutual inductance coil set. SRM 767a is intended to provide fixed points on the 1976 Provisional 0.5 to 30 K Temperature Scale (EPT-76). This SRM should prove particularly valuable to users of ^3He - ^4He dilution refrigerators, in which direct calibrations on the liquid helium vapor pressure-temperature scales are difficult, and to those who wish to determine the temperature reproducibility of physical phenomena or of cryogenic equipment.

| SRM | Type | Material | Nominal Temperature (K) |
|------|---|--|--|
| 767a | Superconductive Thermometric Fixed Point Device | Niobium Lead Indium Aluminum Zinc Cadmium | 9.3 7.2 3.4 1.2 0.9 0.5 |

Freezing Point

SRM's 740a and 741 are defining fixed points for the International Practical Temperature Scale of 1968 (IPTS-68). The secondary reference points are for calibrating thermometers, thermocouples, and other temperature measuring devices. These SRM's are certified per IPTS-68.

Defining Fixed Points

| SRM | Type | Temperature °C | Wt/Unit (grams) |
|-----|------|----------------|-----------------|
| 740 | Zinc | 419.58 | 350 |
| 741 | Tin | 231.9681 | 350 |

Secondary Reference Points

| SRM | Type | Temperature °C | Wt/Unit (grams) |
|-----|----------|----------------|-----------------|
| 43h | Zinc | 419.5 | 350 |
| 44f | Aluminum | 660.3 | 200 |
| 45d | Copper | 1084.8 | 450 |
| 49e | Lead | 327.493 | 600 |
| 743 | Mercury | -38.841 | 680 |

Melting Point

| SRM | Type | Form | Temperature °C | Wt/Unit (grams) |
|------|----------------------------|-------------|----------------|-----------------|
| 742 | Alumina, 99.9 + % | Powder | 2053 | 10 |
| 1968 | Gallium, 99.9999 + % | Sealed Cell | 29.7723 | 25 |
| 1969 | Rubidium, 99.9 + % | Sealed Cell | 39.3 | 154 |
| 1970 | Succinonitrile, 99.999 + % | Sealed Cell | 58.079 | 60 |
| 1971 | Indium, 99.9999 + % | Sealed Cell | 156.635 | 100 |

GM 8000 is issued by NIST in cooperation with the Office of Reference Materials at the National Physical Laboratory (NPL) in Teddington, England. This set of 10 highly purified substances is intended for use in the calibration of thermometry used in determining the melting points of samples in glass capillary tubes. Both the meniscus point and the liquefaction point for each substance are certified by NPL.

| GM | Type | Melting Point | Amount |
|------|---------------------------|---------------|--------|
| 8000 | 4-Nitrotoluene | 52 °C | 1 g |
| | Naphthalene (NPL# M14-11) | 80 | 1 g |
| | Benzil | 95 | 1 g |
| | Acetanilide | 114 | 1 g |
| | Benzoic Acid | 122 | 1 g |
| | Diphenylacetic Acid | 147 | 1 g |
| | Anisic Acid | 183 | 1 g |
| | 2-Chloroanthraquinone | 210 | 1 g |
| | Carbazole | 246 | 1 g |
| | Anthraquinone | 285 | 1 g |

Laboratory Thermometer

This mercury-in-glass thermometer is for use in clinical laboratories. Its main scale extends from 24.00 to 38.00 °C, in 0.05 in °C divisions. It has an auxiliary scale from –0.20 to +0.20 °C.

| SRM | Type | Calibrated Points (°C) | Unit |
|-----|---------------------------------|------------------------|------|
| 934 | Clinical Laboratory Thermometer | 0, 25, 30, 37 | Each |
| | | | |
| | | | |
| | | | |

Thermocouple Material

| SRM | Type | Temperature Range (°C) | Form |
|------|------------------------------------|------------------------|-------------------------------|
| 1967 | Platinum, High-Purity (99.999 + %) | 197 to 1767 | Wire: 0.51 mm D, 1 meter long |
| | | | |
| | | | |
| | | | |

Vapor Pressure

| SRM | Type | Pressure Range (atmosphere) | Temperature Range (K) | Unit Size |
|-----|---------|---------------------------------------|-----------------------|-----------------------|
| 745 | Gold | 10 ⁻⁹ to 10 ⁻³ | 1300–2100 | Wire 1.44 mm × 152 mm |
| 746 | Cadmium | 10 ⁻¹¹ to 10 ⁻⁴ | 350–594 | Rod 6.4 mm × 64 mm |
| 748 | Silver | 10 ⁻¹² to 10 ⁻³ | 800–1600 | Rod 6.4 mm × 64 mm |

Thermal Conductivity

| SRM | Type | Dimension (mm) | Temperature Range (K) | Conductivity at 293 K (W/M·K) |
|------|-------------------|-------------------|-----------------------|-------------------------------|
| 1461 | Stainless Steel | 12.7 D, 50 length | 2-1200 | 14.1 |
| 1462 | Stainless Steel | 34 D, 50 length | 2-1200 | 14.1 |
| 8420 | Electrolytic Iron | 6.4 D, 50 length | 2-1000 | 77.9 |
| 8421 | Electrolytic Iron | 31.7 D, 50 length | 2-1000 | 77.9 |
| 8422 | Sintered Tungsten | 3.2 D, 50 length | 2-3000 | 173 |
| 8423 | Sintered Tungsten | 6.4 D, 50 length | 2-3000 | 173 |
| 8424 | Graphite | 6.4 D, 50 length | 5-2500 | 90.9 |
| 8425 | Graphite | 12.7 D, 50 length | 5-2500 | 90.9 |
| 8426 | Graphite | 25.4 D, 50 length | 5-2500 | 90.9 |

Thermal Expansion

| SRM | Type | Temperature Range (K) | Diameter (mm) | Length (mm) |
|--------|----------------------------|-----------------------|---------------|-------------|
| 731-L1 | Borosilicate Glass | 80-680 | 6.4 | 51 |
| 731-L2 | Borosilicate Glass | 80-680 | 6.4 | 102 |
| 731-L3 | Borosilicate Glass | 80-680 | 6.4 | 152 |
| 738 | Stainless Steel (AISI 446) | 293-780 | 6.4 | 51 |
| 739-L1 | Fused Silica | 80-1000 | 6.4 | 51 |
| 739-L2 | Fused Silica | 80-1000 | 6.4 | 102 |
| 739-L3 | Fused Silica | 80-1000 | 6.4 | 152 |

Thermal Resistance

| SRM | Type | Dimension (mm) | Temperature Range (K) | Thermal Resistance |
|-------|-----------------------|----------------|-----------------------|--------------------------------|
| | | | | at 293 K (m ² ·K/W) |
| 1449 | Fumed Silica Board | 600×600×25.4 | 297.1 | 1.2 |
| 1450b | Fibrous Glass Board | 600×600×25.4 | 100-330 | 0.75 |
| 1451 | Fibrous Glass Blanket | 600×600×25.4 | 100-330 | 0.60 |
| 1452 | Fibrous Glass Blanket | 600×600×25.4 | 297.1 | 0.60 |
| 1459 | Fumed Silica Board | 300×300×25.4 | 297.1 | 1.2 |

Magnetic

Magnetic Susceptibility

| SRM | Type | Gram Susceptibility @297 K | | Form/Unit |
|-----|--------------------|---|--|-----------------------------|
| | | (10 ⁶ χg, cm ³ ·g ⁻¹) | | |
| 763 | Aluminum | 0.604 | | Cylinder 3 mm diameter×3 mm |
| 766 | Manganese Fluoride | 123.3 | | Cube 3×3×3 mm |

Magnetic Moment

| SRM | Type | Magnetic, H (Oe) | Moment, σ | |
|-----|---------------|------------------|--|----------|
| | | | (Oe·cm ³ ·g ⁻¹) | Size |
| 772 | Nickel Sphere | 3500–10,000 | 54.75–54.90 | 2.4 mm D |



Ken Eckerle of the Radiometric Physics Division is shown here making an adjustment to a component of the high accuracy spectrophotometer in preparation for measurements of one of the optical properties SRM's.

OPTICAL

Spectrophotometric

SRM 930D: This SRM consists of three neutral density glass filters. The filters have transmittances of approximately 10, 20, and 30 percent. Each filter is individually certified for transmittance at wavelengths of 440, 465, 546.1, 590, and 635 nm.

SRM 931d: This SRM consists of three sets of four solutions—a blank solution and three concentrations of absorbing liquid. The net absorbances are certified for each concentration at wavelengths of 302, 395, 512, and 678 nm.

SRM 935a: Solutions made with this SRM are certified for apparent specific absorbances at wavelengths of 235, 257, 313, 345, and 350 nm.

SRM 936: A solution made with this SRM is certified for its molecular emission spectrum over the wavelength range of 375 to 675 nm.

SRM 1930: This SRM consists of three neutral density glass filters. The filters have transmittances of approximately 1, 3, and 50 percent. Each filter is individually certified for transmittance at wavelengths of 440, 465, 546.1, 590, and 635 nm.

SRM 1931: This SRM is a set of four luminescent samples and a blank certified for corrected luminescence emission spectra over the wavelength range from 400 to 760 nm. Each sample consists of an inorganic dye in a sintered polytetrafluoroethylene matrix mounted in a cuvette-sized holder.

SRM 2009a: This SRM is for checking the wavelength scale between 400 and 760 nm for bandpasses between 1.5 and 10.5 nm. SRM 2009a is mounted in a standard cuvette-sized holder.

SRM 2031: This SRM consists of three filters mounted in holders and an empty holder; all holders are equipped with shutters. Two of the filters have an evaporated layer of semi-transparent metal sandwiched between two quartz plates assembled by optical contact. The third filter is a single quartz plate. Each filter is individually calibrated at 250, 280, 340, 360, 400, 465, 500, 546.1, 590, and 635 nm.

SRM 2032: Aqueous solutions made with this SRM are certified for specific absorbances from 240 to 280 nm for use as a stray light standard in the ultraviolet region.

SRM 2033: This SRM consists of the same material as SRM 2032 plus a reference beam attenuator for extending the dynamic range of the stray light test.

SRM 2034: This SRM is a solution sealed in a non-fluorescent, fused-silica cuvette for checking the wavelength scale between 240 and 650 nm.

| SRM | Type | Wavelength Range (nm) | Unit |
|-------|---|-----------------------|-----------------------|
| 930D | Glass Filters, Transmittance | 440–635 | 3 filters/4 holders |
| 931d | Liquid Filters, Absorbance | 302–678 | Set: 12 vials |
| 935a | Potassium Dichromate, UV Absorbance | 235–350 | 15 grams |
| 936 | Quinine Sulfate Dihydrate, Fluorescence | 375–675 | 1 gram |
| 1930 | Glass Filters, Transmittance | 440–635 | 3 filters/4 holders |
| 1931 | Fluorescence Corrected Emission Spectra | 400–760 | Set: 4 pieces |
| 2009a | Didymium-oxide Glass, Wavelength | 400–760 | 1 filter/1 holder |
| 2031 | Metal-on-Quartz Filters, Transmittance | 250–635 | 3 filters/4 holders |
| 2032 | Potassium Iodide, Stray Light | 240–280 | 25 grams |
| 2033 | Potassium Iodide with Attenuator | 240–280 | 25 grams w/attenuator |
| 2034 | Holmium-oxide Solution, Wavelength | 240–650 | 1 sealed cuvette |

Reflectance

These SRM's are for calibrating the reflectance scale of integrating sphere reflectometers used to evaluate materials for solar energy collectors and to calibrate reflectometers used in evaluating the appearance of polished metals and metal plated objects.

Specular Spectral Reflectance

| SRM | Type | Wavelength Range (nm) | Size |
|-------|---|-----------------------|-------------|
| 2003 | First Surface, Aluminum on Glass | 250–2500 | 5.1 cm D |
| 2011 | First Surface, Gold on Glass | IN PREP | 5.1 cm D |
| 2023a | Second Surface, Aluminum on Fused Quartz | IN PREP | 5.1 cm D |
| 2025 | Second Surface, Aluminum on Fused Quartz with Wedge | 250–2500 | 2.5×10.2 cm |

Infrared Reflectance

This SRM is for use in establishing the accuracy of the wavelength scale of reflectance spectrophotometers.

| SRM | Type | Wavelength Range (nm) | Size |
|------|--------------------|-----------------------|-----------------|
| 1920 | Near IR Wavelength | 740–2000 | 51 mm D × 12 mm |



Yvonne Barnes prepares to mount a sample on the reference spectral reflectometer to make certification measurements on a reflectance SRM.

Directional-Hemispherical Reflectance

| SRM | Type | Wavelength Range (nm) | Size |
|------|------------------------|-----------------------|-----------------|
| 2015 | Opal Glass | 400–750 | 2.5×5.0×0.64 cm |
| 2016 | Opal Glass | 400–750 | 10×10×0.64 cm |
| 2021 | Black Porcelain Enamel | 280–2500 | 5.1×5.1×0.20 cm |

Refractive Index

SRM's 211c, 2211, and 2213 are certified for refractive index at 20, 25, and 30 °C, from 435.8 to 667.8 nm for seven wavelengths.

SRM 1822 is certified for refractive index at thirteen wavelengths from 404.7 nm to 706.5 nm. This SRM is designed for calibrating refractometers and certifying refractive index immersion liquids. It consists of two rectangular glass slabs: one slab has polished faces and is to be used to check the performance of a refractometer; the second slab is unpolished and can be broken into fragments to certify the refractive index of immersion liquids by microscope methods.

SRM 1823 consists of two silicone liquids that are miscible and span the refractive index range of a variety of glasses and glass fibers. The liquids are suitable for calibrating refractometers and are certified for refractive index at ten wavelengths from 435.8 to 667.8 nm, at temperatures of 20, 40, 60, and 80 °C.

| SRM | Type | n ²⁰ , λ546.1 | Unit Size |
|---------|---|--------------------------|--------------|
| 211c | Toluene | 1.5008 | 5 mL |
| 2211 | Toluene | 1.5008 | 8 mL |
| 2213 | 2,2,4-Trimethylpentane (<i>Isooctane</i>) | 1.3934 | 25 mL |
| 1822 | Glass (Soda-Lime) | 1.5200 | Set: 2 slabs |
| 1823-I | Silicone Liquid (I) | 1.5214 | 60 mL |
| 1823-II | Silicone Liquid (II) | 1.5638 | 60 mL |

Optical Rotation

These SRM's are intended for use in calibrating or checking polarimetric apparatus. In aqueous solution the optical rotation of SRM 17d is certified at three wavelengths, while that of SRM 41c is certified at two wavelengths. SRM 41c is also certified at one wavelength in a dimethyl sulfoxide solution.

| SRM | Type | Optical Rotation In Aqueous Solution (mrad) | | | Unit Size |
|-----|-------------|---|--------|--------|-----------|
| | Wavelengths | 546 | 589 | 633 | |
| 17d | Sucrose | 711.64 | 604.26 | 519.17 | 60 g |
| 41c | Dextrose | 1101.1 | 931.8 | 798.6 | 70 g |

RADIOACTIVITY

These SRM's are shipped express or air freight (shipping charges collect). The amount of a radionuclide in an SRM, at a specified time, is stated as (1) the number of atoms (or the mass, for radium SRM's), (2) the activity, or "decays per second," or (3) the emission rate of a particular radiation, depending on the method of calibration or the intended use. For solution SRM's, the quantity is usually specified per gram of liquid. The active portion of gamma-ray "point-source" standards is usually restricted to the central few millimeters of a low-mass, low-Z support to minimize scattering. Alpha-particle-emitting radionuclides are deposited or plated on metal backings.

The unit for activity has traditionally been the curie (Ci), but simpler relations between activity, emission rate, and counting rate result if the current SI (International System of Units) unit "1 per second" is used. This is symbolized as " s^{-1} " and has been given the special name becquerel (Bq). The relationship between the curie and the becquerel is:


$$1 \text{ Ci} = 3.7 \times 10^{10} \text{ Bq.}$$

Many SRM's are measured and certified in terms of emission rate. In this catalog, αs^{-1} , $\beta^- s^{-1}$, $\beta^+ s^{-1}$, $Kx s^{-1}$, and γs^{-1} are used for the emission rates of alpha particles, negatrons, positrons, K x rays, and gamma rays, respectively.

The SRM's without an asterisk (*) may be ordered singly, without a license, under the general licensing provisions of the Atomic Energy Act of 1954. Those marked by an asterisk are available only under the special licensing provisions of the Atomic Energy Act of 1954.

NOTE: Certain radionuclides are not economical to maintain in stock because of short half lives or low demand. When sufficient demand exists, based on letters of inquiry, these materials are prepared and those who have expressed interest are notified of their availability. If you need any radionuclide not listed, write the Radioactivity Group, Room C114 Radiation Physics Building, National Institute of Standards and Technology, Gaithersburg, MD 20899; or telephone (301) 975-5531.

In addition, chemically stable solutions of most radionuclides may be submitted to NIST for calibration as described in Calibration Services Users Guide, NIST Special Publication 250 (1986-88/Revised ed.). Requests for such tests should be submitted, with full source information, for approval of suitability to the Radioactivity Group.



Bert Coursey and Bill McLaughlin of the Ionizing Radiation Division discuss the preparation and measurement of a radioactivity sample.

Alpha-Particle, Beta-Particle, Gamma-Ray, and Electron-Capture Solutions

| SRM | Radionuclide | Approximate activity, per gram, at time of calibration (month/year) (Bq g ⁻¹) | | Approx. Mass of Solution (g) | Overall Uncertainty (%) |
|--------|--|--|-------|------------------------------------|-------------------------------|
| 4322* | Americium-241 | 38.7 | 11/86 | 5 | 1.0 |
| 4332C* | Americium-243 | 89 | 9/89 | 5 | 1.4 |
| 4251B* | Barium-133 | 5×10^5 | 1/82 | 5 | 1.4 |
| 4233C* | Cesium-137-Barium-137m | 7×10^5 | 8/89 | 5 | 1.4 |
| 4943* | Chlorine-36 | 1×10^4 | 12/84 | 3 | 0.8 |
| 4915D* | Cobalt-60 | 3×10^5 | 3/84 | 5 | 0.8 |
| 4320 | Curium-244 | 57 | 4/89 | 5 | 1.0 |
| 4329* | Curium-243 | 70 | 6/84 | 5 | 1.4 |
| 4370* | Europium-152 | 9×10^4 | 2/87 | 5 | 1.1 |
| 4926D | Hydrogen-3 | 3×10^3 | 5/89 | 18 | 0.6 |
| 4927D | Hydrogen-3 | 7×10^5 | 4/89 | 3 | 0.8 |
| 4947C | Hydrogen-3 | 3×10^5 | 3/87 | 4 | 1.2 |
| 4361B | Hydrogen-3 | 1.1 | 8/87 | 490 | 1.0 |
| 4949B | Iodine-129 | 7×10^3 | 1/82 | 1 | 1.9 |
| 4929D | Iron-55 | 4×10^4 | 8/85 | 5 | 2.6 |
| 4226B* | Nickel-63 | 1×10^6 | 12/84 | 4 | 1.1 |
| 4327* | Polonium-208 | 77 | 6/84 | 1.1 | 1.4 |
| 4326 | Polonium-209 | 370 | 9/89 | 5 | 2.0 |
| 4323* | Plutonium-238 | 33 | 11/86 | 5 | 0.5 |
| 4338* | Plutonium-240 | 18 | 4/80 | 5 | 1.0 |
| 4940C | Promethium-147 | 1×10^5 | 8/85 | 5 | 0.4 |
| 4423 | Strontium-90 | 4×10^6 | 11/85 | 5 | 1.1 |
| 4919F* | Strontium-90 | 4×10^3 | 5/88 | 5 | 1.2 |
| 4288* | Technetium-99 | 4×10^4 | 11/82 | 5 | 1.6 |
| 4328* | Thorium-229 | 884 | 5/84 | 2 | 1.5 |
| 4324* | Uranium-232 | 83 | 2/84 | 5 | 1.5 |
| 4321 | Uranium-238 (Natural) | 263 | 11/86 | 5 | 0.4 |
| 4276C* | Long-Lived Mixed Radionuclide: Antimony-125 Europium-154 Europium-155 | 2×10^4 1×10^4 7×10^3 | 9/88 | 5 | |

*License certification is required by NIST for these radionuclides.

Alpha-Particle Point-Sources

These SRM's consist of a practically weightless deposit of the nuclide on a thin platinum foil cemented to a monel disk.

| SRM | Radionuclide | Approx. α -particle-emission rate into 2π geometry and/or approx. activity at time of calibration (month/year) | Overall uncertainty (%) |
|---------|---------------|---|----------------------------|
| 4904NG* | Americium-241 | 2300 Bq g ⁻¹ | 1.3 |
| 4904SG* | Americium-241 | 2300 Bq g ⁻¹ | 1.0 to 1.3 |
| 4906HC* | Plutonium-238 | 1×10^3 – 5×10^4 | 1.3 |
| 4906C* | Plutonium-238 | 1×10^{-3} – 10^2 | 1.3 |

*License certification is required by NIST for these radionuclides.

Radiocarbon Dating and Ground Water Studies

Contemporary Standard for Carbon-14 Dating Laboratories

| SRM | Material | Description |
|-------|-------------|---|
| 4990C | Oxalic Acid | One-half pound of oxalic acid taken from specially prepared material for use as a common contemporary standard against which world-wide measurements can be compared. |

Low-Level Tritiated-Water Standard

| SRM | Material | Description |
|-------|------------|---|
| 4361B | Hydrogen-3 | Contains 490 grams of ${}^3\text{H}-\text{H}_2\text{O}$ in a flame-sealed bottle. The radioactivity concentration was 1.123 Bq g ⁻¹ , as of the date of the most recent gas-counting measurement—August 1, 1987. The total uncertainty in this value is 1.0. |

Gaseous Materials

| SRM | Radionuclide | Approximate activity or radioactivity concentration at time of calibration (month/year) | Approx. Vol. (cm ³) | Approx. Pressure (atm) | Overall Uncertainty (%) | |
|---------|--------------|---|------------------------------------|------------------------------|-------------------------------|-----|
| 4935C* | Krypton-85 | 5×10^7 Bq mol ⁻¹ | 3/74 | 10 | 1 | 0.9 |
| 4235B* | Krypton-85 | 1×10^7 Bq | 10/86 | 3 | 1 | 1.1 |
| 4308C | Krypton-85 | 2×10^6 Bq | 11/82 | 30 | 0.3 | 3.2 |
| 4415LN* | Xenon-133 | 5×10^8 Bq | 3/90 | 5 | 0.1 | 1.5 |

*License certification is required by NIST for these radionuclides.

Gamma-Ray and X-Ray Point-Sources

These SRM's are usually prepared by depositing the radioactive material and sealing it between two layers of polyester tape, mounted on an aluminum ring. SRM 4206c, Thorium-228, is prepared by depositing and sealing the radionuclide between two layers of gold foil and this sandwich is then sealed between two double layers of polyurethane-film tape.

| SRM | Radionuclide | Principal Photon Energy (MeV) | Approximate activity, Bq, at time of calibration (except MRN) (month/year) | Overall Uncertainty (%) |
|--------------------------------------|---|-------------------------------|--|-------------------------|
| 4202D | Cadmium-109-Silver-109m | 0.088 | 4×10^5 | 10/86 |
| 4241B* | Barium-133 | 0.081-384 | 2×10^5 | 6/81 |
| 4200B | Cesium-137-Barium-137m | 0.662 | 4×10^4 | 9/79 |
| 4207B | Cesium-137-Barium-137m | 0.662 | 3×10^5 | 3/87 |
| 4203D* | Cobalt-60 | 1.173-1.332 | 2×10^4 to 2×10^5 | 3/84 |
| 4218E* | Europium-152 | 0.122 to 1.408 | 5×10^4 to 5×10^5 | 11/82 |
| 4201B | Niobium-94 | 0.702 | 5×10^3 | 4/70 |
| Long-Lived Mixed Radionuclide | | | | |
| 4275C | Antimony-125-Tellurium-125m Europium-154 Europium-155 | 0.027 to 1.596 | 5×10^4 6×10^4 3×10^4 | 9/88 |

*License certification is required by NIST for these radionuclides.

Low-Energy-Photon Point-Sources

SRM's 4260C and 4264B consist of a thin-layer deposit of the radionuclide on a thin stainless steel or platinum foil cemented to a monel disk. SRM 4267 has the same construction as the above gamma-ray point sources.

| SRM | Radionuclide | Principal Photon Energy (MeV) | Approx. emission rate at time of calibration (month/year) | Overall Uncertainty (%) |
|-------|-----------------------|-------------------------------|---|-------------------------|
| 4264B | Tin-121m-Antimony-121 | 0.0372 | $5 \times 10^2 \text{s}^{-1}$ | 11/82 |
| 4267 | Niobium-93m | 0.016 | $8 \times 10^2 \text{ Kx s}^{-1}$ | 11/85 |

Radium-226 Solutions

Radon Analysis

These samples are contained in flame-sealed glass ampoules.

| SRM | Nominal Radium Content (g) | (month/year) | Approx. Mass of Solution (g) | Overall Uncertainty (%) |
|-------|----------------------------|--------------|------------------------------|-------------------------|
| 4952B | 8×10^{-15} | 8/76 | 2 | 68 |
| 4953D | 4×10^{-9} | 6/84 | 5 | 1.2 |
| 4950E | 4×10^{-10} | 6/84 | 5 | 1.3 |

Gamma-Ray Solutions

These samples are contained in flame-sealed glass ampoules.

| SRM | Nominal Radium Content (g) | (month/year) | Approx. Mass of Solution (g) | Overall Uncertainty (%) |
|------|-------------------------------|--------------|---------------------------------|-------------------------------|
| 4956 | 2×10^{-7} | 9/67 | 5 | 4.4 |
| 4957 | 5×10^{-7} | 9/67 | 5 | 1.8 |
| 4958 | 1×10^{-6} | 9/67 | 5 | 1.8 |
| 4959 | 2×10^{-6} | 9/67 | 5 | 1.3 |

Environmental Natural Matrix Materials for Quality Assurance Testing

SRM 4350B—Columbia River Sediment

This material was collected from a river downstream from a nuclear reactor facility. Concentrations of fission and activation products are elevated over typical world-wide levels. $^{239}/^{240}\text{Pu}$ and ^{241}Am are very homogeneously distributed through the sample and are in acid-leachable forms. Inhomogeneity is 3 percent or better for other radionuclides.

SRM 4351—Human Lung

This material contains radioactivity concentrations on the order of 10^{-4} Bq g $^{-1}$. It has been freeze-dried, cryogenically ground, homogenized, and packed in a glass bottle under vacuum. There is significant inhomogeneity in $^{239}/^{240}\text{Pu}$ which is unavoidable because plutonium was taken into the lungs in particulate form. Assessments of accuracy of measurement technique can be improved by averaging over several samples.

SRM 4352—Human Liver

This material contains radioactivity concentrations on the order of 10^{-4} Bq g $^{-1}$. It has been freeze-dried, cryogenically ground, homogenized, and packed in a glass bottle under vacuum.

SRM 4353—Rocky Flats Soil Number 1

This material was collected within 13 centimeters of the soil surface at Rocky Flats, CO. $^{239}/^{240}\text{Pu}$ and ^{240}Am concentrations are about an order of magnitude higher than typical world-wide levels. Approximately 10 percent of the plutonium is in an acid-resistant form. The material also contains "hot" particles and a statistical method is provided for dealing with these. Inhomogeneities, excluding hot particles, are on the order of 3 percent or better.

SRM 4354—Freshwater Lake Sediment

This material (gyttja) contains approximately 25 grams of freeze-dried, pulverized freshwater lake sediment (approximately 50 percent organic by weight) in a polyethylene bottle. The SRM is intended for use in tests of measurements of environmental radioactivity contained in matrices similar to the sample, for evaluating analytical methods, or as a generally available calibrated "real" sample matrix in interlaboratory comparisons.

SRM 4355—Peruvian Soil

This material has non-measurable radioactivity concentrations for many fallout radionuclides and can be used as a blank or for sensitive tests of radioanalytical procedures at low-radioactivity concentrations for other radionuclides. The results of a trace-element study are given for 57 elements.

Radiopharmaceuticals

| SRM | Radionuclide (5 mL solution) | Half Life | Approximate Radioactivity at Time of Dispatch (Bq g ⁻¹) | Overall Uncertainty |
|---------|---------------------------------|--------------|--|---------------------------------------|
| 4400LK* | Chromium-51 | 27.702 | d | 3×10^6 |
| 4416LK* | Gallium-67 | 3.261 | d | 3×10^6 |
| 4421L* | Gold-195 | 183 | d | 5×10^5 |
| 4405LB* | Gold-198 | 2.696 | d | 4×10^6 |
| 4417LI* | Indium-111 | 2.805 | d | 5×10^6 |
| 4414LC* | Iodine-123 | 13.221 | h | 6×10^7 |
| 4407LN* | Iodine-125 | 59.6 | d | 1×10^6 |
| 4401LP* | Iodine-131 | 8.021 | d | 5×10^6 |
| 4411LB* | Iron-59 | 44.51 | d | 8×10^5 |
| 4420LB* | Lead-203 | 51.88 | h | 3×10^6 |
| 4418L* | Mercury-203 | 46.60 | d | 1×10^6 |
| 4412LO* | Molybdenum-99-Technetium-99m | 65.92 | h | 2×10^6 |
| 4406LJ* | Phosphorus-32 | 14.29 | d | 2×10^6 |
| 4409LD* | Selenium-75 | 119.8 | d | 1×10^6 |
| 4403LB* | Strontium-85 | 64.854 | d | 1×10^6 |
| 4410HO* | Technetium-99m | 6.007 | h | 1×10^9 |
| 4404LL* | Thallium-201 | 72.91 | h | 4×10^6 |
| 4402LC* | Tin-113-Indium-113m | 115.08 | d | 1×10^6 |
| 4415LN* | Xenon-133 (5 mL gas) | 5.243 | d | 5×10^8 s ⁻¹ total |
| 4419LC* | Ytterbium-169 | 32.03 | d | 2×10^6 |

*License certification is required by NIST for these radionuclides.



This series of SRM's, ranging from glass composition to spectral wavelength to microhardness, represent some of the diversity of the SRM's now available.

Metallurgical

SRM's 485a, 486, 487, and 488 are for calibrating x-ray diffraction equipment used in determining the amount of retained austenite in ferrous materials. SRM 493 is for calibrating x-ray diffraction and Mössbauer equipment to determine the relative amounts of iron carbide in steel.

| SRM | Type | Form |
|------|--|--|
| 485a | Austenite in Ferrite | 5% |
| 486 | Austenite in Ferrite | 15% |
| 487 | Austenite in Ferrite | 30% |
| 488 | Austenite in Ferrite | 2.5% |
| 489 | Ferrite in Austenite | 10% |
| 493 | Spheroidized Iron Carbide (Fe_3C) in Ferrite | Disk: 21 mm dia. \times 2.4 mm thick Disk: 21 mm dia. \times 2.4 mm thick Disk: 21 mm dia. \times 2.4 mm thick Disk: 21 mm dia. \times 2.4 mm thick IN PREP Wafer: 29 \times 29 \times 2.4 mm |

Abrasive Wear

SRM 1857 is for use in the dry sand/rubber wheel abrasion test per ASTM G65, Procedure A.

| SRM | Type | Form |
|------|----------------|--|
| 1857 | D-2 Tool Steel | 2 blocks: 7.8 \times 25 \times 76 mm |

Corrosion

Electrochemical Potential and Thickness

This SRM is for determining the reliability of step test measurements of electrochemical and thickness of multilayered nickel deposits. It consists of a 50 \times 50 mm plate of copper-plated steel over which a duplex nickel coating has been deposited.

| SRM | Type | Step Test Potential (mV) | Nickel Thicknesses | | |
|------|---------------------------|--------------------------|--------------------|--------|------------|
| | | | Total | Bright | Semibright |
| | | | (micrometers) | | |
| 2350 | Nickel Step Test Standard | 110–150 | 27 | (7) | (20) |

Pitting or Crevice Corrosion

These SRM's are for use in evaluating the pitting or crevice corrosion of surgical implant materials per ASTM F746.

| SRM | Type | Form |
|------|--|--------------------------------|
| 1890 | 316L Stainless Steel Rod and Teflon Collar | 4 sets: 6.4 mm D, 25.4 mm long |
| 1891 | Co-Cr-Mo Alloy Rod and Teflon Collar | 2 sets: 6.4 mm D, 25.4 mm long |

X-ray Fluorescent Emission Target

This SRM is intended for use in determining the detector window absorption in semiconductor x-ray spectrometers according to ANSI-IEEE Standard STO 759. When excited by a ^{55}Fe source this glass target will emit fluorescent x rays in the range 1.0 to 5.2 keV.

| SRM | Type | Form | Unit Size |
|-----|---------------------------|------|--------------------|
| 477 | Glass Fluorescence Source | Disk | 2 \times 25 mm D |

X-ray Diffraction

These SRM's are powdered materials to be used as internal standards for powder diffraction measurements. SRM 674a is a set of five oxides for use in the quantitative analysis (intensity measurement) of materials. See also: SRM's 485a-488, 493 (p. 104), and SRM 1878 (p. 62).

| SRM | Type | Lattice Parameter (25.0 °C) | Unit Size |
|--------|---|---|--------------------------------------|
| 640b | Silicon Powder | 5.430940 Å | 7.5 g |
| 656 | Alpha/Beta Standard | IN PREP | |
| 658 | Tridymite | IN PREP | 5 g |
| 659 | Ceramic Particle Size | IN PREP | |
| 660 | LaB ₆ -2Theta X-Ray Profile | 4.15695 Å | 3 g |
| 674a | Powder Diffraction Intensity Al ₂ O ₃ (α -alumina) CeO ₂ Cr ₂ O ₃ TiO ₂ (Rutile) ZnO | 4.75893 Å 5.41129 Å 4.95916 Å 4.59365 Å 3.24981 Å | 10 g 10 g 10 g 10 g 10 g |
| 675 | Powder Diffraction (Mica) | 9.98104 Å | 5 g |
| 676 | Crystalline Alpha Alumina | IN PREP | |
| 677 | Crystalline X-Ray Line Processing | IN PREP | |
| RM8585 | Kaolin Based Relative XRD | IN PREP | 25 g |

Gas Transmission

SRM 1470 is for use in the measurement of gas transmission rates using a volumetric method (ASTM D1434), manometric method (ASTM D1434), or coulometric method (ASTM D3985) of measurement. The permeances of nitrogen, oxygen, carbon dioxide, and helium through this polyester film at 296.15 K are 0.0421, 0.352, 1.722, and 13.79 pmol·s⁻¹·Pa⁻¹, respectively.

| SRM | Type | Unit Size |
|------|---|-------------------------|
| 1470 | Polyester Plastic Film for Gas Transmission | 15 sheets, 23 cm square |



Jim Cline prepares to insert a X-ray diffraction candidate SRM material into a diffractometer.

Reference Fuel

SRM's 1815a and 1816a are high purity liquids intended for use in maintaining the integrity of the octane rating of motor and aviation fuels as specified in the ASTM Manual for Rating Motor, Diesel and Aviation Fuels.

| SRM | Type | Purity, % | Unit Size |
|-------|------------------------------------|-----------|-----------|
| 1815a | n-Heptane | 99.987 | 100 mL |
| 1816a | Isooctane (2,2,4-Trimethylpentane) | 99.987 | 100 mL |

Electrical Resistivity and Conductivity

Metals

These materials are for evaluating methods of measuring electrical resistance over wide temperature ranges.

| SRM | Type | Temperature Range | Resistivity at 293 K | Form |
|------|-----------------|-------------------|--------------------------------|----------------------------|
| 1461 | Stainless Steel | 5 to 1200 K | 80.5 $\mu\Omega\cdot\text{cm}$ | Rod: 12.7 mm D, 50 mm long |
| 1462 | Stainless Steel | 5 to 1200 K | 80.5 $\mu\Omega\cdot\text{cm}$ | Rod: 34.0 mm D, 50 mm long |
| 8420 | Iron | 6 to 1000 K | 10.1 $\mu\Omega\cdot\text{cm}$ | Rod: 6.4 mm D, 50 mm long |
| 8421 | Iron | 6 to 1000 K | 10.1 $\mu\Omega\cdot\text{cm}$ | Rod: 31.7 mm D, 50 mm long |
| 8422 | Tungsten | 4 to 3000 K | 5.4 $\mu\Omega\cdot\text{cm}$ | Rod: 3.2 mm D, 50 mm long |
| 8423 | Tungsten | 4 to 3000 K | 5.4 $\mu\Omega\cdot\text{cm}$ | Rod: 6.4 mm D, 50 mm long |

Silicon

SRM's 1521, 1522, and 1523 sets of Resistivity Standards will be replaced with SRM 2541-50 series of Resistivity Standards as they become available. This series will be 100 mm in diameter and will have improved accuracy and precision. SRM's 2526, 2527, 2528, and 2529 are mounted on beveling blocks for two-probe test equipment.

| SRM | Type | Resistivity | Form |
|------|--|-------------------------------------|--|
| 2526 | 111 p-Type Silicon, Spreading Resistance | 0.001 to 200 $\Omega\cdot\text{cm}$ | 16 levels, 5 \times 10 \times 0.625 mm |
| 2527 | 111 n-Type Silicon, Spreading Resistance | 0.001 to 200 $\Omega\cdot\text{cm}$ | 16 levels, 5 \times 10 \times 0.625 mm |
| 2528 | 100 p-Type Silicon, Spreading Resistance | 0.001 to 200 $\Omega\cdot\text{cm}$ | 16 levels, 5 \times 10 \times 0.625 mm |
| 2529 | 100 n-Type Silicon, Spreading Resistance | 0.001 to 200 $\Omega\cdot\text{cm}$ | 16 levels, 5 \times 10 \times 0.625 mm |
| 2541 | Silicon Resistivity | | 4" wafers |
| 2542 | Silicon Resistivity | | 4" wafers |
| 2543 | Silicon Resistivity | | 4" wafers |
| 2544 | Silicon Resistivity | | 4" wafers |
| 2546 | Silicon Resistivity | | 4" wafers |
| 2547 | Silicon Resistivity | | 4" wafers |
| 2548 | Silicon Resistivity | | 4" wafers |
| 2549 | Silicon Resistivity | | 4" wafers |
| 2550 | Silicon Resistivity | | 4" wafers |

Residual Resistivity Ratio

This SRM is a set of five aluminum rods that are intended for use in checking four-terminal dc and eddy current decay techniques. The residual resistivity ratio, $\rho(273\text{ K})/\rho(4\text{ K})$, is a sensitive indicator of purity and of the mechanical state of a material.

| SRM | Type | RRR Values | Form |
|-----|----------|----------------------------------|----------------------|
| 769 | Aluminum | 130, 683, 1205, 2650, and 11,000 | 6.4 mm D, 52 mm long |

Eddy Current

These SRM's are intended for use in the calibration of eddy current conductivity meters and of secondary electrical conductivity standards. Eddy current measurements are used in nondestructive inspection of conducting materials and in the sorting of alloys for composition and heat treatment.

| SRM | Type | Conductivity | Form |
|------|--------------------------|--------------|------------------------------|
| 1860 | Aluminum | 60% IACS | $44 \times 44 \times 9.5$ mm |
| 1862 | Aluminum-Magnesium Alloy | 41% IACS | $44 \times 44 \times 9.5$ mm |

Electrolytic Conductance

These SRM's are for calibrating and standardizing conductivity cells and meters used in water purity determinations. They are solutions of high-purity potassium chloride in de-ionized water in equilibrium with atmospheric carbon dioxide.


| SRM | Type | Nominal Conductance (microsiemens per cm, $\mu\text{S}/\text{cm}$) | Unit Size |
|------|--------------------------|--|-----------|
| 3191 | Electrolytic Conductance | 100 | 500 mL |
| 3192 | Electrolytic Conductance | 500 | 500 mL |
| 3193 | Electrolytic Conductance | 1000 | 500 mL |
| 3194 | Electrolytic Conductance | 10000 | 500 mL |
| 3195 | Electrolytic Conductance | 100000 | 500 mL |

Superconducting Critical Current

This SRM is for checking the performance of measurement systems used in superconductor technology. It consists of 2.2 m of a multifilamentary niobium titanium, copper stabilized superconducting wire wound in a single layer onto a spool with a core diameter of 8.7 cm.

| SRM | Type | Magnetic Field (T) | Critical Current (A) |
|------|------------|--------------------|----------------------|
| 1457 | Nb-Ti Wire | 2.000 | 293.30 |
| | | 4.000 | 187.38 |
| | | 6.000 | 124.72 |
| | | 8.000 | 69.72 |

Ambler Thompson is shown here making adjustments on the high accuracy spectrofluorimeter before measuring the relative spectral emission of a fluorescent sample.



Engineering Materials

Standard Rubbers and Rubber-Compounding Materials

These SRM's have been prepared to provide the rubber industry with standard materials for rubber compounding. They are useful for the testing of rubber and rubber-compounding materials in connection with quality control of raw materials and for the standardization of rubber testing.

Each material has been statistically evaluated for uniformity by mixing rubber and rubber compounds, and vulcanizing them in accordance with ASTM Designation D-15 and determining the stress-strain properties of the resulting vulcanizates. Certificates are issued for the rubbers because the properties of different lots are not the same. Replacement lots of rubber-compounding SRM's impart essentially the same characteristics to rubber vulcanizates so that Certificates are not issued for these SRM's.

Rubbers

| SRM | Type | Wt/Unit | Pounds |
|------|------------------------|---------|--------|
| 386k | Styrene-butadiene 1500 | 34 kg | 75 |
| 388p | Butyl | 34 kg | 75 |
| 1495 | Butyl (Low Viscosity) | 34 kg | 75 |

Rubber Compounding Materials

| SRM | Type | Wt/Unit | Pounds |
|------|---|---------|--------|
| 370e | Zinc Oxide | 8 kg | 17.6 |
| 371h | Sulfur | 6 kg | 13.2 |
| 372i | Stearic Acid | 3.2 kg | 7.1 |
| 375g | Channel Black | 28 kg | 61.6 |
| 378b | Oil Furnace Black | 28 kg | 61.6 |
| 382a | Gas Furnace Black | 32 kg | 70.6 |
| 383a | Mercaptobenzothiazole | 3.2 kg | 7.1 |
| 384f | N-tertiary-Butyl-2-benzothiazolesulfenamide | 3.2 kg | 7.1 |

Sizing

Particle Size

SRM's 1003a, 1690, 1691, and 1960 can be used to calibrate various types of particle size measuring instruments including both light and electrical zone flow-through counters. SRM's 1004a, 1017a, 1018a, and 1019a are for calibrating test sieves.

| SRM | Type | Size (μm) | Sieve No. | Wt/Unit |
|-------|---|------------------------|-----------|-----------|
| 1003a | Glass Spheres | 8–58 | — | 25 g |
| 1004a | Glass Spheres | IN PREP | — | |
| 1017a | Glass Spheres | 100–310 | 140–50 | 84 g |
| 1018a | Glass Spheres | 225–780 | 60–25 | 74 g |
| 1019a | Glass Spheres | 760–2160 | 20–10 | 200 g |
| 1690 | Polystyrene Spheres (0.5% wt. concentration in water) | 0.895 | — | 5 mL vial |
| 1691 | Polystyrene Spheres (0.5% wt. concentration in water) | 0.269 | — | 5 mL vial |
| 1960 | Polystyrene Spheres (0.4% wt. concentration in water) | 9.89 | — | 5 mL vial |
| 1961 | Polystyrene Spheres (0.5% wt. concentration in water) | 29.64 | — | 5 mL vial |
| 1962 | Polystyrene Spheres (0.5% wt. concentration in water) | 2.978 | — | 5 mL vial |
| 1965 | Polystyrene Spheres (0.5% wt. concentration in water) | 9.94 | — | 1 slide |

Cement Turbidimetric and Fineness

This SRM is available to calibrate the Blaine fineness meter according to the latest issue of Federal Test Method Standard 158, Method 2101 or ASTM Designation C204; to calibrate the Wagner turbidimeter according to ASTM Designation C115; and to determine sieve residue according to ASTM Designation C430. Each set consists of twenty sealed vials, each containing approximately 10 grams of cement.

| SRM | Type | Properties Certified | Value | Unit |
|------|-----------------|---|----------------------------------|-----------------|
| 114n | Portland Cement | Residue on 45 μm , electroformed sieve wet method | 8.3% | Set of 20 vials |
| | | Surface area (Wagner turbidimeter) | $2020 \text{ cm}^2\text{g}^{-1}$ | |
| | | Surface area (Air-permeability) | $3460 \text{ cm}^2\text{g}^{-1}$ | |

Surface Area of Powders

These materials are for calibrating and checking instruments used to determine the specific surface area of powders by BET. RM's 8005 through 8008 have been certified by the National Physical Laboratory, Teddington, U.K. (and meet the ISO definition for CRM's); RM's 8570, 8571, and 8572 are issued by NBS in cooperation with ASTM, but are not certified.

| RM | Type | Surface Area | Unit Size |
|------|--------------------------------|-------------------------|-----------|
| 8005 | Alpha Alumina (NPL# M11-05/09) | 2.1 m ² /g | 50 g |
| 8006 | Alpha Alumina (NPL# M11-06/10) | 0.3 m ² /g | 50 g |
| 8007 | Alpha Alumina (NPL# M11-07/11) | 0.1 m ² /g | 50 g |
| 8008 | Alpha Alumina (NPL# M11-08/12) | 0.8 m ² /g | 50 g |
| 8570 | Calcined Kaolin | 10.89 m ² /g | 10 g |
| 8571 | Alumina | 158.7 m ² /g | 10 g |
| 8572 | Silica-Alumina | 291.2 m ² /g | 10 g |

PERFORMANCE STANDARDS

Socketed Ball Bar

This SRM is for measuring the performance of coordinate measuring machines (CMM's) as per ASME Standard B89.1.12. It consists of a set of three precision balls pinned and cemented onto threaded shafts, one table-mount magnetic socket, one ram-mount magnetic socket, and 5 partially insulated extension tubes—50, 100, 200, 400, and 800 mm long.

| SRM | Type | Measuring Lengths (50 mm steps) | Unit |
|------|-------------------|------------------------------------|------|
| 2083 | Socketed Ball Bar | 100 to 1650 mm | Set |

Dye Penetrant Test Blocks

These SRM's are for checking the performance of liquid dye penetrants and dye penetrant crack detection techniques. These test blocks have four synthetic cracks, approximately 0.2, 0.5, 1, and 2 µm wide.

| SRM | Type | Surface | Unit Size |
|------|--------------------------|---------------|-----------------------|
| 1850 | Penetrant Test Block | Bright Finish | 5 cm dia., 1 cm thick |
| 1851 | NDE Penetrant Test Block | Matte Finish | 5 cm dia., 1 cm thick |

Surface Roughness

These SRM's are for calibrating stylus instruments that measure surface roughness. These electroless-nickel coated steel blocks have a sinusoidal roughness profile machined on the top surface.

| SRM | Type | Roughness | Wavelength | Unit of Issue |
|------|----------------------|-------------------|-------------------|-----------------|
| 2071 | Sinusoidal Roughness | 0.3 μm | 10 μm | Block, 24×33 mm |
| 2072 | Sinusoidal Roughness | 1.0 μm | 100 μm | Block, 24×33 mm |
| 2073 | Sinusoidal Roughness | 3.0 μm | 100 μm | Block, 24×33 mm |
| 2074 | Sinusoidal Roughness | 1.0 μm | 40 μm | IN PREP |
| 2075 | Sinusoidal Roughness | 1.0 μm | 800 μm | IN PREP |

Charpy V-Notch Test Blocks

These SRM's are test specimens for standardizing the Charpy impact test and for certifying Charpy impact test machines. After the bars are broken in the user's laboratory, they are to be returned to NIST for examination and certification of the test machine.

| SRM | Type | Range | Unit of Issue |
|------|----------------------------|-------------|---------------|
| 2092 | Charpy V-Notch Low Energy | 11–18 ft-lb | 5 bars |
| 2096 | Charpy V-Notch High Energy | 65–80 ft-lb | 5 bars |

Artificial Flaw for Eddy Current NDE

RM 8458 provides a flaw of known size and geometry that closely resembles an actual fatigue crack. It is intended to produce a response suitable for calibrating eddy current nondestructive evaluation (NDE) systems. The flaw size is 3.0 ± 0.1 mm long by 1.0 mm deep in a $7 \times 7 \times 2$ cm block of 7075-T651 aluminum alloy, heat treated to the T6 temper.

Color

These SRM's are available to illustrate a characteristic color for each of the ISCC-NBS color-name blocks in NBS Special Publication 440, COLOR: Universal Language and Dictionary of Names. SRM 2106 consists of 251 color chips on 18 constant-hue centroid color charts, and constitutes a supplement to SP 440. The centroid colors represent a systematic sampling of the whole color solid. Note: The color chips were re-measured in 1984 and are issued with the new data as an addendum. This addendum is available upon request.

| SRM | Type | Unit of Issue |
|------|-----------------------|----------------|
| 2106 | Centroid Color Charts | Set: 18 Charts |

X-ray and Photographic

SRM 1001 is a calibrated x-ray film step tablet of 17 steps that cover the optical density range from 0 to 4; it has a blue tint and emulsion on both sides. SRM 1008 is a calibrated photographic step tablet of 21 steps that cover the optical density range from 0 to 4; it has a black tint and emulsion on a single side.

SRM 1010a, Microcopy Resolution Test Charts, is used to test the resolving power of cameras or of whole microcopying systems. SRM 1010a consists of five charts printed photographically on paper, which have 26 high-contrast five-line patterns ranging in spatial frequency from one cycle per millimeter to 18 cycles per millimeter. Instructions for the use of the charts are supplied with each order.

| SRM | Type | Unit |
|-------|---------------------------------|--------------------|
| 1001 | X-ray Film Step Tablet (0-4) | 1 tablet, 17 steps |
| 1008 | Photographic Step Tablet (0-4) | 1 tablet, 21 steps |
| 1010a | Microcopy Resolution Test Chart | Set of 5 charts |

Magnetic Computer Storage Media

These SRM's are for evaluating the performance of magnetic computer storage media and systems, and for maintaining control over their production. Each SRM is individually calibrated and certified.

| SRM | Description | Unit of Issue |
|------|--|---------------|
| 3200 | Secondary Standard Magnetic Tape—12.7 mm (1/2 in) wide tape, certified for signal amplitude outputs relative to the NIST Standard Reference Amplitudes at 8, 32, and 126 flux transitions per millimeter (200, 800, 3200 flux transitions per inch). | Open Reel |
| 6250 | Secondary Standard High Density Magnetic Tape—12.7 mm (1/2 in) wide tape, certified for signal amplitude output relative to the NIST Standard Reference Amplitude at 356 flux transitions per millimeter (9042 flux transitions per inch). | Open Reel |

Magnetic Computer Storage Media (Continued)

| SRM | Description | Unit of Issue |
|------|---|---------------|
| 1600 | Secondary Standard Magnetic Tape Cassette—3.8 mm (0.15 in) wide tape, certified for signal amplitude output relative to the NIST Standard Reference Amplitude at 63 flux transitions per millimeter (1600 flux transitions per inch). | Cassette |
| 3216 | Secondary Standard Magnetic Tape Cartridge—6.3 mm ($\frac{1}{4}$ in) wide tape, certified for signal amplitude output relative to the NIST Standard Reference Amplitude at 126 flux transitions per millimeter (3200 flux transitions per inch). | Cartridge |
| 3217 | Secondary Standard High Density Magnetic Tape Cartridge—6.3 mm ($\frac{1}{4}$ in) wide tape, certified for signal amplitude outputs relative to the NIST Standard Reference Amplitudes at 252 and 394 flux transitions per millimeter (6400 and 10,000 flux transitions per inch). | Cartridge |

These RM's are certified by the Physikalisch-Technische Bundesanstalt (PTB), Federal Republic of Germany, for signal amplitude, overwrite, and resolution. The RM numbers correspond to the ISO standard number, and the materials conform to relevant ANSI, ISO, and ECMA standards for flexible disk cartridges.

| RM | Description | Unformatted Capacity | Unit/Size |
|------|-------------------------|----------------------|------------------|
| 6596 | Flexible Disk Cartridge | 125 K bytes | 130 mm (5.25 in) |
| 9529 | Flexible Disk Cartridge | 2000 K bytes | 90 mm (3.5 in) |

Centerline Drawings for Optical Character Recognition Style-B Characters

This SRM is an exact copy of the centerline drawings that uniquely define each printed character shape and size used in constant strokewidth Style B Size I Optical Character Recognition (OCR-B) applications in accordance with one or more of the following standards: American National Standard X3.49-1975 (R 1982), Character Set for Optical Character Recognition (OCR-B); Federal Information Processing Standards Publication (OCR), European Computer Manufacturers Association Standard ECMA-11 for the Alphanumeric Character Set OCR-B for Optical Recognition, 3rd Edition, 1976 and International Standard ISO 1073/II-1976, Alphanumeric Character Sets for Optical Recognition Part II: Character Set OCR-B.

This Standard Reference Material contains information on the shape, size, strokewidth, and position relative to the base line of the OCR-B characters.

| SRM | Characters | Sheets | Size | Sheet Size |
|------|------------|--------|---------|---------------|
| 1901 | 118 | 118 | OCR-B I | 32×44×0.01 cm |

NIST Time Software

The NIST automated time service is a telephone time service designed to provide computer with telephone access to NIST time at occurrences approaching one millisecond (0.001 second). Features of the service include automated compensation for telephone-line delay, advanced alert for changes to and from daylight savings time and advanced notice of insertion of leap seconds. The ASC11-character time code should operate with standard modems and most computer systems. While the system can be used to set up computer time-of-day, simple hardware can be developed to set non-computer clock system.

| SRM | Description | Unit of Issue |
|------|--|---------------|
| 8101 | Automated Computer Time Service (Acts) | Diskette |



Cassandra Beck responds to one of the many inquiries or orders she handles for radioactivity SRM's.

FIRE RESEARCH

Surface Flammability

SRM 1002d, Hardboard Sheet, is issued for checking the operation of radiant panel test equipment in accordance with the procedures outlined in ASTM Standard E162-78.

| SRM | Type | Certification | Unit of Issue |
|-------|-----------------|--|------------------------------------|
| 1002d | Hardboard Sheet | Flame Spread Index, I=153 Heat Evolution Factor, Q=36.5 | Set of 4: 6×18× $\frac{1}{4}$ inch |

Smoke Density Chamber

These SRM's are certified for maximum specific optical density and are issued for performing operational checks of smoke density chambers.

| SRM | Type | Maximum Specific Optical Density | Unit of Issue |
|-------|---|-----------------------------------|---------------|
| 1006c | Non-flaming Exposure Condition (α -cellulose) | D _m (corr.)=178 | 6 sheets |
| 1007a | Flaming Exposure Condition (plastic) | D _m (corr.)=421 to 493 | 3 sheets |
| 1048 | Smoke Toxicity-Comp. ABS | | 7 sheets |

Flooring Radiant Panel

This SRM consists of three sheets of kraft paperboard. It is for checking the operation of flooring radiant panel apparatus used to measure critical radiant flux as per ASTM E648.

| SRM | Type | Critical Radiant Flux | Unit Size (cm) |
|------|------------------------|------------------------|------------------|
| 1012 | Flooring Radiant Panel | 0.36 W/cm ² | 104.1×25.4×0.305 |

Tape Adhesion Testing

This material is intended as a uniform source of linerboard for use under ASTM Designation D2860, Procedure A: Adhesion of Pressure Sensitive Tape to Fiberboard at 90 Degree Angle and Constant Stress.

| SRM | Type | Unit |
|-------|--------------------------------------|---------|
| 1810a | Linerboard for Tape Adhesion Testing | IN PREP |

Additional Information

NIST SPECIAL PUBLICATIONS IN THE 260 SERIES

Barber, S. L., Hines, J. K., editors, NIST Standard Reference Materials Catalog (1990-91 edition), NIST Spec. Publ. 260 (January 1990)—SN003-003-00000-0.

Michaelis, R. E., and Wyman, L. L., Standard Reference Materials: Preparation of White Cast Iron Spectrochemical Standards, NBS Misc. Publ. 260-1 (June 1964). COM74-11061**

Michaelis, R. E., Wyman, L. L., and Flitsch, R., Standard Reference Materials: Preparation of NBS Copper-Base Spectrochemical Standards, NBS Misc. Publ. 260-2 (October 1964). COM74-11063**

Michaelis, R. E., Yakowitz, H., and Moore, G. A., Standard Reference Materials: Metallographic Characterization of an NBS Spectrometric Low-Alloy Steel Standard, NBS Misc. Publ. 260-3 (October 1964). COM74-11060**

Alvarez, R., and Flitsch, R., Standard Reference Materials: Accuracy of Solution X-Ray Spectrometric Analysis of Copper-Base Alloys, NBS Misc. Publ. 260-5 (March 1965). PB168068**

Shultz, J. I., Standard Reference Materials: Methods for the Chemical Analysis of White Cast Iron Standards, NBS Misc. Publ. 260-6 (July 1965). COM74-11068**

Bell, R. K., Standard Reference Materials: Methods for the Chemical Analysis of NBS Copper-Base Spectrochemical Standards, NBS Misc. Publ. 260-7 (October 1965). COM74-11067**

* Send order with remittance to Superintendent of Documents, US Government Printing Office, Washington, DC 20402. Remittance from foreign countries should include an additional one-fourth of the purchase price for postage.

** May be ordered from: National Technical Information Services (NTIS), Springfield, Virginia 22161.

Richmond, M.S., Standard Reference Materials: Analysis of Uranium Concentrates at the National Bureau of Standards, NBS Misc. Publ. 260-8 (December 1965). COM74-11066**

Anspach, S. C., Cavallo, L. M., Garfinkel, S. B., Hutchinson, J. M. R., and Smith, C. N., Standard Reference Materials: Half Lives of Materials Used in the Preparation of Standard Reference Materials of Nineteen Radioactive Nuclides Issued by the National Bureau of Standards, NBS Misc. Publ. 260-9 (November 1965). COM74-11065**

Yakowitz, H., Vieth, D. L., Heinrich, K. F. J., and Michaelis, R. E., Standard Reference Materials: Homogeneity Characterization of NBS Spectrometric Standards II: Cartridge Brass and Low-Alloy Steel, NBS Misc. Publ. 260-10 (December 1965). COM74-11064**

Yakowitz, H., Michaelis, R. E., and Vieth, D. L., Standard Reference Materials: Homogeneity Characterization of NBS Spectrometric Standards IV: Preparation and Microprobe Characterization of W-20% Mo Alloy Fabricated by Powder Metallurgical Methods, NBS Spec. Publ. 260-16 (January 1969). COM74-11062**

Paule, R. C., and Mandel, J., Standard Reference Materials: Analysis of Interlaboratory Measurements on the Vapor Pressure of Gold (Certification of Standard Reference Material 745), NBS Spec. Publ. 260-19 (January 1970). PB190071**

Paule, R. C., and Mandel, J., Standard Reference Materials: Analysis of Interlaboratory Measurements on the Vapor Pressures of Cadmium and Silver, NBS Spec. Publ. 260-21 (January 1971). COM74-11359**

Yakowitz, H., Fiori, C. E., and Michaelis, R. E., Standard Reference Materials: Homogeneity Characterization of Fe-3 Si Alloy, NBS Spec. Publ. 260-22 (February 1971). COM74-11357**

Napolitano, A., and Hawkins, E. G., Standard Reference Materials: Viscosity of a Standard Borosilicate Glass, NBS Spec. Publ. 260-23 (December 1970). COM71-00157**

Sappenfield, K. M., Marinenko, G., and Hague, J. L., Standard Reference Materials: Comparison of Redox Standards, NBS Spec. Publ. 260-24 (January 1972). COM72-50058**

Hicho, G. E., Yakowitz, H., Rasberry, S. D., and Michaelis, R. E., Standard Reference Materials: A Standard Reference Material Containing Nominally Four Percent Austenite, NBS Spec. Publ. 260-25 (February 1971). COM74-11356**

Martin, J. F., Standard Reference Materials: National Bureau of Standards-US Steel Corporation Joint Program for Determining Oxygen and Nitrogen in Steel, NBS Spec. Publ. 260-26 (February 1971). 85 cents* PB 81176620

Garner, E. L., Machlan, L. A., and Shields, W. R., Standard Reference Materials: Uranium Isotopic Standard Reference Materials, NBS Spec. Publ. 260-27 (April 1971). COM74-11358**

Heinrich, K. F. J., Myklebust, R. L., Rasberry, S. D., and Michaelis, R. E., Standard Reference Materials: Preparation and Evaluation of SRM's 481 and 482 Gold-Silver and Gold-Copper Alloys for Microanalysis, NBS Spec. Publ. 260-28 (August 1971). COM71-50365**

Geller, S. B., Standard Reference Materials: Calibration of NBS Secondary Standard Magnetic Tape (Computer Amplitude Reference) Using the Reference Tape Amplitude Measurement "Process A-Model 2," NBS Spec. Publ. 260-29 (June 1971). COM71-50282**

Gorozhanina, R. S., Freedman, A. Y., and Shaievitch, A. B. (translated by M. C. Selby), Standard Reference Materials: Standard Samples Issued in the USSR (A Translation from the Russian), NBS Spec. Publ. 260-30 (June 1971). COM71-50283**

Hust, J. G., and Sparks, L. L., Standard Reference Materials: Thermal Conductivity of Electrolytic Iron SRM 734 from 4 to 300 K, NBS Spec. Publ. 260-31 (November 1971). COM71-50563**

Mavrodineanu, R., and Lazar, J. W., Standard Reference Materials: Standard Quartz Cuvettes for High Accuracy Spectrophotometry, NBS Spec. Publ. 260-32 (December 1973). COM74-50018**

Wagner, H. L., Standard Reference Materials: Comparison of Original and Supplemental SRM 705, Narrow Molecular Weight Distribution Polystyrene, NBS Spec. Publ. 260-33 (May 1972). COM72-50526**

Sparks, L. L., and Hust, J. G., Standard Reference Materials: Thermal Conductivity of Austenitic Stainless Steel, SRM 735 from 5 to 280 K, NBS Spec. Publ. 260-35 (April 1972). COM72-50368**

Cali, J. P., Mandel, J., Moore, L. J., and Young, D. S., Standard Reference Materials: A Referee Method for the Determination of Calcium in Serum, NBS SRM 915, NBS Spec. Publ. 260-36 (May 1972). COM72-50527**

Shultz, J. I., Bell, R. K., Rains, T. C., and Menis, O., Standard Reference Materials: Methods of Analysis of NBS Clay Standards, NBS Spec. Publ. 260-37 (June 1972). COM72-50692**

Richmond, J. C., and Hsia, J. J., Standard Reference Materials: Preparation and Calibration of Standards of Spectral Specular Reflectance, NBS Spec. Publ. 260-38 (May 1972). COM72-50528**

Clark, A. F., Denson, V. A., Hust, J. G., and Powell, R. L., Standard Reference Materials: The Eddy Current Decay Method for Resistivity Characterization of High-Purity Metals, NBS Spec. Publ. 260-39 (May 1972). COM72-50529**

McAdie, H. G., Garn, P. D., and Menis, O., Standard Reference Materials: Selection of Thermal Analysis Temperature Standards Through a Cooperative Study (SRM 758, 759, 760), NBS Spec. Publ. 260-40 (August 1972). COM72-50776**

Wagner, H. L., and Verdier, P. H., eds., Standard Reference Materials: The Characterization of Linear Polyethylene, SRM 1475, NBS Spec. Publ. 260-42 (September 1972). COM72-50944**

Yakowitz, H., Ruff, A. W., and Michaelis, R. E., Standard Reference Materials: Preparation and Homogeneity Characterization of an Austenitic Iron-Chromium-Nickel Alloy, NBS Spec. Publ. 260-43 (November 1972). COM73-50760**

Schooley, J. F., Soulen, R. J., Jr., and Evans, G. A., Jr., Standard Reference Materials: Preparation and Use of Superconductive Fixed Point Devices, SRM 767, NBS Spec. Publ. 260-44 (December 1972). COM73-50037**

Greifer, B., Maienthal, E. J., Rains, T. C., and Rasberry, S. D., Standard Reference Materials: Powdered Lead-Based Paint, SRM 1579, NBS Spec. Publ. 260-45 (March 1973). COM73-50226**

Hust, J. G., and Giarratano, P. J., Standard Reference Materials: Thermal Conductivity and Electrical Resistivity Standard Reference Materials: Austenitic Stainless Steel, SRM's 735 and 798, from 4 to 1200 K, NBS Spec. Publ. 260-46 (March 1975). COM75-10339**

Hust, J. G., Standard Reference Materials: Electrical Resistivity of Electrolytic Iron, SRM 797, and Austenitic Stainless Steel, SRM 798, from 5 to 280 K, NBS Spec. Publ. 260-47 (February 1974). COM74-50176**

Mangum, B. W., and Wise, J. A., Standard Reference Materials: Description and Use of Precision Thermometers for the Clinical Laboratory, SRM 933 and SRM 934, NBS Spec. Publ. 260-48 (May 1974). COM74-50533**

Carpenter, B. S., and Reimer, G. M., Standard Reference Materials: Calibrated Glass Standards for Fission Track Use, NBS Spec. Publ. 260-49 (November 1974). COM74-51185**

Hust, J. G., and Giarratano, P. J., Standard Reference Materials: Thermal Conductivity and Electrical Resistivity Standard Reference Materials: Electrolytic Iron, SRM's 734 and 797 from 4 to 1000 K, NBS Spec. Publ. 260-50 (June 1975). COM75-10698**

Mavrodineanu, R., and Baldwin, J. R., Standard Reference Materials: Glass Filters as a Standard Reference Material for Spectrophotometry—Selection, Preparation, Certification, Use—SRM 930, NBS Spec. Publ. 260-51 (November 1975). COM75-10339**

Hust, J. G., and Giarratano, P. J., Standard Reference Materials: Thermal Conductivity and Electrical Resistivity Standard Reference Materials 730 and 799, from 4 to 3000 K, NBS Spec. Publ. 260-52 (September 1975). COM75-11193**

Durst, R. A., Standard Reference Materials: Standardization of pH Measurements, NBS Spec. Publ. 260-53 (December 1975, Revised). PB248127**

Burke, R. W., and Mavrodineanu, R., Standard Reference Materials: Certification and Use of Acidic Potassium Dichromate Solutions as an Ultraviolet Absorbance Standard, NBS Spec. Publ. 260-54 (August 1977). PB272168**

Ditmars, D. A., Cezairliyan, A., Ishihara, S., and Douglas, T. B., Standard Reference Materials: Enthalpy and Heat Capacity; Molybdenum SRM 781, from 273 to 2800 K, NBS Spec. Publ. 260-55 (September 1977). PB272127**

Powell, R. L., Sparks, L. L., and Hust, J. G., Standard Reference Materials: Standard Thermocouple Materials, Pt 67: SRM 1967, NBS Spec. Publ. 260-56 (February 1978). PB277172**

Barnes, J. D., and Martin, G. M., Standard Reference Materials: Polyester Film for Oxygen Gas Transmission Measurements SRM 1470, NBS Spec. Publ. 260-58 (June 1979). PB297098**

Velapoldi, R. A., Paule, R. C., Schaffer, R., Mandel, J., and Moody, J. R., Standard Reference Materials: A Reference Method for the Determination of Sodium in Serum, NBS Spec. Publ. 260-60 (August 1978). PB286944**

Verdier, P. H., and Wagner, H. L., Standard Reference Materials: The Characterization of Linear Polyethylene (SRM 1482, 1483, 1484), NBS Spec. Publ. 260-61 (December 1978). PB289899**

- Soulen, R. J., and Dove, R. B., Standard Reference Materials: Temperature Reference Standard for Use Below 0.5 K (SRM 768), NBS Spec. Publ. 260-62 (April 1979). PB294245****
- Velapoldi, R. A., Paule, R. C., Schaffer, R., Mandel, J., Machlan, L. A., and Gramlich, J. W., Standard Reference Materials: A Reference Method for the Determination of Potassium in Serum, NBS Spec. Publ. 260-63 (May 1979). PB297207****
- Velapoldi, R. A., and Mielenz, K. D., Standard Reference Materials: A Fluorescence Standard Reference Material Quinine Sulfate Dihydrate (SRM 936), NBS Spec. Publ. 260-64 (January 1980). PB-80-132046****
- Marinenko, R. B., Heinrich, K. F. J., and Ruegg, F. C., Standard Reference Materials: Micro-Homogeneity Studies of NBS Standard Reference Materials, NBS Research Materials, and Other Related Samples, NBS Spec. Publ. 260-65 (September 1979). PB300461****
- Venable, W. H., Jr., and Eckerle, K. L., Standard Reference Materials: Didymium Glass Filters for Calibrating the Wavelength Scale of Spectrophotometers—SRM 2009, 2010, 2013 and 2014, NBS Spec. Publ. 260-66 (October 1979). PB-80-104961****
- Velapoldi, R. A., Paule, R. C., Schaffer, R., Mandel, J., Murphy, T. J., and Gramlich, J. W., Standard Reference Materials: A Reference Method for the Determination of Chloride in Serum, NBS Spec. Publ. 260-67 (November 1979). PB-80-110117****
- Mavrodineanu, R., and Baldwin, J. R., Standard Reference Materials: Metal-On-Quartz Filters as a Standard Reference Material for Spectrophotometry—SRM 2031, NBS Spec. Publ. 260-68 (April 1980). PB-80-197486****
- Velapoldi, R. A., Paule, R. C., Schaffer, R., Mandel, J., Machlan, L. A., Garner, E. L., and Rains, T. C., Standard Reference Materials: A Reference Method for the Determination of Lithium in Serum, NBS Spec. Publ. 260-69 (July 1980). PB-80-20917****
- Marinenko, R. B., Biancaniello, F., Boyer, P. A., Ruff, A. W., and DeRobertis, L., Standard Reference Materials: Preparation and Characterization of an Iron-Chromium-Nickel Alloy for Microanalysis, NBS Spec. Publ. 260-70 (May 1981). PB-84-165349****
- Seward, R. W., and Mavrodineanu, R., Standard Reference Materials: Summary of the Clinical Laboratory Standards Issued by the National Bureau of Standards, NBS Spec. Publ. 260-71 (November 1981). PB-82-135161****
- Reeder, D. J., Coxon, B., Enagonio, D., Christensen, R. G., Schaffer, R., Howell, B. F., Paule, R. C., and Mandel, J., Standard Reference Materials: SRM 900, Antiepilepsy Drug Level Assay Standard, NBS Spec. Publ. 260-72 (June 1981). PB-81-220758****
- Interrante, C. G., and Hicho, G. E., Standard Reference Materials: A Standard Reference Material Containing Nominally Fifteen Percent Austenite (SRM 486), NBS Spec. Publ. 260-73 (January 1982). PB-82-215559****
- Marinenko, R. B., Standard Reference Materials: Preparation and Characterization of K-411 and K-414 Mineral Glasses for Microanalysis: SRM 470, NBS Spec. Publ. 260-74 (April 1982). PB-82-221300****
- Weidner, V. R., and Hsia, J. J., Standard Reference Materials: Preparation and Calibration of First Surface Aluminum Mirror Specular Reflectance Standards (SRM 2003a), NBS Spec. Publ. 260-75 (May 1982). PB-82-221367****
- Hicho, G. E., and Eaton, E. E., Standard Reference Materials: A Standard Reference Material Containing Nominally Five Percent Austenite (SRM 485a), NBS Spec. Publ. 260-76 (August 1982). PB-83-115568****
- Furukawa, G. T., Riddle, J. L., Bigge, W. G., and Pfieffer, E. R., Standard Reference Materials: Application of Some Metal SRM's as Thermometric Fixed Points, NBS Spec. Publ. 260-77 (August 1982). PB-83-117325****
- Hicho, G. E., and Eaton, E. E., Standard Reference Materials: Standard Reference Material Containing Nominally Thirty Percent Austenite (SRM 487), NBS Spec. Publ. 260-78 (September 1982). PB-83-115576****
- Richmond, J. C., Hsia, J. J., Weidner, V. R., and Wilmering, D. B., Standard Reference Materials: Second Surface Mirror Standards of Specular Spectral Reflectance (SRM's 2023, 2024, 2025), NBS Spec. Publ. 260-79 (October 1982). PB-84-203447****
- Schaffer, R., Mandel, J., Sun, T., Cohen, A., and Hertz, H. S., Standard Reference Materials: Evaluation by an ID/MS Method of the AACC Reference Method for Serum Glucose, NBS Spec. Publ. 260-80 (October 1982). PB-84-216894****
- Burke, R. W., and Mavrodineanu, R., Standard Reference Materials: Accuracy in Analytical Spectrophotometry, NBS Spec. Publ. 260-81 (April 1983). PB-83-214536****
- Weidner, V. R., Standard Reference Materials: White Opal Glass Diffuse Spectral Reflectance Standards for the Visible Spectrum (SRM's 2015 and 2016), NBS Spec. Publ. 260-82 (April 1983). PB-83-220723****
- Bowers, G. N., Jr., Alvarez, R., Cali, J. P., Eberhardt, K. R., Reeder, D. J., Schaffer, R., and Uriano, G. A., Standard Reference Materials: The Measurement of the Catalytic (Activity) Concentration of Seven Enzymes in NBS Human Serum SRM 909, NBS Spec. Publ. 260-83 (June 1983). PB-83-239509****
- Gills, T. E., Seward, R. W., Collins, R. J., and Webster, W. C., Standard Reference Materials: Sampling, Materials Handling, Processing, and Packaging of NBS Sulfur in Coal Standard Reference Materials 2682, 2683, 2684, and 2685, NBS Spec. Publ. 260-84 (August 1983). PB-84-109552****
- Swyt, D. A., Standard Reference Materials: A Look at Techniques for the Dimensional Calibration of Standard Microscopic Particles, NBS Spec. Publ. 260-85 (September 1983). PB-84-112648****
- Hicho, G. E., and Eaton, E. E., Standard Reference Materials: A Standard Reference Material Containing Two and One-Half Percent Austenite, SRM 488, NBS Spec. Publ. 260-86 (December 1983). PB-84-143296****
- Mangum, B. W., Standard Reference Materials: SRM 1969: Rubidium Triple-Point—A Temperature Reference Standard Near 39.30°C, NBS Spec. Publ. 260-87 (December 1983). PB-84-149996**
- Hust, J. G., Standard Reference Materials: A Fine-Grained, Isotropic Graphite for Use as NBS Thermophysical Property RM's from 5 to 2500 K, NBS Spec. Publ. 260-89 (September 1984). PB-85-112886****

- Hust, J. G., and Lankford, A. B., Standard Reference Materials: Update of Thermal Conductivity and Electrical Resistivity of Electrolytic Iron, Tungsten, and Stainless Steel, NBS Spec. Publ. 260-90 (September 1984). PB-85-115814**
- Goodrich, L. F., Vecchia D. F., Pittman, E. S., Ekin, J. W. and Clark, A. F., Standard Reference Materials: Critical Current Measurements on an NbTi Superconducting Wire Standard Reference Material, NBS Spec. Publ. 260-91 (September 1984). SN003-003-02614-0*
- Carpenter, B. S., Standard Reference Materials: Calibrated Glass Standards for Fission Track Use (Supplement to NBS Spec. Publ. 260-49). NBS Spec. Publ. 260-92 (September 1984). PB-85-113025**
- Ehrstein, J., Standard Reference Materials: Preparation and Certification of Standard Reference Materials for Calibration of Spreading Resistance Probes, NBS Spec. Publ. 260-93 (January 1985). PB-85-177921**
- Gills, T. E., Koch, W. F., Stoltz, J. W., Kelly, W. R., Paulsen, P. J., Colbert, J. C., Kirklin, D. R., Pei, P. T. S., Weeks, S., Lindstrom, R. M., Fleming, R. F., Greenberg, R. R., and Paule, R. C., Standard Reference Materials: Methods and Procedures Used at the National Bureau of Standards to Certify Sulfur in Coal SRM's for Sulfur Content, Calorific Value, Ash Content, NBS Spec. Publ. 260-94 (December 1984). PB-85-165900**
- Mulholland, G. W., Hartman, A. W., Hembree, G. G., Marx, E., and Lettieri, T. R., Standard Reference Materials: Development of a 1 μm Diameter Particle Size Standard, SRM 1690, NBS Spec. Publ. 260-95 (May 1985). SN003-003-02665-4*
- Carpenter, B. S., Gramlich, J. W., Greenberg, R. R., Machlan, L. A., DeBievre, P., Eschbach, H. L., Meyer, H., Van Andenhove, J., Connelly, V. E., Trahey, N. M., and Zook, A. C., Standard Reference Materials: Uranium-235 Isotopic Abundance Standard Reference Materials for Gamma Spectrometry Measurements, NBS Spec. Publ. 260-96 (September 1986). PB-87-108544**
- Mavrodineanu, R., and Gills, T. E., Standard Reference Materials: Summary of the Coal, Ore, Mineral, Rock, and Refractory Standards Issued by the National Bureau of Standards, NBS Spec. Publ. 260-97 (September 1985). SN003-003-02688-3*
- Hust, J. G., Standard Reference Materials: Glass Fiberboard SRM for Thermal Resistance, NBS Spec. Publ. 260-98 (August 1985). SN003-003-02674-3
- Callanan, J. E., Sullivan, S. A., and Vecchia, D. F., Standard Reference Materials: Feasibility Study for the Development of Standards Using Differential Scanning Calorimetry, NBS Spec. Publ. 260-99 (August 1985). SN003-003-02675-1*
- Taylor, J. K., Standard Reference Materials: Handbook for SRM Users, NBS Spec. Publ. 260-100 (September 1985). PB-86-110897**
- Mangum, B. W., Standard Reference Materials: SRM 1970, Succinonitrile Triple-Point Standard: A Temperature Reference Standard Near 58.08°C, NBS Spec. Publ. 260-101 (March 1986). SN003-003-02722-7*
- Weidner, V. R., Mavrodineanu, R., Mielenz, K. D., Velapoldi, R. A., Eckerle, K. L., and Adams, B., Standard Reference Materials: Holmium Oxide Solution Wavelength Standard from 240-650 nm, SRM 2034, NBS Spec. Publ. 260-102 (July 1986). PB-86-245727**
- Hust, J. G., Standard Reference Materials: Glass Fiberblanket SRM for Thermal Resistance, NBS Spec. Publ. 260-103 (September 1985). SN003-003-02687-5*
- Mavrodineanu, R., and Alvarez, R., Standard Reference Materials: Summary of the Biological and Botanical Standards Issued by the National Bureau of Standards, NBS Spec. Publ. 260-104 (November 1985). SN003-003-02704-9*
- Mavrodineanu, R., and Rasberry, S. D., Standard Reference Materials: Summary of the Environmental Research Analysis, and Control Standards Issued by the National Bureau of Standards, NBS Spec. Publ. 260-105 (March 1986). SN003-003-02725-1*
- Koch, W. F., ed., Standard Reference Materials: Methods and Procedures Used at the National Bureau of Standards to Prepare, Analyze, and Certify SRM 2694, Simulated Rainwater, and Recommendations for Use, NBS Spec. Publ. 260-106 (July 1986). PB-86-247483**
- Hartman, A. W., McKenzie, R. L., Standard Reference Materials: SRM 1965, Microsphere Slide (10 μm Polystyrene Spheres), NIST Spec. Publ. 260-107 (November 1988)
- Mavrodineanu, R., and Gills, T. E., Standard Reference Materials: Summary of Gas Cylinder and Permeation Tube Standard Reference Materials Issued by the National Bureau of Standards, NBS Spec. Publ. 260-108 (May 1987).
- Candela, G. A., Chandler-Horowitz, D., Novotny, D. B., Marchiando, J. F., and Belzer, B. J., Standard Reference Materials: Preparation and Certification of an Ellipsometrically Derived Thickness and Refractive Index Standard of a Silicon Dioxide Film (SRM 2530), NIST Spec. Publ. 260-109 (October 1988).
- Kirby, R. K., and Kanare, H. M., Standard Reference Materials: Portland Cement Chemical Composition Standards (Blending, Packaging, and Testing), NBS Spec. Publ. 260-110 (February 1988).
- Gladney, E. S., O'Malley, B. T., Roelandts, I., Gills, T. E. Standard Reference Materials: Compilation of Elemental Concentration Data for NBS Clinical, Biological, Geological, and Environmental, NBS Spec. Publ. 260-111 (November 1987)

NIST Calibration Service Contacts

| Measurement Area | Name | Telephone |
|--|---|----------------|
| General Information | Joe D. Simmons or Measurement Services Staff | (301) 975-2002 |
| Dimensional Measurements | | |
| Angular | Yun H. Wang | (301) 975-3468 |
| API Plug and Ring Gages | Edgar G. Erber | (301) 975-3468 |
| End Standards | Theodore D. Doiron | (301) 975-3468 |
| Gage Blocks | John R. Stoup | (301) 975-3468 |
| Hydrometers | John F. Houser | (301) 975-5956 |
| Length and Diameter; Step Gages | Theodore D. Doiron | (301) 975-3468 |
| Line Standards | William B. Penzes | (301) 975-3468 |
| Measuring Wires | David Stieren | (301) 975-3468 |
| Micrometers; Penetration Needles | Edgar G. Erber | (301) 975-3468 |
| Optical Reference Planes | Jay H. Zimmerman | (301) 975-3468 |
| Plain Conical; Threaded Plug and Ring Gages | Edgar G. Erber | (301) 975-3468 |
| Roundness | Yun H. Wang | (301) 975-3468 |
| Sieves | Theodore D. Doiron | (301) 975-3468 |
| Spherical Diameter; Plug Gages; Ring Gages | Yun H. Wang | (301) 975-3468 |
| Surface Texture | Arie Hartman | (301) 975-3475 |
| Surveying Rods and Tapes; Two Dimensional Gages | Ronald G. Hartsock | (301) 975-3465 |
| Volume and Density | John F. Houser | (301) 975-5956 |
| Electromagnetic Measurements | | |
| AC-DC Voltage/Current Converters (to 1 MHz) | Joseph R. Kinard | (301) 975-4250 |
| AC Resistors | T. Michael Souders | (301) 975-2406 |
| Capacitance Dividers | Robert E. Hebner, Jr. | (301) 975-2403 |
| Coaxial/Waveguide Terminations Reflection Coefficients | John R. Juroshek | (303) 497-5362 |
| Current Transformers | John D. Ramboz | (301) 975-2434 |
| Data Converters | T. Michael Souders | (301) 975-2406 |
| DC Resistance | Ronald F. Dziuba | (301) 975-4239 |
| DC Voltage | June E. Sims | (301) 975-4238 |
| Electromagnetic Field-Strength Parameters | Galen H. Koepke | (303) 497-5766 |
| HF Capacitance/Inductance | George M. Free | (303) 497-3609 |
| High Frequency Resistors | George M. Free | (303) 497-3609 |
| High Frequency Voltage | Robert E. Hebner, Jr. | (301) 975-2403 |
| Inductive Dividers | Bruce F. Field | (301) 975-4230 |
| LF AC Voltmeters and Sources | Nile M. Oldham | (301) 975-2414 |
| LF Capacitance/Inductance | Bruce F. Field | (301) 975-4230 |
| LF Power/Energy | Andrew J. Secula | (301) 975-2416 |
| Microwave Antenna Parameter | Allen C. Newell | (303) 497-3743 |
| Mixed Dividers | Edward F. Kelley | (301) 975-5826 |
| N-Port Scattering | Ronald A. Ginley | (303) 497-3634 |
| Noise Temperature | Sunchana Perera | (303) 497-3546 |
| Phase Angle Meters | Raymond S. Turgel | (301) 975-2420 |
| Power-Frequency Capacitors | William E. Anderson | (301) 975-2423 |
| Pulse Waveform | William L. Gans | (303) 497-3538 |
| Q Standards | George M. Free | (303) 497-3609 |
| Resistive Dividers | Martin Misakian | (301) 975-2426 |
| RF-DC Voltage/Current Converter (100 Hz-1 GHz) | Gregorio Rebuldela | (303) 497-3561 |
| RF/Microwave Attenuators | John R. Juroshek | (303) 497-5362 |
| RF/Microwave Phase Shifters | John L. Workman | (303) 497-3954 |
| RF/Microwave Power Meters | Ronald A. Ginley | (303) 497-3634 |
| VHF Omnidirectional Range | Neil T. Larsen | (303) 497-3711 |
| Voltage Transformers | William E. Anderson | (301) 975-2423 |
| Ionizing Radiation Measurements | | |
| Dosimetry of X rays, Gamma Rays and Electrons | Bert M. Coursey | (301) 975-5584 |
| High-Dose Dosimetry | William L. McLaughlin | (301) 975-5559 |
| Neutron Sources and Dosimeters | E. Dale McGarry | (301) 975-6205 |
| Radioactivity Sources | Jacqueline M. Calhoun | (301) 975-5538 |

| Measurement Area | Name | Telephone |
|--|--------------------------|----------------|
| Mechanical Measurements | | |
| Acoustic | Victor Nedzelnitsky | (301) 975-6638 |
| Acoustic Emission Transducer | Franklin R. Breckenridge | (301) 975-6628 |
| Airspeed | Norman E. Mease | (301) 975-5959 |
| Cryogenic Flow Rate | James A. Brennan | (303) 497-3611 |
| Flow Rate | Kenneth R. Benson | (301) 975-5945 |
| Force | Simone L. Yaniv | (301) 975-6655 |
| Mass | Jerry G. Keller | (301) 975-4218 |
| Ultrasonic Reference Block | Gerald V. Blessing | (301) 975-6627 |
| Ultrasonic Transducer | Steven E. Fick | (301) 975-6629 |
| Vibration | Myroslav R. Serbyn | (301) 975-6646 |
| Optical Radiation Measurements | | |
| Laser Power/Energy | Thomas R. Scott | (303) 497-3651 |
| Photometric | Donald A. McSparron | (301) 975-2321 |
| Radiometric | John K. Jackson | (301) 975-2330 |
| Spectrophotometric | P. Yvonne Barnes | (301) 975-2345 |
| UV Radiometric-Standard Detectors | L. Randall Canfield | (301) 975-3728 |
| UV Radiometric-Standard Sources | J. Mervin Bridges | (301) 975-3228 |
| Thermodynamic Measurements | | |
| Humidity | Gregory E. Scace | (301) 975-2626 |
| Laboratory Thermometers | Jacquelyn A. Wise | (301) 975-4822 |
| Pressure | B. Asoka Ratnam | (301) 975-4857 |
| Radiation Thermometry | Ronald L. Wilkinson | (301) 975-2325 |
| Resistance Thermometry | Gregory F. Strouse | (301) 975-4803 |
| Thermocouples and Pyrometer Indicators | George W. Burns | (301) 975-4817 |
| Vacuum and Low Pressure | Richard W. Hyland | (301) 975-4829 |
| Time and Frequency Measurements | | |
| Frequency Dissemination | George Kamas | (303) 497-3378 |
| Time Dissemination | David W. Allan | (303) 497-5637 |
| Oscillator Characterization | James E. Gray | (303) 497-3209 |



H. M. (Skip) Kingston of the Inorganic Analytical Research Division demonstrates the IR-100 award winning microwave dissolution system he led development of at NIST.

INDICES

Numerical Index

| SRM | SRM Description | Certifi- cate Date | Page | SRM | SRM Description | Certifi- cate Date | Page |
|-----|---|--------------------------|------------|--------|---|--------------------------|-----------|
| 1c | Limestone, Argillaceous | Dec 78 | 79 | 99a | Feldspar, Soda | Aug 81 | 79 |
| 3d | Iron, White | Apr 79 | 27 | 100b | Steel, Mn-2 | Aug 59 | 16 |
| 4k | Iron, Cast | May 76 | 27 | 101g | Steel, Cr18-Ni10 (AISI 304L) | Aug 86 | 19 |
| 5L | Iron, Cast | Nov 70 | 27 | 103a | Chrome Refractory | Sep 62 | 81 |
| 6g | Iron, Cast | Nov 70 | 27 | 105 | Steel, High Sulfur | Aug 81 | 15 |
| 7g | Iron, Cast (High-Phosphorus) | Oct 59 | 27 | 106b | Steel, Cr-Mo-Al | Mar 61 | 16 |
| 8j | Steel, Bessemer (Simulated) 0.1C | Apr 72 | 15 | 107c | Iron, Alloy Cast, Ni-Cr-Mo | May 83 | 27 |
| 11h | Steel, BOH, 0.2C | Apr 74 | 15 | 112b | Carbide, Silicon | Jan 85 | 81 |
| 12h | Steel, BOH, 0.4C | Mar 66 | 15 | 114n | Portland Cement, Fineness | Feb 82 | 122 |
| 13g | Steel, Carbon, 0.6C | Apr 74 | 15 | 115a | Iron, Alloy Cast, Cu-Ni-Cr | Apr 62 | 27 |
| 14f | Steel, Carbon (AISI 1078) | Feb 81 | 15 | 120c | Phosphate Rock (Florida) | Feb 88 | 74, 77 |
| 15h | Steel, BOH, 0.1C | Mar 84 | 15 | 121d | Stainless Steel, Cr17-Ni11-Ti0.3 (AISI 321) | Aug 81 | 19 |
| 16f | Steel, BOH, 1.0C | Mar 83 | 15 | 122h | Iron, Cast, Car Wheel | Apr 83 | 27 |
| 17d | Sucrose (Polarimetric) | Aug 86 | 44, 108 | 123c | Stainless Steel, Cr17-Ni11-Nb0.6 (AISI 348) | Jul 71 | 19 |
| 19h | Steel, Basic Electric, 0.2C | Sep 87 | 15 | 125b | Steel, High-Silicon (Ca-Bearing) | Feb 82 | 16 |
| 20g | Steel, AISI 1045 | Oct 70 | 15 | 126c | Steel, High-Nickel (Ni36) | Dec 77 | 18 |
| 25d | Ore, Manganese | Feb 84 | 77 | 127b | Solder, Sn40-Pb60 | Oct 81 | 34 |
| 27f | Ore, Iron (Sibley) | May 77 | 76 | 129c | Steel, High-Sulfur (SAE 112) | Aug 73 | 16 |
| 30f | Steel, Cr-V (SAE 6150) | Jun 79 | 16 | 131d | Steel, Low-C Si (C & S only) | Dec 86 | 16 |
| 32e | Steel, Ni-Cr | Apr 57 | 16 | 132b | Steel, Tool (AISI M2) | Aug 73 | 19 |
| 33e | Steel, Nickel | Nov 84 | 16 | 133b | Steel, Chromium-Molybdenum | Aug 81 | 19 |
| 36b | Steel, Cr2-Mo1 | Jul 69 | 16 | 134a | Steel, Mo8-W2-Cr4-V1 | May 57 | 19 |
| 39i | Benzoic Acid, Calorimetric | Jul 68 | 99 | 136e | Potassium Dichromate | Jun 89 | 44 |
| 40h | Sodium Oxalate (Reductometric) | May 82 | 44 | 139b | Steel, Cr-Ni-Mo (AISI 8640) | May 78 | 16 |
| 41c | D-Glucose (Dextrose) (Polarimetric) | Nov 84 | 44, 108 | 141c | Acetanilide | Sep 76 | 44 |
| 43h | Zinc Freezing Point | Aug 73 | 102 | 142 | Anisic Acid | Jul 69 | 44 |
| 44f | Aluminum Freezing Point | Apr 73 | 102 | 143c | Cystine | Sep 76 | 44 |
| 45d | Copper Freezing Pt | Dec 71 | 102 | 148 | Nicotinic Acid | Dec 70 | 44 |
| 49e | Lead Freezing Point | Dec 71 | 102 | 152a | Steel, BOH, 0.5C | Oct 65 | 15 |
| 50c | Steel, W18-Cr4-V1 (Tool) | Jun 57 | 19 | 153a | Steel, Co8-Mo9-W2 | Jan 60 | 19 |
| 53e | Bearing Metal, Lead-base | Jan 70 | 34 | 154b | Titanium Dioxide | May 73 | 80 |
| 57a | Silicon Metal | Dec 80 | 26 | 155 | Steel, Cr0.5-W0.5 | Oct 46 | 16 |
| 58a | Ferrosilicon (73-Si, Regular Grade) | Apr 78 | 26 | 158a | Bronze, Silicon | Aug 61 | 31 |
| 59a | Ferrosilicon (48-Si) | Nov 69 | 26 | 160c | Stainless Steel, Cr18-Ni12-Mo2 (AISI 316) | Aug 69 | 19 |
| 64c | Ferrochromium, High-Carbon | Aug 77 | 26 | 163 | Steel, Chromium | Jan 68 | 16 |
| 68c | Ferromanganese, High-Carbon | Aug 79 | 26 | 165a | Sand, Glass | Oct 78 | 79 |
| 69b | Bauxite (Arkansas) | Aug 79 | 77 | 166c | Stainless Steel, Low-C (AISI 316L) | Mar 70 | 19 |
| 70a | Feldspar (Potash) | Aug 81 | 79 | 173b | Ti-Base Alloy 6Al-4V | Dec 84 | 37 |
| 71 | Calcium Molybdate | Feb 29 | 26 | 176 | Ti-Base Alloy 5Al-2.5Sn | Oct 81 | 37 |
| 72g | Steel, Low-Alloy AISI 4130 | Jun 81 | 16 | 178 | Steel, BOF 0.4C | Jul 69 | 15 |
| 73c | Stainless Steel, 13 Cr | Jul 66 | 19 | 180 | Fluorspar, High-Grade | Mar 71 | 75 |
| 76a | Burnt Refractory (Al ₂ O ₃ -39) | May 85 | 80 | 181 | Ore, Lithium (Spodumene) | Oct 81 | 75 |
| 77a | Burnt Refractory (Al ₂ O ₃ -60) | May 85 | 80 | 182 | Ore, Lithium (Petalite) | Oct 81 | 75 |
| 78a | Burnt Refractory (Al ₂ O ₃ -72) | May 85 | 80 | 183 | Ore, Lithium (Lepidolite) | Oct 81 | 75 |
| 79a | Fluorspar (Customs grade) | Jan 80 | 75 | 185f | Potassium Hydrogen Phthalate, pH | Jan 84 | 89 |
| 81a | Sand, Glass (High Iron) | Jan 78 | 79 | 186Id | Potassium Dihydrogen Phosphate, pH | Sep 70 | 89 |
| 82b | Iron, Cast, Ni-Cr | Apr 66 | 27 | 186IId | Disodium Hydrogen Phosphate, pH | Sep 70 | 89 |
| 83d | Arsenic Trioxide, Reductometric | Mar 82 | 44 | 187c | Sodium Tetraborate Decahydrate (Borax), pH | Mar 84 | 89 |
| 84j | Potassium Hydrogen Phthalate | Nov 84 | 44 | 188 | Potassium Hydrogen Tartrate, pH | May 87 | 89 |
| 87a | Silicon-Aluminum Alloy | Aug 81 | 30 | 189a | Potassium Tetroxalate, pH | Apr 86 | 89 |
| 88b | Limestone, Dolomitic | Apr 86 | 79 | 191a | Sodium Bicarbonate, pH | Nov 84 | 89 |
| 89 | Glass, Lead-Barium | Aug 32 | 83 | 192a | Sodium Carbonate, pH | Nov 84 | 89 |
| 90 | Ferrophosphorous | Oct 28 | 26 | 193 | Potassium Nitrate | Mar 74 | 74 |
| 91 | Glass, Opal (Powder) | Oct 82 | 83 | | | | |
| 92 | Glass, Soda-Lime (Powder) | Mar 82 | 83 | | | | |
| 93a | Glass, Borosilicate | Aug 73 | 83 | | | | |
| 94c | Zn-base Die-Casting Alloy | Aug 73 | 38 | | | | |
| 97b | Clay, Flint | Apr 88 | 78 | | | | |
| 98b | Clay, Plastic | Apr 88 | 78 | | | | |

| SRM | SRM Description | Certifi- cate Date | Page | SRM | SRM Description | Certifi- cate Date | Page |
|------|---|--------------------------|------------|------|---|--------------------------|-----------|
| 194 | Ammonium Dihydrogen Phosphate | Jan 74 | 74 | 493 | Iron Carbide in Ferrite | May 85 | 115 |
| 195 | Ferrosilicon (75Si) | Jan 76 | 26 | 494 | Copper I | Jan 78 | 33 |
| 196 | Ferrochromium (Low Carbon) | Nov 70 | 26 | 495 | Unalloyed Copper II | Oct 87 | 33 |
| 198 | Silica Brick (0.16 Al ₂ O ₃) | Jan 60 | 81 | 496 | Unalloyed Copper III | Apr 86 | 33 |
| 199 | Silica Brick (0.48 Al ₂ O ₃) | Jan 60 | 81 | 498 | Copper V | Jan 78 | 33 |
| 200 | Potassium Dihydrogen Phosphate | Aug 74 | 74 | 499 | Unalloyed Copper VI | Mar 86 | 33 |
| 211c | Toluene | Sep 84 | 96, 108 | 500 | Copper VII | Jan 78 | 33 |
| 276a | Carbide, Tungsten | May 80 | 81 | 600 | Bauxite (Australian) | Dec 87 | 77 |
| 277 | Tungsten Concentrate | Oct 78 | 75 | 607 | Potassium Feldspar | May 73 | 86 |
| 278 | Obsidian Rock | Aug 81 | 80 | 610 | Glass, Trace Elements (500 ppm) | Jan 82 | 86 |
| 291 | Steel, Cr-Mo (ASTM A-213) | Oct 75 | 16 | 611 | Glass, Trace Elements (500 ppm) | Jan 82 | 86 |
| 293 | Steel, Cr-Ni-Mo (AISI 8620) | Mar 75 | 16 | 612 | Glass, Trace Elements (500 ppm) | Jan 82 | 86 |
| 330 | Copper, Ore Mill Heads | Jan 77 | 75 | 613 | Glass, Trace Elements (50 ppm) | Jan 82 | 86 |
| 331 | Ore, Copper Mill Tails | Jan 77 | 75 | 614 | Glass, Trace Elements (1 ppm) | Jan 82 | 86 |
| 334 | Iron, Gray Cast | Mar 82 | 27 | 615 | Glass, Trace Elements (1 ppm) | Jan 82 | 86 |
| 335 | Steel, BOH, 0.1C | Apr 66 | 15 | 616 | Glass, Trace Elements (0.02 ppm) | Jan 82 | 86 |
| 337a | Steel, BOH, 1.1C | Apr 85 | 15 | 617 | Glass, Trace Elements (0.02 ppm) | Jan 82 | 86 |
| 338 | Iron, White Cast | Jun 82 | 27 | 620 | Glass, Soda-Lime Flat | Jan 82 | 83 |
| 339 | Steel, Cr17-Ni9-Se0.2 | Jul 65 | 19 | 621 | Glass, Container | Jan 82 | 83 |
| 340 | Ferroniobium | Nov 70 | 26 | 622 | Glass, Soda-Lime-Silicate | Mar 76 | 94 |
| 341 | Iron, Ductile Cast | Mar 62 | 27 | 623 | Glass, Borosilicate | Mar 76 | 94 |
| 342a | Iron, Nodular Cast | Apr 70 | 27 | 624 | Glass, Electrical Resistance | Oct 77 | 94 |
| 343a | Stainless Steel, Cr16-Ni2 (AISI 431) | Jul 85 | 19 | 625 | Zn-Base Alloy A | Apr 64 | 38 |
| 344 | Steel, Cr15-Ni7 (Mo precip harden) | Oct 63 | 18 | 626 | Zn-Base Alloy B | Apr 64 | 38 |
| 345 | Steel, Cr16-Ni4 (Cu precip harden) | Jan 64 | 18 | 627 | Zn-Base Alloy C | Apr 64 | 38 |
| 346a | Steel, Valve (Cr21-Ni3-Mn8) | Oct 85 | 18 | 628 | Zn-Base Alloy D | Apr 64 | 38 |
| 348a | High Temp Alloy A286 (Ni26- Cr15) | Mar 87 | 18 | 629 | Zn-Base Alloy E | Apr 64 | 38 |
| 349a | Waspaloy | Jun 87 | 35 | 630 | Zn-Base Alloy F | Apr 64 | 38 |
| 350a | Benzoic Acid, Acidimetric | Apr 81 | 44 | 631 | Zinc Spelter (mod) | Nov 81 | 38 |
| 352c | Titanium for Hydrogen | * | 39 | 640b | Silicon X-ray Diffraction | Jan 87 | 116 |
| 360b | Zircaloy 2, Zr-Base Alloy | Apr 86 | 38 | 641 | Ti-Base Alloy, 8Mn (A) | Oct 81 | 37 |
| 361 | Steel, AISI 4340 | Feb 81 | 17 | 642 | Ti-Base Alloy, 8Mn (B) | Oct 81 | 37 |
| 362 | Steel, AISI 94B17 (modified) | Feb 81 | 17 | 643 | Ti-Base Alloy, 8Mn (C) | Oct 81 | 37 |
| 363 | Steel, Cr-V (modified) | Feb 81 | 17 | 644 | Ti-Base 2Cr-2Fe-2Mo (A) | Jan 60 | 37 |
| 364 | Steel, High C (modified) | Feb 81 | 17 | 646 | Ti-Base 2Cr-2Fe-2Mo (C) | Jan 60 | 37 |
| 365 | Iron, Electrolytic | Feb 81 | 17 | 647 | Ti-Base 6Al-2Mo-2Sn-4Zr | Aug 86 | 37 |
| 367 | Stainless Steel (AISI 446) | Jul 77 | 19 | 648 | Ti-Base Alloy 5Al-2Sn-2Zr- 4Cr-4Mo | Jun 87 | 37 |
| 368 | Steel, AISI 1211 | Jan 78 | 15 | 650 | Titanium | Nov 85 | 37 |
| 370e | Zinc Oxide | none | 121 | 651 | Titanium | Nov 85 | 37 |
| 371h | Sulfur | none | 121 | 652 | Titanium | Nov 85 | 37 |
| 372i | Stearic Acid | none | 121 | 654a | Titanium Alloy 6A1-4V | Oct 81 | 37 |
| 375g | Channel Black | none | 121 | 656 | Silicon Nitride Alpha/Beta | * | 117 |
| 378b | Oil Furnace Black | none | 121 | 658 | Tridymite Quantitative XRD | * | 117 |
| 382a | Gas Furnace Black | none | 121 | 659 | Ceramic Particle Size | * | 117 |
| 383a | Mercaptobenzothiazole | none | 121 | 660 | LaB ₆ -2 Theta X-Ray Line Profile | * | 117 |
| 384f | n-Tertiary-Butyl-2 | none | 121 | 668 | Steels, Set 661-665 | Sep 81 | 20 |
| 386k | Styrene Butadiene | Jan 85 | 121 | 670 | Ore, Rutile | Jun 85 | 77 |
| 388p | Butyl Rubber | Mar 87 | 121 | 671 | Nickel Oxide 1 | Dec 60 | 36 |
| 393 | Copper "0" | Sep 80 | 33 | 672 | Nickel Oxide 2 | Dec 60 | 36 |
| 394 | Copper I | Jan 78 | 33 | 673 | Nickel Oxide 3 | Dec 60 | 36 |
| 395 | Copper II | Jan 78 | 33 | 674a | Intensity X-ray Diffraction Set | Jun 83 | 117 |
| 396 | Copper III | Jan 78 | 33 | 675 | Mica X-ray Diffraction | Jun 82 | 117 |
| 398 | Copper V | Jan 78 | 33 | 676 | Crystalline Alpha Alumina | * | 117 |
| 399 | Copper VI | Jan 78 | 33 | 677 | Crystalline X-Ray Line Proc- essing | * | 117 |
| 400 | Copper VII | Jan 78 | 33 | 679 | Brick Clay | Jan 87 | 78 |
| 454 | Copper XI | Sep 80 | 33 | 680a | Platinum, High Purity | Mar 77 | 40 |
| 457 | Copper IV | Jan 78 | 33 | 681 | Platinum, Doped | Mar 77 | 40 |
| 470 | Mineral Glasses | Oct 81 | 41, 42 | 682 | Zinc, High Purity | Jul 68 | 40 |
| 473 | Optical Linewidth | * | 91 | 683 | Zinc, Pure | Oct 81 | 40 |
| 475 | AR Cr Optical Linewidth | Apr 81 | 91 | 685 | Gold, High Purity | Oct 81 | 40 |
| 476 | B Cr Optical Linewidth | * | 91 | 688 | Basalt Rock | Aug 81 | 80 |
| 477 | Glass Fluorescence Source | Feb 83 | 116 | 689 | Ferrochromium Silicon | Feb 82 | 26 |
| 480 | Tungsten-Molybdenum | Nov 68 | 41 | 690 | Ore, Iron (Canada) | Oct 78 | 76 |
| 481 | Gold-Silver | Feb 69 | 41 | 691 | Reduced Iron Oxide | Apr 82 | 76 |
| 482 | Gold-Copper | Jun 69 | 41 | 692 | Ore, Iron (Labrador) | Oct 82 | 76 |
| 484f | SEM Magnification | * | 90 | 693 | Ore, Iron (Nimba) | Oct 78 | 76 |
| 485a | 5% Austenite in Ferrite | Oct 81 | 115 | 694 | Phosphate Rock (Western) | Jun 84 | 74, 77 |
| 486 | 15% Austenite in Ferrite | Mar 81 | 115 | 696 | Bauxite (Surinam) | Aug 79 | 77 |
| 487 | 30% Austenite in Ferrite | May 82 | 115 | 697 | Bauxite (Dominican) | Aug 79 | 77 |
| 488 | 2% Austenite in Ferrite | Oct 83 | 115 | 698 | Bauxite (Jamaican) | Aug 79 | 77 |
| 489 | Ferrite in Austenite (10%) | * | 115 | 699 | Alumina, Reduction Grade | Aug 81 | 77 |

*In Prep.

| SRM | SRM Description | Certi- ficate Date | Page | SRM | SRM Description | Certi- ficate Date | Page |
|------|--|--------------------------|-------------------|-------|--|--------------------------|-----------|
| 705 | Polystyrene 179k mol wt | Nov 78 | 97, 98, 100 | 900 | Antiepilepsy Drug Level Assay | Apr 79 | 49 |
| 706 | Polystyrene 258k mol wt | Feb 79 | 97, 98 | 909 | Human Serum | Mar 85 | 49, 50 |
| 708 | Glasses, Stress Optical Coefficient | Sep 73 | 95 | 910 | Sodium Pyruvate | May 81 | 49 |
| 709 | Glass, Extra Dense Lead | Jun 74 | 95 | 911b | Cholesterol | May 88 | 49 |
| 710a | Glass, Soda Lime-Silica | * | 94, 95 | 912a | Urea | Nov 79 | 49 |
| 711 | Glass, Lead-Silica | Jul 64 | 94, 95 | 913 | Uric Acid | Nov 73 | 49 |
| 712 | Glass, Alkali Lead-Silica | Oct 66 | 95 | 914a | Creatinine | Oct 87 | 49 |
| 713 | Glass, Dense Barium Crown | Oct 66 | 95 | 915 | Calcium Carbonate | Nov 73 | 49 |
| 714 | Glass, Alkali Alumina Silica | Oct 66 | 95 | 916a | Bilirubin | * | 49 |
| 715 | Glass, Alkali-free Alumina | Sep 66 | 95 | 917a | d-Glucose (Dextrose) | Sep 73 | 49 |
| 716 | Glass, Neutral | Sep 66 | 95 | 918 | Potassium Chloride | Nov 73 | 49 |
| 717 | Glass, Borosilicate | Nov 69 | 94- | 919a | Sodium Chloride | Nov 73 | 49 |
| | | | 95 | 920 | D-Mannitol | Nov 73 | 49 |
| 723a | Tris(hydroxymethyl) aminomethane, Basimetric | Apr 81 | 44 | 921 | Cortisol (Hydrocortisone) | Dec 73 | 49 |
| 724a | Tris(hydroxymethyl) aminomethane, Calorimetric | Sep 73 | 99 | 922 | Tris(hydroxymethyl) aminomethane, pH | Aug 76 | 89 |
| 726 | Selenium, Inter-Purity | Jan 67 | 40 | 923 | Tris(hydroxymethyl) aminomethane hydrochloride, pH | Aug 76 | 49, 89 |
| 728 | Zinc-Intermediate Purity | Oct 81 | 40 | 924 | Lithium Carbonate | Nov 73 | 49 |
| 731 | Glass, Borosilicate | Jul 72 | 104 | 925 | 4-Hydroxy-3-methoxy-dl-mandelic Acid (VMA) | Dec 73 | 49 |
| 738 | Stainless Steel | Nov 86 | 104 | 926 | Bovine Serum Albumin (Total Protein) | Jul 77 | 49 |
| 739 | Fused Silica | May 71 | 104 | 927a | Bovine Serum Albumin (7% Solution, Total Protein) | Aug 86 | 49 |
| 740a | Zinc Freezing Point | Feb 70 | 102 | 928 | Lead Nitrate | May 76 | 49 |
| 741 | Tin Freezing Point | Jul 72 | 102 | 929 | Magnesium Gluconate Dihydrate | Apr 79 | 49 |
| 742 | Alumina Melting Point | Jul 70 | 102 | 930D | Glass Filters for Spectrophotometry (Visible) | Aug 84 | 106 |
| 743 | Mercury, Triple Point | Apr 76 | 102 | 931d | Liquid Absorbance Filters for UV and Visible Spectrophotometry | Oct 86 | 106 |
| 745 | Gold, Vapor Pressure | May 69 | 102 | 934 | Clinical Laboratory Thermometer | * | 103 |
| 746 | Cadmium, Vapor Pressure | Aug 70 | 102 | 935a | Crystalline Potassium Dichromate for UV Absorbance | Jun 77 | 106 |
| 748 | Silver, Vapor Pressure | Aug 70 | 102 | 936 | Quinine Sulfate Dihydrate | Apr 79 | 106 |
| 763 | Aluminun, Magnetic Susceptibility | Apr 73 | 105 | 937 | Iron Metal | Jun 78 | 49 |
| 766 | Manganese Fluoride, Mag Suscept | Apr 73 | 105 | 938 | 4-Nitrophenol | May 81 | 49 |
| 767a | Thermometric Fix Point Device | Jun 83 | 101 | 951 | Boric Acid | Oct 71 | 44, 88 |
| 769 | Electrical "RRR" Set | Nov 82 | 119 | 952 | Enriched Boric Acid | Oct 71 | 88 |
| 772 | Nickel, Magnetic Moment | Oct 78 | 105 | 953 | Neutron Density Monitor Wire | Mar 69 | 87 |
| 773 | Glass, Liquidus Temperature | Nov 80 | 95 | 956 | Electrolytes in Serum for ISE | * | 49 |
| 774 | Glass, Dielectric Constant | Jul 82 | 94 | 963a | Fission Track Glass (U-1 ppm) | Feb 84 | 87 |
| 853 | Aluminum Alloy 3004 | May 85 | 30 | 968 | Fat Soluble Vitamins in Human Serum | Apr 88 | 49 |
| 854 | Aluminum Alloy 5182 | May 85 | 30 | 975 | Chlorine, Isotopic | Mar 65 | 88 |
| 855 | Aluminum Casting Alloy 356 | Jan 80 | 30 | 976 | Copper, Isotopic | Mar 65 | 88 |
| 856 | Aluminum Casting Alloy 380 | Jan 80 | 30 | 977 | Bromine, Isotopic | Mar 65 | 88 |
| 858 | Aluminum Alloy 6011 (mod) | Jun 80 | 30 | 978a | Silver, Isotopic | Sep 84 | 88 |
| 859 | Aluminum Alloy 7075 | Jun 80 | 30 | 979 | Chromium, Isotopic | May 66 | 88 |
| 864 | Inconel 600 | May 84 | 35 | 980 | Magnesium, Isotopic | Jan 67 | 88 |
| 865 | Inconel 625 | May 84 | 35 | 981 | Lead, Common Isotopic | Apr 73 | 88 |
| 866 | Incoloy 800 | May 84 | 35 | 982 | Lead, Equal-Atom Isotopic | Jun 68 | 88 |
| 867 | Incoloy 825 | May 84 | 35 | 983 | Lead, Radiogenic Isotopic | Jun 68 | 88 |
| 868 | High-Temperature Alloy (Fe-Ni-Co) | May 87 | 18 | 984 | Rubidium Chloride, Assay & Isotopic | Jul 70 | 88 |
| 871 | Phosphor Bronze, CDA 521 | Aug 79 | 31 | 985 | Potassium, Assay & Isotopic | Aug 79 | 88 |
| 872 | Phosphor Bronze, CDA 544 | Aug 79 | 31 | 986 | Nickel, Isotopic | * | 88 |
| 874 | Cupro-Nickel, 10 (CDA 706) (pure) | Jan 78 | 31 | 987 | Strontium, Assay & Isotopic | Oct 82 | 44, 88 |
| 875 | Cupro-Nickel, 10 (CDA 706) (doped) | Jan 78 | 31 | 989 | Rhenium, Assay & Isotopic | Feb 74 | 88 |
| 879 | Nickel Silver, CDA 762 | Jun 79 | 31 | 990 | Silicon, Assay & Isotopic | Aug 75 | 88 |
| 880 | Nickel Silver, CDA 770 | Jun 79 | 31 | 991 | Lead-206 Spike, Assay & Isotopic | Mar 76 | 88 |
| 882 | Ni-Cu Alloy (65Ni-31Cu-3Al) | Aug 79 | 35 | 994 | Gallium, Isotopic | Dec 85 | 88 |
| 887 | Cemented Carbide | Sep 88 | 82 | 997 | Thallium, Isotopic | Jan 86 | 88 |
| 888 | Cemented Carbide | Sep 88 | 82 | 998 | Angiotensin I (Human) | Jan 83 | 49 |
| 889 | Cemented Carbide | Sep 88 | 82 | 999 | Potassium Chloride (Primary) | Sep 72 | 44 |
| 890 | Iron, HA White Cast (HC-250+V) | Apr 82 | 27 | 1001 | X-Ray Film Step Tablet (0-4) | Jun 86 | 125 |
| 891 | Iron, HA White Cast (Ni-Hard I) | Apr 82 | 27 | 1002c | Surface Flammability | Dec 78 | 128 |
| 892 | Iron, HA White Cast (Ni-Hard IV) | Apr 82 | 27 | 1003a | Glass Spheres (8-58 µm) | Sep 84 | 122 |
| 897 | Tracealloy A | Aug 83 | 36 | 1004a | Glass Beads (34-120 µm) | * | 122 |
| 898 | Tracealloy B | Aug 83 | 36 | 1006c | Smoke Density, Nonflame (celulose) | Apr 83 | 128 |
| 899 | Tracealloy C | Aug 83 | 36 | 1007a | Smoke Density, Flame (ABS plastic) | Feb 76 | 128 |
| | | | | 1008 | Photographic Step Tablet (0-4) | Jun 86 | 125 |

*In Prep.

| SRM | SRM Description | | | Certifi- cate Date | Page | SRM | SRM Description | | | Certifi- cate Date | Page |
|--------|---|---------------|------|--------------------------|------|-------|--------------------------------------|--|--|--------------------------|--------|
| 1010a | Microcopy Resolution Charts | Radiant Panel | Test | Jun 82 | 125 | 1157 | Steel, Tool (AISI M2) | | | Aug 73 | 24 |
| 1012 | Glass Beads (100–310 µm) | | | Sep 84 | 128 | 1158 | Steel, High-Nickel (36Ni) | | | Dec 77 | 24 |
| 1017a | Glass Beads (225–780 µm) | | | Sep 71 | 122 | 1159 | Electronic and Magnetic Alloy | | | Aug 81 | 35 |
| 1018a | Glass Spheres (0.76–2.16 mm) | | | May 73 | 122 | 1160 | Electronic and Magnetic Alloy | | | Aug 81 | 35 |
| 1019a | Unalloyed Copper | | | Oct 84 | 122 | 1171 | Stainless Steel (AISI 321) | | | Jul 71 | 24 |
| 1034 | Leaded-Tin Bronze Alloy | | | Feb 82 | 31 | 1172 | Stainless Steel (AISI 348) | | | Jul 71 | 24 |
| 1035 | Smoke Toxicity | | | Feb 82 | 31 | 1173 | Steel, Ni-Cr-Mo-V | | | May 83 | 28 |
| 1048 | Barium Metallo-organic | | | * | 128 | C1173 | Steel, Cast 3 | | | Feb 81 | 28 |
| 1051b | Vanadium Metallo-organic | | | Jul 15 | 72 | 1199 | High Temperature Alloy—L605 | | | Aug 74 | 25 |
| 1052b | Cadmium Metallo-organic | | | Mar 68 | 72 | 1200 | High Temperature Alloy—S816 | | | Aug 74 | 25 |
| 1053a | Tin Metallo-organic | | | Jan 70 | 72 | 1216 | Carbon Modified Silicon | | | Nov 87 | 35 |
| 1057b | Lead Metallo-organic | | | Aug 68 | 72 | 1217 | Steel, Nickel (SAE 4820) | | | Nov 84 | 20 |
| 1059c | Lithium Metallo-organic | | | Sep 87 | 6372 | 1218 | Steel, Silicon, Low C & S | | | Nov 84 | 20 |
| 1060a | Magnesium Metallo-organic | | | Apr 64 | 72 | 1219 | Stainless Steel (AISI 413) | | | Sep 85 | 24 |
| 1061c | Nickel Metallo-organic | | | Oct 81 | 72 | C1221 | Steel, Resulfurized/Rephosphorized | | | Jan 82 | 20 |
| 1065b | Silicon Metallo-organic | | | Nov 67 | 72 | 1222 | Cr, Ni, Mo (AISI 8640) | | | Oct 78 | 20 |
| 1066a | Sodium Metallo-organic | | | Apr 69 | 72 | 1223 | Stainless Steel, High S (AISI 416) | | | Sep 85 | 24 |
| 1069b | Zinc Metallo-organic | | | Feb 69 | 72 | 1224 | Steel, Carbon (AISI 1078) | | | Feb 81 | 20 |
| 1070a | Strontium Metallo-organic | | | Apr 64 | 72 | 1225 | Steel, Low-Alloy (AISI 4130) | | | Mar 83 | 20 |
| 1071b | Phosphorus Metallo-organic | | | Feb 76 | 72 | 1226 | Steel, Low-Alloy (HY 130) | | | Dec 82 | 20 |
| 1073b | Aluminum Metallo-organic | | | Jul 67 | 72 | 1227 | Steel, BOH 1.0 C | | | Mar 83 | 20 |
| 1075a | Silver Metallo-organic | | | Oct 67 | 72 | 1228 | Steel, BOH 0.1 C | | | Sep 82 | 20 |
| 1077a | Chromium Metallo-organic | | | Feb 68 | 72 | 1230 | High Temperature Alloy | | | Jun 87 | 25 |
| 1078b | Iron Metallo-organic | | | Jul 72 | 72 | 1233 | Steel, Valve | | | Feb 86 | 24 |
| 1079b | Copper Metallo-organic | | | Feb 69 | 72 | 1234 | Zirconium A | | | Nov 80 | 38 |
| 1080a | Wear-Metals in Lube Oil (Base Oil) | | | Jul 85 | 73 | 1237 | Zircaloy-4 D | | | Nov 80 | 38 |
| 1085 | Wear-Metals in Lube Oil (300 ppm) | | | Jul 85 | 73 | 1240 | Aluminum Alloy 3004 | | | Jul 85 | 30 |
| 1089 | Gasometric Set (1095—1099) | | | Set | 39 | 1240a | Aluminum Alloy 3004 | | | Jul 85 | 30 |
| 1090 | Oxygen in Ingot Iron | | | Oct 85 | 39 | 1240b | Aluminum Alloy 3004 | | | Jul 85 | 30 |
| 1091a | Oxygen in Stainless Steel (AISI 431) | | | Oct 85 | 39 | 1241a | Aluminum Alloy 5182 | | | Aug 85 | 30 |
| 1093 | Oxygen in Valve Steel | | | Nov 84 | 39 | 1241b | Aluminum Alloy 5182 | | | Aug 85 | 30 |
| 1094 | Oxygen in Maraging Steel | | | Nov 84 | 39 | 1243 | Waspaloy | | | Jan 89 | 35 |
| 1103 | Brass, Free Cutting, A | | | Aug 65 | 32 | 1244 | Inconel 600 | | | May 84 | 25 |
| 1104 | Brass, Free Cutting, B | | | Aug 65 | 32 | 1245 | Inconel 625 | | | May 84 | 25 |
| 1107 | Brass, Naval, B | | | Nov 81 | 32 | 1246 | Inconel 800 | | | May 84 | 25 |
| C1107 | Brass, Naval, B | | | Nov 81 | 32 | 1247 | Inconel 825 | | | May 84 | 25 |
| 1108 | Brass, Naval, C | | | Nov 81 | 32 | C1248 | Nickel Copper Alloy | | | Dec 86 | 35 |
| C1108 | Brass, Naval, C | | | Nov 81 | 32 | 1250 | Pyromet | | | Jun 87 | 25 |
| C1110 | Brass, Red, B | | | Oct 81 | 32 | C1251 | Phosphorized Copper (Cu VIII) | | | Sep 80 | 33 |
| 1111 | Brass, Red, C | | | Oct 81 | 32 | C1252 | Phosphorized Copper (Cu IX) | | | Sep 80 | 33 |
| C1111 | Gilding Metal, A | | | Oct 81 | 32 | C1253 | Phosphorized Copper (Cu X) | | | Sep 80 | 33 |
| C1112 | Gilding Metal, A | | | Oct 81 | 32 | 1254 | Steel, Silicon (Ca only) | | | Apr 82 | 20 |
| 1113 | Gilding Metal, B | | | Oct 81 | 32 | 1255a | Aluminum Casting Alloy 356 | | | Nov 86 | 30 |
| C1113 | Gilding Metal, B | | | Oct 81 | 32 | 1256a | Aluminum Casting Alloy 380 | | | Nov 86 | 30 |
| 1114 | Gilding Metal, C | | | Oct 81 | 32 | C1257 | High Purity Aluminum | | | Jan 87 | 30, 40 |
| C1114 | Gilding Metal, C | | | Oct 81 | 32 | 1258 | Aluminum Alloy 6011 (mod) | | | May 78 | 30 |
| 1115 | Bronze, Commercial, A | | | Nov 81 | 32 | 1259 | Aluminum Alloy 7075 | | | May 78 | 30 |
| C1115 | Bronze, Commercial, A | | | Nov 81 | 32 | 1267 | Stainless Steel (AISI 446) | | | Jan 78 | 24 |
| 1116 | Bronze, Commercial, B | | | Nov 81 | 32 | 1269 | Steel (AISI 1526) Line Pipe (mod) | | | Jun 81 | 20 |
| C1116 | Bronze, Commercial, B | | | Nov 81 | 32 | 1270 | Steel, A336 (F-22) 2.3Cr-1Mo | | | Jun 81 | 20 |
| 1117 | Bronze, Commercial, C | | | Nov 81 | 32 | 1275 | Cupro-Nickel (CDA 706) | | | Mar 80 | 32 |
| C1117 | Bronze, Commercial, C | | | Nov 81 | 32 | 1276 | Cupro-Nickel (CDA 715) | | | Mar 80 | 32 |
| C1119 | Brass, Aluminum, B | | | Jul 82 | 32 | C1285 | Steel, A242 (mod) | | | Jun 82 | 20 |
| 1129 | Solder (63Sn–37Pb) | | | May 89 | 34 | 1286 | Steel, Low Alloy (HY 80) | | | Jun 82 | 20 |
| 1131 | Solder (40Sn–60Pb) | | | Oct 81 | 34 | C1287 | Steel, ACI HK (AISI 310 mod) | | | Jun 81 | 24 |
| 1132 | Bearing Metal, Pb-Base | | | Jan 70 | 34 | C1288 | Steel, ACI CN-7M (A-743) | | | Aug 81 | 24 |
| 1133 | Titanium Base Alloy 5Al–2Sn–2Zr–4Cr–4Mo | | | Jul 87 | 37 | C1289 | Steel, ACI CA-6NM (AISI 414 mod) | | | Jun 81 | 24 |
| 1134 | Steel, High-Silicon | | | Apr 70 | 20 | C1290 | Iron, White Cast (HC-250+V) | | | Jan 85 | 28 |
| 1135 | Steel, High-Silicon | | | Jul 72 | 20 | C1291 | Iron, White Cast (Ni-Hard, Type I) | | | Jan 85 | 28 |
| C1137a | Iron, White Cast | | | Jan 84 | 28 | C1292 | Iron, White Cast (Ni-Hard, Type IV) | | | Jan 85 | 28 |
| 1138a | Steel, Cast, 1 | | | Jan 77 | 28 | 1321 | Cu-Cr Coating (nonmagnetic) on Steel | | | Jul 88 | 92 |
| 1139a | Steel, Cast, 2 | | | Jan 77 | 28 | 1322 | Cu-Cr Coating (nonmagnetic) on Steel | | | Jul 88 | 92 |
| 1144a | Iron, Blast Furnace, 2a | | | Dec 76 | 28 | 1323 | Cu-Cr Coating (nonmagnetic) on Steel | | | Jul 88 | 92 |
| C1145a | Iron, White Cast | | | Jun 87 | 28 | 1357 | Cu-Cr Coating (nonmagnetic) on Steel | | | Aug 84 | 92 |
| C1146a | Iron, White Cast | | | Oct 83 | 28 | 1358 | Cu-Cr Coating (nonmagnetic) on Steel | | | Aug 84 | 92 |
| C1150a | Iron, White Cast | | | Dec 85 | 28 | | | | | | |
| C1151 | Stainless Steel (23Cr-7Ni) | | | Jan 80 | 24 | | | | | | |
| C1152 | Stainless Steel (18Cr-11Ni) | | | Jan 80 | 24 | | | | | | |
| C1153 | Stainless Steel (17Cr-9Ni) | | | Jan 80 | 24 | | | | | | |
| C1153a | Stainless Steel (17Cr-9Ni) | | | Jan 80 | 24 | | | | | | |
| C1154 | Stainless Steel (19Cr-13Ni) | | | Jan 80 | 24 | | | | | | |
| 1155 | Stainless Steel (AISI 316) | | | Aug 69 | 24 | | | | | | |

*In Prep.

| SRM | SRM Description | Certi-ficate Date | Page | SRM | SRM Description | Certi-ficate Date | Page |
|-------|--|-------------------|-------------|-------|--|-------------------|--------|
| 1359 | Cu-Cr Coating (nonmagnetic) on Steel | May 84 | 92 | 1497 | Polyethylene Resin (Pigmented) | Jul 87 | 97 |
| 1360 | Cu-Cr Coating (nonmagnetic) on Steel | May 84 | 92 | 1507 | THC in Freeze-Dried Urine | * | 68 |
| 1361a | Cu-Cr Coating (nonmagnetic) on Steel | May 84 | 92 | 1508 | Cocaine in Urine | * | 68 |
| 1362a | Cu-Cr Coating (nonmagnetic) on Steel | May 84 | 92 | 1514 | Thermal Analysis Purity (DSC) | Jul 84 | 100 |
| 1363a | Cu-Cr Coating (nonmagnetic) on Steel | May 84 | 92 | 1515 | Apple Leaves | * | 53 |
| 1364a | Cu-Cr Coating (nonmagnetic) on Steel | May 84 | 92 | 1543 | GC/MS System Performance | Aug 84 | 67 |
| 1365a | Nickel (magnetic) on Steel | May 84 | 92 | 1547 | Peach Leaves | * | 53 |
| 1366a | Nickel (magnetic) on Steel | May 84 | 92 | 1548 | Total Diet | undated | 51 |
| 1379 | Ultra-thin Gold on Nickel 0.35 mg | May 84 | 93 | 1549 | Non-Fat Milk Powder | Jul 85 | 51 |
| 1380 | Ultra-thin Gold on Nickel 0.55 mg | May 84 | 93 | 1563 | Cholesterol and Fat-Soluble Vitamins in Coconut Oil | Jul 87 | 52 |
| 1387 | Gold Coating on Nickel 2.2 mg | Sep 85 | 93 | 1566a | Oyster Tissue | Oct 89 | 51 |
| 1399b | Gold Coating on Nickel (set) | May 84 | 93 | 1567a | Wheat Flour | Sep 88 | 51 |
| 1411 | Soft Borosilicate Glass | Aug 85 | 83 | 1568a | Rice Flour | Jan 78 | 51 |
| 1412 | Multicomponent Glass | Aug 85 | 83 | 1569 | Brewers Yeast (Cr only) | Sep 76 | 51 |
| 1413 | Glass Sand, High Alumina | Aug 85 | 79 | 1572 | Citrus Leaves | Dec 82 | 53 |
| 1449 | Fumed Silica Board | Jan 89 | 104 | 1573a | Tomato Leaves | ** | 53 |
| 1450b | Thermal Resistance, Fibrous Glass Board | May 85 | 104 | 1575 | Pine Needles | Oct 76 | 53 |
| 1451 | Thermal Resistance, Fibrous Glass Blanket | May 85 | 104 | 1577a | Bovine Liver | Feb 85 | 51 |
| 1452 | Thermal Resistance, Fibrous Glass Batt | Apr 86 | 104 | 1579 | Powdered Lead-Based Paint (Pb only) | Jan 73 | 59 |
| 1457 | Superconducting Critical Current Nb-Ti Wire | Jun 84 | 119 | 1580 | Organics in Shale Oil | Nov 80 | 63, 64 |
| 1459 | Fumed Silica Board | Jan 89 | 104 | 1581 | Polychlorinated Biphenyls in Oil | Jan 82 | 63, 65 |
| 1461 | Thermal Conductivity and Electrical Resistivity, Stainless Steel | May 84 | 104, 118 | 1582 | Petroleum Crude Oil | Jan 84 | 63, 64 |
| 1462 | Thermal Conductivity and Electrical Resistivity, Stainless Steel | May 84 | 104, 118 | 1583 | Chlorinated Pesticides in 2,2,4-Trimethylpentane | Feb 85 | 63, 65 |
| 1470 | Gas Transmission, Polyester Film | Feb 82 | 117 | 1584 | Priority Pollutant Phenols in Methanol | Apr 84 | 63, 65 |
| 1474 | Polyethylene Melt Flow | * | 97 | 1585 | Chlorinated Biphenyls | Jan 86 | 63, 66 |
| 1475 | Linear Polyethylene (52k mol wt) | Dec 78 | 97, 98, 100 | 1586 | Isotopically Labeled and Unlabeled Priority Pollutants in Methanol | Oct 84 | 63, 66 |
| 1476 | Branched Polyethylene (viscosity) | Nov 69 | 97, 98 | 1587 | Nitrated Polycyclic Aromatic Hydrocarbons in Methanol | Jun 85 | 63, 66 |
| 1478 | Polystyrene (37k mol wt) | Jan 79 | 97, 98 | 1588 | Organics in Cod Liver Oil | Jan 89 | 63, 67 |
| 1479 | Polystyrene (1M mol wt) | Mar 81 | 97, 98 | 1589 | PCB's in Human Serum | Nov 85 | 49, 67 |
| 1480 | Polyurethanes Low MW | * | 97 | 1590 | Stabilized Wine | Dec 80 | 52 |
| 1481 | Polyurethanes High MW | * | 97 | 1595 | Tripalmitin | Jul 83 | 49 |
| 1482 | Linear Polyethylene (13k mol wt) | Oct 76 | 97, 98 | 1596 | Dinitropyrene Isomers and I-Nitropyrene in Methylene Chloride | Jul 87 | 63, 67 |
| 1483 | Linear Polyethylene (32k mol wt) | Mar 76 | 97, 98 | 1597 | Complex Mixture of Polycyclic Aromatic Hydrocarbons | Sep 87 | 63, 64 |
| 1484 | Linear Polyethylene (119k mol wt) | Oct 76 | 97, 98 | 1598 | Inorganic Constituents in Bovine Serum | Jan 89 | 49 |
| 1487 | Poly (methylmethacrylate), $M_n \approx 6,000$ | Jun 89 | 97, 98 | 1599 | Anticonvulsant Drug Level Assay | Aug 82 | 49 |
| 1488 | Poly (methylmethacrylate) | Feb 88 | 97, 98 | 1600 | Secondary Standard Magnetic Tape Cassette (Computer Amplitude) | Mar 74 | 126 |
| 1489 | Poly (methylmethacrylate) | Mar 86 | 97, 98 | 1614 | Dioxin in Isooctane | Jul 85 | 63, 67 |
| 1490 | Polyisobutylene Solution in Cetane, Rheology | Dec 77 | 98 | 1616 | Sulfur in Kerosene | Feb 88 | 61 |
| 1491 | Aromatic Hydrocarbons in Hexane | Aug 89 | 63, 64 | 1617 | Sulfur in Kerosene | Feb 88 | 61 |
| 1492 | Chlorinated Pesticides in Hexane | Apr 89 | 63, 65 | 1618 | V and Ni in Residual Fuel Oil | May 85 | 59 |
| 1493 | Polychlorinated Biphenyls in 2,2,4-Trimethylpentane | undated | 63 | 1619 | Sulfur in Residual Fuel Oil (0.7%) | Dec 81 | 61 |
| 1495 | Rubber, Isobutylene-Isoprene (Butyl) (Low Mooney Viscosity) | Mar 81 | 121 | 1620b | Sulfur in Residual Fuel Oil (4.5%) | * | 61 |
| 1496 | Polyethylene Resin (Natural) | Sep 88 | 97 | 1621c | Sulfur in Residual Fuel Oil (0.9%) | Sep 86 | 61 |
| | | | | 1622c | Sulfur in Residual Fuel Oil (1.9%) | Sep 86 | 61 |
| | | | | 1623b | Sulfur in Residual Fuel Oil (0.2%) | * | 61 |
| | | | | 1624b | Sulfur in Distillate (Diesel) Fuel Oil (0.1%) | * | 61 |
| | | | | 1625 | Sulfur Dioxide Permeation Tube, 10 cm | Jan 73 | 57 |
| | | | | 1626 | Sulfur Dioxide Permeation Tube, 5 cm | Aug 71 | 57 |
| | | | | 1627 | Sulfur Dioxide Permeation Tube, 2 cm | Aug 71 | 57 |
| | | | | 1630 | Trace Mercury in Coal | Aug 79 | 59 |

*In Prep.

**Subscription.

| SRM | SRM Description | Certi-ficate Date | Page | SRM | SRM Description | Certi-ficate Date | Page |
|-------|--|-------------------|--------|-------|--|-------------------|--------|
| 1632b | Trace Elements in Coal (Bituminous) | Jun 85 | 62 | 1690 | Polystyrene Spheres, 1 μm | Dec 82 | 122 |
| 1633a | Trace Elements in Coal Fly Ash | Jan 85 | 62 | 1691 | Polystyrene Spheres, 0.3 μm | May 84 | 122 |
| 1634b | Trace Elements in Fuel Oil | Feb 86 | 62 | 1693a | Sulfur Dioxide in Nitrogen, 50 ppm | Jun 89 | 56 |
| 1635 | Trace Elements in Coal (Sub-bituminous) | Aug 79 | 62 | 1694a | Sulfur Dioxide in Nitrogen, 100 ppm | Jun 89 | 56 |
| 1638b | Lead in Reference Fuel | Aug 86 | 59 | 1696 | Sulfur Dioxide in Nitrogen, 3500 ppm | Jul 84 | 56 |
| 1639 | Halocarbons (in methanol) for Water Analysis | Apr 83 | 63, 64 | 1700a | Carbon Dioxide in Nitrogen, 10% Blood Gas | Aug 88 | 49, 54 |
| 1641b | Mercury in Water ($\mu\text{g/mL}$) | Apr 83 | 59 | 1701a | Carbon Dioxide and Oxygen in Nitrogen, 5% and 12%, Blood Gas | Aug 88 | 49, 54 |
| 1643b | Trace Elements in Water | May 84 | 62 | 1702a | Carbon Dioxide and Oxygen in Nitrogen, 5% and 20%, Blood Gas | Aug 88 | 49, 54 |
| 1646 | Estuarine Sediment | Jun 82 | 62 | 1703a | Carbon Dioxide and Oxygen in Nitrogen, 10% and 7%, Blood Gas | Aug 88 | 49, 54 |
| 1647a | Priority Pollutant Polynuclear Aromatic Hydrocarbons | Dec 81 | 63, 64 | 1754 | Steel, Low Alloy (AISI 4320) | Feb 89 | 39 |
| 1648 | Urban Particulate Matter | May 82 | 62 | 1761 | Steel, Low Alloy | Feb 88 | 20 |
| 1649 | Urban Dust/Organics | Apr 82 | 63, 64 | 1762 | Steel, Low Alloy | Feb 88 | 20 |
| 1650 | Diesel Particulate Matter | Feb 85 | 63, 64 | 1763 | Steel, Low Alloy | Feb 88 | 20 |
| 1651 | Heat-Source Powder for Calorimetry Zirconium-Barium Chromate, 1460 | Nov 68 | 99 | 1764 | Steel, Low Alloy | Feb 88 | 20 |
| 1652 | Heat-Source Powder for Calorimetry Zirconium-Barium Chromate, 1632 | Nov 68 | 99 | 1765 | Steel, Low Alloy | Feb 88 | 20 |
| 1653 | Heat Source for Calorimetry Zirconium-Barium Chromate, 1762 | Nov 68 | 99 | 1766 | Steel, Low Alloy | Feb 88 | 20 |
| 1655 | Potassium Chloride for Solution Calorimetry | Mar 81 | 99 | 1767 | Steel, Low Alloy | Feb 88 | 20 |
| 1656 | Thianthrene, Combustion Calorimetry | Jan 85 | 99 | 1804 | Ambient Toxic Organics in Nitrogen | * | 56 |
| 1657 | Synthetic Refuse-Derived Fuel, Combustion Calorimetry | Mar 85 | 99 | 1805 | Benzene in Nitrogen, 0.25 ppm | Dec 82 | 54 |
| 1658a | Methane in Air, 1ppm | Mar 81 | 55 | 1806 | Benzene in Nitrogen, 10 ppm | Dec 82 | 54 |
| 1659a | Methane in Air, 10 ppm | Mar 81 | 55 | 1808 | Tetrachloroethylene in Nitrogen, 0.25 ppm | Jun 83 | 57 |
| 1660a | Methane (4) and Propane (1) in Air | Mar 81 | 55 | 1809 | Tetrachloroethylene in Nitrogen, 10 ppm | Jun 83 | 57 |
| 1661a | Sulfur Dioxide in Nitrogen, 500 ppm | Oct 88 | 56 | 1810a | Linerboard | Dec 83 | 128 |
| 1662a | Sulfur Dioxide in Nitrogen, 1000 ppm | Oct 88 | 56 | 1811 | Aromatic Gases in Nitrogen 0.25 ppm | Nov 85 | 54 |
| 1663a | Sulfur Dioxide in Nitrogen, 1500 ppm | Mar 81 | 56 | 1812 | Aromatic Gases in Nitrogen 10 ppm | Nov 85 | 54 |
| 1664a | Sulfur Dioxide in Nitrogen, 2500 ppm | Mar 81 | 56 | 1813 | Aliphatic Organic Gases in Nitrogen 0.25 ppm | Mar 87 | 54 |
| 1665b | Propane in Air, 3 ppm | Dec 87 | 56 | 1814 | Aliphatic Organic Gases in Nitrogen 10 ppm | Mar 87 | 54 |
| 1666b | Propane in Air, 10 ppm | Dec 87 | 56 | 1815a | n-Heptane, Reference Fuel | Mar 85 | 118 |
| 1667b | Propane in Air, 50 ppm | Dec 87 | 56 | 1816a | Isooctane, Reference Fuel | Mar 85 | 118 |
| 1668b | Propane in Air, 100 ppm | Jan 80 | 56 | 1817b | Catalyst Package for Lubricant Oxidation | Oct 86 | 73 |
| 1669b | Propane in Air, 500 ppm | Jan 80 | 56 | 1818 | Chlorine in Lube Base Oil | Apr 86 | 72 |
| 1670 | Carbon Dioxide in Air, 330-ppm | Dec 82 | 54 | 1819 | Sulfur in Lubricating Base Oil | Jul 85 | 72 |
| 1671 | Carbon Dioxide in Air, 340 ppm | Dec 82 | 54 | 1822 | Refractive Index Glass, Soda-Lime | Nov 84 | 108 |
| 1672 | Carbon Dioxide in Air, 350 ppm | Dec 82 | 54 | 1823 | Refractive Index Silicone Liquids | Dec 76 | 108 |
| 1674b | Carbon Dioxide in Nitrogen, 7% | Jan 80 | 54 | 1825 | Glass Density | * | 96 |
| 1675b | Carbon Dioxide in Nitrogen, 14% | Jan 80 | 54 | 1826 | Glass Density | * | 96 |
| 1677c | Carbon Monoxide in Nitrogen, 10 ppm | Jan 80 | 55 | 1827 | Glass Density | * | 96 |
| 1678c | Carbon Monoxide in Nitrogen, 50 ppm | Jan 80 | 55 | 1828 | Ethanol-Water Solutions | Jun 85 | 52 |
| 1679c | Carbon Monoxide in Nitrogen, 100 ppm | Jan 80 | 55 | 1829 | Alcohols in Reference Fuel | Mar 86 | 60 |
| 1680b | Carbon Monoxide in Nitrogen, 500 ppm | Jan 80 | 55 | 1830 | Soda-Lime Float Glass | Jul 82 | 83 |
| 1681b | Carbon Monoxide in Nitrogen, 1000 ppm | Jan 80 | 55 | 1831 | Soda-Lime Sheet Glass | Jul 82 | 83 |
| 1683b | Nitric Oxide in Nitrogen, 50 ppm | Nov 88 | 55 | 1832 | Thin Glass Film on Polycarbonate for X-ray Fluorescence | May 84 | 69 |
| 1684b | Nitric Oxide in Nitrogen, 100 ppm | Nov 88 | 55 | 1833 | Thin Glass Film on Polycarbonate for X-ray Fluorescence | May 84 | 69 |
| 1685b | Nitric Oxide in Nitrogen, 250 ppm | Jan 80 | 55 | 1834 | Fused Ore Glass for XRF | * | 83 |
| 1686b | Nitric Oxide in Nitrogen, 500 ppm | Jan 80 | 55 | 1835 | Borate Ore | Sep 87 | 77 |
| 1687b | Nitric Oxide in Nitrogen, 1000 ppm | Jan 80 | 55 | 1836 | Nitrogen in Lube Base Oil | * | 72 |

*In Prep.

| SRM | SRM Description | Certi- ficate Date | Page | SRM | SRM Description | Certi- ficate Date | Page |
|-------|--|--------------------------|--------|--------|---|--------------------------|--------|
| 1851 | NDE Penetrant Test Block (Matte) | Apr 84 | 123 | 2011 | Gold Mirror, First Surface, Reflectance | * 107 | |
| 1855 | Ultrasonic Power Transducer | Jan 86 | 97 | 2015 | White Opal Glass Diffuse Spectral Reflectance | May 82 | 108 |
| 1856 | Acoustic Emission Transducer | Jul 85 | 97 | 2016 | White Opal Glass Diffuse Spectral Reflectance | May 82 | 108 |
| 1857 | Tool Steel Abrasive Wear | Mar 83 | 115 | 2021 | Black Porcelain Enamel, Directional-Hemispherical Reflect | Sep 80 | 108 |
| 1860 | Al, Eddy Current 60% IACS | Aug 82 | 119 | 2023a | Aluminum Mirror, Second Surface, Reflectance | * 107 | |
| 1862 | Al, Eddy Current 41% IACS | Aug 82 | 119 | 2025 | Aluminum Mirror with Wedge, Second Surface, Reflectance | Feb 82 | 107 |
| 1866 | Bulk Asbestos, Common | Nov 88 | 71 | 2031 | Metal-on-Quartz Filters for Spectrophotometry | Oct 84 | 106 |
| 1867 | Bulk Asbestos, Uncommon | * | 71 | 2032 | Potassium Iodide Stray Light | Oct 79 | 106 |
| 1871 | Pb-Si Glasses for Microanalysis | May 84 | 41, 42 | 2033 | KI Stray Light with Attenuator | May 80 | 106 |
| 1872 | Pb-Ge Glasses for Microanalysis | May 84 | 41, 42 | 2034 | Holmium Oxide Solution Wavelength | Jun 85 | 106 |
| 1873 | Ba-Zn-Si Glasses for Microanalysis | May 84 | 41, 42 | 2063 | Microanalysis Thin Film Mg-Si-Ca-Fe | Aug 87 | 41, 43 |
| 1874 | Li-Al-Bo Glasses for Microanalysis | Dec 84 | 41, 42 | 2064 | Thin Film Microanalysis | * | 41, 43 |
| 1875 | Al-Mg-P Glasses for Microanalysis | Dec 84 | 41, 42 | 2069a | SEM Performance Standard | Feb 85 | 90 |
| 1876a | Chrysotile Asbestos Fibers | Jun 83 | 71 | 2071 | Sinusoidal Roughness | * | 124 |
| 1878 | Respirable Alpha Quartz | Nov 83 | 70 | 2072 | Sinusoidal Roughness | * | 124 |
| 1879 | Respirable Cristobalite | Jan 88 | 70 | 2073 | Sinusoidal Roughness | Nov 84 | 124 |
| 1880 | Portland Cement, black | Feb 84 | 84 | 2074 | Sinusoidal Roughness | * | 124 |
| 1881 | Portland Cement, white | Feb 84 | 84 | 2075 | Sinusoidal Roughness | * | 124 |
| 1882 | Ca-Al Cement | Jul 86 | 84 | 2083 | Socketed Ball Bar | Aug 85 | 123 |
| 1883 | Ca-Al Cement | Jul 86 | 84 | 2092 | Charpy V-Notch Low Energy | Aug 89 | 124 |
| 1884 | Cement Composition | Sep 89 | 84 | 2096 | Charpy V-Notch High Energy | Aug 89 | 124 |
| 1885 | Cement Composition | Sep 89 | 84 | 2106 | Centroid Color Charts | none | 125 |
| 1886 | Cement Composition | Sep 89 | 84 | 2135c | Ni/Cr Thin-Film Depth Profile | * | 91 |
| 1887 | Cement Composition | Sep 89 | 85 | 2136 | Cr/CrO ₂ Depth Profile | * | 91 |
| 1888 | Cement Composition | Sep 89 | 85 | 2137 | Boron Implant in Silicon Depth Profile | * | 91 |
| 1889 | Cement Composition | Sep 89 | 85 | 2141 | Urea | Aug 70 | 44 |
| 1890 | Stainless Steel for Pitting of Crevice Corrosion | May 83 | 116 | 2142 | o-Bromobenzoic Acid | Sep 70 | 44 |
| 1891 | Co-Cr-Mo Alloy for Pitting of Crevice Corrosion | Sep 85 | 116 | 2143 | p-Fluorobenzoic Acid | Jan 82 | 44 |
| 1893 | Cu Microhardness Knoop | Feb 84 | 96 | 2144 | m-Chlorobenzoic Acid | Apr 73 | 44 |
| 1894 | Cu Microhardness Vickers | Feb 84 | 96 | 2151 | Nicotinic Acid (Calorimetry) | Jan 85 | 99 |
| 1895 | Ni Microhardness Knoop | Feb 84 | 96 | 2152 | Urea (Calorimetry) | Jan 85 | 99 |
| 1896 | Ni Microhardness Vickers | Feb 84 | 96 | 2161 | Low Alloy Steel | * | 17 |
| 1901 | Centerline Drawings for Optical Character Recognition-Type B | Mar 76 | 126 | 2162 | Low Alloy Steel | * | 17 |
| 1905 | Ni Microhardness Knoop-300 | Aug 86 | 96 | 2163 | Low Alloy Steel | * | 17 |
| 1906 | Ni Microhardness Knoop-500 | Sep 86 | 96 | 2164 | Low Alloy Steel | * | 17 |
| 1907 | Ni Microhardness Knoop-1000 | Sep 86 | 96 | 2165 | Low Alloy Steel | * | 17 |
| 1920 | Near Infrared Reflectance Wavelength | Jul 86 | 107 | 2166 | Low Alloy Steel | * | 17 |
| 1923 | Polystyrene Sulfonates | * | 97 | 2167 | Low Alloy Steel | * | 17 |
| 1924 | Polystyrene Sulfonates | * | 97 | 2168 | High Purity Iron | * | 89 |
| 1925 | Polystyrene Sulfonates | * | 97 | 2181 | Heptes | Mar 89 | 89 |
| 1930 | Glass Filters, Transmittance | Mar 87 | 106 | 2182 | Heptes Sodium Salt | Mar 89 | 89 |
| 1931 | Fluorescence Corrected Emission Spectrum | * | 106 | 2185 | Potassium Hydrogen Phthalate, pD | Nov 84 | 89 |
| 1939 | Polychlorinated Biphenyls in Sediment | * | 63 | 2186I | Potassium Dihydrogen Phosphate, pD | May 68 | 89 |
| 1940 | Polychlorinated Biphenyls in Sediment | * | 63 | 2186II | Disodium Hydrogen Phosphate, pD | May 68 | 89 |
| 1941 | Organics in Marine Sediment | * | 63 | 2191a | Sodium Bicarbonate, pD | Nov 84 | 89 |
| 1951a | Cholesterol in Human Serum (Frozen) | * | 49 | 2192a | Sodium Carbonate, pD | Nov 84 | 89 |
| 1952a | Cholesterol in Human Serum (Freeze-dried) | * | 49 | 2201 | Sodium Chloride, pNa & pCl | Mar 84 | 90 |
| 1960 | Polystyrene Spheres, 10 µm | Apr 85 | 122 | 2202 | Potassium Chloride, pK & pCl | Mar 84 | 90 |
| 1961 | Polystyrene Spheres, 30 µm | Jan 87 | 122 | 2203 | Potassium Fluoride, pF | May 73 | 90 |
| 1962 | Polystyrene Spheres, 0.3 µm | * | 122 | 2211 | Toluene 8 mL | Mar 85 | 96, |
| 1965 | Polystyrene Spheres, 10 µm | Jan 87 | 122 | 2213 | 2,2,4-Trimethylpentane 25 mL | Mar 85 | 108 |
| 1967 | High-Purity Platinum Thermoelement | Feb 77 | 103 | 2220 | Tin, Temp and Enthalpy of Fusion | Oct 85 | 96, |
| 1968 | Gallium Melting Point | Jun 77 | 102 | 2221a | Zinc, Temp and Enthalpy of Fusion | Oct 85 | 108 |
| 1969 | Rubidium Triple Point | Jan 84 | 102 | 2222 | Biphenyl, Temp and Enthalpy of Fusion | Sept 87 | 100 |
| 1970 | Succinonitrile Triple Point | Mar 85 | 102 | 2223 | Potassium Nitrate, Temp and Enthalpy of Fusion | * | 100 |
| 1971 | Indium Melting Point | Feb 87 | 102 | 2225 | Temperature and Enthalphy of Fusion-Mercury | Mar 89 | 100 |
| 1974 | Organics in Mussel Tissue | * | 63 | 2321 | Solder Thickness | * | 92 |
| 1975 | Diesel Particulate Bioassay | * | 63 | 2322 | Solder Thickness | * | 92 |
| 2003 | Aluminum Mirror, First Surface, Reflectance | May 85 | 107 | | | | |
| 2009a | Didymium Glass Filter, Wavelength | Jul 84 | 106 | | | | |

*In Prep.

| SRM | SRM Description | Certi-ficate Date | Page | SRM | SRM Description | Certi-ficate Date | Page |
|--------|--------------------------------------|-------------------|------|-------|---|-------------------|--------|
| 2350 | Potential & Thickness Step | Aug 85 | 115 | 2636a | Carbon Monoxide in Nitrogen (250 ppm) | Nov 87 | 55 |
| C2400 | High Alloy Steel (ACI 17/4 PH) | Feb 86 | 25 | 2637a | Carbon Monoxide in Nitrogen (2500 ppm) | Mar 87 | 55 |
| C2401 | High Alloy Steel (ACI-CD-4M-Cu) | Feb 86 | 25 | 2638a | Carbon Monoxide in Nitrogen (5000 ppm) | Mar 87 | 55 |
| C2402 | Hasteloy C | Feb 86 | 25 | 2639a | Carbon Monoxide in Nitrogen (1%) | Mar 87 | 55 |
| C2416 | Bullet Lead | Feb 88 | 34 | 2640a | Carbon Monoxide in Nitrogen (2%) | Jul 79 | 55 |
| C2417 | Lead-Base Alloy | Feb 87 | 34 | 2641a | Carbon Monoxide in Nitrogen (4%) | Jul 79 | 55 |
| C2418 | High-Purity Lead | Feb 87 | 34 | 2642a | Carbon Monoxide in Nitrogen (8%) | Mar 87 | 55 |
| C2423 | Ductile Iron | Nov 85 | 28 | 2645a | Propane in Nitrogen (500ppm) | May 80 | 56 |
| C2423a | Ductile Iron | Nov 85 | 28 | 2646a | Propane in Nitrogen (1000 ppm) | May 80 | 56 |
| C2424 | Ductile Iron | Jul 85 | 28 | 2647a | Propane in Nitrogen (2500 ppm) | May 80 | 56 |
| C2424a | Ductile Iron | Jul 85 | 28 | 2648a | Propane in Nitrogen (5000 ppm) | May 80 | 56 |
| C2425 | Ductile Iron | Jul 85 | 28 | 2649 | Propane in Nitrogen (10,000 ppm) | May 80 | 56 |
| C2425a | Ductile Iron | Jul 85 | 28 | 2650 | Propane in Nitrogen (20,000 ppm) | May 80 | 56 |
| C2430 | Scheelite Ore | Jan 87 | 75 | 2651 | Propane and Oxygen in Nitrogen | Jul 80 | 56 |
| 2526 | 111 p-Type Si, Spreading Resistance | Aug 83 | 118 | 2652 | Propane and Oxygen in Nitrogen | Jul 80 | 56 |
| 2527 | 111 n-Type Si, Spreading Resistance | Aug 83 | 118 | 2654 | Nitrogen Dioxide in Air (500 ppm) | Jun 82 | 55 |
| 2528 | 100 p-Type Si, Spreading Resistance | Jan 84 | 118 | 2655 | Nitrogen Dioxide in Air (1000 ppm) | Jun 82 | 55 |
| 2529 | 100 n-Type Si, Spreading Resistance | May 84 | 118 | 2656 | Nitrogen Dioxide in Air (2500 ppm) | Jun 82 | 55 |
| 2531 | Thin Film Thickness | * | 93 | 2657a | Oxygen in Nitrogen (2%) | Nov 87 | 55 |
| 2532 | Thin Film Thickness | * | 93 | 2658a | Oxygen in Nitrogen (10%) | Nov 87 | 55 |
| 2533 | Thin Film Thickness | * | 93 | 2659a | Oxygen in Nitrogen (20%) | Nov 87 | 55 |
| 2534 | Thin Film Thickness | * | 93 | 2670 | Toxic Metals in Freeze-Dried Urine | Mar 85 | 68 |
| 2541 | Silicon Resistivity | * | 118 | 2671a | Freeze-Dried Urine for Fluorine | Dec 82 | 68 |
| 2542 | Silicon Resistivity | * | 118 | 2672a | Freeze-Dried Urine for Mercury | May 83 | 68 |
| 2543 | Silicon Resistivity | * | 118 | 2676c | Metals on Filter Media (Cd-Mn-Pb-Zn) | Feb 87 | 69 |
| 2544 | Silicon Resistivity | * | 118 | 2677 | Be & As on Filter Media | Oct 85 | 69 |
| 2546 | Silicon Resistivity | * | 118 | 2678 | Membrane Blank Filter | May 88 | 69 |
| 2547 | Silicon Resistivity | * | 118 | 2679a | Quartz on Filter Media | May 84 | 69 |
| 2548 | Silicon Resistivity | * | 118 | 2680 | Membrane Blank Filter | May 88 | 69 |
| 2549 | Silicon Resistivity | * | 118 | 2681 | Ashless Blank Filter | May 88 | 69 |
| 2550 | Silicon Resistivity | * | 118 | 2682 | Sulfur in Coal (0.5) | Feb 85 | 69, 99 |
| 2607 | Carbon Dioxide/Nitrous Oxide in Air | Sep 85 | 54 | 2683a | Sulfur in Coal (1.9) | Feb 85 | 61, 99 |
| 2608 | Carbon Dioxide/Nitrous Oxide in Air | Sep 85 | 54 | 2684a | Sulfur in Coal (3.0) | Feb 85 | 61, 99 |
| 2609 | Carbon Dioxide/Nitrous Oxide in Air | Sep 85 | 54 | 2685 | Sulfur in Coal (4.6) | Feb 85 | 61, 99 |
| 2610 | Carbon Dioxide/Nitrous Oxide in Air | Sep 85 | 54 | 2689 | Coal Fly Ash | Oct 86 | 62 |
| 2612a | Carbon Monoxide in Air (10 ppm) | Jan 80 | 55 | 2690 | Coal Fly Ash | Oct 86 | 62 |
| 2613a | Carbon Monoxide in Air (20 ppm) | Jan 80 | 55 | 2691 | Coal Fly Ash | Oct 86 | 62 |
| 2614a | Carbon Monoxide in Air (45 ppm) | Jan 80 | 55 | 2692 | Sulfur in Coal (0.1) | Nov 88 | 61 |
| 2619a | Carbon Dioxide in Nitrogen (0.5%) | Jan 80 | 54 | 2694 | Simulated Rainwater | * | 59 |
| 2620a | Carbon Dioxide in Nitrogen (1.0%) | Jan 80 | 54 | 2695 | FL Vegetation | * | 53 |
| 2621a | Carbon Dioxide in Nitrogen (1.5%) | Jan 80 | 54 | 2696 | FL Vegetation | * | 53 |
| 2622a | Carbon Dioxide in Nitrogen (2.0%) | Jan 80 | 54 | 2704 | Buffalo River Sediment | Jun 88 | 62 |
| 2623a | Carbon Dioxide in Nitrogen (2.5%) | Jan 80 | 55 | 2709 | Agricultural Soil (Baseline Levels) | * | 79 |
| 2624a | Carbon Dioxide in Nitrogen (3.0%) | Jan 80 | 55 | 2710 | Agricultural Soil (Elevated Levels) | * | 79 |
| 2625a | Carbon Dioxide in Nitrogen (3.5%) | Jan 80 | 55 | 2712 | Lead in Reference Fuel | Sep 88 | 59 |
| 2626a | Carbon Dioxide in Nitrogen (4.0%) | Jan 80 | 55 | 2713 | Lead in Reference Fuel | Sep 88 | 59 |
| 2627a | Nitric Oxide in Nitrogen (5 ppm) | Jun 82 | 55 | 2714 | Lead in Reference Fuel | Sep 88 | 59 |
| 2628a | Nitric Oxide in Nitrogen (10 ppm) | Jun 82 | 55 | 2715 | Lead in Reference Fuel | Sep 88 | 59 |
| 2629a | Nitric Oxide in Nitrogen (20 ppm) | Jun 82 | 55 | 2717 | Sulfur in Residual Fuel Oil | * | 61 |
| 2630 | Nitric Oxide in Nitrogen (1500 ppm) | May 79 | 55 | 2720 | Synthetic Sulfur in Oil | * | 61 |
| 2631 | Nitric Oxide in Nitrogen (3000 ppm) | May 79 | 55 | 2725 | Auto Emission Inspection Maintenance Gas 1.6 mol% | Jan 89 | 55 |
| 2635a | Carbon Monoxide in Nitrogen (25 ppm) | Dec 88 | 55 | 2726 | Auto Emission Inspection Maintenance Gas | Jan 89 | 55 |

*In Prep.

| SRM | SRM Description | Certi- cate Date | Page | SRM | SRM Description | Certi- cate Date | Page |
|------|--|------------------------|------|-------|---|------------------------|------|
| 2727 | Auto Emission Maintenance Gas Inspection | Jan 89 | 55 | 3144 | Rhodium Spectrometric Solution | * | 45 |
| 2728 | Auto Emission Maintenance Gas Inspection | Jan 89 | 55 | 3145 | Rubidium Spectrometric Solution | Nov 86 | 45 |
| 2730 | Hydrogen Sulfide in Nitrogen | May 89 | 56 | 3146 | Ruthenium Spectrometric Solution | * | 46 |
| 2731 | Hydrogen Sulfide in Nitrogen | May 89 | 56 | 3147 | Samarium Spectrometric Solution | Mar 87 | 46 |
| 3101 | Aluminum Spectrometric Solution | Nov 86 | 45 | 3148 | Scandium Spectrometric Solution | Mar 87 | 46 |
| 3102 | Antimony Spectrometric Solution | Nov 86 | 45 | 3149 | Selenium Spectrometric Solution | Dec 86 | 46 |
| 3103 | Arsenic Spectrometric Solution | Nov 86 | 45 | 3150 | Silicon Spectrometric Solution | Dec 86 | 46 |
| 3104 | Barium Spectrometric Solution | Nov 86 | 45 | 3151 | Silver Spectrometric Solution | Dec 86 | 46 |
| 3105 | Beryllium Spectrometric Solution | Nov 86 | 45 | 3152 | Sodium Spectrometric Solution | Nov 86 | 46 |
| 3106 | Bismuth Spectrometric Solution | Dec 86 | 45 | 3153 | Strontium Spectrometric Solution | Nov 86 | 46 |
| 3107 | Boron Spectrometric Solution | Dec 86 | 45 | 3154 | Sulfur Spectrometric Solution | Aug 87 | 46 |
| 3108 | Cadmium Spectrometric Solution | Dec 86 | 45 | 3155 | Tantalum Spectrometric Solution | Jan 88 | 46 |
| 3109 | Calcium Spectrometric Solution | Nov 86 | 45 | 3156 | Tellurium Spectrometric Solution | Aug 87 | 46 |
| 3110 | Cerium Spectrometric Solution | Mar 87 | 45 | 3157 | Terbium Spectrometric Solution | Mar 87 | 46 |
| 3111 | Cesium Spectrometric Solution | Feb 87 | 45 | 3158 | Thallium Spectrometric Solution | Dec 86 | 46 |
| 3112 | Chromium Spectrometric Solution | Nov 86 | 45 | 3159 | Thorium Spectrometric Solution | Jan 88 | 46 |
| 3113 | Cobalt Spectrometric Solution | Dec 86 | 45 | 3160 | Thulium Spectrometric Solution | Mar 87 | 46 |
| 3114 | Copper Spectrometric Solution | Dec 86 | 45 | 3161 | Tin Spectrometric Solution | Nov 86 | 46 |
| 3115 | Dysprosium Spectrometric Solution | Mar 87 | 45 | 3162 | Titanium Spectrometric Solution | Nov 86 | 46 |
| 3116 | Erbium Spectrometric Solution | Mar 87 | 45 | 3163 | Tungsten Spectrometric Solution | Dec 86 | 46 |
| 3117 | Europium Spectrometric Solution | Mar 87 | 45 | 3164 | Uranium Spectrometric Solution | Dec 86 | 46 |
| 3118 | Gadolinium Spectrometric Solution | Mar 87 | 45 | 3165 | Vanadium Spectrometric Solution | Dec 86 | 46 |
| 3119 | Gallium Spectrometric Solution | Mar 87 | 45 | 3166 | Ytterbium Spectrometric Solution | Mar 87 | 46 |
| 3120 | Germanium Spectrometric Solution | Mar 88 | 45 | 3167 | Yttrium Spectrometric Solution | Mar 87 | 46 |
| 3121 | Gold Spectrometric Solution | Nov 86 | 45 | 3168 | Zinc Spectrometric Solution | Nov 86 | 46 |
| 3122 | Hafnium Spectrometric Solution | Apr 88 | 45 | 3169 | Zirconium Spectrometric Solution | Dec 86 | 46 |
| 3123 | Holmium Spectrometric Solution | Mar 87 | 45 | 3171 | Multielement Mix A Solution | Jul 88 | 47 |
| 3124 | Indium Spectrometric Solution | Dec 86 | 45 | 3172 | Multielement Mix B Solution | Jul 88 | 47 |
| 3125 | Iridium Spectrometric Solution | * | 45 | 3173 | Multielement Mix C Solution | * | 47 |
| 3126 | Iron Spectrometric Solution | Nov 86 | 45 | 3174 | Multielement Mix D Solution | Mar 89 | 47 |
| 3127 | Lanthanum Spectrometric Solution | Mar 87 | 45 | 3181 | Sulfate Anion Solution | Jan 87 | 48 |
| 3128 | Lead Spectrometric Solution | Dec 86 | 45 | 3182 | Chloride Anion Solution | Apr 87 | 48 |
| 3129 | Lithium Spectrometric Solution | Nov 86 | 45 | 3183 | Fluoride Anion Solution | Apr 87 | 48 |
| 3130 | Lutetium Spectrometric Solution | Mar 87 | 45 | 3184 | Bromide Anion Solution | * | 48 |
| 3131 | Magnesium Spectrometric Solution | Nov 86 | 45 | 3185 | Nitrate Anion Solution | * | 48 |
| 3132 | Manganese Spectrometric Solution | Dec 86 | 45 | 3186 | Phosphate Anion Solution | * | 48 |
| 3133 | Mercury Spectrometric Solution | Dec 86 | 45 | 3191 | Electrolytic Conductance | May 87 | 119 |
| 3134 | Molybdenum Spectrometric Solution | Nov 86 | 45 | 3192 | Electrolytic Conductance | May 87 | 119 |
| 3135 | Neodymium Spectrometric Solution | Mar 87 | 45 | 3193 | Electrolytic Conductance | May 87 | 119 |
| 3136 | Nickel Spectrometric Solution | Dec 86 | 45 | 3194 | Aqueous Electrolytic Conductance | Mar 89 | 119 |
| 3137 | Niobium Spectrometric Solution | Jan 88 | 45 | 3195 | Aqueous Electrolytic Conductance | Mar 89 | 119 |
| 3138 | Palladium Spectrometric Solution | Nov 86 | 45 | 3200 | Secondary Standard Magnetic Tape (Computer Amplitude Ref) | May 81 | 125 |
| 3139 | Phosphorus Spectrometric Solution | Nov 86 | 45 | 3216 | Secondary Standard Magnetic Tape Cartridge (Computer Amplitude Ref) | Aug 82 | 126 |
| 3140 | Platinum Spectrometric Solution | Nov 86 | 45 | 3217 | Secondary Standard Magnetic Tape Cartridge-High Density (C A Ref) | Jul 87 | 126 |
| 3141 | Potassium Spectrometric Solution | Nov 86 | 45 | 4200B | Cesium-137/Barium-137m Point Source | Dec 79 | 112 |
| 3142 | Praseodymium Spectrometric Solution | Mar 87 | 45 | 4201B | Niobium-94 Gamma-ray | Jun 70 | 112 |
| 3143 | Rhenium Spectrometric Solution | Jul 88 | 45 | 4202D | Cadmium-109-Silver-109m | Oct 86 | 112 |
| | | | | 4203D | Cobalt-60 Point Source | Feb 84 | 112 |
| | | | | 4207B | Cesium-137/Barium-137m Point Source | Mar 87 | 112 |
| | | | | 4218E | Europium-152 Point Source | Nov 82 | 112 |
| | | | | 4226B | Nickel-63 Solution | Dec 84 | 110 |
| | | | | 4233C | Cesium-137 Burn-up Standard | Nov 79 | 110 |

*In Prep.

| SRM | SRM Description | Certi-ficate Date | Page | SRM | SRM Description | Certi-ficate Date | Page |
|--------|--|-------------------|----------|-------|--|-------------------|----------|
| 4235C | Krypton-85 Gaseous | Oct 86 | 111 | 4952B | Radium Standard Blank Solution | Dec 60 | 112 |
| 4241B | Barium-133 Point Source | Apr 82 | 112 | 4953D | Radium-226 Solution | May 84 | 112 |
| 4251B | Barium-133 Solution | Dec 81 | 110 | 4956 | Radium-226 Gamma-ray Solution 0.2 μ g | Mar 68 | 113 |
| 4264B | Tin-121m Point-Source Gamma-ray | Sep 82 | 112 | 4957 | Radium-226 Gamma-ray Solution 0.5 μ g | Mar 68 | 113 |
| 4267 | Niobium-93m Point Source | Oct 85 | 112 | 4958 | Radium-226 Gamma-ray Solution 1.0 μ g | Mar 68 | 113 |
| 4275C | Mixed Radionuclide Point Source | Jul 83 | 112 | 4959 | Radium-226 Gamma-ray Solution 2.0 μ g | Mar 68 | 113 |
| 4276C | Mixed Radionuclide Solution | Jul 83 | 110 | 4990C | Carbon-14 Oxalic Acid | Jul 83 | 111 |
| 4288 | Technetium-99 Solution | Nov 82 | 110 | 6250 | Secondary Standard High Density Magnetic Tape (Comp Amp Ref) | May 82 | 125 |
| 4308C | Krypton-85 Gaseous | Jan 83 | 111 | 6596 | Flexible Disk Cartridge | PTB | 126 |
| 4320 | Curium-244 | Mar 89 | 110 | 8000 | Melting Point Set (NPL CRM M14-11) | NPL | 103 |
| 4321 | Uranium-238 Solution | Nov 86 | 110 | 8005 | Alpha Alumina (Surface Area) | NPL | 103, 123 |
| 4322 | Americium-241 Solution | Nov 86 | 110 | 8006 | Alpha Alumina (Surface Area) | NPL | 123 |
| 4323 | Plutonium-238 Solution | Nov 86 | 110 | 8007 | Alpha Alumina (Surface Area) | NPL | 123 |
| 4324 | Uranium-232 Alpha-particle Solution | May 84 | 110 | 8008 | Alpha Alumina (Surface Area) | NPL | 123 |
| 4326 | Polonium-209 | no date | 110 | 8009 | Automated Computer Time Service | Nov 88 | 127 |
| 4327 | Polonium-208 Alpha-particle Solution | Jan 85 | 110 | 8406 | Mercury in Soil (110 μ g/g) | * | 79 |
| 4328 | Thorium-229 Alpha-particle Solution | May 85 | 110 | 8407 | Mercury in Soil (50 μ g/g) | * | 79 |
| 4329 | Curium-243 Alpha-particle Solution | Mar 85 | 110 | 8408 | Mercury in Soil (0.07 μ g/g) | * | 79 |
| 4332C | Americium-243 Alpha-particle Solution | Feb 84 | 110 | 8410 | Asbestos Research Filter | Apr 84 | 71 |
| 4338 | Plutonium-240 Alpha-particle Solution | Aug 80 | 110 | 8411 | Mixed Asbestos Research Filter | Nov 88 | 71 |
| 4350B | Environmental Radioactivity, River Sediment | Sep 81 | 113 | 8412 | Corn (Zea Mays) Stalk | undated | 53 |
| 4351 | Environmental Radioactivity, Human Lung | Oct 82 | 113 | 8413 | Corn (Zea Mays) Kernel | undated | 53 |
| 4352 | Environmental Radioactivity, Human Liver | Jun 82 | 113 | 8420 | Electrolytic Iron, Thermal Conductivity and Electrical Resistivity | May 84 | 104, 118 |
| 4353 | Environmental Radioactivity, Rocky Flats Soil Number 1 | Dec 80 | 113 | 8421 | Electrolytic Iron, Thermal Conductivity and Electrical Resistivity | May 84 | 104, 118 |
| 4354 | Freshwater Lake Sediment (Gyttja) | Nov 86 | 113 | 8422 | Sintered Tungsten, Thermal Conductivity and Electrical Resistivity | May 84 | 104, 118 |
| 4355 | Environmental Radioactivity, Peruvian Soil | Jun 82 | 113 | 8423 | Sintered Tungsten, Thermal Conductivity and Electrical Resistivity | May 84 | 104, 118 |
| 4361B | Hydrogen-3 Solution | Jan 81 | 110, 111 | 8424 | Graphite, Thermal Conductivity and Electrical Resistivity | May 84 | 104 |
| 4370C | Europium-152 Solution | Mar 87 | 110 | 8425 | Graphite, Thermal Conductivity and Electrical Resistivity | May 84 | 104 |
| 4400L | Chromium-51 Solution | ** | 114 | 8426 | Graphite, Thermal Conductivity and Electrical Resistivity | May 84 | 104 |
| 4401L | Iodine-131 Solution | ** | 114 | 8430 | Aspartate Aminotransferase (AST) Human Erythrocyte Source | Jun 87 | 49 |
| 4404L | Thallium-201 Solution | ** | 114 | 8443 | GC/MS System Performance | Aug 84 | 67 |
| 4405LB | Gold-198 | ** | 114 | 8444 | Cotinine in Freeze-Dried Human Urine | Feb 89 | 68 |
| 4406L | Phosphorus-32 Solution | ** | 114 | 8450 | Polyethylene Piping | Jan 88 | 97 |
| 4407L | Iodine-125 Solution | ** | 114 | 8451 | Polyethylene Piping | Jan 88 | 97 |
| 4410H | Technetium-99m Solution | ** | 114 | 8452 | Polyethylene Piping | Jan 88 | 97 |
| 4412L | Molybdenum-99 Solution | ** | 114 | 8453 | Polyethylene Pipe Socket T Joint | Jan 88 | 97 |
| 4414LC | Iodine-123 | ** | 114 | 8454 | Polyethylene Pipe Butt T | Jan 88 | 97 |
| 4415L | Xenon-133 Gaseous | ** | 111 | 8458 | Artificial Flaw for Eddy Current NDE Calibration | * | 124 |
| 4416L | Gallium-67 Solution | ** | 114 | 8486 | Portland Cement Clinker | May 89 | 85 |
| 4417L | Indium-111 Solution | ** | 114 | 8487 | Portland Cement Clinker | May 89 | 85 |
| 4420LB | Lead-203 | ** | 114 | 8488 | Portland Cement Clinker | May 89 | 85 |
| 4421L | Gold-195 | ** | 114 | 8505 | Vanadium in Crude Oil | — | 59 |
| 4423 | Strontium-90 | ** | 110 | 8531 | Glass Fibers for Microanalysis | Jan 89 | 43 |
| 4904G | Americium-241 Alpha-particle | Apr 82 | 111 | 8570 | Calcined Kaolin (Surface Area) | ASTM | 123 |
| 4906C | Plutonium-238 Point Source | * | 111 | 8571 | Alumina (Surface Area) | ASTM | 123 |
| 4915D | Cobalt-60 Solution | Feb 84 | 110 | 8572 | Silica-Alumina (Surface Area) | ASTM | 123 |
| 4919F | Strontium-99 Solution | May 83 | 110 | 8580 | Nickel on Fresh Alumina | * | 73 |
| 4926D | Hydrogen-3 Tritiated Water | Jan 79 | 110 | 8581 | Nickel on Fresh Alumina | * | 73 |
| 4927D | Hydrogen-3 Tritiated Water | Mar 85 | 110 | 8582 | Nickel on Fresh Alumina | * | 73 |
| 4929D | Iron-55 X-ray Solution | Jul 85 | 110 | 8583 | Platinum on Fresh Alumina | * | 73 |
| 4935C | Krypton-85 Beta-particle Gaseous | Jul 74 | 111 | 8584 | Platinum on Fresh Alumina | * | 73 |
| 4940C | Promethium-147 Beta-particle Solution | Aug 85 | 110 | 8585 | Kaolin Based Relative XRD | * | 116 |
| 4943 | Chlorine-36 Beta-particle Solution | Dec 84 | 110 | 8586 | FCC Catalyst | * | 73 |
| 4947C | Hydrogen-3 Tritiated Toluene | Apr 79 | 110 | 8587 | FCC Catalyst | * | 73 |
| 4949B | Iodine-129 Solution | Feb 82 | 110 | 8588 | Faujasite Y | * | 73 |
| 4950E | Radium-226 Solution | May 84 | 112 | | | | |

*In Prep.

**Subscription.

| SRM | SRM Description | Certifi- cate Date | Page | SRM | SRM Description | Certifi- cate Date | Page |
|--------|--|--------------------------|------|--------|--------------------------------------|--------------------------|------|
| 8589 | FCC Catalyst | * | 73 | GM 758 | ICTA Temp Set 125-435 C (DTA&DSC) | ICTA | 101 |
| 8590 | High Sulfur Gas Oil Feed | * | 61 | GM 759 | ICTA Temp Set 295-675 C (DTA&DSC) | ICTA | 101 |
| 9111 | SRM Compilation | Nov 87 | 132 | GM 760 | ICTA Temp Set 570-940 C (DTA&DSC) | ICTA | 101 |
| 9529 | Flexible Disk Cartridge | PTB | 126 | GM 761 | ICTA Thermogravimetry Set | ICTA | 101 |
| GM 754 | ICTA Polystyrene (DTA & DSC) | ICTA | 101 | RM 1R | Aluminum, Ultra-Purity (rod) | — | 40 |
| GM 757 | ICTA Temp Set Below 350 K (DTA&DSC) | ICTA | 101 | RM 50 | Albacore Tuna | May 77 | 51 |

Alphabetical Index

ABRASIVE WEAR, 115

Tool Steel

ABSORBANCE, see Spectrophotometry

ACETANILIDE, 44

Microchemical

ACIDIMETRIC, 44

(See also, Spectrometry)

Benzoic Acid

Boric Acid

Potassium Hydrogen Phthalate

ACID RAIN, 59

ACOUSTIC EMISSION TRANSDUCER, 97

ADHESION TESTING TAPE, 128

Linerboard

AGRICULTURE, 53

(See also, Food and Beverage)

Ammonium Dihydrogen Phosphate
(Fertilizer), 74

Bovine Liver, 51

Brewers Yeast, 51

Citrus Leaves, 53

Coconut Oil, 52

Corn, 53

Cystine, 44

Dextrose, 44, 108

Mercury in Water, 59

Milk, Non-fat Powdered, 51

Oyster Tissue, 51

Phosphate Rock (Fertilizer), 74, 77

Pine Needles, 53

Potassium Dihydrogen Phosphate
(Fertilizer), 74

Potassium Nitrate (Fertilizer), 74

Rice Flour, 51

Sucrose, 44, 108

Tomato Leaves, 53

Trace Elements in Water, 62

Tuna, Albacore, 51

Wheat Flour, 51

Wine, Stabilized, 52

AIR PARTICULATES, 62

AIR POLLUTION, 54

ALCOHOL

Alcohol in Ref. Fuel, 60

Ethanol in Ref. Fuel, 60

Methanol-Butanol in Ref. Fuel, 60

Methanol in Ref. Fuel, 60

Stabilized Wine, 52

ALLOYS

(See also individual metals)

Ferrous, 15

High Temperature, 30

Nonferrous, 30

Steelmaking, 26

ALPHA QUARTZ, 70

ALUMINA

Melting Point, 102

Reduction Grade, 77

Surface Area, 103, 123

ALUMINUM

Alloys-Composition, 30

Conductivity, 104, 118

First Surface Mirror, 107

Freezing Point, 102

High-Purity, 30, 40

Magnetic Susceptibility, 105

Metallo-Organic, 72

Mirrors, Specular Reflectance, 107

Residual Resistivity Ratio, 119

Spectrometric Solution, 45

Specular Reflectance Mirrors, 107

Ultra-Purity, 40

AMERICIUM, 109-114

Radioactivity

AMMONIUM DIHYDROGEN PHOSPHATE, 74

Fertilizer

ANALYZED GASES, see Gases

ANGIOTENSIN I (Human), 49, 50

Clinical

ANISIC ACID, 44

Microchemical

ANION SOLUTIONS, 48

Bromide

Chloride

Fluoride

Nitrate

Phosphate

Sulfate

ANTICONVULSANT DRUG LEVEL

ASSAY, 49

Clinical

ANTIEPILEPSY DRUG LEVEL ASSAY, 49

Clinical

ANTIMONY,

Radioactivity, 110

Spectrometric Solution, 45

ARGILLACEOUS LIMESTONE, 79

ARSENIC

-on Filter Media, 68

Spectrometric Solution, 45

-Trioxide, Reductometric, 44

ASBESTOS, 71

Bulk-Common, Uncommon

Chrysotile Fibers

Research Filter

ATOMIC ABSORPTION, see Spectrometric Solutions

AUSTENITE, 115

in Ferrite

- BARIUM**
 Metallo-Organic, 72
 Radioactivity, 109–114
 Spectrometric Solution, 45
- BASALT**, 80
 Rock
- BASIMETRIC**, 44
- BAUXITE**, 77
- BEARING METAL**, 34
 Lead Alloy
- BENZOIC ACID**
 Acidimetric, 44
 Calorimetric, 99
- BERYLLIUM**
 Copper Alloys, 32
 -on Filter Media, 69
 Spectrometric Solution, 45
- BET SURFACE AREA**, 123
- BEVERAGE**, see Food and Beverage
- BILIRUBIN**, 49
 Clinical
- BIOLOGICAL**, 51–53
 Agricultural
 Ethanol Solutions
 Food and Beverage
- BIPHENYL**, 100
 Enthalpy
- BISMUTH**
 Spectrometric Solution, 45
- BLACK PORCELAIN ENAMEL**, 108
 Reflectance
- BLOOD GAS**, 49, 54
 Analyzed Gases
- BORATE ORE**, 77
- BORON**
 Boric Acid (Assay & Isotopic), 44, 88
 Boric Acid, Enriched B-10, 88
 Depth Profile, 91
 Implant in Silicon, 91
 Spectrometric Solution, 45
- BOTANICAL**, see Biological
- BOVINE**
 -Liver, 51
 -Serum, Inorganic Constituents, 49
- BOVINE SERUM ALBUMIN**, 49
 Clinical
- BRASS**, see Copper Alloys
- BREWERS YEAST**, 51
- BROMIDE**
 Anion Solution, 48
- BROMINE**, 88
 Isotopic
 o-Bromobenzoic Acid
- BRONZE**, see Copper Alloys
- BURNT REFRactories**, 79
- BUTYL**, 121
 Rubber
- CADMIUM**
 Metallo-Organic, 72
 -on Filter Media, 69
 Spectrometric Solution, 45
 Vapor Pressure, 102
- CALCIUM**
 Carbonate (Clinical), 49
 Molybdate, 26
 Spectrometric Solution, 45
- CALORIMETRY**
 (See also, Fuels & Fossil Fuels)
 Combustion, 99
 Enthalpy, 97, 98, 100
 Heat Capacity, 99, 100
 Heating Value, 99, 100
 Heat Source, 99
 Solution, 99
- CARBIDES**, 81
 Cemented, 82
 Silicon
 Tungsten
- CARBON**
 (See also, Rubber Compounding Materials)
 Carbon Dioxide, 54
 Carbon Monoxide, 55
 Radioactivity, 109–114
 Steels, 15
- CAST IRON**, 27
- CAST STEEL**, 28
- CATALYST CHARACTERIZATION MATERIALS**, 73
- CATALYST PACKAGE FOR LUBRICANT OXIDATION**, 73
- CEMENTS**
 Portland, Composition, 84
 Portland, Fineness, 122
- CENTERLINE DRAWINGS FOR OCR**, 126
- CERAMIC MATERIALS**
 Carbides, 81
 Glasses, 83
 Minerals, 79
 Reflectance, 107
 Refractories, 80
- CERIUM**
 Spectrometric Solution, 45
- CESIUM**
 Radioactivity, 110
 Spectrometric Solution, 45
- CHANNEL BLACK**, 121
 Rubber Compounding Material
- CHARPY V-NOTCH TEST BLOCKS**, 124
- CHEMICAL**
 High Purity Metals, 40
 Microchemical, 44
 Primary, Working, and Secondary, 44
- CHLORIDE**
 Anion Solution, 48
- CHLORINE**
 m-Chlorobenzoic Acid, 44
 in Lube Base Oil, 72
 Isotopic, 88
 Radioactivity, 109–114
- CHOLESTEROL**, 49
 -and Fat Soluble Vitamins in Coconut Oil
 Clinical
 in Egg Powder
 in Human Serum
- CHROMIUM**
 Chrome Refractory, 79
 Cr/CrO₂ Depth Profile, 91
 Isotopic, 88

- Metallo-Organic, 72
 Radioactivity, 109–114
 Spectrometric Solution, 45
 Steel, 16
- CHRYSOTILE**, 71
 Asbestos
- CITRUS LEAVES**, 53
- CLAYS**, 78
 Brick
 Flint
 Plastic
- CLINICAL LABORATORY**, 49
 Angiotensin I (Human)
 Anticonvulsant Drug Level Assay
 Antiepilepsy Drug Level Assay
 AST Human Erythrocyte Source
 Bilirubin
 Blood Gases
 Bovine Serum Albumin
 Bovine Serum, Inorganic Constituents
 Calcium Carbonate
 Cholesterol
 Cortisol (Hydrocortisol)
 Creatinine
 D-Glucose (Dextrose)
 Human Serum
 4-Hydroxy-3-methoxy-D₁-mandelic Acid (VMA)
 Iron Metal
 Lead Nitrate
 Lithium Carbonate
 Magnesium Gluconate Dihydrate
 D-Mannitol
 4-Nitrophenol
 Potassium Chloride
 Sodium Chloride
 Sodium Pyruvate
 Tris(hydroxymethyl)aminomethane, pH
 Tris(hydroxymethyl)aminomethane hydrochloride, pH
 Urea
 Uric Acid
- COAL**
 Calorimetric, 61, 99
 Sulfur in, 61
 Trace Elements in, 62
 Trace Mercury in, 59, 62
- COAL FLY ASH**, 62
 Trace Elements in
- COATING THICKNESS**, 92
 Cu-Cr Coating (nonmagnetic) on Steel
 Nickel (magnetic) on Steel
 Solder Thickness Standards, 92
- COATING WEIGHT**, 93
 Gold on Nickel
 Ultra-Thin Gold on Nickel
- COBALT**
 Metallo-Organic, 72
 Radioactivity, 109–114
 Spectrometric Solution, 45
- COCAINE METABOLYTE**, 68
- COCONUT OIL**, 52
- COD LIVER OIL**, 63
- Co-Cr-Mo ALLOY**, 116
 Pitting/Crevice Corrosion
- COLOR**, 125
 Centroid Color Charts
- COMPUTER, MAGNETIC STORAGE MEDIA**, 125
- CONDUCTANCE (AQUEOUS)**, 119
- CONDUCTIVITY**,
 Aluminum, 119
 Electrical, 118
 Glass, 94
 Graphite, 104
 Iron, 104, 118
 Silicon, 118
 Steel, 104, 118
 Thermal, 104
 Tungsten, 104, 118
- COPPER**
 Alloys, 31
 Benchmarks (Unalloyed Cu), 33
 Brass, 32
 Bronze, 32
 Freezing Point, 102
 Gilding Metal, 32
 Heat Capacity Test Specimen, 100
 Isotopic, 88
 Metallo-Organic, 72
 Microprobe, 41, 42
 Ores, 75
 Spectrometric Solution, 45
- CORN**, 53
- CORROSION**, 115
 Implants
 Pitting or Crevice Corrosion
 Potential and Thickness Step Test
- CORTISOL**, 49
 Clinical
- COTININE IN URINE**, 68
- CREATININE**, 49
 Clinical
- CRUDE OIL**, 59
- CURIUM**, 109–114
 Radioactivity
- CYSTINE**, 44
 Microchemical
- DENSITY**
 Liquids, 96
 Neutral Filters, 106
 Neutron, 87
 Photographic, 125
 Silicon, 96
 Smoke, 128
 X-ray, 125
- DEPTH PROFILE**, 91
 Boron Implant in Silicon
 Cr/CrO₂ Thin-Film
 Ni/Cr Thin-Film
- DEXTROSE**, see D-Glucose
- DIDYMIUM**
 Wavelength, 106

DIFFERENTIAL SCANNING CALORIMETRY, 100

DIFFERENTIAL THERMAL ANALYSIS, 101

DIFFRACTION, X-RAY, 116

DIOXIN IN ISOOCTANE, 63, 67

DISODIUM HYDROGEN PHOSPHATE, 89
pH, pD

DOLOMitic LIMESTONE, 79

DOSIMETRY, 87

DRUG LEVEL ASSAY, see Clinical

DRUGS OF ABUSE, IN URINE, 68

DSC, see Heat

DTA, see Heat

DUST
Urban, 63, 64

DYE PENETRANT CRACK BLOCK, 123

DYSPROSIUM
Spectrometric Solution, 45

EDDY CURRENT
Aluminum, 119
Artificial Flaw, 124
Silicon, 118
Test Block, 123

ELECTRICAL CONDUCTIVITY, 104, 118

ELECTRICAL RESISTIVITY, 104, 118
Graphite
Iron
“RRR” Set
Silicon
Tungsten

ELECTROLYtic CONDUCTANCE, 119

ELECTRONIC AND MAGNETIC ALLOY, 35

ELLIPSOMETRY THIN FILM, 93

ENGINEERING MATERIALS, 121

ENTHALPY, see Heat

ENVIRONMENTAL
Diesel Particulate Matter, 63, 64
Estuarine Sediment, 62
Radioactivity, Natural Matrices, 113
River Sediment, 62
Urban Particulate Matter, 62

ENZYME, 49
Aspartate Aminotransferase Human
Erythrocyte Source

ERBIUM
Spectrometric Solution, 45

ESTUARINE SEDIMENT, 62

ETHANOL SOLUTIONS, 60
Ethanol-Water Solutions

EUROPIUM
Radioactivity, 109–114
Spectrometric Solution, 45

FELDSPAR, 79
Potash

Soda

FERRITE, 115

FERRO-ALLOYS, see Steelmaking Alloys

FERTILIZERS, 74
Ammonium Dihydrogen Phosphate
Phosphate Rock
Potassium Dihydrogen Phosphate
Potassium Nitrate

FILTER MEDIA, 69
Be & As on
Metals on (Cd-Mn-Pb-Zn)
Quartz on
Toxic Metals on

FILTERS, 106

FINENESS, 122
Portland Cement

FIRE RESEARCH, 128

FISSION TRACK, 87

FLAMMABILITY, 128

FLEXIBLE DISK CARTRIDGE
Magnetic Computer Storage Media, 125

FLOORING RADIANT PANEL, 128

FLOUR
Rice, 51
Wheat, 51

FLUORESCENCE
Glass, 116
Quinine Sulfate Dihydrate, 106
X-ray, Glass Target, 116

FLUORINE, 68
in Freeze-Dried Urine
p-Fluorobenzoic Acid

FLUORSPAR, 75
Customs Grade
High Grade

FLY ASH, COAL, 62
Lime
Trace Elements

FOOD AND BEVERAGE
(See also, Agriculture)
Coconut Oil, 52
Cod Liver Oil, 52
Corn, 53
Liver, Bovine, 51
Non-Fat Milk Powder, 51
Oyster Tissue, 51
Rice Flour, 51
Tuna, Albacore, 51
Wheat Flour, 51
Wine, Stabilized, 52
Yeast, Brewers, 51

FOSSIL FUELS
(See also, Fuels)
Chlorine in Lube Base Oil, 72
Heating Values, 99
Lime in Coal Fly Ash, 62
Nitrogen in Lube Base Oil, 72
Sulfur in Coal, 61
Sulfur in Distillate Fuel Oil, 61
Sulfur in Kerosine, 61
Sulfur in Residual Fuel Oil, 61
Trace Elements in Coal, 62
Trace Elements in Coal Fly Ash, 62
Trace Elements in Fuel Oil, 62
Trace Mercury in Coal, 59
V and Ni in Residual Fuel Oil, 59

- FREEZING POINT**, 102
 (See also, Heat & Melting Point)
- Aluminum
 - Copper
 - Lead
 - Mercury
 - Tin
 - Zinc
- FUELS**
 (See also, Calorimetry)
- Alcohol in Reference Fuel (Gasoline), 60
 - Ethanol in Reference Fuel (Gasoline), 60
 - Heating Values, 99
 - n-Heptane, 118
 - Isooctane, 118
 - Lead in Reference Fuel (Gasoline), 59
 - Lime in Coal Fly Ash, 62
 - Methanol-Butanol in Reference Fuel (Gasoline), 60
 - Methanol in Reference Fuel (Gasoline), 60
 - Reference Fuels, 118
 - Synthetic Refuse-Derived Fuel (Gasoline), 99
 - Sulfur in Coal, 61
 - Sulfur in Distillate Fuel Oil, 61
 - Sulfur in Kerosene, 61
 - Sulfur in Residual Fuel Oil, 61
 - Synthetic Sulfur in Oil, 61
 - Trace Elements in Coal, 62
 - Trace Elements in Coal Fly Ash, 62
 - Trace Elements in Fuel Oil, 62
 - Trace Mercury in Coal, 59
 - V and Ni in Residual Fuel Oil, 59
- FUMED SILICA BOARD**, 104
- GADOLINIUM**
 Spectrometric Solution, 45
- GALLIUM**
 Isotopic, 44
 Melting Point, 102
 Radioactivity, 109–114
 Spectrometric Solution, 45
- GAS CHROMATOGRAPHY**, 67
 GC/MS System Performance
- GAS FURNACE BLACK**, 121
- GASES**, 54
- Aliphatic Organic
 - Aromatic Organic
 - Benzene in Nitrogen
 - Blood Gases
 - Carbon Dioxide in Air
 - Carbon Dioxide in Nitrogen
 - Carbon Dioxide in Nitrogen, Blood Gas
 - Carbon Dioxide and Oxygen in Nitrogen, Blood Gas
 - Carbon Monoxide in Air
 - Carbon Monoxide in Nitrogen
 - Methane and Propane in Air
 - Methane in Air
 - Nitric Oxide in Nitrogen
 - Nitrogen Dioxide in Air
 - Nitrogen Dioxide Permeation Device
- Organic
 Oxygen in Nitrogen
 Propane and Oxygen in Nitrogen
 Propane in Air
 Propane in Nitrogen
 Sulfur Dioxide in Nitrogen
 Sulfur Dioxide Permeation Tubes
 Tetrachloroethylene in Nitrogen
 Tetrachloroethylene Permeation Device
- GASES IN METALS**, 39
- Hydrogen
 - Oxygen
- GASOLINE**, see Fuels
- GAS TRANSMISSION**, 117
- Polyester Film
- GEOLOGICAL**, 79
- GERMANIUM**
 Spectrometric Solution, 45
- GLASS BEADS**, see Sizing
- GLASSES, CHEMICAL COMPOSITION**
- Glass, 83
 - Microanalytical, 41
- GLASSES, PHYSICAL PROPERTIES**
- Dielectric Constant, 94
 - Electrical Resistance, 94
 - Fixed Points, 95
 - Fluorescence Source, 116
 - Leaching Resistance, 94
 - Liquidus Temperature, 95
 - Stress Optical Coefficient, 95
 - Thermal Expansion, 104
 - Viscosity, 94
- GLASS SAND**, 79
- GLASS SPHERES**, see Sizing
- D-GLUCOSE (Dextrose)**
 Clinical, 49
 Primary Chemical, 44
- GOLD**
- Coating Weight, 93
 - High Purity, 40
 - Microprobe, 41
 - Mirror (Reflectance), 107
 - Radioactivity, 109–114
 - Reflectance (Specular), 107
 - Spectrometric Solution, 45
 - Vapor Pressure, 102
- GRAPHITE**, 104
 Thermal Conductivity
- HAFNIUM**, 38
- in Zircaloy
 - in Zirconium
 - Spectrometric Solution, 45
- HARDNESS**, 96
- Microhardness
- HASTELOY**, 25
- HEAT**
 (See also, Fuels)
- Calorimetric, 99
 - Differential Scanning Calorimetry, 100
 - Differential Thermal Analysis, 101
 - Enthalpy, 100
 - Freezing Points, 102

Heat Capacity, 100
Heat Source, 99
Heating Values, 99
Melting Points, 102
Superconducting Fixed Points, 101

Thermal Conductivity, 104
Thermal Expansion, 104
Thermocouple Material, 103
Thermogravimetry, 101
Vapor Pressure, 102

HEPES (and salts), 89

n-HEPTANE, 118

HIGH ALLOY STEEL, 18

HIGH PURITY METALS, 18
(See also, Spectrometric Solutions)
Aluminum (See Al Alloys)
Cadmium (See Vapor Pressure)
Copper (See Cu Benchmarks)
Gallium (See Stable Isotopes)
Gold
Iron (See Iron, Electrolytic)
Lead (See Freezing Points)
Magnesium (See Stable Isotopes)
Mercury (See Freezing Points)
Nickel (See Magnetic Moment)
Platinum
Rhenium (See Stable Isotopes)
Selenium
Silicon (See X-ray Powder Diffraction)
Silver (See Vapor Pressure)
Tin (See Freezing Points)
Titanium (See Ti Alloys)
Tungsten (See Thermal Conductivity)
Zinc
Zirconium (See Zr Alloys)

HIGH TEMPERATURE ALLOYS, 25

HOLMIUM

Oxide, Wavelength, 106
Wavelength, 106
Spectrometric Solution, 45

HUMAN

Liver (Radioactivity), 113
Lung (Radioactivity), 113
Serum (Clinical), 49, 50

HUMAN SERUM, 49, 50

AST Enzyme
Cholesterol in
Clinical
PCB's in

HYDROGEN

Gases in Metals, 39
Radioactivity, 109-114
4-HYDROXY-3-METHOXY-DL-
MANDELIC ACID (VMA), 49

IMAGE QUALITY INDICATOR

Radiographic, 125

IMPLANTS, see Corrosion

INDIUM

Melting Point, 102
Radioactivity, 109-114
Spectrometric Solution, 45

INDUSTRIAL HYGIENE, 68

INFRARED, NEAR, 107

Reflectance

Wavelength

INSTRUMENT PERFORMANCE

Atomic Absorption Spectrophotometer,
45

Beta-Backscatter Gages, 93

Calorimetry, 99

Conductivity Meters, 119

Coordinate Measuring Machines, 123

Differential Scanning Calorimetry, 100

Differential Thermal Analysis, 101

Dilatometers, 104

Electron Microprobe, 41, 42

Ellipsometers, 93

Gas Analysis, 54

GC/MS, 67

Hardness, Micro-, 96

Optical Emission Spectrometers, 15-40

Optical Microscopes, 91

Particle Counters, 122

pH Meters, 89

Polarimetry, 108

Profilometers, 124

Radioactivity Counting Systems, 109-114

Reflectometers, 107

Sacchrimeters, 108

SEM's, 90

Spectrophotometers, 106

Surface Analyzers, 123

Thermometers, 103

Turbidimeters, 122

Viscometers, 94, 95

X-ray Diffraction, 116

X-ray Fluorescence, 116

IODINE, 109-114

Radioactivity

ION ACTIVITY, 90

Potassium Chloride

Potassium Fluoride

Sodium Chloride

IRIDIUM

Spectrometric Solution, 45

IRON

Alloy Cast, 27, 28

Blast Furnace, 28

Cast, 27, 28

Clinical, 49, 50

Ductile Cast, 27, 28

Electrical Resistivity, 104, 118

Electrolytic, 17

Gray Cast, 27

High Alloy White Cast, 28

Metallo-Organic, 72

Nodular Cast, 27

Radioactivity, 109-114

Spectrometric Solution, 45

Thermal Conductivity, 104

White, 27, 28

ISOBUTYLENE-ISOPRENE, 121

Rubber

ISOTOPIC (STABLE), 88

Boron

Bromine

Chlorine

Chromium

Copper
Gallium
Lead, Common
Lead, Equal-Atom
Lead, Radiogenic
Magnesium
Nickel
Potassium
Rhenium
Rubidium
Silicon
Silver
Strontium
Thallium

KAOLIN, 123
 Surface Area
KEROSINE, 61
 Sulfur in
KNOOP, 96
 Microhardness
KRYPTON, 109-114
 Radioactivity

LANTHANUM
 Spectrometric Solution, 45

LEAD
 Alloys, 34
 Base Material, 34
 -Based Paint, Powdered, 59
 Freezing Point, 102
 -in Reference Fuels, 59
 Isotopic, 88
 -on Filter Media, 69
 Metallo-Organic, 72
 Nitrate (Clinical), 49
 Radioactivity, 109-114
 Spectrometric Solution, 45

LEAVES, 53
 Citrus
 Pine Needles
 Tomato

LIMESTONE, 79
 Argillaceous
 Dolomitic

LINERBOARD, 128
 Tape Adhesion Testing

LINEWIDTH MEASUREMENT, 91

LITHIUM
 Carbonate (Clinical), 49
 Metallo-Organic, 72
 Ores, 75
 Spectrometric Solution, 45

LIVER
 Bovine, Biological, 51
 Human, Radioactivity, 113

LUBRICANT OXIDATION PACKAGE, 73

LUBRICATING OIL

 Chlorine in, 72
 Ni and V in, 59
 Nitrogen in, 72
 Sulfur in, 72
 Wear-Metals in, 73

LUNG, 113

 Human, Radioactivity

LUTETIUM

 Spectrometric Solution, 45

MAGNESIUM

 Clinical, 45
 Isotopic, 88
 Metallo-Organic, 72
 Spectrometric Solution, 45

MAGNETIC COMPUTER STORAGE

MEDIA, 125

 Cartridge: High-Density Tape
 Cartridge: Tape
 Cassette: Tape
 Flexible Disk
 Reel: High-Density Tape
 Reel: Tape

MAGNETIC MOMENT, 105

 Nickel

MAGNETIC SUSCEPTIBILITY, 105

 Aluminum
 Manganese Fluoride
 Palladium

MANGANESE

 -on Filter Media, 69
 Spectrometric Solution, 45

D-MANNITOL, 49

 Clinical

MARIAJUANA METABOLITE (THC IN FREEZE DRIED URINE), 68

MARINE CHEMISTRY

 (See also, Environmental, Sediment, & Water Analysis)
 Albacore Tuna, 51
 Anion Solution, 48
 Estuarine Sediment, 62
 Limestone, 79
 Oyster Tissue, 51
 PCB's in Sediment, 63
 Spectrometric Solution, 45

MATERIALS ON FILTER MEDIA, 69

 Arsenic
 Beryllium
 Cadmium
 Lead
 Manganese
 Quartz
 Zinc

MELTING POINT, 102

 (See also, Heat & Freezing Point)

 Alumina
 Gallium
 Indium
 Rubidium
 Succinonitrile

- MERCAPTOBENZOTHIAZOLE**, 121
- MERCURY**
- Freezing Point, 102
 - in Coal, 59
 - in Freeze-Dried Urine, 68
 - in Soil, 79
 - in Water, 59
 - Radioactivity, 109–114
 - Spectrometric Solution, 45
- METALLO-ORGANIC COMPOUNDS**, 72
- Aluminum
 - Barium
 - Cadmium
 - Calcium
 - Chromium
 - Copper
 - Iron
 - Lead
 - Lithium
 - Magnesium
 - Nickel
 - Phosphorus
 - Silicon
 - Silver
 - Sodium
 - Strontium
 - Tin
 - Vanadium
 - Zinc
- METALLURGICAL**, 115
- Austenite in Ferrite
 - Iron Carbide in Ferrite
- METALS ON FILTER MEDIA**, 69
- Arsenic
 - Beryllium
 - Cadmium
 - Lead
 - Manganese
 - Zinc
- METHANE**, 55
- METROLOGY**
- Depth Profile, 91
 - Optical Linewidth, 91
 - Particle Size, 122
 - SEM Magnification, 90
 - Socketed Ball Bar, 123
- MICROANALYTICAL**, 41, 42
- Fe-Cr-Ni Alloy
 - Gold-Copper
 - Gold-Silver
 - Iron-Silicon
 - Mg-Si-Ca-Fe
 - Mineral Glasses
 - Thin Film
 - Tungsten-Molybdenum
- MICROCHEMICAL**, 44
- Acetanilide
 - Anisic Acid
 - o-Bromobenzoic Acid
 - m-Chlorobenzoic Acid
 - Cystine
 - p-Fluorobenzoic Acid
 - Nicotinic Acid
 - Urea
- MICROHARDNESS**, 96
- MICROSCOPY**, 91, 122
- Optical
- Scanning Electron
- Knoop
- Vickers
- MICROPROBE**, see Microanalytical
- MICROSPHERE**, 122
- MILK, NON-FAT POWDERED**, 51
- MINERALS**, 79
- MOLECULAR WEIGHT-MELT FLOW**, 97, 98
- Polymers
- MOLYBDENUM**
- Concentrate, 75
 - Heat Capacity, 97, 98, 100
 - Radioactivity, 109–114
 - Spectrometric Solution, 45
- NEODYMIUM**
- Spectrometric Solution, 45
- NICKEL**
- Alloys, 35
 - Isotopic, 88
 - Magnetic Moment, 105
 - Metallo-Organic, 72
 - Oxides, 36
 - Radioactivity, 109–114
 - Spectrometric Solution, 45
- NICOTINIC ACID**, 44
- Microchemical
- NIOBIUM**, 109–114
- Radioactivity
 - Spectrometric Solution, 45
- NIST TIME SOFTWARE**, 127
- NITRATE**, 48
- Anion Solution
- NITRIC OXIDE**, 55
- NITROGEN**
- in Lubricating Base Oil, 72
 - Nitrogen Dioxide, 55
- NITROUS OXIDE**, 55
- 4-NITROPHENOL**, 49
- Clinical
- NONDESTRUCTIVE EVALUATION (NDE)**, 123
- Dye Penetrant Test Blocks
- NONFERROUS ALLOYS**, 30
- N-TERTIARY-BUTYL-2-**
- BENZOTHIAZOLESULFENAMIDE**, 121
 - Rubber Compounding Material
- NUCLEAR MATERIALS**, 87
- Fission Track Glass
 - Neutron Density Monitor Wire
- NUTRITION**, see Food and Beverage, and Diet
- OBSIDIAN ROCK**, 80
- OIL**
- Base Oil, Wear Metals in, 73
 - Chlorine in Lube Base Oil, 72

- Distillate (Diesel) Fuel Oil (Sulfur),** 61
Nitrogen in Lube Base Oil, 72
Petroleum Crude Oil (Organics), 63, 64
Residual Fuel Oil (Sulfur), 61
Residual Fuel Oil (Trace Elements), 62
Residual Fuel Oil (V and Ni), 61
Sulfur in Lube Base Oil, 72
- OIL FURNACE BLACK,** 121
- OPAL GLASS**
- Composition, 83
 - Reflectance, 107
- OPTICAL**
- (See also, Reflectance and Spectrophotometry)
 - Centerline Drawings for OCR, 126
 - Linewidth Measurement, 91
- ORES, 75**
- Alumina
 - Bauxites
 - Copper
 - Fluorspar
 - Iron
 - Lithium
 - Manganese
 - Molybdenum
 - Phosphate Rock
 - Reduced Iron Oxide
 - Rutile
 - Scheelite
 - Tungsten
- ORGANIC CONSTITUENTS, 66, 67**
- Chlorinated Biphenyls
 - Chlorinated Pesticides
 - Dinitropyrene Isomers
 - Dioxin in *Isooctane*
 - Generator Columns for PAH's
 - Halocarbons for Water Analysis
 - in Cod Liver Oil
 - in Marine Sediment
 - I-Nitropyrene in Methylene Chloride
 - Isotopically Labeled and Unlabeled Priority Pollutants
 - Nitrated Polycyclic Aromatic Hydrocarbons
 - Organics in Shale Oil
 - PAH's (from Coal Tar)
 - Petroleum Crude Oil
 - Polychlorinated Biphenyls in Human Serum
 - Polychlorinated Biphenyls in Oil
 - Polychlorinated Biphenyls in Sediments
 - Priority Pollutant Phenols in Methanol
 - Priority Pollutant PAH's
 - Urban Dust
- OXALIC ACID, 111**
- Radioactivity
- OXYGEN, 55**
- Analyzed Gas
- OYSTER TISSUE, 51**
- PACKAGE**
- Catalyst Package for Lubricant Oxidation, 73
- PAINT**
- Lead-Based, 59
- PALLADIUM**
- Spectrometric Solution, 45
- PARTICLE SIZE, 122**
- PARTICULATES**
- Metals on Filter Media, 69
 - Urban Dust/Organics, 63, 64
 - Urban, Trace Elements, 62
- pD, 89**
- Disodium Hydrogen Phosphate
 - Hepes
 - Potassium Dihydrogen Phosphate
 - Potassium Hydrogen Phthalate
 - Sodium Bicarbonate
 - Sodium Carbonate
- PERMEATION DEVICES, 57**
- PERUVIAN SOIL, 113**
- Radioactivity
- PESTICIDES, 65**
- pH, 89**
- Disodium Hydrogen Phosphate
 - Potassium Tetroxalate
 - Sodium Bicarbonate
 - Sodium Carbonate
 - Sodium Tetraborate Decahydrate (Borax)
 - Tris(hydroxymethyl)aminomethane
 - Tris(hydroxymethyl)aminomethane hydrochloride
- PHOSPHATE**
- Anion Solution, 48
- PHOSPHORUS**
- Metallo-Organic, 72
 - Phosphate Rock, 74
 - Radioactivity, 109–114
 - Spectrometric Solution, 45
- PHOTOGRAPHIC, 125**
- Photographic Step Tablet (0–4)
 - Microcopy Resolution Test Charts
- PINE NEEDLES, 53**
- PLASTIC, see Polymer**
- PLASTIC PIPE, 97**
- PLATINUM**
- Doped, 40
 - High Purity, 40
 - Spectrometric Solution, 45
 - Thermoelement, 103
- PLUTONIUM**
- Radioactivity, 109–114
- POLLUTANTS**
- Atmospheric Gases, 54–56
 - Dust, 63, 64
 - Inorganics, 62
 - Organics, 68
 - Rainwater, Simulated, 59
 - Sediments, 62
 - Spectrometric Solution, 45
- POLONIUM, 109–114**
- Radioactivity
- POLYCHLORINATED BIPHENYLS (PCB)**
- in Human Serum, 49, 67
 - in Oil, 63, 65
 - in Sediment, 63
- POLYESTER FILM, 117**
- POLYETHYLENE, 97**
- Molecular Weight

| | |
|---|--|
| POLYISOBUTYLENE , 98 | |
| Rheology | |
| POLYMER | |
| Particle Size, 122 | |
| Polyethylene (Molecular Weight), 97 | |
| Polyisobutylene (Rheology), 98 | |
| Poly(methyl methacrylate), 97, 98 | |
| Polystyrene (Molecular Weight), 97 | |
| POLY(METHYL METHACRYLATE) , 97, 98 | |
| POLYSTYRENE | |
| Differential Thermal Analysis, 101 | |
| Molecular Weight, 97, 98, 100 | |
| Particle Size, 122 | |
| POTASSIUM | |
| -Chloride (Calorimetric), 99 | |
| -Chloride (Clinical), 49 | |
| -Chloride (Electrolytic Conductance), 119 | |
| -Chloride (Isotopic), 88 | |
| -Chloride (Primary), 44 | |
| -Dichromate (Primary), 44 | |
| -Dichromate (Spectrophotometry), 106 | |
| -Dihydrogen Phosphate (Fertilizer), 74 | |
| -Dihydrogen Phosphate (pH), 89 | |
| -Hydrogen Tartrate (pH), 89 | |
| -Hydrogen Phthalate (pH), 89 | |
| -Hydrogen Phthalate (Primary), 44 | |
| -Iodide (Stray Light), 106 | |
| -Nitrate (Enthalpy), 100 | |
| -Nitrate (Fertilizer), 74 | |
| Spectrometric Solution, 45 | |
| -Tetroxalate (pH), 89 | |
| POTENTIAL AND THICKNESS STEP TEST , 115 | |
| POWDERED LEAD-BASED PAINT , 59 | |
| PRAESODYMIUM | |
| Spectrometric Solution, 45 | |
| PRIMARY CHEMICALS , 44 | |
| Arsenic Trioxide | |
| Benzoic Acid | |
| Boric Acid | |
| Dextrose (D-glucose) | |
| Potassium Chloride | |
| Potassium Dichromate | |
| Potassium Hydrogen Phthalate | |
| Sodium Oxalate | |
| Strontium Carbonate | |
| Sucrose | |
| Tris(hydroxymethyl)aminomethane | |
| PRIORITY POLLUTANTS , see Organics | |
| PROMETHIUM , 109–114 | |
| Radioactivity | |
| QUARTZ , 70 | |
| Respirable Alpha | |
| QUININE SULFATE DIHYDRATE , 106 | |
| Fluorescence | |
| RADIOACTIVITY , 109–114 | |
| Americium | |
| Barium | |
| Carbon-14 | |
| Cesium-137/Barium-137m Point Source | |
| Chlorine | |
| Chromium | |
| Cobalt | |
| Curium | |
| Environmental | |
| Europium | |
| Gallium | |
| Gold | |
| Human Liver | |
| Human Lung | |
| Hydrogen | |
| Indium | |
| Iodine | |
| Iron | |
| Krypton | |
| Lead | |
| Mercury | |
| Mixed Radionuclides | |
| Molybdenum | |
| Natural Matrix | |
| Nickel | |
| Niobium | |
| Peruvian Soil | |
| Phosphorus | |
| Plutonium | |
| Polonium | |
| Promethium | |
| Radium | |
| River Sediment | |
| Rocky Flats Soil Number 1 | |
| Selenium | |
| Strontium | |
| Sulfur | |
| Technetium | |
| Thallium | |
| Thorium | |
| Tin | |
| Uranium | |
| Xenon | |
| Ytterbium | |
| RADIOGRAPHIC | |
| X-ray Film Step Tablet, 125 | |
| RADIUM , 109–114 | |
| Radioactivity | |
| RAINWATER , 59 | |
| Simulated | |
| REFERENCE FUELS , see Fuels | |
| REFLECTANCE , 107 | |
| (See also, Spectrophotometry) | |
| Aluminum Mirrors | |
| Aluminum Mirror with Wedge | |
| Black Porcelain Enamel | |
| Directional-Hemispherical | |
| First Surface Mirror | |
| Gold Mirror | |
| Mirrors | |
| Near IR | |
| Second Surface Mirror | |
| White Ceramic Tile | |
| White Opal Glass | |
| REFRACTIVE INDEX , 108 | |
| Ellipsometrically Derived, SiO ₂ on Si | |

- Isooctane*
 Silicone Liquids
 Soda-Lime Glass
 Toluene
 2,2,4-Trimethylpentane
- REFRACTORIES**, 80
- RESIDUAL RESISTIVITY RATIO**, 119
- RESISTIVITY**
- Electrical, 118
 - RRR, 119
 - Silicon, Eddy Current Testers, 118
 - Silicon, Spreading Resistance, 118
 - Thermal, 104
- RHENIUM**, 109-114
- Isotopic
 - Spectrometric Solution, 45
- RHODIUM**
- Spectrometric Solution, 45
- RICE FLOUR**, 51
- RIVER SEDIMENT**
- Environmental, 62
 - PCB's in, 63
 - Radioactivity, 113
- ROCKS**
- Basalt, 80
 - Obsidian, 80
 - Phosphate, 74, 77
- ROCKY FLATS SOIL**, 113
- Radioactivity
- RRR**, see **Residual Resistivity Ratio**
- RUBBER**, 121
- Butyl
 - Isobutylene-Isoprene (Butyl)
 - Styrene Butadiene
- RUBBER COMPOUNDING MATERIALS**, 121
- Channel Black
 - Gas Furnace Black
 - Mercaptobenzothiazole
 - n-Tertiary-Butyl-2
 - Oil Furnace Black
 - Stearic Acid
 - Sulfur
 - Zinc Oxide
- RUBIDIUM**
- Isotopic, 88
 - Melting Point, 102
 - Spectrometric Solution, 45
- RUTHENIUM**
- Spectrometric Solution, 46
- SAMARIUM**
- Spectrometric Solution, 46
- SAND, GLASS**, 79
- SCANDIUM**
- Spectrometric Solution, 46
- SCANNING ELECTRON MICROSCOPE (SEM)**, 90
- SEM Magnification
 - SEM Performance
- SCHEELITE ORE**, 75
- SEDIMENT**
- Estuarine, 62
- Marine, Organics in, 63
- PCB's in, 63
- River, 62
- River (Radioactivity), 113
- SELENIUM**
- Intermediate Purity, 40
 - Radioactivity, 109-114
 - Spectrometric Solution, 46
- SEM**, see **Scanning Electron Microscope**
- SEMICONDUCTOR MANUFACTURING**
- SILICA**
- Silica Refractory, 81
 - Thermal Expansion, 104
- SILICON**
- Carbide, 81
 - Density, 96
 - Eddy Current, 118
 - Isotopic, 88
 - Metallo-Organic, 72
 - Resistivity, 118
 - Silicon Metal, 26
 - Spectrometric Solution, 46
 - X-ray Powder Diffraction, 117
- SILICON DIOXIDE**
- Thickness and Refractive Thickness, 93
- SILICONE**, 108
- Refractive Index Liquids
- SILVER**
- Isotopic, 88
 - Metallo-Organic, 72
 - Microprobe, 41, 42
 - Spectrometric Solution, 46
 - Vapor Pressure, 102
- SIZING**, 122
- Glass Beads
 - Glass Spheres
 - Polystyrene Spheres
 - Portland Cement
- SMOKE DENSITY**, 128
- Flaming (ABS plastic)
 - Nonflaming (cellulose)
- SOCKETED BALL BAR**, 123
- SODA LIME GLASS**, 83
- Float
 - Sheet
- SODIUM**
- Bicarbonate (pH), 89
 - Carbonate (pH), 89
 - Chloride (Clinical), 49
 - Metallo-Organic, 72
 - Oxalate (Reductometric), 44
 - Pyruvate (Clinical), 49
 - Spectrometric Solution, 46
 - Tetraborate Decahydrate (pH) [Borax], 89
- SOLDER**, 34
- Lead-Base Alloy
- SPECTRAL REFLECTANCE**, see **Reflectance**
- SPECTROMETRIC SOLUTIONS**, 45
- Aluminum
 - Antimony
 - Arsenic
 - Barium
 - Beryllium
 - Bismuth
 - Boron

| | |
|--|--|
| Cadmium | Potassium Dichromate (UV Absorbance) |
| Calcium | Potassium Iodide Stray Light |
| Cerium | Quinine Sulfate Dihydrate (Fluorescence) |
| Cesium | Stray Light |
| Chromium | Ultra Violet |
| Cobalt | Visible |
| Copper | SPECULAR REFLECTANCE , see Reflectance |
| Dysprosium | STAINLESS STEEL |
| Erbium | Chip Form, 19 |
| Europium | Pitting/Crevice Corrosion, 116 |
| Gadolinium | Solid Form, 24 |
| Gallium | Thermal Expansion, 104 |
| Germanium | STEARIC ACID , 121 |
| Gold | STEELMAKING ALLOYS , 26 |
| Hafnium | Ferrochromium, High-Carbon |
| Holmium | Ferrochromium, Low Carbon |
| Indium | Ferrochromium Silicon |
| Iridium | Ferromanganese, High-Carbon |
| Iron | Ferroniobium |
| Lanthanum | Ferrophosphorous |
| Lead | Ferrosilicon (48-Si) |
| Lithium | Ferrosilicon (73-Si) |
| Lutetium | Ferrosilicon (75-Si) |
| Magnesium | STEEL, CHEMICAL COMPOSITION |
| Manganese | Chip form, 15 |
| Mercury | Gasimetric, 39 |
| Molybdenum | Solid form, 20 |
| Neodymium | STEP TABLET , 125 |
| Nickel | Photographic |
| Niobium | Radiographic |
| Palladium | X-ray |
| Phosphorus | STEP TEST , 115 |
| Platinum | Potential & Thickness Step |
| Potassium | STRONTIUM |
| Rubidium | Carbonate, 44, 88 |
| Ruthenium | Isotopic, 44, 88 |
| Samarium | Metallo-Organic, 72 |
| Scandium | Radioactivity, 109–114 |
| Selenium | Spectrometric Solution, 46 |
| Silicon | STYRENE BUTADIENE , 121 |
| Silver | Rubber |
| Sodium | SUCCINONITRILE , 102 |
| Strontium | Triple Point Cell |
| Sulfur | SUCROSE , 44, 108 |
| Tantalum | Primary Chemical |
| Tellurium | SULFATE |
| Terbium | Anion Solution, 48 |
| Thallium | SULFUR |
| Thorium | -in Base Oil, 72 |
| Thulium | -in Coal, 61 |
| Tin | -in Fossil Fuels, 61 |
| Titanium | -in Kerosene, 61 |
| Tungsten | -in Lubricating Oil, 72 |
| Uranium | Rubber Compounding Material, 121 |
| Vanadium | Spectrometric Solution, 46 |
| Ytterbium | SULFUR DIOXIDE |
| Yttrium | Analyzed Gas, 54 |
| Zinc | Permeation Tubes, 57 |
| Zirconium | SUPERCONDUCTING |
| SPECTROPHOTOMETRY , 106 | Thermometric Fixed Point Device, 101 |
| (See also, Reflectance) | SURFACE AREA , 123 |
| Fluorescence Corrected Emission | SURFACE FINISH , 124 |
| Spectrum | Sinusoidal Roughness |
| Glass Filters (Visible) | |
| Liquid Absorbance Filters (UV and | |
| Visible) | |
| Metal-on-Quartz Filters (UV and Visible) | |

- SURFACE FLAMMABILITY**, 128
- TANTALUM**
Spectrometric Solution, 46
- TAPE ADHESION TESTING**
Linerboard, 128
- TECHNETIUM**, 109–114
Radioactivity
- TELLURIUM**
Spectrometric Solution, 46
- TERBIUM**
Spectrometric Solution, 46
- TETRACHLOROETHYLENE**
Analyzed Gas, 54
Permeation Device, 57
- TETRAHYDROCANNABINOL IN FREEZE-DRIED URINE**, 68
- THALLIUM**
Isotopic, 88
Radioactivity, 109–114
Spectrometric Solution, 46
- THC IN URINE**, 68
- THERMAL ANALYSIS**
Differential Scanning Calorimetry, 100
Differential Thermal Analysis, 101
- THERMAL CONDUCTIVITY**, 104
Fibrous Glass
Graphite
Iron
Stainless Steel
Tungsten
- THERMAL EXPANSION**, 104
Borosilicate Glass
Fused Silica
Stainless Steel
Tungsten
- THERMAL RESISTANCE**, 104
Fumed Silica Board
Glass Batt
Glass Blanket
Glass Board
- THERMOGRAVIMETRY**, 101
- THERMOMETER**, 103
Clinical Laboratory
- THERMOMETRIC FIXED POINTS**, 103
- THIANTHRENE**, 99
Calorimetric
- THICKNESS**
Coating, 92
Ellipsometrically Derived, SiO₂ on Si, 93
Step Test, Potential, 115
- THORIUM**, 109–114
Radioactivity
Spectrometric Solution, 46
- THULIUM**
Spectrometric Solution, 46
- TIME SERVICE DISKETTE**, 127
- TIN**
Enthalpy, 97, 98, 100
Freezing Point, 102
Metallo-Organic, 72
Radioactivity, 109–114
Spectrometric Solution, 46
- TITANIUM**
Alloys, 37
-Dioxide, 80
Spectrometric Solution, 46
Unalloyed Titanium, 37
Unalloyed Titanium for Gases, 39
- TOLUENE**
Density, 96
Refractive Index, 108
- TOMATO LEAVES**, 53
- TOXIC METALS IN FREEZE-DRIED URINE**, 68
- TRACE ELEMENTS**
(See also, Environmental and Organic Constituents)
-in Coal, 62
-in Estuarine Sediment, 62
-in Glass, 86
-in Oil, 62
-in Potassium Feldspar, 86
-in River Sediment, 62
-in Urban Particulate, 62
-in Water, 62
- TRANSMISSION, GAS**, 117
- 2,2,4-TRIMETHYLPENTANE (*Isooctane*)**
Calorimetric, 99
Density, 96
Reference Fuel, 118
Refractive Index, 108
- TRIDYMITE**
Quantitative XRD, 117
- TRIPALMITIN**, 49
Clinical
- TRIS(HYDROXYMETHYL)AMINOMETHANE**
Basimetric, 44
Calorimetric, 99
Clinical pH, 99
pH, 89
- TRIS(HYDROXYMETHYL)AMINOMETHANE HYDROCHLORIDE**
Clinical pH, 49
pH, 89
- TUNA, ALBACORE**, 51
- TUNGSTEN**
Carbide, 81
Concentrate, 75
Spectrometric Solution, 46
Thermal Conductivity, 104
Thermal Expansion, 104
- TURBIDIMETRIC**, 122
- ULTRASONICS**, 97
Acoustic Emission Transducer
Power Standard
- URANIUM**
Radioactivity, 109–114
Spectrometric Solution, 46
- UREA**
Calorimetric, 99
Clinical, 49, 50
Microchemical, 44

- URIC ACID**, 49
 Clinical
- URINE, FREEZE-DRIED**
 Cocaine Metabolite, 68
 Cotinine, 68
 Fluorine, 68
 Industrial Hygiene, 68
 Mercury, 68
 THC in, 68
 Toxic Metals, 68
- VANADIUM**
 -in Crude Oil, 59
 -in Residual Fuel Oil, 59
 Metallo-Organic, 72
 Spectrometric Solution, 46
- VAPOR PRESSURE**, 102
 Cadmium
 Gold
 Silver
- VICKERS**
 Microhardness, 96
- VISCOOSITY**, 94, 95
- VITAMINS**
 -in Coconut Oil, 52
- VMA** (see Clinical)
- WASPALOY**, 35
- WATER ANALYSIS**
 (See also, Electrolytic Conductance)
 Acid Rain, 59
 Anion Solutions, 48
 Halocarbons, 63, 64
 Mercury, 59
 Rainwater, Simulated, 59
 Spectrometric Solution, 46
 Trace Elements, 62
- WAVELENGTH**, 106
 Didymium Glass Filter
 Holmium Oxide Solution
 Near IR Reflectance
- WEAR, ABRASIVE**
 Tool Steel, 115
- WEAR-METALS IN OIL**, 73
- WHEAT FLOUR**, 51
- WHITE CERAMIC TILE**
 Reflectance, 122
- WHITE OPAL GLASS**
 Reflectance, 122
- WINE, STABILIZED**, 52
- XENON**, 109–114
 Radioactivity
- X-RAY**
 Diffraction, 116
 Fluorescence, Glass Target, 116
 Step Tablet, 125
 Thin-Glass Film on Polycarbonate Filter, 69
- YEAST, BREWERS**, 51
- YTTERBIUM**
 Radioactivity, 109–114
 Spectrometric Solution, 46
- YTTRIUM**
 Spectrometric Solution, 46
- ZINC**
 Alloys, 38
 Enthalpy, 97, 98, 100
 Freezing Point, 102
 High Purity, 40
 Metallo-Organic, 72
 -on Filter Media, 71
 Oxide, Rubber Compounding Material, 121
 Spectrometric Solution, 46
 Spelter, 38
- ZIRCONIUM**
 Alloys, 38
 Spectrometric Solution, 46



Stephen Wise, Reenie Parris and Stephen Chesler (left to right), of the Organic Analytical Research Division, evaluate chromatographic data for measurement of organic contaminants in an environmental SRM.

U.S. DEPARTMENT OF COMMERCE
National Institute of Standards and Technology
Standard Reference Materials Program
Bldg. 202, Room 204
Gaithersburg, MD. 20899

Return and Forwarding Postage Guaranteed

Address Correction Requested

Official Business

Penalty for Private Use \$300

**Bulk Rate
Bound Printed Matter
Postage and Fees Paid
NIST
Permit No. G195**

